

**Flex Integral /QF - Flex Integral /5
BC/FC/FCV/NTC**

**Reservedele - Spare parts
Ersatzteile - Pièces détachées
Pezzi di ricambio**

List of versions YA-001

Go to
page:

	Reserve dele Indhold	Spare parts Contents	Ersatzteile Inhalt	Pièces détachées Sommaire	Pezzi di ricambio Indice	Side - Seite Page-Pagina
8	Stander:	Unit Stand:	Standkörper:	Borne :	Colonna:	
10	Dæksler, type 89	Covers, type 89	Abdeckungen Typ 89	Carters, type 89	Carter, tipo 89	 0.6.00
14	Dæksler, type 92	Covers, type 92	Abdeckungen Typ 92	Carters, type 92	Carter, tipo 92	 0.6.01
18	Mekaniske dele	Mechanical Parts	Mechanik	Pièces mécaniques	Pezzi meccanici	 0.6.02
	Elektriske dele	Electrical Parts	Elektronik	Pièces électriques	Pezzi elettrici	 0.6.03
20	Select Vac	Select Vac	Select Vac	Select Vac	Select Vac	 0.6.04
22	HygiFlex Vac	HygiFlex Vac	HygiFlex Vac	HygiFlex Vac	HygiFlex Vac	 0.6.05
24	HygiFlex Thermo/Rinse	HygiFlex Thermo/Rinse	HygiFlex Thermo/Rinse	HygiFlex Thermo/Rinse	HygiFlex Thermo/Rinse	 0.6.06
26	Fontæne	Cuspidor	Speifontäne	Crachoir	Sputacchiera	 0.6.07
30	Armsystem	Arm System	Armsystem	Système de bras	Braccio	 0.6.08
32	Instrumentbro:	Instrument Bridge:	Instrumenten- träger:	Porte- instrument:	Faretra:	
36	Ventiler	Valves	Ventile	Vannes	Valvole	 0.6.09
	Slanger	Tubes	Schläuche	Tuyaux	Tubi Circuiti	
	Print	Printed Circuits	Platinen	Cartes	Stampati	
	Flowmeter	Flowmeter	Flowmeter	Flowmeter	Flussometro	 0.6.10
	Selectventil	Select Valve	Select Ventil	Vanne select	Valvola Select	
38	Kirurgianlæg	Surgery Plant	Chirurgie- vorrichtung	Module chirurgical	Impianto chirurgico	 0.6.11
40	Amalgam- separator:	Amalgam Separator:	Amalgam- separator:	Séparateur d'amalgame:	Separatore d'amalgama:	 0.6.12
	Metasys	Metasys	Metasys	Pompe d'écoulement	Metasys	
	Afløbspumpe	Drain Pump	Abfluspumpe	Metasys	Pompa di drenaggio	
42	Fodkontakt, type 89 BC/FC/FCV/NTC	Foot Control, type 89 BC/FC/FCV/NTC	Fussanlasser, Typ 89 BC/FC/FCV/NTC	Pédale, type 89 BC/FC/FCV/NTC	Pedale, tipo 89 BC/FC/FCV/NTC	 0.6.13
44	Instrumenter:	Instruments:	Instrumente:	Instruments :	Strumentazione:	 0.6.14
	QF-kobling	QF Coupling	QF-Kupplung	Raccord QF	Raccordo QF	
46	Motor type 89, MC2	Motor Type 89, MC2	Motor Typ 89, MC2	Moteur type 89, MC2	Motore tipo 89, MC2	 0.6.15
	Tandrenser type L	Scaler Type L	Zahnreiniger Typ L	Détartreur type L	Ablatore tipo L	
	LC-lampe	Composite Lamp	Kompositlicht	Lampe composite	Lampada compositi	
48	Sprøjte, type M	Syringe, Type M	Spritze, Typ M	Seringue, type M	Siringa, tipo M	 0.6.16
	Kold/varm	Cold/Warm	Kalt/warm	Froide/chaude	Fredda/calda	
50	Lampe:	Lamp:	OP-Lampe:	Eclairage :	Lampada:	 0.6.17
54	Unitmonteret	On the Unit	Unitmontage	... sur l'unit	... sull'unità	 0.6.18
	Loftmonteret	On the Ceiling	Deckenmontage	... au plafond	... sul soffitto	
	Patientstol, type 89:	Patient chair, type 89:	Patientenstuhl Typ 89:	Fauteuil, type 89:	Poltrona, tipo 89:	 0.6.19
56	Mekaniske dele	Mechanical Parts	Mechanik	Pièces mécaniques	Pezzi meccanici	
	Elektriske dele	Electrical Parts	Elektronik	Pièces électroniques	elettrici	
58	Polster	Upholstery	Polster	Rembourrage	Tappezzeria.	 0.6.20
	Kunststofdele	Plast Parts	Kunststoffteile	Pièces en matière synthétique	Pezzi in materiale sintetico	

	Reservedele Indhold	Spare parts Contents	Ersatzteile Inhalt	Pièces détachées Sommaire	Pezzi di ricambio Indice	Side - Seite Page-Pagina
60 62	VarioFlex: Mekaniske dele Polster Kunststofdele	VarioFlex: Mechanical Parts Upholstery Plast Parts	VarioFlex: Mechanik Polster Kunststoffteile	VarioFlex : Pièces mécaniques Rembourrage Pièces en matière synthétique	VarioFlex: Pezzi meccanici Tappezzeria Pezzi in materiale sintetico	 0.6.21 0.6.22
64	ErgoFlex	ErgoFlex	ErgoFlex	ErgoFlex	ErgoFlex	 0.6.23
66	Plejemidler. Specielt værktøj Sikringer	Care & Cleaning Special Tools Fuses	Pflegemittel Spezialwerkzeug Sicherungen	Produits d'entretien Outillage spécial Fusibles	Prodotti di pulizia e manutenzione Attrezzi speciali	 0.6.24
70	Filmbetrakter FlexView	X-ray Viewer FlexView	Filmbetrachter FlexView	Négatoscope FlexView	Negativoscopio FlexView	 0.6.25
72	Tandrenser: Supprasson	Scaler: Supprasson	Zahnreiniger: Supprasson	Détartreur: Supprasson	Ablatore: Supprasson	 0.6.26
74	Amalgam- separator: Dürr 7110	Amalgam Separator: Dürr 7110	Amalgam- separator: Dürr 7110	Séparateur d'amalgame : Dürr 7110	Separatore d'amalgama: Dürr 7110	 0.6.27
76	Berøringsfri funktionsskærm: NTC/FCV	No touch control monitor: NTC/FCV	Berührungsfrei Monitor: NTC/FCV	Ecran de commande: NTC/FCV	Schermo de comando: NTC/FCV	 0.6.28
78	Patientstol, type 93: Polster Kunststofdele	Patient chair, Type 93: Upholstery Plast Parts	Patientenstuhl Typ 93: Polster Kunststoffteile	Fauteuil, type 93: Rembourrage Pièces en matière synthétique	Poltrona, tipo 93: Tappezzeria Pezzi in materiale sintetico	 0.6.29
82	Mekaniske dele Elektriske dele	Mechanical parts Electrical parts	Mechanik Elektronik	Pièces mécaniq Pièces électroniq	Pezzi meccanici Pezzi elettrici	 0.6.30
84	Tandrenser: Type O, Odontoson	Scaler: Type O, Odontoson	Zahnreiniger: Typ O, Odontoson	Détartreur: Type O, Odontoson	Ablatore: Tipo O, Odontoson	 0.6.31
86	HygiFlex Vac Ultra	HygiFlex Vac Ultra	HygiFlex Vac Ultra	HygiFlex Vac Ultra	HygiFlex Vac Ultra	 0.6.32
88	Turbine: Type B, type P	Turbine: Type B, Type P	Turbine: Typ B, Typ P	Turbine: Type B, type P	Turbina: Tipo B, tipo P	 0.6.33
90	Håndinstrument- bord, Holdere	Hand Instrument Table, Adapter	Handinstrumenten- tische, Adapter	Plateaux Adaptateur	Tavoletta per gli strumenti a mano, Supporto	 0.6.34
92	Sugeslangeholder, fontæne. Lynkoblinger	Suction hose holder, cupidor. Quick Connections	Saugschlauchhalter, Fontäne. Schnellkupplungen	Porte-tuyaux d'aspi- ration, crachoir. Raccords rapides	Porto-tubo d'aspira- zione, sputacchiera. Raccordi rapidi	 0.6.35
94	Sprøjte, type S: Kold/varm	Syringe, type S: Cold/warm	Spritze, Typ S: Kalt/warm	Seringue, type S: Froide/chaude	Siringa, tipo S: Fredda/calda	 0.6.36
96	Vandvarmer, vand- glasfylder	Water heater, cup filler	Wasserwärmer, Glasfüller	Réchauffeur d'eau, jet de verre	Riscaldatore per l'acqua de bicchiere, il dispositivo per riempimento del bicchiere	 0.6.37

	Reservedele Indhold	Spare parts Contents	Ersatzteile Inhalt	Pièces détachées Sommaire	Pezzi di ricambio Indice	Side - Seite Page-Pagina
98	Ventiler, BC	Valves, BC	Ventilen, BC	Vannes, BC	Valvole, BC	 0.6.38
100	Sugeslængeholder, patientstol.	Suction hose holder, patient chair.	Saugschlauchhalter, Patientenstuhl.	Porte-tuyaux d'aspi- ration, fauteuil den- taire.	Porto-tubo d'aspira- zione, poltrona.	 0.6.39
102	Fontæneventil: Dürr	Cuspidor valve: Dürr	Speifontæneventil: Dürr	Vanne crachoir: Dürr	Valvola sputacchiera Dürr	 0.6.40
104	Spuleenhed: Type 98, Dürr	Flushing unit: Type 98, Dürr	Spüleinheit: Typ 98, Dürr	Unité de rinçage: Type 98, Dürr	Sistema di lavaggio: Tipo 98, Dürr	 0.6.41
106	Integral /5: Elektriske dele EPROM Mekaniske dele	Integral /5: Electrical Parts EPROM mechanical Parts	Integral /5: Elektronik EPROM Mechanik	Integral /5: Pièces électroniques EPROM Pièces mécaniques	Integral /5: Pezzi elettrici EPROM Pezzi meccanici	 0.6.42
108	Tandrenser: Amdent	Scaler: Amdent	Zahnreiniger: Amdent	Détartreur: Amdent	Ablatore: Amdent	 0.6.43
110	LC-lampe: Demetron	Composite Lamp: Demetron	Kompositlicht: Demetron	Lampe à composite: Demetron	Lampada compositi: Demetron	 0.6.44
112	Motor: Type 98, MC3	Motor: Type 98, MC3	Motor: Typ 98, MC3	Moteur: Type 98, MC3	Motore: Tipo 98, MC3	 0.6.45
114	Fodkontakt, type 98 BC/FC/FCV/NTC	Foot Control, type 98 BC/FC/FCV/NTC	Fussanlasser, Typ 98 BC/FC/FCV/NTC	Pédale, type 98 BC/FC/FCV/NTC	Pedale, tipo 98 BC/FC/FCV/NTC	 0.6.46
116	Tandrenser: Odontoson	Scaler: Odontoson	Zahnreiniger: Odontoson	Détartreur: Odontoson	Ablatore: Odontoson	 0.6.47
118	Tandrenser: EMS	Scaler: EMS	Zahnreiniger: EMS	Détartreur: EMS	Ablatore: EMS	 0.6.48
120	Videoskærm, 12"	Monitor, 12" video	Monitor, 12" Video	Moniteur, vidéo 12"	Monitor, 12" video	 0.6.49
	Ombytningsdele	Exchange Parts	Austauschteile	Pièces d'échange	Parti in sostituzione	 0.6.99
126	Indeks over alle Flex reservedels kataloger. Sidetal, der begynder med 0.6., finder du i dette Integral katalog. Ellers henvises der til Flex cc kataloget.	Index for all Flex spare part catalogues. Pages starting with 0.6. are in this catalogue. Other Pages are in the Flex cc catalogue.	Index für alle Flex Ersatzteilkataloge. Seitenzahlen, die mit 0.6. beginnen, beziehen sich auf diesen Integral Katalog. Alle anderen sind aus dem Flex-cc Katalog.	Index des pièces détachées de tous les catalogues Flex. Les pages dont le numéro commence par 0.6 sont dans le présent catalogue Intégral. Se reporter sinon au catalogue Flex cc.	Indice per tutti i cataloghi Flex dei pezzi di ricambio (troverete le pagine che cominciano per 0.6. all'interno di questo catalogo. Le altre appartengono al catalogo Flex cc).	

Flex ombytningsordning

Alle varenumre, som er mærket med *, findes også som ombytningsdel (oversigt på side 0.6.99.B). Ombytningsdele har samme varenummer som nye dele, dog er første bogstav erstattet med Z, fx AC-103 flowmeter * = ZC-103 flowmeter.

Betingelserne for FLEX' ombytningsordning fremgår af Bilag 2 i Leveringsbetingelsernes punkt 8-11 som følger:

8. Bestiller FORHANDLEREN en af de dele, der er omfattet af FLEX' ombytningsordning, skal FLEX senest 8 uger efter fakturadatoen for ombytningsdelen modtage en tilsvarende defekt del, der kan reparerer.
9. Dersom den defekte del ikke kan reparerer, fakturerer FLEX FORHANDLEREN differencen mellem ombytningsdelen og en tilsvarende ny del.
10. Dersom FLEX ikke modtager den defekte del senest 8 uger efter fakturadatoen for ombytningsdelen, debiteres FORHANDLEREN differencen mellem ombytningsdelen og en tilsvarende ny del tillagt et gebyr, på 10% af nyprisen, dog min. DKK 100
11. Modtager FLEX den defekte del senere end 8 uger efter fakturadatoen for ombytningsdelen, kan den defekte del ikke krediteres.

Flex Exchange System

All article numbers marked with an * are found as exchange parts as well (survey on page 0.6.99.B). The article number of the exchange parts is the same as the number for new parts, however the first letter has been replaced by Z, e.g. AC-103 flowmeter * = ZC-103 flowmeter.

The conditions for the Flex exchange system are given in appendix 2 in the Delivery Conditions, points 8 - 11, as follows:

8. If the DEALER orders one of the parts belonging to the Flex Exchange System FLEX must receive a corresponding, defective part which can be repaired not later than 8 weeks after the invoice date of the exchange part.
9. If the defective part can not be repaired, FLEX will invoice the difference between the exchange part and a corresponding new part.
10. If FLEX does not receive the defective part by 8 weeks after the invoice date of the exchange part, the DEALER will be invoiced the difference between the exchange part and a corresponding, new part plus a fee of 10% of the list price, however, at least DEM 25.
11. If FLEX receives the defective part later than 8 weeks after the invoice date of the exchange part, the defective part can not be credited to the DEALER.

Not valid

Flex Umtauschbestimmungen

Alle mit einem * gekennzeichneten Waren sind auch Umtauschartikel (Übersicht auf Seite 0.6.99.B). Umtauschartikel haben bis auf das Z als ersten Buchstaben die gleichen Warennummern wie entsprechende Neuteile. Z. B. AC-103 Flowmeter * = ZC-103 Flowmeter.

Die Bedingungen der Flex Umtauschregelung ergeben sich wie folgt aus Anlage 2 der Lieferbedingungen, Pkt. 8-11:

8. Bestellt der HÄNDLER ein Teil, das von FLEX' Umtauschbestimmungen umfaßt ist, soll FLEX spätestens 8 Wochen nach Fakturadatum des Umtauschteils ein entsprechendes, mangelhaftes Teil zur Reparatur bekommen.
9. Wenn das mangelhafte Teil nicht repariert werden kann, berechnet FLEX dem HÄNDLER die Differenz zwischen dem Umtauschteil und einem entsprechenden, neuen Teil.
10. Wenn FLEX das mangelhafte Teil nicht spätestens 8 Wochen nach Fakturadatum des Umtauschteils erhält, berechnet FLEX dem HÄNDLER die Differenz zwischen dem Umtauschteil und einem entsprechenden, neuen Teil, zuzüglich einer Gebühr von 10% des Neupreises, jedoch mindestens DEM 25.
11. Erhält FLEX das mangelhafte Teil später als 8 Wochen nach dem Fakturadatum des Umtauschteils, kann das mangelhafte Teil nicht gutgeschrieben werden.

Règlement pour l'échange standard

Tous les numéros de référence marqués d'un * sont aussi pièces d'échange (se reporter éventuellement à la page 0.6.99.B). Le numéro de référence des pièces d'échange est le même que celui des pièces neuves, sauf que la première lettre est remplacée Z, par exemple AC-103 flowmètre * = ZC-103 flowmètre.

Les conditions du règlement pour l'échange Flex résultent de l'annexe 2 des conditions de livraison, points 8 - 11, comme suit :

8. Si le concessionnaire commande une des pièces soumise au règlement Flex pour l'échange, au plus tard 8 semaines après la date de facturation de la pièce d'échange Flex doit avoir en main une pièce défectueuse équivalente et réparable.
9. Si la pièce défectueuse n'est pas réparable, Flex facture au concessionnaire la différence entre la pièce d'échange et une pièce neuve équivalente.
10. Si Flex n'a pas en main la pièce défectueuse au plus tard 8 semaines après la date de facturation de la pièce d'échange, le concessionnaire est débité de la différence entre la pièce d'échange et la pièce neuve équivalente avec une majoration de 10% du prix catalogue et de FRF 100 au minimum.
11. Si Flex a en main la pièce défectueuse plus tard que 8 semaines après la date de facturation de la pièce d'échange, la pièce défectueuse ne peut pas être portée à crédit.


Regole Flex per la sostituzione

Tutti i codici Flex contrassegnati * anche sono parti in sostituzione (inoltre, puoi consultare la pagina 0.6.99.B). Il codice delle parti in sostituzione è lo stesso di quello delle parti nuove, solamente la prima lettera viene sostituita da uno Z, per esempio AC-103 flussometro * = ZC-103 flussometro.

Le condizioni delle regole Flex per la sostituzione risultano dall'allegato 2 delle condizioni di consegna, punti 8 - 11 come segue:

8. Se il concessionario ordina una delle parte soggetta alle regole Flex per la sostituzione, al più tardi 8 settimane dopo la data della fattura della parte in sostituzione, la parte difettosa equivalente e riparabile deve pervenire alla Flex.
9. Se la parte difettosa non è riparabile, la Flex emette una fattura al concessionario della differenza fra la parte in sostituzione e una nuova parte equivalente.
10. Se la parte difettosa non perviene alla Flex al più tardi 8 settimane dopo la data della fattura della parte in sostituzione, il concessionario verrà addebitato della differenza fra la parte in sostituzione e una nuova parte equivalente maggiorata di una tassa del 10% del prezzo nuovo, comunque LIT 23.000 al minimo.
11. Se la parte difettosa perviene alla Flex più tardi di 8 settimane dopo la data della fattura della parte di sostituzione, la parte difettosa non può essere portata al credito.

Type 89



Stander: Dæksler, type 89

Unit Stand: Covers, type 89

YA-001


Strandkørrer: Abdeckung Typ 89

Borne: Carters, type 89

0.6.00.B

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZICHNUNG	DESCRIPTION
1	AC-429	DÆKKNAP, RAL 9018	COVER PLUG, 9018, AC-336, AC-590	DECKSTÖPSEL, 9018, AC-336, AC-590	BOUCHON, 9018, AC-336, AC-590
2	AE-610	DÆKSEL, STANDERFOD, CLASSIC	COVER, FOOT, INTEGRAL II CLASSIC	ABDECKUNG, FUSS, INTE. II CLASSIC	CARTER, EMBRAS, INTE. II CLASS
2	AC-763	DÆKSEL, STANDERFOD, INTEGR. II	COVER, FOOT, INTEGRAL II	ABDECKUNG, FUSS, INTE. II	CARTER, EMBASE, INTE. II
3	AC-467	DÆKSEL, STANDERFOD, EL-SIDE	COVER, FOOT, E-S, INTEGRAL	ABDECKUNG, FUSS, E-S, INTEGRAL	CARTER, EMBASE, COTE, E-S, INTE
4	AC-468	DÆKSEL, STANDERFOD, MEK-SIDE	COVER, FOOT, M-S INTEGRAL	ABDECKUNG, FUSS, M-S INTEGRAL	CARTER, EMBASE, M-S INTEGRAL
5	MC-063	DÆKSEL, STANDERF.	COVER, FOOT	ABDECKUNG, FUSS	CARTER, EMBASE
6	MC-590	DÆKSEL, STANDERF., M/KABELUDG.	COVER, FOOT, WITH CABLE OUTLET	ABDECKUNG, FUSS, KABELAUSGANG	CARTER, EMBASE, SORTIE CABLE
7	AC-469	DÆKSEL, AMALGAMOPSAMLER	COVER, AMALGAM COLLECTOR	ABDECKUNG, AMALGAMSAMMLER	CARTER, RECEPTACLE D'AMALGAM
8	AE-441	BLENDPLADE, FODDÆKSEL, CLASSIC	COVER PLATE, FOOT, CLASSIC	VERBLENDPLATTE, FUSS, CLASSIC	PLAQUE, EMBRASE, CLASSIC
8	MC-144	BLENDPLADE, FODDÆKSEL, INTEG.	COVER PLATE, FOOT	VERBLENDPLATTE, FUSS	PLAQUE, EMBASE
9	AE-461	LÅGE, STANDER, BA-202	DOOR, UNIT STAND, BA-202,	TÜR, STANDKÖRPER, BA-202,	PORTE, BORNE, BA-202,
9	AE-462	LÅGE, STANDER, BA-203	DOOR, UNIT, STAND, BA-203	TÜR, STANDKÖRPER, BA-203	PORTE, BORNE, BA-203
9	AE-463	LÅGE, STANDER, BA-204	DOOR, UNIT, STAND, BA-204	TÜR, STANDKÖRPER, BA-204	PORTE, BORNE, BA-204
9	AE-464	LÅGE, STANDER, BA-205	DOOR, UNIT, STAND, BA-205	TÜR, STANDPÖRPER, BA-205	PORTE, BORNE, BA-205
9	AE-465	LÅGE, STANDER, BA-206	DOOR, UNIT, STAND, BA-206	TÜR, STANDPÖRPER, BA-206	PORTE, BORNE, BA-206
9	AE-466	LÅGE, STANDER, BA-207	DOOR, UNIT, STAND, BA-207	TÜR, STANDPÖRPER, BA-207	PORTE, BORNE, BA-207
9	AC-112	LÅGE, STANDER, BA-015	DOOR, UNIT STAND, BA-015	TÜR, STANDKÖRPER, BA-015	PORTE, BORNE, BA-015
9	AC-236	LÅGE, STANDER, BA-049	DOOR, UNIT STAND, BA-049	TÜR, STANDKÖRPER, BA-049	PORTE, BORNE, BA-049
9	AC-237	LÅGE, STANDER, BA-026	DOOR, UNIT STAND, BA-026	TÜR, STANDKÖRPER, BA-026	PORTE, BORNE, BA-026
9	AC-238	LÅGE, STANDER, BA-048	DOOR, UNIT STAND, BA-048	TÜR, STANDKÖRPER, BA-048	PORTE, BORNE, BA-048
9	AC-938	LÅGE, STANDER, BA-102	DOOR, UNIT STAND, BA-102	TÜR, STANDKÖRPER, BA-102	PORTE, BORNE, BA-102
10	MC-195	BESLAG, HÆNGSEL, LÅGE	HINGE FITTINGS, DOOR	ANGELBESCHLAG, TÜR	CHARNIÈRE, PORTE
11	UA-118	SKRUE M 2,9 X 9,5 PH PLASTTITE	SCREW, M 2,9 X 9,5	SCHRAUBE, M 2,9 X 9,5	VIS, M 2,9 X 9,5
12	UC-208	GUMMIBØSNING, LÅS SERIE 1000	RUBBER BUSH, LOCK SERIES 1000	GUMMIBUCHSE, VERSCHLUSS S. 1000	DOUILLE, VERROUILLAGE S. 100
13	UC-210	LÅSEKNAP SERIE 1000, RAL 9018	LOCK BUTTON	EINRASTKNOPF	BOUTON VERROUILLAGE.
13	AE-616	LÅSEKNAP SERIE 1000, CLASSIC	LOOK BUTTON, CLASSIC	EINRASTKNOPF, CLASSIC	BOUNTON VERROUILLAGE, CLASSIC
14	AC-332	DÆKSEL STAND. E-S, BA-015	COVER, U. STAND. E-S, BA-015	ABDECK, STANDK. E-S, BA-015	CARTER, BORNE, COTE E-S, BA-0
14	AC-242	DÆKSEL STAND. E-S, BA-049	COVER, U. STAND E-S, BA-049	ABDECK. STANDK. E-S, BA-049	CARTER, BORNE, COTE E-S, BA-0
14	AC-243	DÆKSEL STAND. E-S, BA-026	COVER, U. STAND. E-S, BA-026	ABDECK., STANDK. E-S, BA-026	CARTER, BORNE, COTE E-S, BA-0
14	AC-244	DÆKSEL STAND, E-S, BA-048	COVER, U. STAND E-S, BA-048	ABDECK., STANDK. E-S, BA-048	CARTER, BORNE, COTE E-S, BA-04
14	AC-903	DÆKSEL STAND. E-S, BA-102	COVER, STAND. E-S, BA-102	ABDECK., STANDK. E-S, BA-102	CARTER, BORNE, COTE E-S, BA-102
14	AE-656	DÆKSEL STAND. E-S, BA-202	COVER, U. STAND. E-S, BA-202	ABDECK., STANDK. E-S, BA-202	CARTER, BORNE, COTE E-S, BA-20
14	AE-657	DÆKSEL STAND. E-S, BA-203	CIVERMU. STAND, E-S, BA-203	ABDECK., STANDK. E-S, BA-203	CARTER, BORNE, COTE E-S, BA-20
14	AE-658	DÆKSEL STAND. E-S, BA-204	COVER, U. STAND. E-S, BA-204	ABDECK., STANDK., E-S, BA-204	CARTER, BORNE, COTE E-S, BA-2
14	AE-659	DÆKSEL STAND. E-S, BA-205	COVER, U. STAND, E-S, BA-205	ABDECK., STANDK, E-S, BA-205	CARTER, BORNE COTE E-S, BA-205
14	AE-660	DÆKSEL STAND. E-S, BA-206	COVER, U. STAND. E-S, BA-206	ABDECK., STANDK. E-S, BA-206	CARTER, BORNE, COTE E-S, BA-2
14	AE-661	DÆKSEL STAND. E-S, BA-207	COVER U. STAND, E-S, BA-207	ABDECK., STANDK, E-S, BA-207	CARTER, BORNE, COTE E-S, BA-20
15	AC-240	DÆKSEL STAND. M-S, BA-026	COVER, U. STAND M-S, BA-026	ABDECK. STANDK. M-S, BA-026	CARTER, BORNE, COTE M-S, BA-0
15	AC-241	DÆKSEL STAND. M-S, BA-048	COVER, U. STAND M-S, BA-048	ABDECK. STANDK. M-S, BA-048	CARTER, BORNE, COTE M-S, BA-0
15	AC-333	DÆKSEL, STAND. M-S BA-015	COVER, U. STAND M-S, BA-015	ABDECK. STANDK. M-S, BA-015	CARTER, BORNE. COTE M-S, BA-0
15	AC-239	DÆKSEL, STAND. M-S BA-049	COVER, U. STAND M-S, BA-049	ABDECK. STANDK. M-S, BA-049	CARTER, BORNE. COTE M-S, BA-0
15	AE-662	DÆKSEL, STAND. M-S, BA-202	COVER, U. STAND M-S, BA-202	ABDECK. STANDK. M-S, BA-202	CARTER, BORNE. COTE M-S, BA-2
15	AE-663	DÆKSEL, STAND. M-S, BA-203	COVER, U. STAND M-S, BA-203	ABDECK. STANDK. M-S, BA-203	CARTER, BORNE. COTE M-S, BA-2
15	AE-664	DÆKSEL, STAND. M-S, BA-204	COVER, U. STAND M-S, BA-204	ABDECK. STANDK. M-S, BA-204	CARTER, BORNE. COTE M-S, BA-2
15	AE-665	DÆKSEL, STAND. M-S, BA-205	COVER, U. STAND M-S BA-205	ABDECK. STANDK. M-S BA-205	CARTER BORNE COTE M-S BA-205
15	AE-666	DÆKSEL, STAND. M-S, BA-206	COVER, U. STAND M-S, BA-206	ABDECK. STANDK. M-S, BA-206	CARTER, BORNE. COTE M-S, BA-2
15	AE-667	DÆKSEL, STAND. M-S, BA-207	COVER, U. STAND M-S, BA-207	ABDECK. STANDK. M-S, BA-207	CARTER, BORNE. COTE M-S, BA-2
15	AC-904	DÆKSEL, STAND. M-S BA-102	COVER, U. STAND M-S, BA-207	ABDECK. STANDK. M-S, BA-207	CARTER, BORNE. COTE M-S, BA-2
16	AE-702	DÆKSEL KABELH. M/AFLAST. CLASSIC	CABLE HOLDER, RIGHT, CLASSIC	KABELHALTER, RECHTS, CLASSIC	PORTE-CABLE, DROIT, CLASSIC
16	AC-336	DÆKSEL KABELH. M/AFLASTN. HØJRE	COVER, CABLE HOLDER, RIGHT	ABDECKUNG, KABELHALTER, RECHTS	CARTER, PORTE-CABLE, DROIT
16	AC-590	DÆKSEL KABELH. M/AFLASTN. VENST	COVER, CABLE HOLDER, LEFT	ABDECKUNG, KABELHALTER, LINKS	CARTER, PORTE-CABLE, GAUCHE
17	UC-355	CELLEGUMMIL. 8 X 8 SORT M. KLÆB	ADHESIVE RUBBER 8 X 8 BLACK	KLEBGUMMI 8 X 8 SCHWARZ	CAOUTCHOUC ADHESIF 8 X 8 NOI
18	MC-274	BESLAG FOR FODDÆKSEL	BRACKET FOR COVERS	BESCHLAG FÜR FUSSABDECKUNG	GARNITURE, CARTER, EMBASE
9	AF-131	Låge, stander, BA-244	Door, unit stand, BA-244	Tür, Standkörper, BA-244	
14	AF-119	Dæksel, stand E-S, BA-244	Cover, Unit stand E-S, BA-244	Abdeck. Standk. E-S, BA-244	
15	AF-120	Dæksel, stand M-S, BA-244	Cover, Unit stand M-S, BA-244	Abdeck. Standk. M-S, BA-244	

Type 92



Flex Dental A/S	Unit Stand: Covers, Typ 92	1999-06-28
Stander: Dæksler, type 92	Borne : Carters, type 92	YA-001
Standkôrper: Abdeckungen, type 92		0.6.01.B


No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	MC-540	BLÆNDPROP, SUG, INTEGRAL UNIT	COVER PLUG, SUCTION	DECKSTÖPSEL, ABSAUGER	OBTURATEUR, ASPIRATION
2	MC-748	AFBLÆNDING, GULDFANG/LAMPE INT.	STOPPER, GOLD TRAP, LAMP/INT.	STÖPSEL, GOLDFÅNGER, LAMPE/INT	OPERCULE, TRAPPE OR/LAMPE INT
3	AE-610	DÆKSEL, STANDERFOD, CLASSIC	COVER, FOOT, INTEGRAL II CLASSIC	ABDECKUNG, FUSS, INTE. II CLASSIC	CARTER, EMBRAS, INTE. II CLASS
3	AC-763	DÆKSEL, STANDERFOD, INTEGR. II	COVER, FOOT, INTEGRAL II	ABDECKUNG, FUSS, INTE. II	CARTER, EMBASE, INTE. II
4	AB-379	FRONTPLADE, SELECT VAC, BA-202	FRONT PLATE, SELECT, BA-202	FRONTPLATTE, SELECT, BA-202	PLAQUE, FR., SELECT, BA-202
4	AB-380	FRONTPLADE, SELECT VAC, BA-203	FRONT PLATE, SELECT, BA-203	FRONTPLATTE, SELECT, BA-203	PLAQUE, FR., SELECT, BA-203
4	AB-381	FRONTPLADE, SELECT VAC, BA-204	FRONT PLATE, SELECT, BA-204	FRONTPLATTE, SELECT, BA-204	PLAQUE, FR., SELECT, BA-204
4	AB-382	FRONTPLADE, SELECT VAC, BA-205	FRONT PLATE, SELECT, BA-205	FRONTPLATTE, SELECT, BA-205	PLAQUE, FR., SELECT, BA-205
4	AB-383	FRONTPLADE, SELECT VAC, BA-206	FRONT PLATE, SELECT, BA-206	FRONTPLATTE, SELECT, BA-206	PLAQUE, FR., SELECT, BA-206
4	AB-384	FRONTPLADE, SELECT VAC, BA-207	FRONT PLATE, SELECT, BA-207	FRONTPLATTE, SELECT, BA-207	PLAQUE, FR., SELECT, BA-207
4	AC-902	FRONTPLADE, SELECT VAC, BA-102	FRONT PLATE, SELECT VAC, BA-102	FRONTPLATTE, SELECT VAC, BA-102	PLAQUE FR., SELECT VAC, BA-1
4	AC-265	FRONTPLADE, SELECT VAC, BA-015	FRONT PLATE, SELECT VAC, BA-015	FRONTPLATTE, SELECT VAC, BA-015	PLAQUE FR., SELECT VAC, BA-0
4	AC-266	FRONTPLADE, SELECT VAC, BA-049	FRONT PLATE, SELECT VAC, BA-049	FRONTPLATTE, SELECT VAC, BA-049	PLAQUE FR., SELECT VAC, BA-0
4	AC-267	FRONTPLADE, SELECT VAC, BA-026	FRONT PLATE, SELECT VAC, BA-026	FRONTPLATTE, SELECT VAC, BA-026	PLAQUE FR., SELECT VAC, BA-0
4	AC-268	FRONTPLADE, SELECT VAC, BA-048	FRONT PLATE, SELECT VAC, BA-048	FRONTPLATTE, SELECT VAC, BA-048	PLAQUE FR., SELECT VAC, BA-0
5	AB-442	DÆKPLADE, STANDERDÆKSEL, CLASSIC	COVER PLATE, SIDE, CLASSIC	BLLENDE, STANDKÔRPERABD. CLASSIC	PLAQUE, CARTER, BORNE. CLASSIC
5	MC-140	DÆKPLADE, STANDERDÆKSEL	COVER PLATE, SIDE	BLLENDE, STANDKÔRPERABD.	PLAQUE, CARTER, BORNE.
6	MC-778	TAP DÆKSEL STANDER INTEGRAL II	PIVOT TUBE COVER STAND II	ZAPPEN ABDECK. STANDK II	TENON CARTER BORNE II
7	MC-859	CLIPS, STANDERDÆKSEL TY II	CLIPS, COVER UNIT STAND	CLIPS, STANDKÔRPERABDECKUNG	CLIPS, CARTER, BORNE
8	UC-209	GUMMIBØSNING, LÅS SERIE 600	RUBBER BUSH, LOCK SERIES 600	GUMMIBUCHSE, VERSCHLUSS S. 600	DOUILLE, VERROUILLAGE S. 600
9	UC-214	LÅSEKNAP SERIE 600	LOCK BUTTON SERIES 600	BINRASTKNOPF SERIE 600	BOUTON DE VERROUILLAGE S. 60
10	MC-264	BLADFJEDER, STANDERDÆKSEL TY I	LAMINATED SPRING, UNIT ST. COVER	BLATTFEDER, STANDKÔRPERDECKEL	RESSORT A FEUILLE, CARTER BOR
11	MC-265	CLIPS, STANDERDÆKSEL TY I	CLIPS, COVER UNIT STAND	CLIPS, STANDKÔRPERABDECKUNG	CLIPS, CARTER, BORNE
12	AB-443	DÆKPLADE, STANDER, CLASSIC	COVER PLATE, UNIT STAND CLASSIC	ABDECKUNG, STANDKÔRPER, CLASSIC	PLAQUE, BORNE, CLASSIC
12	MC-141	DÆKPLADE, STANDER, 9018	COVER PLATE, UNIT STAND, 9018	ABDECKUNG, STANDKÔRPER, 9018	PLAQUE, BORNE, 9018
13	ME-010	DÆKPL. STANDERDÆKSEL, 100V STIK	COVER PLADE, SIDE 100V	NUR FUER JAPAN	ONLY FOR JAPAN
14	UC-206	CLIPS, LÅS SERIE 1000	CLIPS, LOCK SERIES 1000	CLIPS, VERSCHLUSS SERIE 1000	COLLIERS, VERROUILLAGE S. 100
15	UC-207	CLIPS, LÅS SERIE 600	CLIP, LOCK SERIES 600	CLIP, VERSCHLUSS SERIE 600	COLLIER, VERROUILLAGE S. 600
16	AB-674	DÆKPL. STANDERDÆK. AMALG. CLASSIC	COVER PLATE, COMBIS, CLASSIC	BLLENDE STANDK. COMBIS, CLASSIC	PLAQUE, CARTER, COMBIS, CLASSI
16	MC-779	DÆKPL. STANDERDÆK, AMALG.	COVER PLATE, SIDE COVER, COMBIS.	BLLENDE, STANDKÔRPERABD., COMBIS.	PLAQUE, CARTER DE LA BORNE, CO
17	AB-337	DÆKSEL, ARMBASIS, BA-202	COVER ARM BASE, BA-202	ABDECKUNG, ARMBASIS, BA-202	CARTER, BASE DU BRAS, BA-202
17	AB-338	DÆKSEL, ARMBASIS, BA-203	COVER ARM BASE, BA-203	ABDECKUNG, ARMBASIS, BA-203	CARTER, BASE DU BRAS, BA-203
17	AB-339	DÆKSEL, ARMBASIS, BA-204	COVER ARM BASE, BA-204	ABDECKUNG, ARMBASIS, BA-204	CARTER, BASE DU BRAS, BA-204
17	AB-340	DÆKSEL, ARMBASIS, BA-205	COVER ARM BASE, BA-205	ABDECKUNG, ARMBASIS, BA-205	CARTER, BASE DU BRAS, BA-205
17	AB-341	DÆKSEL, ARMBASIS, BA-206	COVER ARM BASE, BA-206	ABDECKUNG, ARMBASIS, BA-206	CARTER, BASE DU BRAS, BA-206
17	AB-342	DÆKSEL, ARMBASIS, BA-207	COVER ARM BASE, BA-207	ABDECKUNG, ARMBASIS, BA-207	CARTER, BASE DU BRAS, BA-207
17	AC-895	DÆKSEL, ARMBASIS, BA-102	COVER, ARM BASE, BA-102	ABDECKUNG, ARMBASIS, BA-102	CARTER, BASE DU BRAS, BA-102
17	MC-134	DÆKSEL, ARMBASIS, BA-015	COVER, ARM BASE, BA-015	ABDECKUNG, ARMBASIS, BA-015	CARTER, BASE DU BRAS, BA-015
17	AC-216	DÆKSEL, ARMBASIS, BA-049	COVER, ARM BASE, BA-049	ABDECKUNG, ARMBASIS, BA-049	CARTER, BASE DU BRAS, BA-049
17	AC-217	DÆKSEL, ARMBASIS, BA-026	COVER, ARM BASE, BA-026	ABDECKUNG, ARMBASIS, BA-026	CARTER, BASE DU BRAS, BA-026
17	AC-218	DÆKSEL, ARMBASIS, BA-048	COVER, ARM BASE, BA-048	ABDECKUNG, ARMBASIS, BA-048	CARTER, BASE DU BRAS, BA-048
18	AC-894	DÆKSEL, FONTÆNEBASIS, BA-102	COVER, CUSPIDOR BASE, BA-102	ABDECK. SPÛLSCHALENBASIS, BA-102	CARTER, BASE DU CRACHOIR, BA-1
18	MC-133	DÆKSEL, FONTÆNEBASIS, BA-015	COVER, CUSPIDOR BASE, BA-015	ABDECK. SPÛLSCHALENBAS. BA-015	CARTER, BASE CRACHOIR, BA-015
18	AC-214	DÆKSEL, FONTÆNEBASIS, BA-026	COVER, CUSPIDOR BASE, BA-026	ABDECK., SPÛLSCHALENBAS. BA-026	CARTER, BASE DU CRACHOIR, BA-0
18	AC-213	DÆKSEL, FONTÆNEBASIS, BA-049	COVER, CUSPIDOR BASE, BA-049	ABDECK. SPÛLSCHALENBASIS, BA-049	CARTER, BASE DU CRACHOIS, BA-0
18	AC-215	DÆKSEL, FONTÆNEBASIS, BA-048	COVER, CUSPIDOR BASE, BA-048	ABDECK. SPÛLSCHALENBASIS, BA-048	CARTER, BASE DU CRACHOIR, BA-0
18	AB-331	DÆKSEL, FONTÆNEBASIS, BA-202	COVER, CUSPIDOR BASE, BA-202	ABDECK. SPÛLSCHALENBAS. BA-202	CARTER, BASE CRACHOIR, BA-202
18	AB-332	DÆKSEL, FONTÆNEBASIS, BA-203	COVER, CUSPIDOR BASE, BA-203	ABDECK. SPÛLSCHALENBAS. BA-203	CARTER, BASE CRACHOIR, BA-203
18	AB-333	DÆKSEL, FONTÆNEBASIS, BA-204	COVER, CUSPIDOR BASE, BA-204	ABDECK. SPÛLSCHALENBAS. BA-204	CARTER, BASE CRACHOIR, BA-204
18	AB-334	DÆKSEL, FONTÆNEBASIS, BA-205	COVER, CUSPIDOR BASE, BA-205	ABDECK. SPÛLSCHALENBAS. BA-205	CARTER, BASE CRACHOIR, BA-205
18	AB-335	DÆKSEL, FONTÆNEBASIS, BA-206	COVER, CUSPIDOR BASE, BA-206	ABDECK. SPÛLSCHALENBAS. BA-206	CARTER, BASE CRACHOIR, BA-206
18	AB-336	DÆKSEL, FONTÆNEBASIS, BA-207	COVER, CUSPIDOR BASE, BA-207	ABDECK. SPÛLSCHALENBAS. BA-207	CARTER, BASE CRACHOIR, BA-207
19	AB-939	BASIS, ARM, BA-202	ARM BASE, BA-202	ARMBASIS, BA-202	BASE DU BRAS, BA-202
19	AB-949	BASIS, ARM, BA-203	ARM BASE, BA-203	ARMBASIS, BA-203	BASE DU BRAS, BA-203
19	AB-940	BASIS, ARM, BA-204	ARM BASE, BA-204	ARMBASIS, BA-204	BASE DU BRAS, BA-204

Flex Dental A/S
 Stander: Dæksler, type 92
 Standkørper: Abdeckungen, type 92

Unit Stand: Covers, Typ 92
 Borne : Carters, type 92


1999-06-28
 YA-001
 0.6.01.C

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
19	AE-941	BASIS, ARM, BA-206	ARM BASE, BA-206	ARMBASIS, BA-206	BASE DU BRAS, BA-206
19	AE-942	BASIS, ARM, BA-205	ARM BASE, BA-205	ARMBASIS, BA-205	BASE DU BRAS, BA-205
19	AE-943	BASIS, ARM, BA-207	ARM BASE, BA-207	ARMBASIS, BA-207	BASE DU BRAS, BA-207
19	AC-479	BASIS, ARM, BA-015	ARM BASE, BA-015	ARMBASIS, BA-015	BASE DU BRAS, BA-015
19	AC-480	BASIS, ARM, BA-026	ARM BASE, BA-026	ARMBASIS, BA-026	BASE DU BRAS, BA-026
19	AC-481	BASIS, ARM, BA-049	ARM BASE, BA-049	ARMBASIS, BA-049	BASE DU BRAS, BA-049
19	AC-482	BASIS, ARM, BA-048	ARM BASE, BA-048	ARMBASIS, BA-048	BASE DU BRAS, BA-048
19	AE-105	BASIS, ARM, BA-102	ARM BASE, BA-102	ARMBASIS, BA-102	BASE DU BRAS, BA-102
20	AC-483	STOPSTIFT SÆT, BASIS, ARM	STOP PIN, SET, BASE ARM	ANSCHLAGSTIFT, BASIS ARM	GOUPILLE D'ARRÊT, BASE BRAS
21	AE-361	BASIS FOR FONTÈNE, BA-202	CUSPIDOR BASE, BA-202	SPÜLSCHALENBASIS, BA-202	BASE DU CRACHOIR, BA-202
21	AE-362	BASIS FOR FONTÈNE, BA-203	CUSPIDOR BASE, BA-203	SPÜLSCHALENBASIS, BA-203	BASE DU CRACHOIR, BA-203
21	AE-363	BASIS FOR FONTÈNE, BA-204	CUSPIDOR BASE, BA-204	SPÜLSCHALENBASIS, BA-204	BASE DU CRACHOIR, BA-204
21	AE-364	BASIS FOR FONTÈNE, BA-205	CUSPIDOR BASE, BA-205	SPÜLSCHALENBASIS, BA-205	BASE DU CRACHOIR, BA-205
21	AE-365	BASIS FOR FONTÈNE, BA-206	CUSPIDOR BASE, BA-206	SPÜLSCHALENBASIS, BA-206	BASE DU CRACHOIR, BA-206
21	AE-366	BASIS FOR FONTÈNE, BA-207	CUSPIDOR BASE, BA-207	SPÜLSCHALENBASIS, BA-207	BASE DU CRACHOIR, BA-207
21	AC-622	BASIS FOR FONTÈNE, BA-015	CUSPIDOR BASE, BA-015	SPÜLSCHALENBASIS, BA-015	BASE DU CRACHOIR, BA-015
21	AC-219	BASIS FOR FONTÈNE, BA-049	CUSPIDOR BASE, BA-049	SPÜLSCHALENBASIS, BA-049	BASE DU CRACHOIR, BA-049
21	AC-899	BASIS FOR FONTÈNE, BA-102	CUSPIDOR BASE, BA-102	SPÜLSCHALENBASIS, BA-102	BASE DU CRACHOIR, BA-102
21	AC-220	BASIS FOR FONTÈNE, BA-026	CUSPIDOR BASE, BA-026	SPÜLSCHALENBASIS, BA-026	BASE DU CRACHOIR, BA-026
21	AC-209	BASIS FOR FONTÈNE, BA-048	CUSPIDOR BASE, BA-048	SPÜLSCHALENBASIS, BA-048	BASE DU CRACHOIR, BA-048
22	MC-200	CLIPS, FODDÆKSEL	CLIPS, FOOT COVER	CLIPS, FUSSABDECKUNG	CLIPS, CARTER, EMBASE
23	MC-274	BESLAG FOR FODDÆKSEL	BRACKET FOR COVERS	BESCHLAG FÜR FUSSABDECKUNG	GARNITURE, CARTER, EMBASE
3	AF-147	Dæksel, standerfor AB-386	Cover, foot, BA-386	Abdeck. Fuss, BA-386	Plaque, Select Vac 9018
4	AF-118	Frontplade, select v Ral9018	Front plate, select vac 9018	Frontplatte, Select Vac 9018	
4	UC-129	Fingerskrue 4mm sort	Finger screw 4mm black	Fingerschraube 4mm schwarz	
4	UA-249	Skrue 4x16 for UC-129	Screw 4x16 for UC-129	Schraube 4x16 für UC-129	
4	UA-082	Skrue 4x8 for Select Vac	Screw 4x8 for Select Vac	Schraube 4x8 für Select Vac	Clips, Borne, COMBISEP.
9	AF-131	Låge, stander, BA-244	Door, unit stand, BA-244	Tür, Standkørper, BA-244	
14	AF-119	Dæksel, stand E-S, BA-244	Cover, Unit stand E-S, BA-244	Abdeck. Standk. E-S, BA-244	
11	MC-899	Clips, Standerdksl amalg	Clips, cover unit combi sep.	Klamme, Standkørper Amalg	
15	AF-120	Dæksel, stand M-S, BA-244	Cover, Unit stand M-S, BA-244	Abdeck. Standk. M-S, BA-244	



No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
	AE-103	MONT.S�T STANDER,MONT�R	MOUNTING KIT,STAND,TECHNICIAN	MONTAGEKIT, ST�NDER, TECHNIKER	JEU DE MONTAGE, BORNE, TECHN
	AE-194	MONT.S�T STANDER,MONT�R, JAPAN	MOUNTING KIT, STAND JAPAN	NUR F�R JAPAN	ONLY FOR JAPAN
	AE-790	INDSATS FOR FILTER, AC-168	BUSH FOR FILTER AC-168	EINSATZ F�R FILTER AC-168	MANCHON POUR FILTRE AC-168
1	AC-160	VENTILATOR	FAN	VENTILATOR	VENTILATEUR
2	AC-345	AFL�SSAMLER, KOMPLET	DRAIN JUNCTION, COMPLETE	ABFLUSSAMMLER, KOMPLETT	JONCTION DES �COULEMENTS,COM
3	UC-159	BLINDH�TTE, W 110	BLIND CAP, W 110	BLINDKAPPE, W 110	CAPUCHON D'OBTURATION, W 11
4	AC-166	VANDTRYKREGULATOR, SAMLET	PRESSURE GOVERNOR, WATER, INT.	DRUCKREGLER, WASSER, INT.	REGULATEUR DE PRESSION A EAU
5	AC-167	LUFTTRYKREGULATOR, SAMLET	PRESSURE GOVERNOR, AIR, INT.	DRUCKREGLER, LUFT, INT.	REGULATEUR DE PRESSION A AIR
6	AC-168	HOVEDVANDFILTER, SAMLET INT.	MAIN WATER FILTER	HAUPTWASSERFILTER	FILTRE � EAU PRINCIPAL
7	SD-255	FILTERSKRUE INCL. O-RING	FILTER SCREW + O-RING	FILTERSCHRAUBE + O-RING	VIS DU FILTRE + JOINT TORIQUE
8	SD-256	O-RING, FILTERSKRUE	O-RING, FILTER SCREW	O-RING, FILTERSCHRAUBE	JOINT TORIQUE, VIS DU FILTRE
9	SD-257	FILTERGLAS	FILTER BOWL	FILTERGLAS	CLOCHE DE FILTRE
10	SD-258	OML�BERM�TRIK, FILTERHUS	UNION, FILTER HOUSING	�BERWURFMUTTER, FILTERGEH�USE	�CROU, CARTER DU FILTRE
11	SD-259	O-RING, FILTERHUS	O-RING, FILTER HOUSING	O-RING, FILTERGEH�USE	JOINT TORIQUE,CARTER DU FILT
12	SD-232	FILTER, VAND	FILTER, WATER	FILTER, WASSER	FILTRE, EAU
13	SD-245	SKRUE, FILTER	SCREW, FILTER	SCHRAUBE, FILTER	VIS, FILTRE
14	UC-032	MANOMETER �40-1/8"-0-6 (LUFT)	MANOMETER �40-1/8"-0-6 (AIR)	MANOMETER �40-1/8"-0-6 (LUFT)	MANOMETRE �40-1/8"-0-6 (AIR)
15	SD-233	MEMBRANS�T, LUFT, INTEGRAL	MEMBRANE SET, AIR	MEMBRANSATZ, LUFT	JEU DE MEMBRANES, AIR
16	SD-247	DEFLEKTOR, FILTER	DEFLECTOR, FILTER	DEFLECTOR, FILTER	DEFLECTEUR, FILTRE
17	SD-492	FILTERGLAS INCL.L�S DELE, BC	FILTER BOWL /BC	FILTERGLAS/BC	CLOCHE DE FILTRE/BC
18	SD-237	FILTER, LUFT	FILTER, AIR	FILTER, LUFT	FILTRE, AIR
19	UC-033	MANOMETER �40-1/8"-0-2,5(VAND)	MANOMETER �40-1/8"0-2.5(WATER)	MANOMETER �40-1/8"0-2,5 WASSER	MANOMETRE �40-1/8"-0-2,5(EAU
20	SD-249	MEMBRANS�T,VAND	MEMBRANE SET, WATER	MEMBRANSATZ, WASSER	JEU DE MEMBRANES, EAU
21	SD-236	O-RING, FILTERGLAS	O-RING, FILTER BOWL	O-RING, FILTERGLAS	JOINT TORIQUE,CLOCHE DE FILT
22	AC-484	VENTILBLOK INKL. FITTINGS	VALVE BLOCK + FITTINGS	VENTILBLOCK + FITTINGS	BLOC VANNE + RACCORDS
23	AC-138	KABEL, SEPARATIONSTANK	CABLE, SEPARATION AUTOMAT	KABEL, SEPARIERAUTOMAT	CABLE, SEPARATEUR AUTOMATIQUE
24	AC-249	SLANGE,TILSL�TNING AF SUG	TUBE, CONNECTION OF SUCTION	SCHLAUCH, ABSAUGVORRICHTUNG	TUYAU,RACCORDEMENT ASPIRATIO
25	AD-419	PRINT, SEP. TANK CONNECTIONS *	PRINT.CIRC."SEPARATION AUTM" *	PLATINE "SEPARATION AUTOMAT" *	CARTE "SEPARATION AUTOMAT"
26	HA-160	SEPARATIONSTANK	* SEPARATION AUTOMAT, INTEGRAL *	* SEPARIERAUTOMAT, INTEGRAL *	* SEPARATEUR AUTOMATIQUE,INT.
27	ML-183	MAGNETVENTIL, STANDER, FLANGE	SOLENOID VALVE, UNIT STAND	MAGNETVENTIL, STANDK�RPER	ELECTROVANNE, BORNE
28	AC-524	MAGNETVENTIL INKL,FITTING,VAND	SOLENOID VALVE +FITTINGS,WATER	MAGNETVENTIL + FITTINGS,WASSER	�LECTROVANNE + RACCORDS, EAU
29	AC-525	MAGNETVENTIL INKL.FITTING,LUFT	SOLENOID VALVE + FITTINGS, AIR	MAGNETVENTIL + FITTINGS, LUFT	�LECTROVANNE + RACCORDS, AIR
30	AC-775	S�T, OMBYGNING AF AFL�SSPUMPE	SET, REPLACEMENT OF DRAIN PUMP	SATZ, AUSTAUSCH ABFLUSSPUMPE	JEU, ECHANGE DE POMPE D'�COU
31	SD-292	ANKER, VAND, ML-183,AC-524	PLUNGER, WATER, ML-183,AC-524	PLUNGER, WASSER, ML-183,AC-524	PLONGEUR, EAU, ML-183,AC-524
32	SD-293	FJEDER, SD-292	SPRING, SD-292	FEDER, SD-292	RESSORT, SD-292
33	SD-294	ANKER, LUFT, FOR AC-525	PLUNGER, AIR, AC-525	PLUNGER, LUFT, AC-525	PLONGEUR, AIR, AC-525
34	SD-297	MEMBRAN,SEP.TANK HA-160	MEMBRANE, SEPARAT. AUT. HA-160	MEMBRAN,SEPARIERAUTOMAT HA-160	MEMBRANE, SEP. AUTOMAT.HA-16
35	SD-298	MAGNETVENTIL,SEP.TANK HA-160	SOLENOID VALVE,SEP.AUT. HA-160	MAGNETVENTIL,SEPARIERA. HA-160	ELECTROVANNE,SEPARAT. HA-160
36	SD-299	MAGNETVENTIL,SEP.TANK HA-160	SOLENOID VALVE,SEP.AUT. HA-160	MAGNETVENTIL,SEPARIERA. HA-160	ELECTROVANNE,SEPARAT. HA-160
37	SD-300	OML�BER, SEP. TANK HA-160	UNION FOR SEP. TANK HA-160	�BERWURF F�R SEP. TANK HA-160	�CROU/SEP. TANK HA-160
38	SD-301	STUDS, SEP.TANK HA-160	CONNECTOR, SEPARAT. AUT.HA-160	STUTZEN, SEPARIERAUT. HA-160	RACCORD, SEPARATEUR HA-160
39	SD-302	SP�NDRING, SEP. TANK HA-160	LOCK RING, SEPARAT.AUT. HA-160	SPANNRING,SEPARIERAUT. HA-160	BAGUE DE SERRAGE,SEPAR.HA-16
40	UC-526	AFD�KNINGSPROP K35	PLUG, DRAIN JUNCTION	ST�PSEL, ABFLUSSSAMLER	BOUCHON,BAC CENTRAL.�COULEME
41	UC-602	MONT.BOLT M10X70+TR� GEVIND	ATTACHMENT BOLT, M 10 X 70	MONTAGEBOLZEN, M 10 X 70	BOULON DE FIXATION, M 10 X 7
42	UA-025	M�TRIK M 10	NUT M 10	SCHRAUBENMUTTER M 10	�CROU M 10
43	AC-317	VANDL�S MED KAPSLER , INTEGRAL	TRAP WITH CAPS	SIPHON MIT KAPSELN	SIPHON AVEC CAPSULES
44	AC-601	EKVIPOTENTIALKLEMME M/SKRUE	EQUIPOTENTIAL TERMINAL + SCREW	�QUIPOTENTIALKLEMME + SCHRAUBE	TERMINAL EQUIPOTENTIEL + VIS
45	MH-242	HOVEDAFBRYDER, STOL/UNIT 220V	MAIN SWITCH, UNIT/CHAIR 220V	HAUPTSCHALTER, UNIT/STUHL 220V	INTERRUPTEUR,UNIT/FAUTEUIL22
45	MH-745	HOVEDAFBRYDER, STOL/UNIT 100V	MAIN SWITCH, UNIT/CHAIR 100V	ONLY FOR JAPAN	ONLY FOR JAPAN
46	SD-382	SERVICE KIT VENTILBLOK AC-484	STANDARD KIT,VALVE BL. AC-484	STANDARDKIT, VENTILBL. AC-484	KIT STANDARD,BLOCK VANNE AC4
47	AC-264	BL�NDH�TTE+PLANDPLADE,FLANGEV.	BLIND CAP	BLINDKAPPE	CAPUCHON D'OBTURATION
48	SD-136	SAMLER�R INCL. O-RINGE	CONNECTION PIPE 32/32	SAMMELROHR 32/32	RACCORD TUYAUX 32/32
49	MR-041	90 GRADER VINKEL �32 PVC	ANGLE 90 DEGREES � 32 PVC	WINKEL 90 GRAD � 32 PVC	COUDE 90 DEGREES � 32 PVC
50	MB-284	SLANGESTUDS M 4 X 3/16"	HOSE CONNECTOR M 4 X 3/16"	SCHLAUCHSTUTZEN M 4 X 3/16"	RACCORD POUR TUYAU M 4 X 3/1
51	MC-286	SLANGESTUDS F. SEPARATIONSTANK	HOSE CONNECTOR, SEP. AUTOMAT	SCHLAUCHSTUTZEN, SEPARIERAUT.	RACCORD POUR TUYAU,S�PARATEU

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
52	UC-434	SPULENHED. INTEGRAL	FLUSHING UNIT, INTEGRAL	SP�ULEINHEIT, INTEGRAL	UNITE DE RINCAGE INTEGRAL
53	MC-381	MELLEMR�R FOR SUG UDEN SEPTANK	ADAPTER, SUCTION	ZWISCHENROHR, ABSAUGANLAGE	TUBE INTERMEDIAIRE, ASPIRATIO
54	MC-584	BESLAG, SPULENHED	FITTINGS, FLUSHING UNIT	FITTINGS, SP�ULEINHEIT	RACCORD, UNITE DE RINCAGE
55	SG-233	PAKNING, FOR ANKER S�DE SD-292	GASKET, SEAT, VALVE SD-292	DICHTUNG, SITZ PLUNGER SD-292	JOINT, SIEGE, VANNE SD-292
56	SD-488	FITTINGS, BANJO, AC-484	FITTINGS, BANJO, AC-484	FITTINGS, BANJO, AC-484	RACCORD, BANJO, AC-484
57	SD-291	O-RING, AC-524+AC-525+ML-183	O-RING, AC-524+AC-525+ML-183	O-RING, AC-524+AC-525+ML-183	JOINT TOR., AC-524+AC-525+ML1
58	UC-126	O-RING � 1,78 X 2,9	O-RING � 2.90 X 1.78	O-RING � 2,90 X 1,78	JOINT TORIQUE � 2,90 X 1,78
59	AC-589	KABEL, SEPARAT SPULENHED	CABLE, SEPARATE FLUSHING UNIT	KABEL, SEPARATE SP�ULEINHEIT	CABLE, UNITE DE RINCAGE SEPAR
60	UC-120	D�RR SEPAMATIC STYRINGSPRINT	D�RR SEPAMATIC CONTROL PCB	D�RR SEP. STEUERUNG (PLATINE)	CARTE D�RR, SEPAMATIC
61	SD-235	FILTERGLAS / NTC - FC	FILTER BOWL / NTC - FC	FILTERGLAS / NTC - FC	CLOCHE DE FILTRE / NTC - FC
62	AB-034	LUFTTRYKREGULATOR, SAMLET, BC	PRESSURE GOVERNOR, AIR, BC	DRUCKREGLER, LUFT, BC	REGULATEUR DE PRESSION A AIR
63	AC-250	SLANGE, TILSLUTNING AF AFL�B	TUBE, CONNECTION OF OUTLET	SCHLAUCH, ANSCHLUSS VON ABFLUSS	TUYAU, RACCORDEMENT ECOULEMEN
64	AC-041	HOVEDMAGNETVENTILER V+L, INT	M. SOLENOID VALVES W+L, INTEGRAL	H. MAGNETVENTILE W+L, INTEGRAL	ELECTROVANNES PR. E+A, INTEGR
65	AC-137	KABEL HOVEDMAGNETVENTILER INT	CABLE, MAIN SOLENOID VALVES	KABEL, HAUPTMAGNETVENTILE	CABLE, ELECTROVANNES PRINC.
66	UC-762	T�TNINGSRING TYPE OL-1/8"	WASHER TYPE OL-1/8"	DICHTUNGSRING TYP OL-1/8"	RONDELLE TYPE OL-1/8"
67	AD-448	VALVE CONTROL	VALVE CONTROL	VALVE CONTROL	VALVE CONTROL
19	UC-903	Fjeder, manometer	Spring, mano meter	Feder, Manometer	
41	AE-103	Montagekit, komplet	Mounting kit, complete	Montagekit, Komplett	
22	AC-145	Kabels�t Ventilblok	Cable set for valve block	Kabelsatz Ventilblock	
MS-356		Slange �24, sort	Hose, black	Schlauch, schwarz	



Flex Dental A/S
 Stander: Elektriske dele
 Standk rper: Elektronik


Unit Stand: Electric Parts
 Borne : Pi ces  lectriques

1999-06-28
 YA-001
 0.6.03.B

QF

QF


No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZICHNUNG	DESCRIPTION
1	MC-454	ISOLATORPLADE, HOVEDPRINT	INSULATING PLATE,MAIN PR.CIRC.	ISOLIERPLATTE, HAUPTPLATINE	PLAQUE ISOLANTE,CARTE PRINCI
2	AE-075	TRAFO, 100V UNIT INTEGRAL	TRANSFORMER, 100V UNIT	NUR FUER JAPAN	ONLY FOR JAPAN
2	AE-287	TRAFO, UNIT	TRANSFORMER, UNIT INTEGRAL	TRANSFORMATOR, UNIT INTEGRAL	TRANSFORMATEUR, UNIT INTEGRA
3	MH-439	KLEMRE�KKE WEIDMULLER BL5	TERMINAL BLOCK WEIDMULLER BL5	KLEMMENLEISTE WEIDMULLER BL5	BARRETTE DE RACCORDEMENT BL5
4	AD-447	PRINT, MAIN NTC/BC/FC	* PR.CIRCUIT,MAINCONT,NTC/BC/FC*	PLATINE,MAINCONTROL,NTC/BC/FC*	CARTE,MAIN CONTROL,NTC/BC/FC
5	MH-427	SIKRING TR�G 200 MA	FUSE T 200 MA	SICHERUNG T 200 MA	FUSIBLE T 200 MA
5	MH-429	SIKRING TR�G 800 MA	FUSE T 800 MA	SICHERUNG T 800 MA	FUSIBLE T 800 MA
5	MH-435	SIKRING TR�G 315 MA	FUSE T 315 MA	SICHERUNG T 315 MA	FUSIBLE T 315 MA
5	ML-065	SIKRING TR�G 5 A	FUSE T 5 A	SICHERUNG T 5 A	FUSIBLE T 5 A
5	ML-066	SIKRING TR�G 2,5 A	FUSE T 2.5 A	SICHERUNG T 2,5 A	FUSIBLE T 2,5 A
5	ML-068	SIKRING TR�G 4 A	FUSE T 4 A	SICHERUNG T 4 A	FUSIBLE T 4 A
5	ML-069	SIKRING TR�G 6,3 A	FUSE T 6.3 A	SICHERUNG T 6,3 A	FUSIBLE T 6,3 A
6	AD-412	PRINT, MAINS CONNECTIONS/2	* PR. CIRCUIT, MAINS CONNECT/2 *	PLATINE, MAINS CONNECTIONS/2 *	CARTE, MAINS CONNECTIONS/2
6	AD-413	PRINT, MAINS CONNEX./2 NORGE	* PCB, MAINS CONNEX./2 NORWAY *	NUR FUER NORWEGEN	ONLY FOR NORWAY
6	AD-439	PRINT,MAINS CONNECTIONS 100V	PCB, MAINS CONNECTIONS, 100 V	NUR FUER JAPAN	ONLY FOR JAPAN
7	MH-438	KLEMRE�KKE BL3	TERMINAL BLOCK BL3	KLEMMENLEISTE BL3	BARRETTE DE RACCORDEMENT BL3
8	MH-327	KL�RK 2P RIACON 31049102	CONNECTOR 2P RIACON 0201.4901B	STECKER 2P RIACON 0201.4901B	PRISE 2P RIACON 0201.4901B
9	AE-107	E-PROM PROGRAMMERET, BC/FC 400	E-PROM PROGRAMMED, BC/FC 400	E-PROM PROGRAMMIERT, BC/FC 400	E-PROM PROGRAMME BC/FC 400
9	AE-203	E-PROM PROGRAMM. BC/FC AD-440	E-PROM PROGRAMM. BC/FC AD-440	E-PROM PROGRAMM. BC/FC AD-440	E-PROM PROGRAMME BC/FC AD-44
9	AE-108	E-PROM PROGRAMMERET, NTC, 440	E-PROM PROGRAMM. NTC/AD-440	E-PROM PROGRAMM. NTC/AD-440	E-PROM PROGRAMMEMME NTC/AD-4
9	AE-237	E-PROM PROG. JAPAN NTC,AD- 446	E-PROM PROG.JAPAN, NTC/AD-446	NUR FUER JAPAN	ONLY FOR JAPAN
9	AE-908	E-PROM PROGRAMM. FC-V AD-440	EPROM FC/V-440	EPROM FC/V-440	EPROM FC/V 440
9	AE-909	E-PROM PROG.JAPAN FC-V AD-446	EPROM FC/V-446	ONLY JAPAN	ONLY JAPAN
9	AE-681	E-PROM PROG. JAPAN FC, AD-446	EPROM FC-446	ONLY JAPAN	ONLY JAPAN
10	AC-429	D�KKNAP, RAL 9018	COVER PLUG,9018,AC-336+AC-590	DECKST�PSEL,9018,AC-336+AC-590	BOUCHON, 9018, AC-336+AC-590
11	MH-242	HOVEDAFBRYDER, STOL/UNIT 220V	MAIN SWITCH, UNIT/CHAIR 220V	HAUPTSCHALTER, UNIT/STUHL 220V	INTERRUPTEUR,UNIT/FAUTEUIL22
11	MH-745	HOVEDAFBRYDER, STOL/UNIT 100V	MAIN SWITCH, UNIT/CHAIR 100V	ONLY FOR JAPAN	ONLY FOR JAPAN
12	AC-738	TRAFOSK�RM KOMPLET INTEGRAL I	TRANSFORMATOR COVER I	ABDECK. DES TRANSFORMATORS I	PROTECTION DU TRANSFORMA. I
12	AE-071	TRAFOSK�RM KOMPLET INTEGRAL II	TRANSFORMATOR COVER COMPL.,II	ABD.,TRANSFORMATOR,KOMPLETT,II	PROTECT.,TRANSFORM,COMPLET I
13	MH-215	NVRAM THERMO/RINSE BC/FC	NVRAM THERMO/RINSE BC/FC	NVRAM THERMO/RINSE BC/FC	NVRAM THERMO/RINSE BC/FC
13	AE-239	NVRAM NTC/FC-V	NVRAM NTC/FC-V	NVRAM NTC/FC-V	NVRAM NTC/FC-V
13	MH-535	NVRAM BC/FC	NVRAM BC/FC	NVRAM BC/FC	NVRAM BC/FC
14	AE-273	V-LISTE, STANDER	DRAIN LIST, UNIT STAND	ABLUSSRINNE, STANDK�RPER	PROT�GE-CARTE, BORNE
15	AD-448	VALVE CONTROL	VALVE CONTROL	VALVE CONTROL	VALVE CONTROL
16	ME-014	ISOLATORPLADE, MAINS CONNECT.	INSULATING PLATE,MAINS CONNEX.	ISOLIERPLATTE, MAINS CONNECT.	PLAQUE ISOLANTE, MAINS CONN.
17	AD-447	PRINT, MAIN NTC/BC/FC	* PR.CIRCUIT,MAINCONT,NTC/BC/FC*	PLATINE,MAINCONTROL,NTC/BC/FC*	CARTE,MAIN CONTROL,NTC/BC/FC
17	AD-446	PRINT,MAIN CONTROL 100V,	PCB, MAIN CONTROL, 100V NTC	NUR FUER JAPAN	ONLY FOR JAPAN
18	AC-117	TRAFO, OPERATIONSLAMPE INTEG.	TRANSFORMER, OPERATING LAMP	TRANSFORMATOR, OPERATIONSLAMPE	TRANSFORMATEUR,ECLAIRAGE BUC
18	AE-073	TRAFO, 100V OPERATIONSLAMPE	TRANSFORMER, OP. LAMP, 100 V	NUR FUER JAPAN	ONLY FOR JAPAN



Flex Dental A/S		1999-06-28
Select Vac	Select Vac	YA-001
Select Vac	Select Vac	0.6.04.B

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	AC-107	SLANGE Ø 17,M/NIPLER,SELECT V.	TUBE Ø 16, NIPPLES, SELECT VAC	SCHLAUCH Ø 16,ANSÄTZE,SEL.VAC	CORDON Ø16,EMBOUTS,SELECT VA
1	AC-108	SLANGE Ø 12,M/NIPLER,SELECT V.	TUBE Ø 11, NIPPLES, SELECT VAC	SCHLAUCH Ø 11,ANSÄTZE,SEL.VAC	CORDON Ø11,EMBOUTS,SELECT VA
2	AC-106	FILTERENHED SELECT VAC, INT.	FILTER COMPLETE SELECT VAC	FILTEREINHEIT SELECT VAC	FILTRE COMPLET SELECT VAC
3	MR-075	FILTERRØR	FILTER TUBE	FILTERROHR	TUYAU DE FILTRE
3	AC-370	FILTERRØR, 12 X MR-075	FILTER TUBE, 12 X MR-075	FILTERROHR, 12 X MR-075	TUYAU DE FILTRE, 12 X MR-075
3	SD-320	FILTERRØR, 20 X MR-075	FILTER TUBE, 20 X MR-075	FILTERROHR, 20 X MR-075	TUYAU DE FILTRE, 20 X MR-075
4	AC-153	FILTERHOLDER,SUG,INKL.O-RINGE	FILTER HOLDER,SUCTION,+O-RINGS	FILTERHALTER,ABSAUGER,+O-RINGE	PORT.FILTRE,ASP.,+JOINTS TOR
5	AC-154	FILTERPATRON,SUG,INKL.O-RINGE	FILTER CARTR.,SUCTION,+O-RINGS	FILTERPATRONE,ABSAUG.,+O-RINGE	CARTOUCHE, ASPIR.,+JOINTS TO
6	UC-040	O-RING Ø1,78 X 4,48 EPDM/70SH.	O-RING Ø 4.48 X 1.78	O-RING Ø 4,48 X 1,78	JOINT TORIQUE Ø 4,48 X 1,78
7	UC-063	O-RING Ø 1,78 X 17,17	O-RING Ø 17.17 X 1.78	O-RING Ø 17,17 X 1,78	JOINT TORIQUE Ø 17,17 X 1,78
8	UC-064	O-RING Ø 1,78 X 28,3	O-RING Ø 28.30 X 1.78	O-RING Ø 28,30 X 1,78	JOINT TORIQUE Ø 28,30 X 1,78
9	AC-186	SUGENIPPEL, Ø 11	SUCTION TUBE NIPPLE Ø 11	SAUGSCHLAUCHANSATZ Ø 11	EMBOUT DE TUYAU D'ASPIR. Ø 1
10	AC-185	SUGENIPPEL, Ø 16	SUCTION TUBE NIPPLE Ø 16	SAUGSCHLAUCHANSATZ Ø 16	EMBOUT DE TUYAU D'ASPIR. Ø 1
11	MC-188	KAPPE, SUGENIPPEL Ø 16	JACKET,SUCTION TUBE NIPPLE Ø16	KAPPE, SAUGSCHLAUCHANSATZ Ø 16	MANCHON, EMBOUT ASPIR. Ø 16
12	MC-190	KAPPE, SUGENIPPEL, Ø 11	JACKET,SUCTION TUBE NIPPLE Ø11	KAPPE, SAUGSCHLAUCHANSATZ Ø 11	MANCHON, EMBOU ASPIR. Ø 11
13	MC-196	NIPPEL, SUGESLANGE Ø 16	NIPPLE, SUCTION TUBE Ø 16	ANSATZ, SAUGSCHLAUCH Ø 16	EMBOUT,TUYAU D'ASPIRATION Ø1
14	MC-197	NIPPEL, SUGESLANGE Ø 11	NIPPLE, SUCTION TUBE Ø 11	ANSATZ, SAUGSCHLAUCH Ø 11	EMBOUT,TUYAU D'ASPIRATION Ø1
15	AC-280	SUGESLANGE, Ø 12 GRÅ, 1650MM	SUCTION HOSE 11 MM, GREY	SAUGSCHLAUCH 11 MM, GRAU	TUYAU D'ASPIRATION 11 MM,GRI
15	MS-341	SUGESL. Ø 12 GRÅ, RAL7042	HOSE, SUCTION, Ø 11, GREY	SCHLAUCH, ABSAUGER, Ø 11, GRAU	TUYAU, ASPIRATION, Ø 11, GRI
15	AC-369	SUGESLANGE, Ø 12 GRÅ 50 METER	SUCTION TUBE, Ø 11, 50 M, 9018	SAUGSCHLAUCH, Ø 11, 50 M, 9018	TUYAU ASPIR., Ø 11, 50 M,901
16	AC-279	SUGESLANGE, Ø 17 GRÅ 1650MM	SUCTION HOSE Ø 16 MM, GREY	SAUGSCHLAUCH Ø 16 MM, GRAU	TUYAU D'ASPIRATION Ø16MM,GRI
16	MS-348	SUGESL. Ø 17 GRÅ, RAL7042	HOSE, SUCTION, Ø 16, GREY	SCHLAUCH, ABSAUGER, Ø 16, GRAU	TUYAU, ASPIRATION, Ø 16, GRI
16	AC-368	SUGESLANGE, Ø 17 GRÅ 50 METER	SUCTION HOSE, Ø 16, 50 M, 9018	SAUGSCHLAUCH, Ø 16, 50 M, 9018	TUYAU ASPIR., Ø 16, 50 M,901
17	AC-105	SELECT VAC, ENKELT, INTEGRAL	SELECT VAC	SELECT VAC	SELECT VAC
18	AB-124	HOLDER M. BALG, SELECTVAC	HOLDER WITH BELLOWS, SELECTVAC	HALTER MIT BALG, SELECTVAC	SUPPORT AVEC SOUFFLET, SELEC
19	MC-105	LÅSEFJEDER, SELECT VAC, INT.	RETAINING SPRING, SELECT VAC	HALTEFEDER, SELECT VAC	RESSORT DE SERRURE, SELECT V
20	UC-582	MAGNETVENTIL, 3-VEJS, M 5	SOLENOID VALVA	MAGNETVENTIL	ELECTROVANNE
21	MS-001	SLANGE Ø 5/Ø 3 KLAR PVC	TUBE Ø 5/Ø 3 TRANSPARENT PVC	SCHLAUCH Ø5/Ø3 DURCHSICHT. PVC	TUYAU Ø 5/Ø 3 PVC TRANSPAREN
22	AC-033	SLANGESTUDS, DRIVLUFT, Ø 4,7	TUBE CONNECTOR, AIR, QF TUBE	SCHLAUCHSTUTZEN, LUFT, QF.SCHL	RACCORD TUYAU, AIR, TUYAU QF
23	UC-126	O-RING Ø 1,78 X 2,9	O-RING Ø 2.90 X 1.78	O-RING Ø 2,90 X 1,78	JOINT TORIQUE Ø 2,90 X 1,78
24	AC-152	UNDERPART SELECT V.INKL.O-RING	LOWER PART SELECT V., + O-RING	UNTERTEIL SELECT V., + O-RING	PARTIE INF.SELECT V.,+JOINT
25	MC-080	FILTERHUS, SUG, INTEGRAL	FILTER HOUSING, SUCTION	FILTERGEHÄUSE,ABSAUGEVORRICHT.	CARTER DU FILTRE DE L'ASPIR.
26	MC-263	REDUKTION, SUG Ø 11-Ø 7, INT.	REDUCTION, SUCTION Ø 11 - 7	REDUKTION, ABSAUGER Ø 11 - 7	ADAPTATEUR, ASPIRAT. Ø 11 -7
27	AC-159	MANIFOLD SELECT VAC MONT.,INT.	MANIFOLD, SELECT VAC	ANSCHLUSSBLOCK, SELECT VAC	REPARTITEUR SELECT VAC
28	UC-037	BLÆNDPLADE,FLANGEVENTIL ML-183	COVER, VALVE ML-183	ABDECKUNG, VENTIL ML-183	CARTER, VANNE ML-183
29	MR-210	O-RING 45 X 2	O-RING 45 X 2	O-RING 45 X 2	JOINT TORIQUE 45 X 2
30	MC-540	BLÆNDPROP, SUG, INTEGRAL UNIT	COVER PLUG, SUCTION	DECKSTÖPSEL, ABSAUGER	OBTURATEUR, ASPIRATION
31	UC-127	Y-FORBINDELSSESSTUDS, NORMA YS4	Y-SHAPED CONNECTION PIECE	Y-FÖRMIGES ANSCHLUSSSTÜCK	RACCORD EN Y
32	AC-712	ENGANGSFILTER SUG, 1 STK.	FILTER+HOLDER,SUCTION,INTEGRAL	FILTER,HALTER,ABSAUGER,INTEGR.	FILTRE,SUPPORT,ASPIR.,INTEGR
32	SD-400	ENGANGSFILTER SUG, 6 STK.	FILTERHOLDER+6 FILTERS,SUCTION	FILTERHALTER+6 FILTER,ABSAUGER	PORTE-FILTRE+6 FILTRES,ASPIR
32	SD-347	ENGANGSFILTER SUG, 20 STK.	FILTERS (20), SUCTION	FILTER (20), ABSAUGVORRICHTUNG	FILTRES (20), ASPIRATION
33	SD-348	SUGEFILTERHOLDER, FLEX CC/1/2	SUCTION FILTER HOLDER,F.CC/1/2	SAUGFILTERHALTER, FLEX CC/1/2	PORTE-FILTRE, FLEX CC/1/2
34	AC-367	KAPPER, SUGENIPPEL, 6 SÆT	JACKET,SUC.TUBE NIPPLE(6 SETS)	KAPPE, SAUGSCHL.ANSATZ(6 SATZ)	MANCHON, EMBOUT ASPIR. 6 JEU
35	MC-248	SLANGESTUDS Ø4, SELECTVAC	CONNECTOR Ø 4, SELECT VAC	SCHLAUCHSTUTZEN Ø 4,SELECT VAC	RACCORD POUR TUYAU Ø 4,SEL.V
36	UC-062	O-RING Ø 18,0 X 2,0	O-RING Ø 18.0 X 2.0	O-RING Ø 18,0 X 2,0	JOINT TORIQUE Ø 18,0 X 2,0
37	SD-446	MANIFOLD INCL. MAGNETVENTIL	MANIFOLD INCL. SOLENOID VALVE	SAMMELROHR MIT MAGNETVENTIL	REPARTITEUR AVEC ELECTROVANN
38	MC-223	BESLAG FOR MANIFOLD	FITTING FOR MANIFOLD	BESCHLAG FÜR ANSCHLUSSBLOCK	PATTE / REPARTITEUR
39	AC-423	SL.STUDS F.BJEKTOR M. O-RINGE	TUBE CONNECTOR,EJECTOR,O-RING	SCHLAUCHSTUTZEN,EJEKTOR,O-RING	RACCORD TUYAU,POMPE SAL.,JOI
40	AC-040	EJEKTOR, INTEGRAL	EJECTOR/INTEGRAL	EJEKTOR/INTEGRAL	POMPE A SALIVE/INTEGRAL
41	AC-146	KABEL SELECT VAC VENTILER INT	CABLE FOR SELECT VAC VALVE INT	KABEL FÜR SELECT VAC VENTIL	CABLE,SELECT VAC SOQPAPE INT
42	AC-420	SUGNIPPEL F.EJEKTOR M. O-RING	SUCTION TUBE NIPPLE + O-RING	SAUGSCHLAUCHANSATZ + O-RING	EMBOUT T. D'ASPIRATION +JOIN

1	AE-979	Slange Ø17, 2 meter	Tube Ø17, 2 meters	Schlauch Ø17, 2 Meter
	MS-360	Inderslange Ø17	Internal hose Ø17,	Innenschlauch Ø17




Flex Dental A/S
HygiFlex Vac
HygiFlex Vac

HygiFlex Vac
HygiFlex Vac

1999-06-28
YA-001
0.6.05.B

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	AC-143	KABEL, HYGIFLEX VAC PCB	CABLE, HYGIFLEX VAC	KABEL, HYGIFLEX VAC	CABLE, HYGIFLEX VAC
2	AC-182	VAC CLEAN, PUMPE	VAC CLEAN, PUMP	VAC CLEAN, PUMPE	VAC CLEAN, POMPE
3	UC-289	O-RING Ø29,82 X 2,62 VITON 75S	O-RING Ø 29.82 X 2.62	O-RING Ø 29,82 X 2,62	JOINT TORIQUE Ø 29,82 X 2,62
4	UC-398	PAKNING, OMLØBER	GASKET, UNION	DICHTUNG, ÜBERWURFMUTTER	JOINT D'ECROU
5	MC-231	OMLØBER Ø 43	UNION Ø 43	ÜBERWURFMUTTER Ø 43	ECROU Ø 43
6	MC-112	BLANDEBEHOLDER, VAC CLEAN	MIXING TANK, VAC CLEAN	MISCHBEHÄLTER, VAC CLEAN	RECIPIENT A MELANGER, VACCLEA
7	AD-416	PRINT, HYGIFLEX VAC	PRINTED CIRCUIT "HYGIFLEX VAC"	PLATINE "HYGIFLEX VAC"	CARTE "HYGIFLEX VAC"
8	ML-183	MAGNETVENTIL, STANDER, FLANGE	SOLENOID VALVE, UNIT STAND	MAGNETVENTIL, STANDKÖRPER	ELECTROVANNE, BORNE
9	MS-009	SLANGE Ø4 X 8	TUBE Ø4X8 , SEPARATION AUTOMAT	SCHLAUCH Ø4X8 ,SEPARIERAUTOMAT	CORDON Ø4X8,SEPARATEUR AUTOM
10	AC-521	BESLAG, HYGIFLEXVAC	* FITTINGS, HYGIFLEX VAC, MEC.	* FITTINGS, HYGIFLEX VAC,MECH.	* GARNITURE, HYGIFLEX VAC,MEC.
11	AC-183	KONTRAVENTIL, VAC CLEAN	NON-RETURN VALVE, VAC CLEAN	RÜCKSCHLAGVENTIL, VAC CLEAN	SOUPAPE RETENUE, VAC CLEAN
12	MB-280	FITTINGS, BANJO, 6/4 X M 5	FITTINGS, BANJO, 6/4 X M 5	FITTINGS, BANJO, 6/4 X M 5	RACCORD, BANJO, 6/4 X M 5
13	AC-259	GENNEMFØRINGSDÆKSEL MED STUDBE	FEED-THROUGH COVER+CONNECTORS	DURCHFÜHRUNGSDECKEL + STUTZEN	CARTER PASSE-TUYAUX + RACCOR
14	MC-377	NIPPEL FOR KONTRAV. OROTOLPUMPE	NIPPLE, NON-RETURN VALVE, OROTOL	NIPPEL, RÜCKSCHLAGVENTIL, OROTOL	EMBOUT, SOUPAPE DE RETENUE
15	UC-707	NUTRING VI N 25-035 X 8 VITON	GASKET, PUMP	DICHTUNG, PUMPE	JOINT, POMPE
16	AC-260	BLANDETANK MED MONTERET DÆKSEL	MIXING TANK WITH COVER	MISCHBEHÄLTER MIT DECKEL	RECIPIENT A MELANGER+COUVERC
17	UC-126	O-RING Ø 1,78 X 2,9	O-RING Ø 2.90 X 1.78	O-RING Ø 2,90 X 1,78	JOINT TORIQUE Ø 2,90 X 1,78
18	SD-291	O-RING, AC-524+AC-525+ML-183	O-RING, AC-524+AC-525+ML-183	O-RING, AC-524+AC-525+ML-183	JOINT TOR., AC-524+AC-525+ML1




Flex Dental A/S		1999-06-28
HygiFlex Thermo/Rinse	HygiFlex Thermo/Rinse	YA-001
HygiFlex Thermo/Rinse	HygiFlex Thermo/Rinse	0.6.06.B

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	MC-161	FLANGE MED VARMELEGME	HEATING ELEMENT	HEIZSPIRALE	CARTOUCHE CHAUFFANTE
1	ME-015	FLANGE MED VARMELEGME 100V	HEATING ELEMENT, 750W, 100V	ONLY FOR JAPAN	ONLY FOR JAPAN
2	AC-352	SIKKERHEDSVENTIL 45 PSI	SAPETY VALVE	SICHERHEITSVENTIL	VANNE DE SECURITE
3	UC-049	UDLUFTNINGSVENTIL, VARMTV.BEH	RELIEF VALVE, BOILER	ENTLÜFTUNGSVENTIL, BOILER	VIS DE PURGE, CHAUFFE-BAU
4	UC-105	O-RING Ø 1,78 X 47,35, VITON	O-RING Ø 47.37 X 1.78	O-RING Ø 47,37 X 1,78	JOINT TORIQUE Ø 47,37 X 1,78
5	AC-628	INSTRUMENTH. HF THERMO/RINSE	THERMO SUPPORT, THERMO/RINSE	THERMOHALTER, THERMO/RINSE	BAC-THERMO THERMO/RINSE
6	MB-127	FITTINGS 6/4 X 1/8 VINKEL	FITTINGS 6/4 X 1/8 ANGLE	FITTINGS 6/4 X 1/8 WINKEL	RACCORD 6/4 X 1/8 COUDE
7	MB-128	FITTINGS 6/4 X 1/4	FITTINGS 6/4 X 1/4	FITTINGS 6/4 X 1/4	RACCORD 6/4 X 1/4
8	AC-591	SILICONESL. Ø8XØ12 (200MM)	TUBE, THERMO/RINSE Ø 8 X Ø 12	SCHLAUCH,THERMO/RINSE Ø8 X Ø12	CORDON,THERMO/RINSE Ø8 X Ø 1
8	UC-124	SILICONESLANGE Ø8 X Ø12	SILICONE TUBE Ø 8 X Ø 12	SILIKONSCHLAUCH Ø 8 X Ø 12	CORDON EN SILICONE Ø 8 X Ø 1
9	AC-522	NTC-FØLER M/KABEL, VARMTV.BEH.	NTC DETECTOR + CABLE, BOILER	NTC FÜHLER + KABEL, BOILER	CAPTEURS NTC.CABLE,CHAUFFE-B
10	SD-232	FILTER, VAND	FILTER, WATER	FILTER, WASSER	FILTRE, EAU
11	MC-165	VARMVANDSBEHOLDER, SVEJSET	BOILER	BOILER	CHAUFFE-BAU
12	MH-215	NVRAM THERMO/RINSE BC/FC	NVRAM THERMO/RINSE BC/FC	NVRAM THERMO/RINSE BC/FC	NVRAM THERMO/RINSE BC/FC
12	AE-239	NVRAM NTC/FC-V	NVRAM NTC/FC-V	NVRAM NTC/FC-V	NVRAM NTC/FC-V
13	UC-762	TÆTNINGSRING TYPE OL-1/8"	WASHER TYPE OL-1/8"	DICHTUNGSRING TYP OL-1/8"	RONDELLE TYPE OL-1/8"
14	AC-630	FILTER, HYGIFLEX THERMO	FILTER, HYGIFLEX THERMO	FILTER, HYGIFLEX THERMO	FILTRE HYGIFLEX THERMO
15	AE-013	OVERLØBSSIKRING, HF-THERMO	OVERFLOW PROTECTION, HF-THERMO	ÜBERLAUFSCHUTZ, HYGIFLEX THERMO	SECURITE ANTIDEBORD, HF-THERM
16	BA-087	ADAPTOR, HYGIFLEXT.SUPRASSON	ADAPTOR, HYGIFLEXT, SUPRASSON	ADAPTOR, HYGIFLEXT, SUPRASSON	ADAPTATEUR, HYGIFLEXT, SUPRASS
17	BA-085	ADAPTOR, HYGIFLEXT.4+4	ADAPTOR, HYGIFLEX THERMO 4+4	ADAPTER, HYGIFLEX THERMO 4+4	ADAPTATEUR, HYGIFLEX THERMO 4
18	BA-086	ADAPTOR, HYGIFLEX, MOTOR	ADAPTOR, HYGIFLEX, MOTOR	ADAPTER, HYGIFLEX, MOTOR	ADAPTATEUR, HYGIFLEX, MOTEUR
19	BA-083	ADAPTOR, HYGIFLEXT.SJØDING SP91	ADAPTOR, HYGIFLEXT, SJØDING SP91	ASAPTOR, HYGIFLEXT, SJØDING SP91	ADAPTATEUR, HYGIFL.SJØDING SP
20	BA-082	ADAPTOR, HYGIFLEX, SPRØJTE	ADAPTOR, HYGIFLEX, SERINGE	ADAPTOR, HYGIFLEX, SPRITZE	ADAPTEUR, HYGIFLEX, SERINQUE
21	AC-043	VARMTVANDSBEHOLDER, INTEGRAL	BOILER, INTEGRAL	HEISSWASSERBOILER, INTEGRAL	RECHAUFFEUR D'EAU, INTEGRAL
21	AE-243	VARMTVANDSBEHOLDER, INT. 100 V	BOILER, HYGIFLEX THERMO, 100 V	NUR FÜR JAPAN	ONLY FOR JAPAN
22	UA-707	SKRUE M 3 X 4 RUSTFRI	SCREW M 3 X 4	SCHRAUBE M 3 X 4	VIS M 3 X 4
5	MR-435	O-RING/Holder	O-RING for Inst. holder	O-RING für Inst.Halter	
20	BA-341	ADAPTOR, SPRØJTE S	ADAPTOR SYRINGE S	ADAPTOR SPRITZE S	
20	AE-811	ADAPTOR, SPRØJTE O	ADAPTOR, SYRINGE O	ADAPTOR, SPRITZE O	

QF + Turb.

Typ M



No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	MR-218	O-RING Ø 19,5 X 3	O-RING Ø 19.5 X 3	O-RING Ø 19,5 X 3	JOINT TORIQUE Ø 19,5 X 3
2	MR-226	O-RING Ø 4 X 1	O-RING Ø 4 X 1	O-RING Ø 4 X 1	O-RING Ø 4 X 1
3	AC-025	VANDGLASSKÅL KOMPLET	WATER GLASS BOWL	GLASSCHALE, KOMPLETT	REPOSE-VERRE COMPLET
4	MR-077	O-RING Ø 23,4 X 3,53	O-RING 23.4 X 3.53	O-RING 23,4 X 3,53	JOINT TORIQUE 23,4 X 3,53
5	MC-436	FJEDER, RØRSTUDS VENTILB./FONT	SPRING,TUBE CONNECTOR,CUSPIDOR	FEDER, ROHRSTUTZEN,SPEIFONTÅNE	RESSORT, TUYAU RAC., CRACHOI
6	AC-026	FONTÆNESKÅL KOMPLET	CUSPIDOR BOWL, COMPLETE	SPEISCHALE, KOMPLETT	CUVETTE DU CRACHOIR COMPLETE
7	MR-214	O-RING Ø 36,1 X 3,53	O-RING Ø 36.1 X 3.53	O-RING Ø 36,1 X 3,53	JOINT TORIQUE Ø 36,1 X 3,53
8	AC-130	FONTÆNESKYLLERØR INKL O-RINGE	FLUSHING PIPE,CUSPIDOR,O-RINGS	SPÜLROHR, FONTÅNE, + O-RINGE	TUYAU RINC. CRACHOIR,J.TORIQ
8	AC-795	SKYLLERØR, FONTÆNE FC-BASIC	FLUSHING PIPE, BC	SPÜLROHR, BC	TUYAU DE RINCAGE, BC
9	UC-051	O-RING Ø 6 X 1	O-RING Ø 6 X 1	O-RING Ø 6 X 1	JOINT TORIQUE Ø 6 X 1
10	MC-465	DYSE, SKYLLERØR	NOZZLE, FLUSHING TUBE	DÜSE, SPÜLROHR	BUSE, TUYAU D'ARRIVEE
11	MC-245	TRYKFJEDER, VANDGLASFYLDER	COMPRESSION SPRING, CUP FILLER	DRUCKFEDER, GLASFÜLLER	RESSORT DE PRESSION, JET VER
12	UC-053	SEEGERRING DIN 471-8X0,8	SEEGER-A-RING DIN 471-8 X 0.8	SEEGER-A-RING DIN 471-8 X 0,8	BAGUE SEEGER-A DIN 471-8X0,
13	UC-058	STØTTESKIVE SS 8 X 14 X 1,2	SUPPORTING DISC SS 8X14X1.2	STÜTZSCHEIBE SS 8 X 14 X 1,2	RONDELLE SUPPORT SS 8X14X1,2
14	UC-059	PINOLSKRUE M 3 X 4 RUSTFRI	POINTED SCREW M 3 X 4	REITSTOCKSCHRAUBE M 3 X 4	VIS POINTUE M 3 X 4
15	AC-131	VANDGLASFYLD.RØR INKL. O-RINGE	CUP FILLER TUBE WITH O-RINGS	GLASFÜLLROHR MIT O-RINGEN	TUYAU JET VERRE/JOINTS TORIQ
15	AC-796	SKYLLERØR, VANDGLASFYLDER BC	FLUSHING PIPE,CUP FILLER, BC	SPÜLROHR, GLASFÜLLER, BC	TUYAU RINCAGE,JET DE VERRE,B
16	AC-132	MODTAGER, VANDGLASFYLDERØR	HOLDER, CUP FILLING TUBE	ANSCHLUSSROHR, GLASFÜLLROHR	TUYAU SUPPORT,TUYAU JET VERR
17	MR-209	O-RING 17,1 X 1,6	O-RING 17.1 X 1.6	O-RING 17,1 X 1,6	JOINT TORIQUE 17,1 X 1,6
18	AB-373	KAPPE, NIVEAUINDIKATOR, BA-202	COVER, LEVEL INDICATOR, BA-202	ABDECKUNG,NIVEAUINDIK. BA-202	CAPOR,INDICATEUR NIVEAU BA-2
18	AB-374	KAPPE, NIVEAUINDIKATOR, BA-203	COVER, LEVEL INDICATOR, BA-203	ABDECKUNG,NIVEAUINDIK. BA-203	CAPOR,INDICATEUR NIVEAU BA-2
18	AB-375	KAPPE, NIVEAUINDIKAT. BA-204	COVER, LEVEL INDICATOR, BA-204	ABDECKUNG,NIVEAUINDIK. BA-204	CAPOR,INDICATEUR NIVEAU BA-2
18	AB-376	KAPPE, NIVEAUINDIKAT. BA-205	COVER, LEVEL INDICATOR, BA-205	ABDECKUNG,NIVEAUINDIK. BA-205	CAPOR,INDICATEUR NIVEAU BA-2
18	AB-377	KAPPE, NIVEAUINDIKAT. BA-206	COVER, LEVEL INDICATOR, BA-206	ABDECKUNG,NIVEAUINDIK. BA-206	CAPOR,INDICATEUR NIVEAU BA-2
18	AB-378	KAPPE, NIVEAUINDIKAT. BA-207	COVER, LEVEL INDICATOR, BA-207	ABDECKUNG,NIVEAUINDIK. BA-207	CAPOR,INDICATEUR NIVEAU BA-2
18	AC-901	KAPPE, NIVEAUINDIKATOR, BA-102	COVER, LEVEL INDICATOR, BA-102	ABDECKUNG,NIVEAUINDIKAT.BA-102	CAPOT,INDICATEUR NIVEAU,BA-1
18	AC-253	KAPPE, NIVEAUINDIKATOR, BA-015	COVER, LEVEL INDICATOR, BA-015	ABDECKUNG,NIVEAUINDIKAT.BA-015	CAPOT,INDICATEUR NIVEAU,BA-0
18	AC-224	KAPPE, NIVEAUINDIKATOR, BA-049	COVER, LEVEL INDICATOR, BA-049	ABDECKUNG,NIVEAUINDIK. BA-049	CAPOR,INDICATEUR NIVEAU BA-0
18	AC-225	KAPPE, NIVEAUINDIKATOR, BA-026	COVER, LEVEL INDICATOR, BA-026	ABDECKUNG,NIVEAUINDIK. BA-026	CAPOR,INDICATEUR NIVEAU BA-0
18	AC-226	KAPPE, NIVEAUINDIKATOR, BA-048	COVER, LEVEL INDICATOR, BA-048	ABDECKUNG, NIVEAUINDIK. BA-048	CAPOT,INDICATEUR NIVEAU,BA-0
18	MC-666	KAPPE, SKYLLERØR, BC	COVER, FLUSHING PIPE, BC	ABDECKUNG, SPÜLROHR, BC	CAPOT, TUYAU DE RINCAGE, BC
19	MC-101	AFLØBSSTUDS, FONTÆNE	DRAIN CONNECTOR, CUSPIDOR	ABFLUSSSTUTZEN, SPEIFONTÅNE	RACCORD ECOULEMENT, CRACHOIR
20	MC-102	AFLØBSSTUDS, VANDGLASSKÅL	DRAIN CONNECTOR, W. GLASS BOWL	ABFLUSSSTUTZEN, GLASSCHALE	RACCORD ECOULEM.,REPOSE-VERR
21	MR-076	O-RING Ø 47,22 X 3,53	O-RING 47.22 X 3.53	O-RING 47,22 X 3,53	JOINT TORIQUE 47,22 X 3,53
22	MR-217	O-RING Ø 61,91 X 3,53 NBR 70SH	O-RING Ø 63.09 X 3.53	O-RING Ø 63,09 X 3,53	JOINT TORIQUE Ø 63,09 X 3,53
23	AB-349	FONTÆNEARM, BA-202	ARM, CUSPIDOR, BA-202	ARM, SPEIFONTÅNE, BA-202	BRAS, CRACHOIR, BA-202
23	AB-350	FONTÆNEARM, BA-203	ARM, CUSPIDOR, BA-203	ARM, SPEIFONTÅNE, BA-203	BRAS, CRACHOIR, BA-203
23	AB-351	FONTÆNEARM, BA-204	ARM, CUSPIDOR, BA-204	ARM, SPEIFONTÅNE, BA-204	BRAS, CRACHOIR, BA-204
23	AB-352	FONTÆNEARM, BA-205	ARM, CUSPIDOR, BA-205	ARM, SPEIFONTÅNE, BA-205	BRAS, CRACHOIR, BA-205
23	AB-353	FONTÆNEARM, BA-206	ARM, CUSPIDOR, BA-206	ARM, SPEIFONTÅNE, BA-206	BRAS, CRACHOIR, BA-206
23	AB-354	FONTÆNEARM, BA-207	ARM, CUSPIDOR, BA-207	ARM, SPEIFONTÅNE, BA-207	BRAS, CRACHOIR, BA-207
23	AC-897	FONTÆNEARM, BA-102	ARM, CUSPIDOR, BA-102	ARM, SPEIFONTÅNE, BA-102	BRAS, CRACHOIR, BA-102
23	AC-331	FONTÆNEARM BA-015	ARM, CUSPIDOR, BA-015	ARM, SPEIFONTÅNE, BA-015	BRAS, CRACHOIR, BA-015
23	AC-229	FONTÆNEARM, BA-048	ARM, CUSPIDOR, BA-048	ARM, SPEIFONTÅNE, BA-048	BRAS, CRACHOIR, BA-048
23	AC-227	FONTÆNEARM, BA-049	ARM, CUSPIDOR, BA-049	ARM, SPEIFONTÅNE, BA-049	BRAS, CRACHOIR, BA-049
23	AC-228	FONTÆNEARM, BA-026	ARM, CUSPIDOR, BA-026	ARM, SPEIFONTÅNE, BA-026	BRAS, CRACHOIR, BA-026
24	MC-108	DÆKSEL, FONTÆNEARM	COVER, CUSPIDOR ARM	ABDECKUNG, SPEIFONTÅNE	CARTER, BRAS CRACHOIR
24	ME-037	DÆKPLADE, FONTÆNE M. SUG	PLATE, SUCTION AT THE CUSPIDOR	PLATTE, ABSAUGER AN FONTÅNE	PLAQUE,ASPIR. SUR LE CRACHOI
25	AD-409	PRINT, CUPFILLER SENSOR	* PRINTED CIRC."CUPFILL.SENSOR *	PLATINE "CUPFILLER SENSOR" *	CARTE "CUPFILLER SENSOR"
26	FC-005	GULDFANG	GOLD TRAP	GOLDFÄNGER	TRAPPE À OR
27	UC-059	PINOLSKRUE M 3 X 4 RUSTFRI	POINTED SCREW M 3 X 4	REITSTOCKSCHRAUBE M 3 X 4	VIS POINTUE M 3 X 4
28	MR-075	FILTERRØR	FILTER TUBE	FILTERROHR	TUYAU DE FILTRE
29	AC-517	MODTAGER, FONTÆNESKYLLERØR	HOLDER, CUSPIDOR FLUSH	EINSTECKROHR, SPÜLROHR	TUYAU SUPPORT,TUYAU DE RINCA
30	AC-518	VENTILBLOK, FONTÆNE M/FITTING	VALVE BLOCK,CUSPIROD,+FITTINGS	VENTILBLOCK, FONTÅNE, +FITTINGS	BLOC VANNE,CRACHOIR,+ RACCOR
31	SD-401	FILTERINDSATS, GULDFANG	FILTER, GOLD TRAP	FILTER, GOLDFÄNGER	FILTRE, TRAPPE À OR

23 MS-357 Slange Ø25 Fontænearm Hose Ø25, cuspidor arm Schlauch Ø25 Speifontænearm

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
32	AC-520	AFLØBSRØR, FONT/VAND. M/O-RING	DRAINTUBE, CUSPIDOR, WASSER	ABFLUSSROHR, SPEIFONTÅNE, WASSER	TUYAU ÉCOULEMENT, CRACHOIR, EA
33	UC-538	O-RING Ø 1,78 X 17,17. NITRIL	O-RING Ø 1,78 X 17,17 NITRIL	O-RING Ø 1,78 X 17,17 NITRIL	O-RING Ø 1,78 X 17,17 NITRIL
34	AC-141	KABEL, CUPFILLER SENSE PCB	CABLE, CUP FILLER, PHOTO CELL	KABEL, GLASFÜLLER, FÜHLER	CABLE, JET VERRE, C. PHOTO-SENS
35	AC-415	DÆKSEL, Ø 116 PONTB. M/O-R. HVID	COVER Ø 116, CUSPIDOR BASE	ABDECKUNG Ø 116, FONTÅNEBASIS	CARTER Ø 116, EMBASE CRACHOI
36	UC-068	O-RING Ø 94,5 X Ø 3	O-RING Ø 94.3 X Ø 3	O-RING Ø 94,3 X Ø 3	JOINT TORIQUE Ø 94,3 X Ø 3
37	MC-131	GLIDESKIVE, FONTÆNEARM	SLIDING DISC, CUSPIDOR ARM	GLEITSCHIEBE, SPEIFONTÄNENARM	RONDELLE DE GLISSEMENT, BRAS
38	MC-938	HOLDER, GULDFANG	HOLDER, GOLD TRAP	HALTER, GOLDFÄNGER	SUPPORT, TRAPPE À OR
39	UA-707	SKRUE M 3 X 4 RUSTFRI	SCREW M 3 X 4	SCHRAUBE M 3 X 4	VIS M 3 X 4

AE-620 VANDGLASSKÅL KOMPLET, KERAMIK CUPFILLERBOWL, COMPL, PORCELAIN GLASSCHALE, KOMPL, PORZELAN REPOSE-VERRE, COMPKL. PORCELAIN

AE-621 FONTÆNESKÅL KOMPLET, KERAMIK CUSPIDOR BOWL, COMPL, PORCELAIN SPEISCHALE, KOMPL PORZELAN CUVETTE DU CRACHOIR, PORCELAIN


MC-809 GULDFANG, TYPE 99, INTEGRAL GOLD TRAP, TYPE 99, INTEGRAL GOLDFANGER, TYPE 99, INTEGRAL TRAPPE A OR, TYPE 99, INTEGRAL

ME-036 KAPPE, GULDFANG, TY 99, INTEGRAL CAP, GOLD TRAP, TUPE 99 KAPPE, GOLDFANGER, TYPE 99 CAPOUCHON, TRAPPE A OR, TYPE 99

3	HA-467	Rist for vandglasskål,	Grid for cupfiller bowl,	Rost für Glasschale
3	ME-225	Holder for AE-620, kunststof	Support for AE-620, plastics	Ausser Stützung für AE-620, Kunststoff
6	ME-222	Holder for AE-621, kunststof	Support for AE-621, plastics	Ausser stützung für AE-621, Kunststoff
7	MR-060	O-Ring for keramik vandglasskål Ø44x5	O-Ring for ceramic cupfil.bowl Ø44x5	O-Ring für Porcelain Glasschale Ø44x5
18	AF-117	Kappe, niveauindikat. BA-244,	Cover, level ind. BA-244	Abdeckung, Niveua. BA-244
22	MR-062	O-Ring for keramik fontæneskål Ø54x5	O-Ring for ceramic bowl Ø54x5	O-Ring für Porcelain Speischale Ø54x5
23	MS-357	Slange Ø25 Fontænearm	Hose Ø25, cuspidor arm	Schlauch Ø25 Speifontænearm
23	AF-849	Fontænearm BA-386 sølv	Arm, cuspidor BA-386 silver	Arm, Speifontäne BA-386 Silber



No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
	UC-298	PRODUKTFOD SJ 5302 Ø7.6 X 2.2	PRODUCT FOOT, ROUND Ø 8 X 1.6	PRODUKTFUSS, RUND Ø 8 X 1,6	PIED CIRCULAIRE Ø 8 X 1,6
	AC-756	JUSTERSKRUE+AKSEL, INT., ARM	ADJUSTM. SCREW,AXLE, INT., ARM	JUSTIERSCHR. + ACHSE, INT.,ARM	VIS REGLAGE + AXE BRAS INT.
	SD-532	ARMSYSTEM, U/DÆKSLER CLASSIC *	ARM SYSTEM,WITHOUT COV.CLASSIC	ARMSYSTEM,OHNE DECKPLA.CLASSIC	BRAS SANS CARTERS, CLASSIC
	SD-414	ARMSYSTEM, U/DÆKSLER	* ARM SYSTEM,INT,WITHOUT COVERS*	ARMSYSTEM,INT,OHNE DECKPLATTE*	BRAS INTEGRAL SANS CARTERS
	MC-033	GLIDESKIVE FOR ARM	WASHER FOR ARM	SCHEIBE FÜR ARM	RONDELLE POUR BRAS
1	AC-275	DÆKSEL, BAGARM, ØVERST ,M.LÅS	COVER, REAR ARM, TOP	ABDECKUNG, HINTERARM, OBEN	CARTER, BRAS ARRIERE,AU-DESS
1	AE-944	DÆKSEL,BAGARM,ØVERST ,CLASSIC	COVER, REAR ARM TOP, CLASSIC	ABDECK.HINTERARM,OBEN CLASSIC	CART.BRAS ARRIE.AU-DESSU.CLA
2	AC-276	DÆKSEL, BAGARM,UNDERST,M.LÅSET	COVER, REAR ARM, BOTTOM	ABDECKUNG, HINTERARM, UNTEN	CARTER, BRAS ARRIERE,AU-DESSO
2	AE-945	DÆKSEL, BAGARM,UNDERST,CLASSIC	COVER,REAR ARM,BOTTOM,CLASSIC	ABDECK.HINTERARM,UNTEN,CLASSIC	CARTER, BRAS ARRIERE,CLASSIC
3	AC-175	BREMSE, ARM, KOMPLET	BRAKE, ARM, COMPLETE	BREMSE, ARM, KOMPLETT	FREIN, BRAS, COMPLET
4	HA-147	GASFJEDER, ARM	GAS SPRING, ARM	GASFEDER, ARM	RESSORT A GAZ, BRAS
5	MC-047	PYNTESKIVE, ARMLED	DECORATION DISC, ARM JOINT	DEKORATIONSSCHEIBE, GELENK	RONDELLE DECORATIVE, JOINT
6	UA-060	PINOLSKRUE M 3 X 5	POINTED SCREW M 3 X 5	REITSTOCKSCHRAUBE M 3 X 5	VIS POINTUE M 3 X 5
7	UA-067	SEEGERRING 12 X 1,00 (A12)	SEEGER CIRCLIP 12 X 1.00(A 12)	SEEGERRING 12 X 1,00 (A 12)	CLIP CIRCULAIRE 12 X 1,00(A1
8	MC-028	DÆKSEL, FORARM	COVER, FRONT ARM	ABDECKUNG, VORDERARM	CARTER, BRAS AVANT
8	AE-322	DÆKSEL, FORARM CLASSIC	COVER,FRONT ARM, CLASSIC	ABDECKUNG VORDERARM, CLASSIC	CARTER, BRAS AVANT, CLASSIC
9	AC-527	MODTAGER, LÅSETAP M/FJEDER	RECEIVER, DOWEL + SPRING	EMPFÄNGER, SPERRBOLZEN + FEDER	RECEPT.,GOUJON BLOCAGE+RESSO
10	MC-397	LÅSETAP, ARMDÆKSEL	DOWEL, ARM COVER	SPERRBOLZEN, ARMADECKUNG	GOUJON DE BLOCAGE, CARTER
11	AC-316	KABELSÆT, ARMSYSTEM	CABLE SET, ARM SYSTEM	KABELSATZ, ARMSYSTEM	JEU DE CABLES, SYSTEME BRAS
11	AC-121	MAIN CABLE POWER	MAIN CABLE, POWER	HAUPTKABEL, POWER	CABLE PRINCIPAL PUISSANCE
11	AE-246	KABEL, INT.SIGNAL(FERRIT) APFR	MAIN CABLE,SIGNAL,NTC.	HAUPTKABEL,SIGNAL,NTC.	CABLE PRINCIPAL,SIGNAL,NTC.
12	AC-599	MØTRIK OG SKIVE, FASTG. INSTRB	NUT AND DISC, FIX. INSTR.BRIDGE	SCHRAUBENM.UND SCHEIBE,BEFEST.	ECROU ET RONDELLE, FIXATION
13	UC-505	KABELBØJLE RS1. 7/9 W1.	CABLE CLIPS RS1. 7/9 W1	SCHELLE RS1. 7/9 W1	SERRE-CABLE RS1. 7/9 W1
14	MC-073	DÆKSEL, ARM	COVER, ARM	ABDECKUNG, ARM	CARTER, BRAS
14	AE-445	DÆKSEL, ARM CLASSIC	COVER, ARM CLASSIC	ABDECK., ARM CLASSIC	CARTER, BRAS, CLASS
14	UA-011	SKRUE FOR DÆKSEL	SCREWS FOR COVER	SCHRAUBEN F. ABDECK.	



QF


QF

QF


QF

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
		AC-884 MONT. KIT FOR AC-081 + AC-967	MOUNTING KIT FOR AC-081+AC-967	MONTAGEKIT FÜR AC-081 + AC-967	JEU DE MONTAGE, AC-081+AC-96
1	AC-293	PROGRAMERINGSTAST, INTEGRAL	PROGRAMMING BUTTON	PROGRAMMIERUNGSKNOPF	BOUTON DE PROGRAMMATION
2	AD-402	PRINT, INDICATOR LAMPS, INTBGR.*	PRINTED CIRC.INDICATOR LAMPS*	PLATINE "INDICATOR LAMPS"	* CARTE "INDICATOR LAMPS"
3	AD-423	PRINT, INDICATOR LAMPS, INTBGR.*	PRINTED CIRC.INDICATOR LAMPS *	PLATINE "INDICATOR LAMPS"	* CARTE "INDICATOR LAMPS"
4	MC-111	HÅNDTAG, BRO/LAMPE, 1 STK.	HANDLE, INSTRUMENT BRIDGE	HANDGRIFF, INSTRUMENTENTRÄGER	POIGNEE, PORTE-INSTRUMENT
4	SD-399	HÅNDTAG, BRO/LAMPE, 4 STK.	HANDLES (4) BRIDGE/LAMP	HANDGRIFFE (4) TRÄGER/LAMPE	POIGNEES (4) PORTE-INSTR./ECL
4	AC-371	HÅNDTAG, BRO/LAMPE, 24 STK.	HANDLES (24), BRIDGE/LAMP	HANDGRIFFE (24), TRÄGER/LAMPE	POIGNEES (24)PORTE-INSTR./ECL
5	MC-186	MOTAGER, HÅNDTAG, INSTRUMENTBRO	HOLDER, HANDLE, INSTR. BRIDGE	EINSTECKROHR, GRIFF, INST. TRÄGER	TUYAU SUPP., POIGNEE, PORTE-IN
6	AC-530	DÆKSEL. INSTRUMENTBRO 9018	COVER, INSTRUMENT BRIDGE,	ABDECKUNG, INSTR.TRÄGER,	CARTER, PORTE-INSTRUMENTS,
6	AB-946	DÆKSEL. INSTRUMENTBRO CLASSIC	COVER, INSTR.BRIDGE, LAMPS,	ABDECK., INSTR.TRÄGER, LAMP,	CARTER, PORTE-INSTR, LAMPES,
7	AD-408	PRINT INSTRUMENT SELECTION	* PRINTED CIRC."INSTR.SELECT."	* PLATINE "INSTRUMENT SELECT."	* CARTE "INSTRUMENT SELECTION"
8	AD-404	PRINT, INSTR. CONTROL, BC/FC	* PRINTED CIRC.INSTRUMENT CONT *	* PLATINE "INSTRUMENT CONTROL"	* CARTE "INSTRUMENT CONTROL"
8	AD-444	PRINT, INSTR. CONTROL NYC	* PR.CIRCUIT, INSTR.CONTROL, NTC *	* PLATINE, INSTR.CONTROL, NTC *	* CARTE, INSTRUMENT CONTROL, NTC
9	AC-081	SUGESLANGEOPHÆNG, INT. HØJRE	SUCTION HOSE HOLDER, RIGHT	SCHLAUCHHALTER, RECHTER	PORTE-TUYAUX DROIT
9	AC-967	SUGESLANGEOPHÆNG, INT. VENSTRE	SUCTION HOSE HOLDER, LEFT	SCHLAUCHHALTER, LINKER	PORTE-TUYAU GAUCHE
10	AC-533	GENNEMFØRING, SUGETELESKOP	LEAD-IN BUSH, ASPIR. TELESCOPE	FÜHRUNG, SAUGTELESKOP	TUBE DE LIAISON, TELESCOPE AS
11	AC-147	KABEL, SUGESLANGEHOLDER, INT.	CABLE, SUCTION TUBE HOLDER	KABEL, SAUGSCHLAUCHHALTER	CABLE, PORTE-TUYAUX D'ASPIR.
12	AD-403	PRINT, SUGESLANGEHOLDER	* PRINTED CIRC."SUCTION HOLDER"*	* PLATINE "SUCTION HOLDER"	* CARTE "SUCTION HOLDER"
13	MC-128	HOLDER, SUGENIPLER	SUCTION HOSE HOLDER	SAUGSCHLAUCHHALTER	PORTE-TUYAUX D'ASPIRATION
14	AB-232	BESLAG, SUGESLANGEHOLD. KOMPL.	FITTING, SUCTION HOLDER, COMP.	BESCHLAG, SCHLAUCHHALTER, KOM.	GARNITURE COMPL., PORTE-TUYAU
15	MC-458	GEVINDTAP, INSTR. RULLE FJEDER	THREADED PIN, ROCKER PIECE	GEWINDEZAPPEN, KIPPTEILFEDER	TENON TARAUDE, PIECE BASCULE
16	ML-008	TRYKKONTAKT, FONTÆNE	PUSH BUTTON SWITCH, CUSPIDOR	DRUCKSCHALTER, SPEIFONTÂNE	INTERRUPT.A PRESSION, CRACHOI
17	MC-804	MØTRIK FOR FORLED, BRO	NUT FOR INSTRUMENTBRIDGE	MUTTER FÜR VORDERGLIED	ECROU,
18	MC-805	SKIVE FOR FORLED, BRO	WASHER, INSTRUMENTBRIDGE	SCHIBE, VORDERGLIED, BRÜCKE	RONDELLE, TABLET
19	UC-843	CH SKRUE 6 KANT M3 X 6 RUSTFRI	SCREW, M3 X 6, STAINLES	SCHRAUBE, M3 X 6, ROSTFREI	VIS, M3 X 6, INOX
20	AC-595	SLANGEOPHÆNG + SJØDINGSPRØJTE	SUSPENSION + SJØDING SYRINGE	NUR FUER SCHWEDEN	ONLY FOR SWEDEN
20	AC-893	OPHÆNG FOR SJØDINGSPRØJTE	SUSPENSION FOR SJØDING SYRINGE	NUR FUER SCHWEDEN	ONLY FOR SWEDEN
20	AC-035	BALANCEOPHÆNG, 4+4, INTEGRAL	* BALANCED SUSPENSION, 4+4	* BALANCEAUFHÄNGUNG, 4+4	* FOUET, 4+4
20	AB-271	MOTORS LANGE, FAST, INTEGRAL	* MOTOR TUBE, FLEX MOTOR	* MOTORSCHLAUCH, FLEX MOTOR	* CORDON MOTEUR FLEX INTEGRAL
20	AC-771	BALANCEOPHÆNG SUPPRAS. SCALER	BALANCED SUSPENSION, SUPRASSON	BALANCEAUFHÄNGUNG, SUPRASSON	FOUET SUPRASSON
20	SD-228	BALANCEOPHÆNG, 3 FKT.SPRØJTE	* BALANCED SUSPENSION, 3-WAY SYR*	BALANCEAUFHÄNG., 3-FUNKT SPR.*	* FOUET, SERINGUE, 3 FONCTIONS
20	SD-229	BALANCEOPHÆNG, 6 FKT.SPRØJTE	* BALANCED SUSPENSION, 6-WAY SYR*	BALANCEAUFHÄNG., 6-FUNKT.SPR *	* FOUET, SERINGUE, 6 FONCTIONS
20	AB-753	BALANCEOPHÆNG, TURBINE MIDWEST*	BALANCED SUSPENSION, MID-WEST	BALANCEAUFHÄNGUNG, MID-WEST	FOUET, MID-WEST
21	AC-187	INSTRUMENTANLÆG, INTEGRAL	INSTRUMENT PAD WITH MAGNETS	INSTRUMENTENABLAGE + MAGNETE	REPOSE-INSTRUMENT AVEC AIMAN
21	AC-387	INSTRUMENTANLÆG, INTEGRAL, 6 STK	6 INSTRUMENT PADS	6 INSTRUMENTENABLAGEN	6 REPOSE-INSTRUMENTS
22	MC-352	DÆKSEL, SUGETELESKOP	COVER, SUCTION TELESCOPE	ABDECKUNG, SAUGTELESKOP	CARTER, TELESCOPE ASPIRATION
23	AC-599	MØTRIK OG SKIVE, FASTG. INSTRB	NUT AND DISC, FIX. INSTR.BRIDGE	SCHRAUBENM.UND SCHIBE, BEFEST.	ECROU ET RONDELLE, FIXATION
24	UC-844	SEBGER PÅSKIVE Ø12 X Ø18 X 1	SEBGER CIRCULAR CLIP Ø12XØ18X1	SEBGERRING Ø12 X Ø18 X 1	CLIP CIRCULAIRE Ø12 X Ø18 X
25	MC-457	FJEDER DT=0,8 D=6 L0=30,6	SPRING	FEDER	RESSORT
26	AB-245	RULLE M.MAGN.AC-035/SD-229/28	ROCKER PIECE, AC-035/SD-29/28	KIPPTEIL, AC-03/SD-229/28	PIECE BASCUL.AC-035/SD-229/2
26	AB-264	RULLE M.MAGNET MOTOR (AB-271)	ROCKER PIECE, FOR AB-271	KIPPTEIL FÜR AB-271	PIECE A BASCULER, AB-271
27	MC-109	GUMMILISTE, INSTRUMENTBORD	RUBBER, INSTRUMENT TABLE	GUMMILISTE, INSTRUMENTENTISCH	LISTBAU EN CAOUTCHOUC, TABLET
28	AB-276	KIT, SUGESLANGEHOLDER, BRO	KIT, SUCTION HOSE HOLDER, BRIDGE	KIT, SAUGSCHLAUCHHALTER, TRÄGER	JEU, PORTE-TUYAUX, PORTE-INSTR
29	ML-061	SIKRING TRÆG 1,6A	FUSE T 1.6 A	SICHERUNG T 1,6 A	FUSIBLE T 1,6 A
30	MH-430	SIKRING FLINK 6,3 A	FUSE F 6.3 A	SICHERUNG F 6,3 A	FUSIBLE F 6,3 A
31	AB-263	CODE CHIP, TANDRENSER TYPE O	CODE CHIP, SCALER TYPE O	CODE CHIP, ZEG TYP O	BLOC-CODE, DETARTEUR TYPE O
32	MC-176	DÆKSEL, OLIESUMP, INST. BRO	COVER, OIL WELL, INSTR.BRIDGE	ABDECKUNG, INSTRUMENTENTRÄGER	CARTER, PORTE-INSTRUMENTS
33	SD-495	KIT, AFDÆKNING AF AD-408	KIT. COVER AD-408	KIT. ABDECKUNG AD-408	KIT. CARTER AD-408
34	UC-216	FJEDER, SUGOPHÆNG, INTEGRAL	SPRING, SUCTION HOSE, INTEGRAL	FEDER SAUGSCHLAUCH, INTEGRAL	RESSORT TOYAUX D'ASPIRATION
35	ME-023	CLIPS, INSTRUMENTOPHÆNG, BRO	CLIPS, BALANCESUSPENSION, BRIDGE	CLIPS, BALANCEAUFHANG, BRÜCKE	CLIPS, FOUET, PORTE-INSTRUMENT
36	UC-532	SKRUE Ø 4 X 16 PH, PT-SKRUE	SCREW Ø 4 X 16 PH, PT-SCREW	SCHRAUBE Ø4 X 16 PH, PT-SCHRAUB	VIS Ø 4 X 16 PH, PT-VIS
37	UC-533	SKRUE Ø 4 X 12 PH, PT-SKRUE	SCREW Ø 4 X 12 PH, PT-SCREW	SCHRAUBE Ø4X12 PH, PT-SCHRAUBE	VIS Ø 4 X 12 PH, PT-VIS
38	AB-277	KIT, KABEL, SUGESLANGEOPHÆNG	KIT, CABLE, SUCTION TUBE HOLDER	KIT, KABEL, SAUGSCHLAUCHHALTER	KIT, CABLE, PORTE-TUYAUX D'ASP
39	MC-459	SVAMP, OLIESUMP, INTEGRAL	SPONGE, OIL WELL	SCHWAMM, ÖLWANNE	EPONGE, CHAMBRE A HUILE


No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
40	AE-915	KIT, GEVINDTAP OG M-CLIPS M.M.	KIT, BALANCED SUSPENSION	KIT, BALANCEAUFHÄNGUNGEN	KIT DE FOUET
41	AE-281	CODE CHIP, FLEX SCALER	FLEX SCALER CODE CHIP	FLEX SCALER/ZEG CODE CHIP	FLEX SCALER CODE CHIP
42	MC-351	DÆKSEL, ARM TIL SKÆRM	COVER, SCREEN ARM	ABDECKUNG, SCHIRM ARM	OPERCULE, BRAS ECRAN
42	AE-452	DÆKSEL, ARM TIL SKÆRM, CLASSIC	COVER, SCREEN ARM	ABDECKUNG, SCHIRM ARM	OPERCULE, BASE ECRAN
43	MC-948	BLÆNDDÆKSEL, OPHÆNG	COVER PLATE, INSTR. SUSPENSION	VERBLENDPLATTE, INST. AUFHÄNGUNG	CACHE FOUET
43	AE-615	BLÆNDDÆKSEL, OPHÆNG, CLASSIC	COVER PLATE, INST. SUSP. CLASSIC	VERBLENDPLATTE, AUFHÄNGUNG. CLASS	CACHET, FOUET, CLASSIC
44	MC-949	LÅSEPLADE FOR BLÆNDDÆKSEL	LOCKING PLATE, COVER PLATE	HALTEBLECH, VERBLENDPLATTE	PLAQUETTE VERROUILLAGE, CACH
45	MC-261	CLIPS FOR TRISSE, SUGEOPHÆNG	CLIPS, PULLEY, SUCTION HOSE	CLIPS, ROLLE, ABSAUGSSCHLAUCH	CLIPS, POULIE CABLE SUPP. ASP
46	MC-266	LEJEBØSNING FOR TRISSE TIL SUG	BEARING BUSH, SUCTION PULLEY	LAGERBUCHSE, ABSAUGSROLLE	DOUILLE-ROULEMENT, POULIE ASP



No	NUMBER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	SA-555	ANKER	PLUNGER	PLUNGER	PLONGEUR
2	SA-556	FJEDER	SPRING	FEDER	RESSORT
3	AC-103	TRYKMÅLER KOMPLET	* PRESSURE GAUGE, COMPLETE	* DRUCKMESSER, KOMPLETT	* INDICATEUR PRESSION COMPL.
4	AC-114	DYSE, FLOWMÅLER, MED O-RINGE	NOZZLE, FLOW METER, + O-RINGS	DÜSE, MENGENMESSER, + O-RINGE	BUSE, DEBITMETRE, + JOINTS TOR
5	UC-162	O-RING Ø 4,6 X 2,0	O-RING Ø 4.6 X 2.0	O-RING Ø 4,6 X 2,0	JOINT TORIQUE Ø 4,6 X 2,0
6	AC-115	REGULERINGSVENTIL VAND/LUFT	REGULATING VALVE WATER/AIR	REGULIERVENTIL WASSER/LUFT	VANNE DE REGLAGE EAU/AIR
7	MR-227	O-RING Ø 3 X 1	O-RING Ø 3 X 1	O-RING Ø 3 X 1	JOINT TORIQUE Ø 3 X 1
8	AB-972	PROPORTIONAL VAND VENTIL, R	PROPORTIONAL VALVE - WATER	PROPORTIONAL VENTIL, WASSER	ELECTROVANNE, REGLABLE; EAU
9	AC-164	MAGNETVENTIL, 2 VEJS, M/STIK	SOLENOID VALVE, 2-WAY/PLUG	MAGNETVENTIL, 2-WEG/STECKER	ELEVTVANNE, 2-VOIES/PRISE
10	SD-303	O-RING, 12,5 X 1,1 N -O B75	O-RING	O-RING	JOINT TORIQUE
11	SD-307	O-RING, 12 X 1,5 B	O-RING	O-RING	JOINT TORIQUE
12	SD-308	O-RING	O-RING	O-RING	JOINT TORIQUE
13	MB-382	FITTINGS, BANJO, 6/4 X M 5	FITTINGS, BANJO, 6/4 X M 5	FITTINGS, BANJO, 6/4 X M 5	RACCORD, BANJO, 6/4 X M 5
14	AB-973	PROPORTIONAL LUFT VENTIL, R	PROPORTIONAL VALVE, AIR	PROPORTIONALVENTIL, LUFT	ELECTROVANNE, REGLABLE, AIR
15	AB-265	FILTERBLOK, INSTRUMENTBRO	FILTER BLOCK, INSTRU.BRIDGE	FILTERBLOCK, INSTRUM.TRÄGER	BLOC FILTRE, PORTE-INSTRUMENT
16	MR-203	O-RING 11,5 X 1,5	O-RING 11.5 X 1.5	O-RING 11,5 X 1.5	JOINT TORIQUE 11,5 X 1,5
17	AC-011	SELECTVENTIL KOMPLET	SELECT VALVE	SELECT VENTIL	VANNE SELECT
18	AC-568	SELECTVENTIL, U/MAGNETVENTIL	* VALVE BLOCK WITHOUT SOLENOID	* VENTILBLOCK OHNE MAGNETVENT.	* BLOC VANNE SANS ELECTROVANNE
19	AC-165	MAGNETVENTIL, 3 VEJS, M/STIK	SOLENOID VALVE, 3-WAY/CONNECTOR	MAGNETVENTIL, 3-WEG/STECKER	ELECTROVANNE, 3-VOIES/PRISE
20	MR-150	PAKNING, SÆDE, SELECT VENTIL	GASKET, SEAT, SELECT VALVE	DICHTUNG, SITZ, SELECT VENTIL	JOINT, SIEGE, VALVE SELECT
21	AC-076	VENTILBLOK INSTRUMENTBRO	VALVE BLOCK, INSTRUMENT BRIDGE	VENTILBLOCK, INSTRUMENTENTRÄGER	BLOCK VANNE, PORTE-INSTRUMEN
22	MR-135	FILTER ELEMENT, VAND	FILTER, WATER	FILTER, WASSER	FILTRE A EAU
23	MC-605	FILTERHOLDER, BRO, INTEGRAL	FILTER HOLDER, BRIDGE,	FILTERHALTER, INSTR.TRÄGER,	SUPPORT DU FILTRE, PORTE-INST
24	MB-329	FITTING M5 6/4 BANJO	BANJO FITTING M5 6/4	BANJOFITTING M5 6/4	RACCORD BANJO M5 6/4
25	UC-807	GUMMITÆTNING Ø8 X Ø5 X 1.5	RUBBER SEAL Ø8 X Ø5 X 1.5	GUMMIABDICHTUNG Ø8 X Ø5 X 1,5	JOINT CAOUTCHOUC, Ø8 X Ø5 X 1
26	UC-433	BANJO BOLT M5 X 11	FITTINGS, BANJO, M5 X 11	FITTINGS, BANJO, M5 X 11	RACCORD, BANJO, M5 X 11
27	AC-745	VENTILBLOK M. KUGLER OG NIPPEL	VALVEBLOCK	VENTILBLOCK	BLOC, VANNE
28	MC-176	DÆKSEL, OLIESUMP, INST. BRO	COVER, OIL WELL, INSTR.BRIDGE	ABDECKUNG, INSTRUMENTENTRÄGER	CARTER, PORTE-INSTRUMENTS
29	MC-459	SVAMP, OLIESUMP, INTEGRAL	SPONGE, OIL WELL	SCHWAMM, ÖLWANNE	EPONGE, CHAMBRE A HUILE
30	AC-632	FILTER M.KRAVE, HYGIFLEX THERMO	FILTER WITH RIM, HYGIFL. THERMO	FILTER MIT FLANSCH, HYG.THERMO	FILTRE AVEC COLLET, HYG.THERM
30	UC-523	FILTERHOLDER	SUPPORT FOR FILTER	HALTER FÜR FILTER	



No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
	MA-354	BLIK FOR STIK INCL. EL.ZINK	LOCK FOR EL. CONNECTOR	VERSCHLUSS FÜR EL. STECKER	SERRURE POUR PRISE ELECTRIQUE
	BA-062	SALTVANDSPOSER 8 X 1 LITER	8XBAGS(1/1 L) STERILE SALINE S	8 X 1/1 L STERILE NACL-LÖSUNG	8XPOCHES(1/1 L)SERUM PHYSIOL
	BA-061	SPRAYSLANGER 50 STK KIRURGI	50 TUBE SETS FOR SURGERY PLANT	50 SCHLAUCHSETS, CHIRURGIE	50 TUYAUX IRRIGATION,BISTOUR
	AC-432	BASIS BØJLE/HÅNDTAG KIRURGI	BASE SUPPORT/HANDLE, SURGERY	BASISBÜGEL/GRIFF, CHIRURGIE	BASE ETRIER/POIGNEE,CHIRURGI
1	AC-542	KLEMMVENTIL, SALTVAND M/KABEL	PINCH VALVE, SURGERY, + CABLE	ABKLEMMUNG, CHIRURGIE, + KABEL	SYST.PINCEMENT,CHIRUR.,+CABL
2	UC-217	TRYKMANCHET, KIRURGI, GRÅ	SLEEVE, SURGERY, GREY	MANSCHETTE, CHIRURGIE, GRAU	MANCHON, CHIRURGIE, GRIS
3	MC-551	CLIPS, SALTVANDSSLANGE	CLIPS, SALINE SOLUTION TUBE	CLIPS, NACL-LÖSUNG SCHLAUCH	CLIPS, CORDON SERUM PHYSIOL.
4	MH-107	KABEL, FLEXVIEW/KIRURGI VENTIL	CABEL, COILED, FLEXVIEW	SPIRALLEITUNG, FLEXVIEW	CORDON TÉLÉPHONE, FLEXVIEW
5	AD-422	SURGICAL SALINE CONNECTOR	PLUG AND CABLE, SURGERY	STECKER UND KABEL, CHIRURGIE	PRISE ET CABLE, CHIRURGIE
6	AC-433	BØJLE, KIRURGIANLÆG	SUPPORT, SURGERY	BÜGEL, CHIRURGIE	ETRIER, CHIRURGIE
7	UC-218	KUGLELÅS M 4 X 10	BALL BLOCKING M 4 X 10	KUGELSCHLOSS M 4 X 10	OBTURATION A BILLES M 4 X 10
8	AC-566	STYREPLADE M/GUMMISTOP,KIRURGI	GUIDE PLATE W. STOP,SURGERY	FÜHRUNGSPLATTE M.STOPP,KIRURG.	PLAQUE GUIDAGE A/ARRET,CHIRU
8	AC-405	BASIS, SALTVANDSANLÆG, INT.	BASE, SALINE SOLUTION,INTEGRAL	BASIS, SALZLÖSUNG, INTEGRAL	EMBASE,SERUM PHYSIOLOGIQUE,I
9	SD-470	PUMPE FOR TRYKMANCHET, KIRURGI	PUMP, SURGERY PLANT	PUMPE, CHIRURGIEVORRICHTUNG	POMPE DU MANCHON, CHIRURGIE



Amalgamseparator: Metasys. Afløbspumpe

Amalgam Separator: Metasys. Drain Pump


YA-001

Amalgamseparator: Metasys. Abflusspumpe

Séparateur d'amalgame : Metasys. Pompe d'écoulement


0.6.12.B

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	AC-447	AFLØBSPUMPE Ø 25-Ø 25	OUTLET PUMP Ø 25 - Ø 25	ABFLUSSPUMPE Ø 25 - Ø 25	POMPE D'ÉCOULEMENT Ø 25 - Ø 25
2	AC-586	KABINET, AMALGAMUDSK. METASYS	CABINET, AMALGAM SEP., METASYS	GEHÄUSE, AMALGAMABSCHIEDER, MET.	COFFRET, SEP. AMALGAME, METASYS
2	AC-451	KABINET, AMALGAMUDSK. MERCURY M.	CABINET, AMALGAM SEP. MERCURY	GEHÄUSE, AMALGAMSEP. MERCURY	COFFRET, SEP. AMALGAME MERCURY
3	MC-590	DÆKSEL, STANDERF., M/KABELUDG.	COVER, FOOT, WITH CABLE OUTLET	ABDECKUNG, FUSS, KABELAUSGANG	CARTER, EMBASE, SORTIE CABLE
4	AC-401	SLANGESTUDS Ø 25/Ø 25, M. O-R.	TUBE CONNECTOR Ø25/Ø25+O-RINGS	SCHLAUCHSTUTZEN Ø25/Ø25+O-RING	RACCORD TUYAU Ø25/Ø25+JOINT
5	AC-248	SLANGESTUDS Ø32/25 M.O-RINGE	TUBE CONNECTOR Ø 25 + O-RINGS	SCHLAUCHSTUTZEN Ø 25 + O-RINGS	RACCORD TUYAU Ø 25+JOINTS TO
6	AC-402	MELLEMSTYKKE, AFLØBSSAMLER/Ø 25	INTERMED. PIECE, DRAIN JUNCTION	ZWISCHENSTÜCK, ABFLUSSSAMMLER	PIECE INTERMED., JONCTION ECO
7	AC-384	LÅSESKRUE PUSH-TURN KOMPLET	LOCK SCREW PUSH-TURN	VERSCHLUSSSCHRAUBE	VIS DE BLOCAGE
8	SC-688	SPENDBÅND 13-20, RUSTFRI	TIGHTEN. STRAP 13-20, RUSTPROOF	SPANNBAND 13-20, ROSTFREI	RUBAN DE SERRAGE 13-20 INOX.
9	UC-279	BETJENINGSPANEL, METASYS	CONTROL PANEL, METASYS	KONTROLLPANEL, METASYS	PANNEAU DE CONTRÔLE METASYS
10	UC-281	CENTRIFUGE, METASYS	CENTRIFUGE, METASYS	ZENTRIFUGE, METASYS	CONTRIFUGE METASYS
11	UC-282	SEPARATIONSDEL, METASYS	PART, SEPARATION AUT., METASYS	TEIL, SEPARIERAUTOMAT, METASYS	PIECE, SEPARATEUR AUT. METASYS
12	UC-283	BASIS, PRINT, METASYS	BASE, PRINTED CIRCUIT, METASYS	BASIS, PLATINE, METASYS	BASE, CARTE, METASYS
13	UC-500	ALMAGAM AFFALDS BEHOLD. METASYS	AMALGAM COLLECTOR	AMALGAMSAMMLER	RECEPTACLE D'AMALGAME
14	MH-539	TRAFO METASYS	TRANSFORMER METASYS	TRANSFORMATOR METASYS	TRANSFORMATEUR METASYS
15	SD-402	SI FOR METASYS (10 STK.)	STRAINER, METASYS (10)	SIEB, METASYS (10)	FILTRE METASYS (10)
16	MC-525	BESLAG FOR AFLØBSPUMPE	FITTINGS, OUTLET PUMP	BESCHLAG, ABFLUSSPUMPE	RACCORD, POMPE D'ÉCOULEMENT
17	MC-563	BESLAG, METASYS AMALGAMOPSAMLER	FITTINGS, METASYS AMALGAM COLL.	BESCHLAG, METASYS AMALGAMSAM.	RACCORD, RECEPTACLE AMALGAME
18	AC-469	DÆKSEL, AMALGAMOPSAMLER	COVER, AMALGAM COLLECTOR	ABDECKUNG, AMALGAMSAMMLER	CARTER, RECEPTACLE D'AMALGAM
19	MC-063	DÆKSEL, STANDERF.	COVER, FOOT	ABDECKUNG, FUSS	CARTER, EMBASE
20	AC-436	SLANGESÆT, METASYS AMALGAMOPS.	TUBE SET, METASYS AMALGAM COL.	SCHLAUCHSATZ, METASYS AMALGAMA.	JEU DE TUYAUX, RECEPT. AMALGA
21	AC-965	SAMLESTYKKE KOMPLET, II	CONNECTION COMPLETE, II	VERBINDUNGSTÜCK KOMPLETT, II	RACCORD COMPLET, II
22	ML-035	SIKRING, PLUG-IN 3,15AT Ø10	FUSE, PLUG-IN 3,15AT Ø 10	SICHERUNG	FUSIBLE



No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
	AB-166	KIT, SORT GUMMI FODKONT./STOLII	KIT, BLACK RUBBER, FOOT CONTROL KIT, FUSSSCHALTER, SCHWARZ, GUMMI	KIT, CAOUTCHOUC NOIR, PEDALE	
1	AC-176	KUGLEPEDAL, FODKONTAKT, KOMP, INT	BALL PEDAL, FOOT CONTR, COMPLETE	KUGLEPEDAL, FUSSANLASSER, KOMP.	PEDAL BOULE, PEDALE, COMPLET
2	AC-142	KABEL FODKONTAKT	CABLE, FOOT CONTROL	KABEL, FUSSANLASSER	CABLE, PEDALE
3	MC-177	FJEDER, VIPPEDEL, FODKONT	SPRING, ROCKER PIECE, FOOT CONT.	FEDER, PEDALMECHANIK, FUSSANLAS.	RESSORT, PIECE A BASCUL., PED.
4	AB-129	SÆT, PRINT, FODKONTAKT	* PR. CIRCUIT, SET, FOOT CONTROL*	PLATINENSATZ, FUSSANLASSER	* CARTE, JEU, PEDALE
5	MC-370	DÆKSEL FODKONTAKT INTEGRAL U/B	FOOT PEDAL COVER, INTEGRAL	FUSSPEDALABDECKUNG, INTEG.	CAPOT PEDALE SANS ETRIER, INT
5	AC-181	DÆKSEL, FODKONTAKT, FC/NTC M/B	FOOT PEDAL COVER, FC/NTC	FUSSPEDALABDECKUNG, FC/NTC	CAPOT PEDALE, FC/NTC
5	AB-111	DÆKSEL, FODKONTAKT M/BØJLE, BC	FOOT PEDAL COVER WITH HOOP, BC	FUSSANLASSERABD. MIT BÜGEL, BC	CAPOT PEDALE AVEC ETRIER, BC
6	AC-457	DÆKSEL FODK. UNDERST M/SKRIDG.	LOWER FOOT COVER, NON-SKID RUB.	UNTERE FUSSDECKEL, FRIKT. GUMMI	CAPOT PEDALE BAS, AVEC CAOUTC
7	MR-101	SKRIDGUMMI, FIRKANTET	NON-SKID RUBBER, SQUARE	FRIKTIONSGUMMI, VIERECKIG	CAOUTCHOUC SOUS PEDALE, CARR
8	MC-346	DÆKKAPPE, KUGLEPEDAL	RUBBER COVER, BALL PEDAL	GUMMIABDECKUNG, KUGLEPEDAL	CAPOT EN CAOUTC., PEDALE BOUL
9	MC-347	DÆKKAPPE, FODPEDAL	COVER, FOOT PEDAL	ABDECKUNG, FUSSANLASSER	CARTER, PEDALE
10	AC-528	LÅSEMEKANISME, FODKONTAKTBØJLE	LOCKING DEVICE, FOOT PEDAL HOOP	SPERRE, TRAGEBÜGEL, FUSSANLASSER	FERMETURE, ETRIER DE LA PEDAL
11	AC-529	FODRULLE M/SKRUE OG SKIVE	ROLL + SCREW AND DISC	ROLLE + SCHRAUBE UND SCHEIBE	ROULETTE + VIS ET RONDELLE
12	MC-359	KNOP FOR FODKONTAKTBØJLE HØJRE	KNOB, HOOP ON PEDAL, RIGHT	KNOPF, FUSSANLASSERBÜGEL, RECHTS	BOUTON, ETRIER DROIT
13	MC-363	KNOP F. FODKONTAKTBØJLE, VENSTRE	KNOB, HOOP ON PEDAL, LEFT	KNOPF, FUSSANLASSERBÜGEL, LINKS	BOUTON, ETRIER GAUCHE
14	MC-361	FODKONTAKTBØJLE, NAT. ELOX.	HOOP, FOOT CONTROL, NAT. ANO.	FUSSANLASSERBÜGEL, ELOX. NATUR	ETRIER ANODISE TEINTE NATUR.
15	AB-064	ADAPTOR, FODKONTAKT, BC	ADAPTOR, FOOT CONTROL, BC	ADAPTER, FUSSANLASSER, BC	ADAPTATEUR, PEDALE, BC
16	MC-367	FJEDER F. VIPPEARM, FODK. BØJLE	SPRING, LEVEL ARM, FOOTCON. STRAP	FEDER, KIPPHEBEL, FUSSANL. BÜGEL	RESSORT, LEVIER BASC, ETRIER P
17	AC-097	FODKONTAKT, FC/NTC	FOOT CONTROL, FC/NTC	FUSSANLASSER, FC/NTC	PEDALE FC/NTC
17	AB-110	FODKONTAKT M. BØJLE, BC	FOOT CONTROL WITH HOOP, BC	FUSSANLASSER MIT BÜGEL, BC	PEDALE AVEC ETRIER, BC
17	AC-743	FODKONTAKT, BC	FOOT CONTROL, BC	FUSSANLASSER, BC	PEDALE BC
18	SD-494	KIT, MODIFIKATION AF FODKONT.	KIT, MODIFICATION FOOT CONTR.	KIT, ÄNDERUNG VON FUSSANLASSER	KIT, MODIFICATION DE PEDALE
19	AC-179	VIPPEDEL, FODKONTAKT, INTEGRAL	ROCKER PIECE, FOOT CONTROL	PEDALMECHANIK, FUSSANLASSER	PIECE A BASCUL., PEDALE
20	AC-378	AKSEL FOR FK BØJLE HØJRE	SPINDLE FOR HOOP, FOOT C. RIGHT	ACHSE FÜR BÜGEL, FUSSANL. RECHT	AXE ETRIER PEDALE, DROIT
21	AC-379	AKSEL FOR FK BØJLE VENSTRE	SPINDLE FOR HOOP, FOOT C., LEFT	ACHSE FÜR BÜGEL, FUSSANL. LINKS	AXE ETRIER PEDALE, GAUCHE
22	MC-678	DÆKSEL, FODKONTAKT U/BØJLE, BC	FOOT PEDAL COVER, NO HOOP, BC	FUSSPEDALABD. OHNE BÜGEL, BC	CAPOT PEDALE SANS ETRIER, BC
23	MC-984	BØSNING Ø5 X Ø4 X 11,7	BUSH, Ø5 X Ø4 X 11,7	BUCHSE, Ø5 X Ø4 X 11,7	DUILLE, Ø5 X Ø4 X 11,7
24	UA-012	UNDERLAGSSKIVE Ø 4 Z	WASCHER Ø 4Z	SCHEIBE Ø 4Z	RONDELLE Ø 4Z
25	UA-249	SKRUE M 4 X 16 CH. INDV. 6KT. Z	SCREW M 4 X 16	SCHRAUBE M 4 X 16	VIS M 4 X 16
26	UC-149	FIBERSKIVE Ø10 X Ø5,1 X 1	FIBRE DISC Ø 10 X Ø 5,1 X 1	FIBERSCHEIBE Ø 10 X Ø 5,1 X 1	RONDELLE FIBRE Ø 10 X Ø 5,1 X
27	AB-165	AFDÆKNING, KUGLEPEDAL, LAKERET	COVER JOYSTICK FOOT PEDAL	ABDECKUNG, KUGLEPEDAL, FUSSANLA.	CAPOT PEDALE JOYSTICK
28	MC-988	TRAVERS, AFDÆKNING, FODK. BC	BC FOOT-CONT. COVER CROSS-BAR	QUERTRÄGER, BC FUSSANLASSER	TRAVERSE CAPOT PEDALE BC

UC-054 Stift for AC-378/AC-379 Ø3x6 Pin for AC-378/AC-379 Ø3x6 Stift für AC-378/AC-379 Ø3x6mm



Instrumenter: QF-kobling. Turbine type 89, MC2

Instruments: QF Coupling. Turbine Type 89, MC2


YA-001

Instrumente: QF-Kupplung. Turbine Typ 89, MC2r

Instruments : Raccord QF. Turbine type 89,MC2

0.6.14.B

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	SC-740	O-RING, UNIFIX/QF, STOR	O-RING, UNIFIX COUPLING, BIG	O-RING, UNIFIX KUPPLUNG, ,GROSS	JOINT TORIQUE,RAC.UNIFIX,GRA
2	SD-276	KAPPE, MOTOR INTRA	MOTOR COVER, INTRA	MOTORGEHÄUSE, INTRA	CAPOT DU MOTEUR, INTRA
3	SA-024	O-RING, BIEN AIR/FLEX MOTOR	O-RING, BIEN AIR/FLEX MOTOR	O-RING, BIEN AIR/FLEX MOTOR	JOINT TOR.,MOTEUR BIEN A./FL
4	HB-005	PÆRE, LUX MOTOR/TURBINE	BULB, LUX MOTOR/TURBINE	BIRNE, LUX MOTOR/TURBINE	AMPOULE, MOTEUR LUX/TURBINE
5	SD-005	KUL MED FJEDER, FLEX MOTOR	CARBON+SPRING, FLEX MOTOR	KOHLEBÜRSTE+FEDER, FLEX MOTOR	CHARBONS+RESSORT, MOTEUR FLE
6	MC-191	MELLEMKAPPE,QF 4+4 MOTOR HT065	INTERMEDIATE COVER, MOTOR	ZWISCHENKAPPE, MOTOR	CARTER INTERMEDIARE, MOTEUR
7	SA-042	VÆRKTØJ, KULSKIFT	TOOL, CHANGE OF CARBON	WERKZEUG, KOHLEN AUSWECHSELN	TIGE, CHANGEMENT DE CARBON
8	HT-065	KOBLING 4+4, QF MOTOR	* COUPLING 4+4,QUICKFLEX MOTOR *	* KUPPLUNG 4+4,QUICKFLEX MOTOR *	* RACCORD 4+4,MOTEUR QUICKFLEX
9	SD-339	PAKNING, MOTORKOBLING	GASKET, MOTOR COUPLING	DICHTUNG, MOTORKUPPLUNG	JOINT, RACCORD MOTEUR
10	SD-216	KAPPE, MOTOR ISO	MOTOR COVER, ISO	MOTORGEHÄUSE, ISO	CAPOT DE MOTOR, ISO
11	SD-439	O-RING FOR FLEX MOTOR	O-RING FOR FLEX MOTOR	O-RING FÜR FLEX MOTOR	JOINT TORIQUE P. MOTEUR FLEX
12	SD-210	PAKNING, QF KOBLING	GASKET, QUICKFLEX CONNECTION	DICHTUNG, QUICKFLEX KUPPLUNG	ROND. D'ETANC. RACC.QUICKFLE
13	SC-748	SPRAYRØR, FLEX LUX KOBLING	SPRAY TUBE, FLEX LUX COUPLING	SPRAYROHR, FLEX LUX KUPPLUNG	TUYAU SPRAY, RACCORD FLEX LU
15	SC-998	ANKER, FLEX MOTOR	ROTOR, FLEX MOTOR	SCHAUFELRAD, FLEX MOTOR	ROUE MOBILE, MOTEUR FLEX
17	SA-036	KUGLELEJE	BALL BEARING	KUGELLAGER	ROULEMENT A BILLES
18	SA-035	KUGLELEJE	BALL BEARING	KUGELLAGER	ROULEMENT A BILLES
20	HT-070	QF KOBLING FOR 4+4 SLANGE	QUICKFLEX COUPLING, 4+4 TUBE	QUICKFLEX KUPPLUNG, 4+4 SCHL.	RACCORD QUICKFLEX, CORDON 4+
22	HT-097	INT. MOTOR, 40.000 INTRA,LUX *	* MOTOR,INT.,40,000 INTRA, LUX *	* MOTOR,INT.,40.000 INTRA, LUX *	* MOTEUR,INT,40.000 INTRA,LUX
23	HT-092	INT. MOTOR, 40.000, ISO, LUX *	* MOTOR, INT., 40,000, ISO,LUX *	* MOTOR, INT., 40.000, ISO,LUX *	* MOTEUR,INT, 40.000, ISO,LUX
24	SC-821	SPRAYRØR, ISOLUXKOBLING	SPRAY TUBE, ISOLUXKOBLING	SPRAYROHR, ISOLUXKUBLUNG	TUYAU SPRAY, RACCORD ISOLUX
25	SC-750	INTRAMATIC KOBLING,FLEX LUX M.	INTRAMATIC COUPLING, FLEX LUX	INTRAMATIC KUPPLUNG, FLEX LUX	RACCORD INTRAMATIC, FLEX LUX
26	SC-817	FLEX ISOLUX KOBLING	FLEX ISOLUX COUPLING	FLEX ISOLUX KUPPLUNG	RACCORD FLEX ISOLUX
27	SD-312	RING,FOR MOTOR	RING, FOR MOTOR	RING, FÜR MOTOR	JOINT POUR MOTEUR
28	SD-313	O-RING FOR SD-312	O-RING FOR SD-312	O-RING FÜR SD-312	JOINT TORIQUE POUR SD-312



Instrumenter: Tandrenser.type L, LC lampe

Instruments: Scaler Type 1, Composite Lamp

YA-001

Instrumente: Zahnreiniger Typ L, Kompositlicht

Instruments : Détartreur type 1, Lampe composite

0.6.15.B

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZICHTUNG	DESCRIPTION
1	SD-227	SERVICESÆT INTEGRAL TANDRENSER	SERVICE KIT, SCALER	SERVICESATZ, ZAHNREINIGER	SERVICE KIT, DETARTREUR
2	SD-217	SPIDS, TANDRENSER	TIP, SCALER	SPITZE, ZAHNREINIGER	INSERT, DETARTREUR
3	SD-218	VÆRKTØJ, INCL.SPIDSMÅLER TANDR.	TOOL, CHANGE OF SCALER TIP	WERKZEUG, WECHSEL, Z.E.G-SPITZE	OUTIL, CHANGEMENT INSERT DETA
4	MC-987	RULLESIKRING, TANDRENSER	ANTI-ROLL DEVICE	ANTIROLLANORDNUNG, ZEG	DISPOSITIF ANTI-ROULEMENT
5	UC-321	BLYPLADE, TEST AF TANDRENSER	LEAD PLATE FOR TESTING SCALER	BLIPLATTE F. PRÜFUNG DES ZEG	PLAQUE EN PLOMB, TEST DETARTR
6	WH-004	PÆRE, KOMPOSITLAMPE	BULB, COMPOSITE LAMP	BIRNE, KOMPOSITLAMPE	AMPOULE, LAMPE À COMPOSITE
7	SD-220	FIBERLYSSTAV, STANDARD, KOMPOSIT	FIBRE ROD, STANDARD, COMPOSITE	STANDARD FIBERLICHTSTAB, KOMPO.	TIGE EN FIBRE, STANDARD, COMP
8	SD-221	FIBERLYSSTAV 8 MM, KOMPOSIT	FIBRE ROD 8 MM, COMPOSITE	FIBERLICHTSTAB, KOMPOSIT	TIGE EN FIBRE 8 MM, KOMPOSIT
9	SD-222	FIBERLYSSTAV 13 MM, KOMPOSIT	FIBRE ROD 13 MM, COMPOSITE	FIBERLICHTSTAB 13 MM, KOMPOSIT	TIGE EN FIBRE 13MM, KOMPOSIT
10	SD-223	BESKYTTELSSESBRILLER, (KOMPOSIT)	PROTECTIVE GLASSES (COMPOSITE)	SCHUTZBRILLEN (KOMPOSIT)	LUNETTES DE PROTECTION (COMP
11	UC-278	LYSFILTER, KOMPOSITLAMPE	LIGHT FILTER, COMPOSITE LAMP	LICHTFILTER, KOMPOSITLICHT	FILTRE A LUMIERE, LAMPE COMP
12	UC-429	LYSTRAGT, KOMPOSITLAMPE	LIGHT FUNNEL, COMPOSITE LAMP	LICHTTRICHTER, KOMPOSITLICHT	ENTONNOIR A LUMIERE, LAMPE C
13	UC-430	Q-FLEX HUNDEL, LC.M/ELEKTRONIK	QF SOCKET, COMPOSITE + ELECT.	QF BUCHSE, KOMPOSIT. + ELECTRO.	RACCORD QF, COMPOSITE, ELECTRO
14	SD-205	BESKYTTELSSESKJOLD, KOMPOSITL.	PROTECTION SHIELD, COMPOSITE	SCHUTZSCHILD, KOMPOSITLICHT	ECRAN DE PROTECTION, COMPOSIT
15	UC-732	SKRUE, FOR MC-987 RUSTFRI	SCREW, STAINLESS, FOR MC-987	SCHRAUBE, ROSTFREI, FÜR MC-987	VIS INOX POUR MC-987
16	SD-426	TEST ADAPTER, LC-LAMPE	ADAPTOR, MEASUREMENT, LC-VOLTAGE	ADAPTER, MESSUNG, LC-SPANNUNG	ADAPTAT., MESURE, TENSION, LAMP
17	SD-410	Q-FLEX HUNDEL, SCALER.M DIODE	QF SOCKET, SCALER + DIODE	QF BUCHSE, ZEG + DIODE	RACCORD QF, DETARTREUR + DIOD
18	UC-512	NÆSEMØTRIK, KOMPOSITLAMPE INT.	UNION NUT, COMPOSITE LAMP, INT.	ÜBERWURFMUTTER, KOMPOSITLICHT	ECROU DE SERRAGE, L. COMPOSIT
19	SD-411	TANDRENSER MED QF	* SCALER WITH QUICKFLEX	* ZEG MIT QUICKFLEX	* DETARTREUR AVEC QUICKFLEX
20	SD-415	KOMPOSIT, U/LYSSTAV, INTOR.	* COMPOSIT CURING LIGHT	* KOMPOSITLAMPE	* LAMPE A PHOTOPOLYMERISER
21	UC-665	KOMPOSITLAMPE TESTER	TESTER, COMPOSITE CURING LIGHT	TESTER, KOMPOSITLAMPE	TESTEUR, DURCISSEMENT


NEW LAMP (BLACK & WHITE)

NEUE LAMPE (SCHWARZ & WEISS)

6 WH-006 Pære, Nye type sort/hvid
 7 FA-064 Fiberlysstav, nye type sort/hvid
 20 HT-125 Kompositlampe QF, ny type

Bulb, new type black/white
 Fibre rod, new type black/white
 Composite lamp QF, new type

Birne, neue Typ schwarz/weiss
 Fiberlichtstab, neue Typ, schwarz/weiss
 Kompositlampe QF, neue Typ



No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
	AE-915	KIT, GEVINDTAP OG M-CLIPS M.M.	KIT, BALANCED SUSPENSION	KIT, BALANCEAUFHÄNGUNGEN	KIT DE FOUET
1	SD-213	KAPPE, 3 FUNKTIONSPRØJTE	* JACKET, 3-WAY SYRINGE	* SCHUTZKAPPE, 3-FUNKT.-SPRITZE	* MANCHON, SERINGUE 3 FONCTIONS
2	SD-226	KAPPE, 6 FUNKTIONSPRØJTE	* JACKET, 6-WAY SYRINGE	* SCHUTZKAPPE, 6-FUNKT.SPRITZE	* MANCHON, SERINGUE 6 FONCTIONS
3	SD-214	SPIDS, INTEGRAL SPRØJTE	TIP, SYRINGE	SPITZE, SPRITZE	EMBOUT, SERINGUE
4	SD-215	VENTILBLOK, SPRØJTE VARM	VALVE BLOCK, SYRINGE, WARM	VENTILBLOCK, SPRITZE, WARM	BLOC VANNE, SERINGUE, CHAUDE
4	SD-489	VENTILBLOK, SPRØJTE KOLD	VALVE BLOCK, SYRINGE COLD	VENTILBLOCK, SPRITZE, KALT	BLOC VANNE, SERINGUE, FROIDE
5	SD-322	SKRUE M 2, SPRØJTE	SCREW M 2, SYRINGE	SCHRAUBE M 2, SPRITZE	VIS M 2, SERINGUE
6	SD-325	FORBINDELSBESTYKKE, SPRØJTE	CONNECTION PIECE, SYRINGE	VERBINDUNGSTÜCK, SPRITZE	RACCORD, SERINGUE
7	SD-326	KONTAKTFJEDER, SPRØJTE	CONTACT SPRING, SYRINGE	KONTAKTFEDER, SPRITZE	RESSORT DE CONTACT, SERINGUE
8	SD-327	O-RING, SPRØJTE	O-RING, SYRINGE	O-RING, SPRITZE	JOINT TORIQUE, SERINGUE
9	SD-328	STUDS, LUFT, SPRØJTE	CONNECTOR, AIR, SYRINGE	STUTZEN, LUFT, SPRITZE	RACCORD, AIR, SERINGUE
10	SD-329	STUDS, VAND, SPRØJTE	CONNECTOR, WATER, SYRINGE	STUTZEN, WASSER, SPRITZE	RACORD, EAU, SERINGUE
11	SD-330	FJEDER, SD-332	SPRING, SD-332	FEDER, SD-332	RESSORT, SD-332
12	SD-331	FJEDER, VENTIL	SPRING, VALVE	FEDER, VENTIL	RESSORT, VANNE
13	SD-332	KONTAKTSTYKKE, SPRØJTE	CONTACT PIECE, SYRINGE	KONTACKSTÜCK, SPRITZE	PIECE DE CONTACT, SERINGUE
14	SD-333	VENTILKEGLE, SPRØJTE	VALVE CONE, SYRINGE	VENTILKONUS, SPRITZE	CONE VANNE, SERINGUE
15	SD-334	TRYKSTYKKE, SPRØJTE	PRESSURE PIECE, SYRINGE	DRUCKSTÜCK, SYRINGE	PIECE DE PRESSION, SYRINGE
16	SD-335	OMLØBER, SPIDS, SPRØJTE	UNION, SYRINGE	ÜBERWURFMUTTER, SPRITZE	ECROU, SERINGUE
17	SD-336	NYLON KAPPE, INTEGRAL, SPRØJTE	NYLON JACKET, SYRINGE	NYLONDECKKAPPE, SPRITZE	MANCHON EN NYLON, SERINGUE
18	SD-337	OMSKIFTER, KOLD-VARM SPRØJTE	SWITCH, COLD/WARM SYRINGE	SCHALTER, KALTE/WARME SPRITZE	INVERSEUR, SERINGUE FR./CHAU
19	SD-228	BALANCEOPHÆNG, 3 FKT.SPRØJTE	* BALANCED SUSPENSION, 3-WAY SYR*	BALANCEAUFHÄNG., 3-FUNKT SPR.*	FOUET, SERINGUE, 3 FONCTIONS
19	SD-229	BALANCEOPHÆNG, 6 FKT.SPRØJTE	* BALANCED SUSPENSION, 6-WAY SYR*	BALANCEAUFHÄNG., 6-FUNKT.SPR*	FOUET, SERINGUE, 6 FONCTIONS
21	AE-245	RULLE M.MAGN.AC-035/SD-229/28	ROCKER PIECE, AC-035/SD-29/28	KIPPTEIL, AC-03/SD-229/28	PIECE BASCUL.AC-035/SD-229/2
21	AE-264	RULLE M.MAGNET MOTOR (AE-271)	ROCKER PIECE, FOR AE-271	KIPPTEIL FÜR AE-271	PIECE A BASCULER, AE-271
22	MC-457	FJEDER DT=0,8 D=6 L0=30,6	SPRING	FEDER	RESSORT
23	UC-726	O-RING Ø3X1, SILICONE, 60 SH.	O-RING Ø 3 X 1, SILICONE	O-RING Ø 3 X 1, SILIKON	JOINT TORIQUE Ø3X1, SILICONE
24	MC-458	GEVINDTAP, INSTR.RULLE FJEDER	THREADED PIN,ROCKER PIECE	GEWINDEZAPFEN, KIPPTEILFEDER	TENON TARAUDE, PIECE BASCULE

Rubber/Gummi ME-052



No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	MC-111	HÅNDTAG, BRO/LAMPE, 1 STK.	HANDLE, INSTRUMENT BRIDGE	HANDGRIFF, INSTRUMENTENTRÄGER	POIGNEE, PORTE-INSTRUMENT
1	SD-399	HÅNDTAG, BRO/LAMPE, 4 STK.	HANDLES (4) BRIDGE/LAMP	HANDGRIFFE (4) TRÄGER/LAMPE	POIGNEES(4) PORTE-INSTR./ECL
1	AC-371	HÅNDTAG, BRO/LAMPE, 24 STK.	HANDLES (24), BRIDGE/LAMP	HANDGRIFFE (24), TRÄGER/LAMPE	POIGNEES(24)PORTE-INSTR./ECL
2	MC-186	MODTAGER, HÅNDTAG, INSTRUMENTBRO	HOLDER, HANDLE, INSTR. BRIDGE	EINSTECKROHR, GRIFF, INST. TRÄGER	TUYAU SUPP., POIGNEE, PORTE-IN
3	UC-331	PINOLSKRUE M8 X 20 RUSTFRI	POINTED SCREW M8X20, STAINLESS	REITSTOCKSCHRAUBESX20, ROSTFREI	VIS POINTUE M8X20, INOXYDABL
4	HA-304	GASFJEDER, ARM, LAMPE, M.1 ØJE	GAS SPRING, LAMP ARM	GASFEDER, LAMPENARM	RESSORT A GAZ, BRAS DE ECLAI
5	AC-329	RØNGTENRØR, INTEG. 9018 KORT	X-RAY TUPE, 9018 SHORT	RÖNTGENADAPTER, 9018 KURZ	TYAU RADIO, 9018 COUR
5	AB-583	RØNGTENRØR, INTEG. CLASSIC KORT	X-RAY TUPE, CLASSIC SHORT	RÖNTGENADAPTER, CLASSIC KURZ	TYAU RADIO, CLASSIC COUR
5	AB-981	RØNGTENRØR, INTEGRAL LANG	X-RAY TUBE, 9018	RÖNTGEN ROHR, 9018	TUYAU RADIO, 9018
5	AB-982	RØNGTENRØR, INTEG CLASSIC LANG	X-RAY TUBE, CLASSIC	RÖNTGENROHR, CLASSIC	TUYAU RADIO, CLASSIC
6	AK-049	BØSNING TYPE F (Ø 32,1)	BUSH TYPE F (Ø 32.1 MM)	BUCHSE TYP F (Ø 32,1 MM)	DOUILLE TYPE F (Ø 32,1 MM)
6	AK-050	BØSNING TYPE B (Ø 35,1)	BUSH TYPE B (Ø 35.1 MM)	BUCHSE TYP B (Ø 35,1 MM)	DOUILLE TYPE B (Ø 35,1 MM)
6	AK-052	BØSNING TYPE A (Ø 34,95)	BUSH TYPE A (Ø 34.95 MM)	BUCHSE TYP A (Ø 34,95 MM)	DOUILLE TYPE A (Ø 34,95 MM)
6	AK-053	BØSNING TYPE C (Ø 38,1)	BUSH TYPE C (Ø 38.1 MM)	BUCHSE TYP C (Ø 38,1 MM)	DOUILLE TYPE C (Ø 38,1 MM)
6	AK-054	BØSNING TYPE D (Ø 41,1)	BUSH TYPE D (Ø 41.1 MM)	BUCHSE TYP D (Ø 41,1 MM)	DOUILLE TYPE D (Ø 41,1 MM)
6	AK-055	BØSNING TYPE E (Ø 34,85)	BUSH TYPE E (Ø 34.85 MM)	BUCHSE TYP E (Ø 34,85 MM)	DOUILLE TYPE E (Ø 34,85 MM)
6	AK-056	BØSNING TYPE G (Ø 43,1)	BUSH TYPE G (Ø 43.1 MM)	BUCHSE TYP G (Ø 43,1 MM)	DOUILLE TYPE G (Ø 43,1 MM)
6	AK-057	BØSNING TYPE H (Ø 41,4)	BUSH TYPE H (Ø 41.4 MM)	BUCHSE TYP H (Ø 41,4 MM)	DOUILLE TYPE H (Ø 41,4 MM)
6	AK-058	BØSNING TYPE K (Ø 39,75)	BUSH TYPE K (Ø 39,75)	BUCHSE TYP K (Ø 39,75)	DOUILLE TYPE K (Ø 39,75)
6	AK-059	BØSNING TYPE I (Ø 30,1)	BUSH TYPE I (Ø 30.1 MM)	BUCHSE TYP I (Ø 30,1 MM)	DOUILLE TYPE I (Ø 30,1 MM)
6	AK-060	BØSNING TYPE L (Ø 22,1)	BUSH TYPE L (Ø 22.1 MM)	BUCHSE TYP L (Ø 22,1 MM)	DOUILLE TYPE L (Ø 22,1 MM)
6	AK-061	BØSNING TYPE M (Ø 20,1)	BUSH TYPE M (Ø 20.1 MM)	BUCHSE TYP M (Ø 20,1 MM)	DOUILLE TYPE M (Ø 20,1 MM)
6	AK-063	BØSNING TYPE X (Ø 24,1)	BUSH TYPE X (Ø 24.1 MM)	BUCHSE TYP X (Ø 24,1 MM)	DOUILLE TYPE X (Ø 24,1 MM)
6	AK-064	BØSNING TYPE Z (Ø 29,1)	BUSH TYPE Z (Ø 29.1 MM)	BUCHSE TYP Z (Ø 29,1 MM)	DOUILLE TYPE Z (Ø 29,1 MM)
6	AK-065	BØSNING TYPE V (Ø 25,1)	BUSH TYPE V (Ø 25.1 MM)	BUCHSE TYP V (Ø 25,1 MM)	DOUILLE TYPE V (Ø 25,1 MM)
7	AB-947	DÆKSEL, LAMPE, UNDERST CLASSIC	COVER, LAMP, LOWER, CLASSIC	ABDECKUNG, LAMPE UNTEN, CLASSIC	CARTER, BRAS, DESSOUS, CLASSI
7	AB-948	DÆKSEL, LAMPE ØVERST CLASSIC	COVER, LAMP, UPPER, CLASSIC	ABDECKUNG LAMPEN OBEN, CLASSIC	CARTER, AVANT, DESSOUS, CLASSI
7	AC-252	DÆKSEL, LAMPEFORARM, UNDERST	COVER, LAMP FRONT ARM, LOWER	ABDECKUNG, VORDERARM, UNTEN	CARTER, BRAS AVANT, DESSOUS
7	AC-985	DÆKSEL, LAMPEFORARM, ØVERST	COVER, LAMP FRONT ARM, UPPER	ABDECKUNG, VORDERARM, OBEN	CARTER, BRAS AVANT, DESSUS
8	SD-397	PLASTHÅNDGREB, INTEG.U/GEVIND	HANDLE, 9018, WITHOUT THREAD	HANDGRIFF, 9018, OHNE GEWINDE	POIGNEE, 9018, SANS FILETAGE
8	SD-125	PLASTHÅNDGREB, INTEG. M/GEVIND	HANDLE, FLEX LAMP, 9018	HANDGRIFF, FLEX LAMPE, 9018	POIGNEE, ECLAIRAGE FLEX, 901
9	SD-127	PLASTKAPPE, INTEGRAL 9018	CAP, FLEX LAMP, 9018	KAPSEL, FLEX LAMPE, 9018	CARTER, ECLAIRAGE FLEX, 9018
10	SD-130	KONTAKT, INTEGRAL, 9018 TY.S-308	SWITCH, INTEGRAL, 9018 TY S-308	SCHALTER, INTEGRAL, 9018 TY S308	INTERRUPTEUR, ECLAIRAGE, 901
11	SD-126	BAGDÆKSEL, INTEGRAL, 9018	REAR COVER, FLEX LAMP, 9018	HINTERE ABDECKUNG, LAMPE 9018	CARTER ARRIERE, ECLAIRAGE, 901
12	WH-001	PÆRE, MULTILUX, INTEGRAL	BULB, MULTILUX	BIRNE, MULTILUX	AMPOULE, MULTILUX
13	MC-205	BØSNING, LAMPEADAPTER	BUSH, LAMP ADAPTOR	BUCHSE, LAMPENADAPTER	DOUILLE, COLONNE DE LAMPE
14	AB-454	ADAPTOR, FREMMED LAMPE, CLASSIC	ADAPTOR, OTHER MAKE, CLASSIC	ADAPTER, FREMDHERST., CLASSIC	ADAPTATEUR, AUTRE FABR. CLASSI
14	AC-388	ADAPTOR, FREMMED, LAMPE 9018	ADAPTOR, OTHER MAKE, 9018	ADAPTER, FREMDHERSTELLER, 9018	ADAPTATEUR, AUTRE FABRIC. 901
15	AB-453	ADAPTOR, LAMPE, INT. CLASSIC	ADAPTOR, LAMP. INT. CLASSIC	ADAPTER, LAMPE, INT. CLASSIC	ADAPTATEUR, INTEGRAL CLASSIC
15	AC-257	ADAPTOR, LAMPE, INTEGRAL 9018	ADAPTOR, INTEGRAL, 9018	ADAPTER, INTEGRAL, 9018	ADAPTATEUR, INTEGRAL, 9018
16	AC-117	TRAFO, OPERATIONS LAMPE INTEG.	TRANSFORMER, OPERATING LAMP	TRANSFORMATOR, OPERATIONS LAMPE	TRANSFORMATEUR, ECLAIRAGE BUC
16	AB-073	TRAFO, 100V OPERATIONS LAMPE	TRANSFORMER, OP. LAMP, 100 V	NUR FUER JAPAN	ONLY FOR JAPAN
17	AC-541	BREMSE, LAMPE	BRAKE, LAMP	BREMSE, LAMPE	FREIN, ECLAIRAGE BUCCAL
18	MC-562	JUST. SKRUE, BREMSE, LAMPE	ADJUSTING SCREW, BRAKE, LAMP	STELLSCHRAUBE, BREMSE, LAMPE	VIS DE REGLAGE, FREIN, ECLAIR.
19	MC-748	AFBLANDING, GULDFANG/LAMPE INT.	STOPPER, GOLD TRAP, LAMP/INT.	STÖPSEL, GOLDFÄNGER, LAMPE/INT	OPERCULE, TRAPPE OR/LAMPE INT
20	WR-001	FRONTDÆKSEL, MULTILUX/INTEGRAL	FRONT SHIELD, MULTILUX	FRONTDECKPLATTE, MULTILUX	CACHE AVANT, MULTILUX
21	WB-009	DÆKSEL, PÆRE, MULTILUX/INTEGRAL	COVER, BULB	DECKEL, LAMPE	CACHE, AMPOULE
22	MR-310	O-RING 45 X 2 RTG./LAMPE BØSN.	O-RING 45 X 2	O-RING 45 X 2	JOINT TORIQUE 45 X 2
23	UC-579	SKRUE, MULTILUX / INTEGRAL	FLAT HEAD SCREW 4MM	SCHRAUBE	VIS
24	AC-162	KABEL FOR OPERATIONS LAMPE	CABLE, OPERATING LAMP	KABEL, OPERATIONS LAMPE	CABLE, ECLAIRAGE BUCCAL
25	MC-070	PYNTESKIVE, LAMPELED	DECORATION DISC, ARTICULATION	ZIERSCHIBE, GELENK	RONDELLE ENJOLIVEUR, ARTICUL
26	WB-100	KABELSÆT LAMPEHOVED	CABLE SET, LAMP HEAD	KABELSATZ, LAMPENKOPF	CABLES JEU, TÊTE DE LAMPE
27	MC-249	LÅSEPÅL FOR LAMPERØR	LOCKING PAWL, LAMP TUBE	SPERRHAKEN, LAMPENROHR	CLIQUET DE FERMETURE, SUPPOR
28	MC-758	SLIDSKIVE FOR LAMPERØR	WEARING WASHER, LAMP TUBE	SCHUTZSCHIBE, LAMPENROHR	RONDELLE D'USURE, SUPPORT
29	SC-917	BØLGEFJEDER MULTILUX/INTEGRAL	CORRUGATED SPRING, MULTILUX	WELLFEDER, MULTILUX	RESSORT ONDULE MULTILUX

Flex Dental A/S
Lampe: Unitmonteret
Lampe: Unitmontage

Lamp: On the Unit
Eclairage: sur l'unit

1999-06-28
YA-001
0.6.17.C

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
30	SC-946	SKIVE FOR MULTILUX/INTEGRAL	CORRUGATED WASHER, MULTILUX	WELLSCHLEIBE, MULTILUX	RONDELLE ONDULEE MULTILUX
31	AB-585	DEKSEL, VINKEL, LAMPE, CLASSIC	COVER, ANGLE, LAMP, CLASSIC	ABDECK. WINKEL, LAMPE, CLASSIC	CARTER, COUDE, CLASSIC
31	AC-626	DEKSEL, VINKEL, LAMPE, RAL 9018	COVER, ANGLE, LAMP, RAL 9018	ABDECKUNG, WINKEL, LAMPE, RAL 9018	CARTER, COUDE, RAL 9018
32	SD-496	MELLEMLÆD, LAMPE INT. RAL 9018	INTERMEDIATE LAMP, RAL 9018	ZWISCHENSTÜCK, LAMP, RAL 9018	PIÈCE INTERME. ECLAIRAGE 901
32	SD-533	MELLEMLÆD LAMPE INT. CLASSIC	INTERMEDIATE LAMP, CLASSIC	ZWISCHENSTÜCK LAMPE, CLASSIC	PIÈCE INTERME, ECLAIRAGE, CLAS
33	SD-508	FORLED, LAMPE, INT. CLASSIC	FRONT JOINT, LAMP, CLASSIC	VORDERGLIED, LAMPE, CLASSIC	CHARNIERE AVANT ECLAIR, CLASS
33	SD-422	FORLED, LAMPE, INT. RAL 9018	FRONT JOINT, LAMP, INT. RAL 9018	VORDERGLIED, LAMPE, INT. RAL 9018	CHARNIERE AVANTECLAIR, RAL 90
34	AC-664	SKIVE, FORLED/LAMPE, RAL 9018	WASHER, FR. JOINT/LAMP, RAL 9018	SCHLEIBE, VORDERGL/LAMPE, RAL 9018	RONDELLE, CHARN/ECLAIR, RAL 90
34	AB-694	SKIVE, FORLED/LAMPE, CLASSIC	WASHER, FR. JOINT/LAMP, CLASSIC	SCHLEIBE, VORDERGL/LAMPE, CLASSIC	RONDELLE, CHARN/ECLAIR, CLAISS
35	AC-741	LAMPERØR, INTEGRAL, RAL 9018	LAMP TUBE, INTEGRAL, RAL 9018	LAMPENROHR, INTEGRAL, RAL 9018	TUYAU ECLAIRAGE, RAL 9018
35	AB-586	LAMPERØR, INTEGRAL, CLASSIC	LAMP TUBE, CLASSIC	LAMPENROHR, CLASSIC	TUYAU ECLAIRAGE, CLASSIC
36	AC-627	VINKEL, LAMPELED, LAKERET, 9018	ANGLE, LAMP, RAL 9018	WINKEL, LAMPE, RAL 9018	COUDE, ECLAIRAGE, RAL 9018
36	AB-584	VINKEL, LAMPELED, CLASSIC	ANGLE, LAMP, CLASSIC	WINKEL, LAMPE, CLASSIC	COUDE, ECLAIRAGE, CLASSIC
37	WB-014	GAFFEL INTEGRAL 9018	HOLDER (FORK) INTEGRAL LAMP 9018	GABEL INTEGRALLAMPE 9018	SUPPORT (FOURCHE) ECLAIRAGE 90
38	SC-895	GUMMI PAKNING, MULTILUX/INTEG.	RUBBER GASKET, MULTILUX	GUMMIPACKUNG, MULTILUX	JOINT EN CAOUTCHOUC, MULTILU
39	WR-005	REFLEKTOR, MULTILUX, INTEGRAL	REFLECTOR, MULTILUX	REFLEKTOR, MULTILUX	RÉFLECTEUR, MULTILUX
40	AC-527	MODTAGER, LÅSETAP M/FJEDER	RECEIVER, DOWEL + SPRING	EMPFÄNGER, SPERRBOLZEN + FEDER	RECEPT., GOUJON BLOCAGE+RESSO
41	SC-890	SAMLEBOLT FOR MULTILUX/INTEG.	JOINT BOLT	JOINT BOLT	JOINT BOLT
42	SC-919	MØTRIK FOR MULTILUX/INTEGRAL	NUT, MULTILUX	MUTTER, MULTILUX	ECROU, MULTILUX
43	SC-920	NYLON SKIVE, MULTILUX/INTEGRAL	WASHER, NYLON, FOR MULTILUX	UNTERLAGSSCHLEIBE, NYLON, MULTIL.	RONDELLE NYLON, MULTILUX
44	WH-003	FATNING, MULTILUX / INTEGRAL	SOCKET, MULTILUX	FASSUNG, MULTILUX	DOUILLE, MULTILUX
45	WB-004	FJEDER, MULTILUX / INTEGRAL	SPRING CLIP	FEDER	RESSORT
46	WB-016	SIDESKÆRM, INTEGRAL/MULTILUX	SIDE COVER, INTEGRAL/MULTILUX	ABDECKUNG, INTEGRAL/MULTILUX	CACHE, INTEGRAL/MULTILUX
47	UC-365	PINOLSKRUE M 6 X 8 BSLOK	POINTED SCREW M 6 X 8 BSLOK	PINOLSCHRAUBE M 6 X 8 BSLOK	VIS POINTUE M 6 X 8 BSLOK



No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	MC-111	HÅNDTAG, BRO/LAMPE, 1 STK.	HANDLE, INSTRUMENT BRIDGE	HANDGRIFF, INSTRUMENTENTRÄGER	POIGNEE, PORTE-INSTRUMENT
1	SD-399	HÅNDTAG, BRO/LAMPE, 4 STK.	HANDLES (4) BRIDGE/LAMP	HANDGRIFFE (4) TRÄGER/LAMPE	POIGNEES (4) PORTE-INSTR./ECL
1	AC-371	HÅNDTAG, BRO/LAMPE, 24 STK.	HANDLES (24), BRIDGE/LAMP	HANDGRIFFE (24), TRÄGER/LAMPE	POIGNEES (24) PORTE-INSTR./ECL
2	UC-331	PINOLSKRUE M8 X 20 RUSTFRI	POINTED SCREW M8X20, STAINLESS	REITSTOCKSCHRAUBESX20, ROSTFREI	VIS POINTUE M8X20, INOXYDABL
3	AC-170	LOFTOPHÆNG, U/SKÆRM	CEILING SUSPENS. WITHOUT COVER	DECKENAUFHÄNG. OHNE ABDECKUNG	FIXATION PLAFONIERE SANS CAC
4	AE-600	LOFTSKÆRM, CLASSIC	CEILING COVER, CLASSIC	ABDECKUNG, CLASSIC	CACHE, CLASSIC
4	AC-404	LOFTSKÆRM, RAL 9018	CEILING COVER, 9018	ABDECKUNG, 9018	CACHE, 9018
5	MR-115	O-RING 50, 16 X 5,33	O-RING 50, 16 X 5,33	O-RING 50, 16 X 5,33	JOINT TORIQUE 50, 16 X 5,33
6	MC-205	BØSNING, LAMPEADAPTER	BUSH, LAMP ADAPTOR	BUCHSE, LAMPENADAPTER	DOUILLE, COLONNE DE LAMPE
7	MH-432	KABEL 2 X 2,5 MM2 FLEKSIB. GRÅ	CABLE 2 X 2,5 MM2, GREY	KABEL 2 X 2,5 MM2, GRAU	CABLE 2 X 2,5 MM2, GRIS
7	MH-433	KABEL 3 X 0,25 MM2 SKÆRMET GRÅ	CABLE 3X0.25MM2, SHIELDED, GREY	KABEL 3X0,25MM2, GESCHÜTZT, GRAU	CABLE BLINDE 3X0,25 MM2, GRI
8	UC-101	PASSKIVE PS 20 X 28 X 1,2	ADAPTER DISC	ADAPTERSCHEIBE	RONDELLE ADAPTATEUR
9	UC-133	SEEGERRING DIN 471-20 X 1,2	SEEGER-A-RING DIN 471-20 X 1,2	SEEGER-A-RING DIN 471-20 X 1,2	BAGUE SEEGER-A DIN 471-20X2,
10	AE-582	LAMPERØR, LOFT 2000 MM, CLASSIC	LAMP CEILING, 2000MM, CLASSIC	LAMPENROHR, DECKE, 2000MM, CLAS.	LAMPE, PLAFOND, 2000MM, CLAS
10	AE-581	LAMPERØR, LOFT 925 MM, CLASSIC	LAMP CEILING, 925MM, CLASSIC	LAMPENROHR, DECKE, 925MM, CLASSI	LAMPE, PLAFOND, 925MM, CLASS
10	AC-416	LAMPERØR, LOFT 925 MM 9018	LAMP CEILING, 925MM,	LAMPENROHR, DECKE, 925MM,	LAMPE, PLAFOND, 925MM,
10	AC-421	LAMPERØR, LOFT 2000 MM 9018	LAMP, CEILING, 2,000 MM	LAMPENROHR, DECKE, 2.000 MM	LAMPE, PLAFOND, 2.000 MM
11	AC-117	TRAFO, OPERATIONS-LAMPE INTEG.	TRANSFORMER, OPERATING LAMP	TRANSFORMATOR, OPERATIONS-LAMPE	TRANSFORMATEUR, ECLAIRAGE BUC
11	AE-073	TRAFO, 100V OPERATIONS-LAMPE	TRANSFORMER, OP. LAMP, 100 V	NUR FUER JAPAN	ONLY FOR JAPAN
12	MC-748	AFBLÆNDING, GULDFANG/LAMPE INT.	STOPPER, GOLD TRAP, LAMP/INT.	STÖPSEL, GOLDFÄNGER, LAMPE/INT	OPERCULE, TRAPPE OR/LAMPE INT
13	UA-210	MONTAGEBOLT M 10 X 95 MM	ATTACHMENT BOLT M 10 X 95 MM	MONTAGEBOLZEN M 10 X 95 MM	BOULON DE FIXATION M 10X95MM
14	UB-070	SKIVE M 10,5	WASHER M 10,5	SCHIBE M 10,5	RONDELLE M 10,5
15	UC-132	MELLEMBØSNING 20 X 10,5 X 22	DISTANCE PIECE 20 X 10,5 X 22	ABSTANDSTÜCK 20 X 10,5 X 22	PIECE D'ECARTEMENT 20X10,5X2
16	AC-541	BREMSE, LAMPE	BRAKE, LAMP	BREMSE, LAMPE	FREIN, ECLAIRAGE BUCCAL
17	MC-562	JUST. SKRUE, BREMSE, LAMPE	ADJUSTING SCREW, BRAKE, LAMP	STELLSCHRAUBE, BREMSE, LAMPE	VIS DE REGLAGE, FREIN, ECLAIR.
18	UA-025	MØTRIK M 10	NUT M 10	SCHRAUBENMUTTER M 10	ECROU M 10
19	MC-186	MODTAGER, HÅNDTAG, INSTRUMENTBRO	HOLDER, HANDLE, INSTR. BRIDGE	EINSTECKROHR, GRIFF, INST. TRÄGER	TUYAU SUPP., POIGNEE, PORTE-IN
20	MH-341	STIK 6P MOLEX 5557-NR	PLUG 6P MOLEX 5557-NR	STECKER 6P MOLEX 5557-NR	PRISE 6P MOLEX 5557-NR
21	MH-310	TERMINAL HUN AMP 163557-2	TERMINAL, FEMALE, AMP 163557-2	TERMINAL, WEIB., AMP 163557-2	TERMINAL, FEM., AMP 163557-2
22	UC-190	PINOLSKRUE M 5 X 8 RUSTFRI	POINTED SCREW M 5X8, RUSTPROOF	REITSTOCKSCHRAUBE 5X8, ROSTFREI	VIS POINTUE M 5X8 INOXYDABLE

This exploded view diagram illustrates the components of a mechanical assembly, numbered 1 through 40. The parts are arranged to show their relative positions and assembly sequence. Key components include:


- 1:** A small rectangular component, possibly a sensor or actuator.
- 2:** A printed circuit board (PCB) with various electronic components.
- 3:** A cylindrical motor or actuator.
- 4:** A larger motor or actuator with a threaded rod attached.
- 5:** A black cylindrical component, possibly a connector or plug.
- 6:** A small black component, possibly a switch or button.
- 7:** A small black component, possibly a switch or button.
- 8:** A small black component, possibly a switch or button.
- 9:** A small black component, possibly a switch or button.
- 10:** A black cylindrical component, possibly a connector or plug.
- 11:** A white cylindrical component, possibly a connector or plug.
- 12:** A small black component, possibly a switch or button.
- 13:** A black cable or wire.
- 14:** A small black component, possibly a switch or button.
- 15:** A black cylindrical component, possibly a connector or plug.
- 16:** A long white vertical component, possibly a frame or support.
- 17:** A small black component, possibly a switch or button.
- 18:** A small black component, possibly a switch or button.
- 19:** A white vertical component, possibly a frame or support.
- 20:** A small black component, possibly a switch or button.
- 21:** A small black component, possibly a switch or button.
- 22:** A long black vertical component, possibly a frame or support.
- 23:** A small black component, possibly a switch or button.
- 24:** A small black component, possibly a switch or button.
- 25:** A small black component, possibly a switch or button.
- 26:** A small black component, possibly a switch or button.
- 27:** A long black horizontal component, possibly a frame or support.
- 28:** A white cylindrical component, possibly a connector or plug.
- 29:** A small black component, possibly a switch or button.
- 30:** A small black component, possibly a switch or button.
- 31:** A small black component, possibly a switch or button.
- 32:** A black cylindrical component, possibly a connector or plug.
- 33:** A small black component, possibly a switch or button.
- 34:** A black cylindrical component, possibly a connector or plug.
- 35:** A black cylindrical component, possibly a connector or plug.
- 36:** A small black component, possibly a switch or button.
- 37:** A small black component, possibly a switch or button.
- 38:** A small black component, possibly a switch or button.
- 39:** A black cable or wire.
- 40:** A small black component, possibly a switch or button.

Patientstol, type 89: Mekaniske dele. Elektriske dele
 Patientenstuhl, type 89: Mechanik. Elektronik

Patient Chair, type 89: Mechanical Parts. Electrical Parts YA-001
 Fauteuil, type 89 : Pièces mécaniques. Pièces électroniques 0.6.19.B

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
	UC-558	SÆDE, CHASSIS, INTEGRAL STOL	PAINTED SEAT FRAME INTEGRAL I	SITZ, CHASSIS, INTEGRAL-STUHL	ASSISE, CHASSIS, FAUTEUIL II
	UC-255	DÆKSEL, LØFTEARM, STOL I	COVER, LIFTING ARM, CHAIR	ABDECKUNG, HUBARM, STUHL	CARTER, BRAS DE LEVAGE, FAUT.
	UC-564	DÆKPLADE, UC-255	COVER PLATE, UC-255	DECKPLATTE, UC-255	PLAQUE, UC-255
	SD-374	BASIS, SÆDERAMME, I + ALPHAFLEX	BASE, SEAT FRAME, ALPHAFLEX	BASIS, SITZRAHMEN, ALPHAFLEX	BASE, CHASSIS, ALPHAFLEX
1	AC-366	E-PROM PROGRAMMERET AD-417	EPROM PROGRAMMED AD-417	EPROM PROGRAMMIERT AD-417	EPROM PROGRAMME AD-417
2	AD-417	PRINT, STOL I	* PRINTED CIRC."CHAIR I"	* PLATINE "CHAIR I"	* CARTE "CHAIR I"
3	SC-842	LØFTEMOTOR, FLEX CHAIR	* LIFTING MOTOR, FLEX CHAIR	* HUBMOTOR, FLEX STUHL	* MOTEUR DE LEVAGE, FAUT.FLEX
4	SC-843	RYGMOTOR, FLEX CHAIR	* BACKREST MOTOR, FLEX CHAIR	* RÜCKLEHNENMOTOR, FLEX STUHL	* MOTEUR DU DOSSIER, FAUT.FLEX
5	SC-944	PROGRAMVÆLGER, STOL I + ALPHA.	PROGRAM SELECTOR, ALPHAFLEX LUX	PROGRAMMWÄHLER, ALPHAFLEX LUX	SÉLECT. PROGRAMME, ALPHAFLEX
6	MC-979	NAKKESTØTTE, INTEGRAL, HØJRE	* HEAD REST, INTEGRAL, RIGHT-HAND	* KOPFSTÜTZE, INTEGRAL, RECHTSH.	* TÊTIÈRE, INTEGRAL, DROITIER
6	MC-791	NAKKESTØTTE, INTEGRAL VENSTRE	* HEADREST INTEGRAL, LEFT-HAND	* KOPFSTÜTZE INTEGRAL, LINKSHÄN.	* TÊTIÈRE, INTEGRAL, GAUCHER
7	SC-841	SIKRINGSHOLDER, STOL I + ALPHA	FUSE HOLDER, 220 V ALPHAFLEX	SICHERUNGSHALTER, 220 V ALPHA	SUPPORT FUSIBLE 220 V ALPHA
9	SC-849	GUMMI, FODK. STOL I + ALPHA.	MICROSWITCH, FLEX CHAIR	MIKROSCHALTER, FLEX STUHL	MICRO-INTERRUPT. FAUTEUIL FLE
10	SC-845	GUMMI, FODK. STOL I + ALPHA.	RUBBER, ALPHAFLEX PEDAL	GUMMI, ALPHAFLEX PEDAL	CAOUTCHOUC, LEVIER ALPHAFLEX
11	UC-271	SPINDELMØTRIK, RYGMOTOR INT. I	SPINDLE NUT, PATIENT CHAIR	SPINDELSCHRAUBENMUTTER, STUHL	ECROU DE TRANSLATION, FAUTEUI
12	SD-380	UDLØSERKNAP, NAKKESTØTTE	RELEASE NUT, HEAD REST	AUSLÖSEMUTTER, KOPFSTÜTZE	ECROU DE DECLENCHEMENT, TETIE
13	SD-423	PROGRAMKNAP FOR STOL I, II	PROGRAM BOTTON, CHAIR II	PROGRAMMIERKNOPF, STUHL II	BOUTON PROGRAMMES, FAUTEUIL I
14	MH-242	HOVEDAFBRYDER, STOL/UNIT 220V	MAIN SWITCH, UNIT/CHAIR 220V	HAUPTSCHALTER, UNIT/STUHL 220V	INTERRUPTEUR, UNIT/FAUTEUIL 22
14	MH-745	HOVEDAFBRYDER, STOL/UNIT 100V	MAIN SWITCH, UNIT/CHAIR 100V	ONLY FOR JAPAN	ONLY FOR JAPAN
15	SD-388	TAPNØGLE, JUST, NAKKESTØTTE	KEY, ADJUSTMENT HEADREST	SCHLÜSSEL, REGULIER. KOPFSTÜTZE	CLEF, REGLAGE TÊTIÈRE
16	SD-394	TRÆKSTANG STOL I	LEVER FOOTREST MOVEMENT INT. I	PLEUELSTANGE, STUHL I	BARRE DE COUPLAGE, FAUTEUIL
17	SD-395	AKSEL, STOL I	PIN, CHAIR I	ACHSE, STUHL I	ARBRE, FAUTEUIL I
18	SD-396	AKSEL, TRÆKSTANG INTEGRALSTOL	ROD-PIN FOR LEGREST MOVEMENT	ACHSE, PLEUELSTANGE, STUHL	ARBRE, BARRE COUPLAGE, FAUTEUI
19	SD-436	MEKANIK, FODSTØTTE, INTEGRAL I	MECHANICS, FOOTREST, CHAIR I	MECHANIK, FUSSLEHNE, STUHL I	MÉCANIQUE, APPUIE-PIEDS, FAUT.
20	SD-403	RULLE FOR BENSTØTTE, INTEGRAL	ROLL, LEGREST, CHAIR I	ROLLE, BEINSTÜTZE, STUHL I	ROULEAU, APPUI-JAMBES, FAUT.
21	SD-416	UDLØSERARM, NAKKESTØTTE	RELEASE ARM, HEAD REST	AUSLÖSE ARM, KOPFSTÜTZE	DECLENCHEMENT BRAS, TÊTIÈRE
22	SD-417	TRYKSTANG FOR LÅSEPAL STOLE	PRESSURE BAR, LOCKING PAWL	DRUCKSTANGE, SPERRHAKEN	BATON DE PRESSION, CLIQUET
23	SD-418	FJEDER FOR LÅSEPAL STOL	SPRING, LOCKING PAWL, INTEGRAL	FEDER, SPERRHAKEN, INTEGRAL	RESSORT, CLIQUET FERMETURE
24	SD-419	RØRSTIFT 3MM FOR LÅSEPAL STOLE	PIN 3 MM, LOCKING PAWL	STIFT 3 MM, SPERRHAKEN	GOUPILLE 3MM, CLIQUET FERMETU
25	SD-420	LÅSEPAL NAKKESTØTTE STOL	LOCKING PAWL, INTEGRAL CHAIR	SPERRHAKEN, INTEGRAL STUHL	CLIQUET FERMETURE, FAUTEUIL
26	SD-421	UNBRÆKOSKRUE FOR NAKKEST. STOL	SCREW FOR HEAD REST	SCHRAUBE FÜR KOPFSTÜTZE	VIS, TÊTIÈRE
27	SD-437	AKSEL FOR MEK. FODSTØT. STOL I	HINGE-PIN F/INTE.LEGREST G/W	SCHARNIERACHSE, FUSSLEHNE, ST. I	AXE MEC., APPUIE-PIEDS, FAUT. I
28	SD-485	SPINDELMØTRIK, PLAST INT I	SPINDLE NUT, PLASTIC, INT. I	SPINDELSCHR.MUTTER, PLAST INT I	ECROU TRANSLAT., PLAST, INT I
29	MH-428	SIKRING TRÆG 500 MA	FUSE T 500 MA	SICHERUNG T 500 MA	FUSIBLE T 500 MA
30	ML-069	SIKRING TRÆG 6,3 A	FUSE T 6.3 A	SICHERUNG T 6,3 A	FUSIBLE T 6,3 A
31	SD-480	PLASTKNAP FOR SD-423	PLASTIC BUTTON FOR SD-423	PLASTKNOPF FÜR SD-423	BOUTON PLASTIQUE, POUR SD-423
32	SD-493	KNOP, PROGRAMVÆLGER STOL I-ALPA	KNOP, PROG.SELECTOR CHAIR	KNOPF, PROGRAMMWÄHLER, STUHL	BOUTON, SELECTEUR PROG. FAUTEU
33	UC-563	SKRUE, INTEGRALSTOL	SCREW, INTEGRAL PATIENT CHAIR	SCHRAUBE, INTEGRAL STUHL	VIS, FAUTEUIL INTEGRAL
34	SD-451	AKSEL FOR INTEGRAL STOL I	HINGE-PIN	ACHSE	AXE
35	SD-452	AKSEL FOR INTEGRAL STOL I	HINGE-PIN	ACHSE	AXE
36	SD-449	VINKELBESLAG INTEGRAL STOL	ANGLE PLATE, CHAIR, INTEGRAL	WINKELBESCHLAG, STUHL, INTEGRAL	PATTE DE FIXATION, INTEGRAL
37	SD-471	GEVINDSKRUE FOR UC-271	ATTACHMENT ALLEN PIN MOTOR	GEWINDESCHRAUBE FÜR UC-271	VIS TARAUE POUR UC-271
38	SD-481	STILLESKRUE MED MØTRIK STOL I	SCREW/NUT	SCHRAUBENMUTTER	VIS/ECROU
39	SC-886	PLASTGREB, INTEG.NAKKESTØT. I, II	LOCKING HANDLE FOR HEADREST	VERSCHLUSSGRIFF FÜR KOPFSTÜTZE	POIGNÉE FERMETURE POUR TÊTIÈ
40	SD-381	NYLON STYR, RYG BASIS STOL I	NYLON GUIDEWAY, BACK BASE	NYLON STEUERUNG, RÜCKLEHN BASIS	GUIDAGE NYLON, DOSSIER
	SD-413	Kondensator, Rygmotor	Condensator, bagrest motor	Kondensator, Rückenmotor	
	SD-431	Kondensator, Løftmotor	Condensator, lifting motor	Kondensator, Hubmotor	

type I



Patientstol, type 89 Polster Kunststofdele


Patient chair, type 89 Upholstery. Plast Parts

YA-001

Patientenstuhl, Typ 89 Polster. Kunststoffteile

Fauteuil type 89, Rembourrage. Pièces en matière synthétique 0.6.20.B

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
	AC-144	KABEL STOL KOMM INT	CABLE CHAIR CONTROL	KABEL STUHLSTEUERUNG	CABLE COMMANDE FAUTEUIL
	UC-347	SILICONESLANGE Ø 21 X 16,8	SILICONE TUBE Ø 21 X 16,8	SILIKONSCHLAUCH Ø 21 X 16,8	TUYAU SILICONE Ø 21 X 16,8
	AC-886	KIT FOR BENSTØTTE STOL I	KIT, LEGREST	KIT, BEINSTÜTZE	KIT, APPUIE-JAMBES
	BA-137	TILLÆG, PLAST, BENSTØTTE	PROTECTION COVER OF PLASTIC	SCHUTZBEZUG AUS KUNSTSTOFF	PROTECTION EN PLASTIC
1	UC-242	BASIS, RYG, STOL I	BASE, BACKREST, INT. I CHAIR	BASIS, RÜCKENLEHNE, INT. I STUHL	BASE, DOSSIER, FAUTEUIL INT.
2	UC-257	BASIS, BENSTØTTE, STOL I-II	BASE, LEGREST, PATIENT CHAIR II	BASIS, BEINSTÜTZE, STUHL II	BASE, APPUIE-JAMBES, FAUTEUIL
3	UC-294	ARMLÆN, STOL I	ARMREST, PATIENT CHAIR I	ARMLEHNE, PATIENTENSTUHL	ACCOUDOIR, FAUTEUIL
4	UC-291	BASIS, ARMLÆN VENSTRE STOL I	BASE, ARMREST, PAT. CHAIR, LEFT	ARMLEHNENBASIS, STUHL, LINKS	BASE, ACCOUDOIR GAUCHE, FAUTEUIL
5	UC-293	BASIS, ARMLÆN HØJRE STOL I	BASE, ARMREST, PAT. CHAIR, RIGHT	ARMLEHNENBASIS, STUHL, RECHTS	BASE, ACCOUDOIR DROIT, FAUTEUIL
6	UC-639	POLST. NAKKE. STOL I+II WB-265	UPHOLSTERY, HEADREST, I+IIWB-265	POLSTER, KOPFSTÜTZE, I+II WB-265	REMBOURRAGE, TETIERE, I+IIWB-2
6	UC-227	POLST. NAKKE. STOL I+II WB-230	UPHOLSTERY, HEADREST, I+I WB-230	POLSTER, KOPFSTÜTZE, I+II WB-230	REMBOURRAGE, TETIERE, I+IIWB-2
6	UC-631	POLST. NAKKE. STOL I+II WB-264	UPHOLSTERY, HEADREST, I+IIWB-264	POLSTER, KOPFSTÜTZE, I+II WB-264	REMBOURRAGE, TETIERE, I+IIWB-2
6	UC-635	POLST. NAKKE. STOL I+II WB-263	UPHOLSTERY, HEADREST, I+IIWB-263	POLSTER, KOPFSTÜTZE, I+II WB-263	REMBOURRAGE, TETIERE, I+IIWB-2
6	UC-229	POLSTER, NAKKE. I+II WB-233	UPHOLSTERY, HEADREST I+IIWB-233	POLSTER, KOPFSTÜTZE, I+IIWB-233	REMBOURRAGE, TETIERE, I+IIWB-2
6	UC-230	POLST. NAKKE. STOL I+II WB-234	UPHOLSTERY, HEADREST, I+IIWB-234	POLSTER, KOPFSTÜTZE, I+II WB-234	REMBOURRAGE, TETIERE, I+IIWB-2
6	UC-231	POLSTER, NAKKE, I+II WB-235	UPHOLSTERY, HEADREST I+IIWB-235	POLSTER, KOPFSTÜTZE I+IIWB-235	REMBOURRAGE, TETIERE I+IIWB-2
7	UC-640	POLSTER, RYG, STOL I WB-265	UPHOLSTERY, BACKREST I WB-265	POLSTER, RÜCKLEHNE I WB-265	REMBOURRAGE, DOSSIER, I WB-2
7	UC-232	POLSTER, RYG, STOL I WB-230	UPHOLSTERY, BACKREST I WB-230	POLSTER, RÜCKLEHNE, I WB-230	REMBOURRAGE, DOSSIER, I WB-23
7	UC-632	POLSTER, RYG STOL I WB-264	UPHOLSTERY, BACK, I WB-264	POLSTER, RÜCKLEHNE, I WB-264	REMBOURRAGE, DOSSIER, I WB-2
7	UC-636	POLSTER, RYG, INT. STOL I WB-263	UPHOLSTERY, BACKREST, I WB-263	POLSTER, RÜCKLEHNE, I WB-263	REMBOURRAGE, DOSSIER, I WB-26
7	UC-234	POLSTER, RYG, I WB-233	UPHOLSTERY, BCKREST, I WB-233	POLSTER, RÜCKLEHNE, I WB-233	REMBOURRAGE, DOSSIER, I WB-2
7	UC-235	POLSTER, RYG, STOL I WB-234	UPHOLSTERY, BACK, I WB-234	POLSTER, RÜCKLEHNE, I WB-234	REMBOURRAGE, DOSSIER, I WB-2
7	UC-236	POLSTER, RYG, I WB-235	UPHOLSTERY, BACKREST, I WB-235	POLSTER, RÜCKLEHNE, I WB-235	REMBOURRAGE, DOSSIER, I WB-23
8	UC-642	POLSTER, SÆDE, STOL I WB-265	UPHOLSTERY, SEAT, I WB-265	POLSTER, SITZ, I WB-265	REMBOURRAGE, SIEGE, I WB-265
8	UC-237	POLSTER, SÆDE, STOL I WB-230	UPHOLSTERY, SEAT, I WB-230	POLSTER, SITZ, I WB-230	REMBOURRAGE, SIEGE, I WB-230
8	UC-634	POLSTER, SÆDE, STOL I WB-264	UPHOLSTERY, SEAT, I WB-264	POLSTER, SITZ, I WB-264	REMBOURRAGE, SIEGE, I WB-264
8	UC-638	POLSTERSÆDE, STOL I WB-263	UPHOLSTERY, SEAT, I WB-263	POLSTER, SITZ, I WB-263	REMBOURRAGE, SIEGE, I WB-263
8	UC-239	POLSTER, SÆDE, STOL I WB-233	UPHOLSTERY, SEAT, I WB-233	POLSTER, SITZ, I WB-233	REMBOURRAGE, SIEGE, I WB-233
8	UC-240	POLSTER, SÆDE, STOL I WB-234	UPHOLSTERY, SEAT, I WB-234	POLSTER, SITZ, I WB-234	REMBOURRAGE, SIEGE, I WB-234
8	UC-241	POLSTER, SÆDE, STOL I WB-235	UPHOLSTERY, SEAT, I WB-235	POLSTER, SITZ, I WB-235	REMBOURRAGE, SIEGE, I WB-235
10	UC-253	DÆKSEL, VENSTRE SIDE, STOL I	COVER, LEFT, PATIENT CHAIR	ABDECKUNG, LINKS, PATIENTENSTUHL	CARTER GAUCHE, FAUTEUIL
11	UC-252	DÆKSEL, HØJRE SIDE, STOL I	COVER, RIGHT, PATIENT CHAIR	ABDECKUNG, RECHTS, PATIENT. STUHL	CARTER DROIT, FAUTEUIL
12	UC-254	DÆKSEL, BAGSTE, STOL I	COVER, REAR, PATIENT CHAIR	ABDECKUNG, HINTEN, PATIENT. STUHL	CARTER ARRIERE, FAUTEUIL
13	UC-255	DÆKSEL, LØFTEARM, STOL I	COVER, LIFTING ARM, CHAIR	ABDECKUNG, HUBARM, STUHL	CARTER, BRAS DE LEVAGE, FAUT.
14	UC-256	BASIS, SÆDE, STOL I	BASE, SEAT, PATIENT CHAIR	BASIS, SITZ, PATIENTENSTUHL	BASE, SIEGE, FAUTEUIL
15	UC-269	DÆKSEL, FORREST, STOL I	COVER, FRONT, PATIENT CHAIR	ABDECKUNG, VORNE, STUHL	CARTER AVANT, FAUTEUIL
16	UC-270	BASIS, SÆDERAMME, STOL I	BASE, PAINTED " TRIANGLE "	BASIS, SITZRAHMEN, PAT. STUHL	BASE, CHASSIS, FAUTEUIL
17	UC-272	DÆKSEL, INGEN ARMLÆN, STOL I	COVER, NO ARMREST, CHAIR	ABDECKUNG, KEINE ARMLEHNE, STUHL	CARTER, AUCUN ACCOUDOIR, FAU
18	MC-979	NAKKESTØTTE, INTEGRAL, HØJRE *	HEAD REST, INTEGRAL, RIGHT-HAND *	KOPFSTÜTZE, INTEGRAL, RECHTSH. *	TETIERE, INTEGRAL, DROITIER
18	MC-791	NAKKESTØTTE, INTEGRAL VENSTRE *	HEADREST INTEGRAL, LEFT-HAND *	KOPFSTÜTZE INTEGRAL, LINKSHÄN. *	TETIERE, INTEGRAL, GAUCHER
19	UC-277	POLSTER, BEN, STOL I+II WB-235	UPHOLSTERY, LEGREST, I+II WB-235	POLSTER, BEINSTÜTZE, I+II WB-235	REMBOURRAGE, APPUIE, I+II WB-2
19	UC-641	POLSTER, BEN, STOL I+II WB-265	UPHOLSTERY, LEGREST, I+II WB-265	POLSTER, BEINSTÜTZE, I+II WB-265	REMBOURRAGE, APPUIE. I+II WB-2
19	UC-273	POLSTER, BEN, STOL I+II WB-230	UPHOLSTERY, LEGREST, I+IIWB-230	POLSTER, BEINSTÜTZE, I+IIWB-230	REMBOURRAGE, APPUIE. I+IIWB-2
19	UC-633	POLSTER, BEN. STOL I+II WB-264	UPHOLSTERY, LEGREST, I+IIWB-264	POLSTER, BEINSTÜTZE, I+II WB-264	REMBOURRAGE, APPUIE. I+II WB-2
19	UC-637	POLSTER, BEN STOL I+II WB-263	UPHOLSTERY, LEGREST, I+IIWB-263	POLSTER, BEINSTÜTZE, I+IIWB-263	REMBOURRAGE, APPUIE. I+II WB-2
19	UC-275	POLST. BENSTØT. STOL I WB-233	UPHOLSTERY, LEGREST, BLUE	POLSTER, BEINSTÜTZE, BLAU	REMBOURRAGE, APPUIE-JAMB. , BLE
19	UC-276	POLSTER, BEN, STOL I+II WB-234	UPHOLSTERY, LEGREST, I+IIWB-234	POLSTER, BEINSTÜTZE, I+IIWB-234	REMBOURRAGE, APPUIE. I+IIWB-2
20	SD-427	BASIS PLADE INTEGR. I STOL	BASE BOARD, CHAIR I	BASTSPLATTE, STUHL I	PLAQUE DI MONTAGE, FAUTEUIL I
21	AC-424	DÆKKNAP STOLEBASIS INTEG. I+II	COVER, CHAIR BASE	ABDECKUNG, STUHLBASIS	CAPOT, BASE DU FAUTEUIL




No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	UC-516	FODKRYDS, VARIOFLEX IV	FOOT CROSS, OPERATOR'S STOOL	FUSSKREUZ, ARBEITSSTUHL	CROIX DU PIED DU SIÈGE
2	UC-519	CENTERPART, VARIO IV, HØJ	CENTER PART, VARIOFLEX IV, LONG	MITTELTEIL, VARIOFLEX IV, LANG	PARTIE CENTR., VARIO IV, LONGU
3	UC-606	HJUL. U/BREMSE, VARIOFLEX IV	WHEEL WITHOUT BRAKE, VARIOF. IV	ROLLE OHNE BREMSE, VARIOFLEX IV	ROULETTE SANS FREIN, VARIOF. I
4	FA-367	UDLØSERPLADE, VARIOF. FORSTÆRKET	RELEASE PLATE, VARIOFLEX	AUSLÖSEPLATTE, VARIOFLEX	PLAQUE DECLenchement, VARIO
5	AC-718	HÅNDTAG SÆDE/VIP. NCS4502-G G30	HANDLE, SEAT, INCLINATION	HANDGRIFF, SITZ, NEIGUNG	POIGNEE, SIEGE, INCLINAISON
6	MC-496	UDLØSERPLADE SÆDEREG., VARIOF.	RELEASE PLATE, SEAT, VARIOFLEX	AUSLÖSEPLATTE, SITZ, VARIOFLEX	PLAQUE DECLench., SIEGE, VARI
7	MC-673	KNOP, HÅNDTAG, NCS 4502-G	KNOB, HANDLE, NCS 4502-G	KNOPF, GRIFF, NCS 4502-G	BOUTON, MANETTE, NCS 4502-G
8	AC-720	HÅNDTAG HØJDE. NCS 4502-G. G30	HANDLE, HEIGHT, NCS 4502-G. G30	GRIFF, HÖHE, NCS 4502-G. G30	MANETTE, MONTEE, NCS 4502-G G3
9	AE-274	KIT, GASFJEDER, SÆDE, VARIOFL.	KIT, GAS SPRING, SEAT, VARIOFLEX	KIT, GASFEDER, SITZ, VARIOFLEX	KIT, RESSORT A GAZ, SIEGE, VARI
10	UC-404	SEEGERRING DIN 471-14 X 1	SEEGER-A-RING DIN 471-14 X 1	SEEGER-A-RING DIN 471-14 X 1	BAGUE SEEGER-A DIN 471-14 X
11	UC-406	SEEGERRING DIN 471-8 X 0,8	SEEGER-A-RING DIN 471-8 X 0.8	SEEGER-A-RING DIN 471-8 X 0,8	BAGUE SEEGER-A DIN 471-8 X0,
12	AB-275	KIT, GASFJEDER, RYG, VARIOFL.	KIT, GAS SPRING, BACK REST, VARIO	KIT, GASFEDER, RÜCKLEHNE, VARIOF.	KIT, RESSORT A GAZ, DOSSIER, VF
13	UC-517	CENTERPART, VARIO IV, LAV	CENTRE PART, VARIO IV, LOW	MITTELTEIL, VARIO IV, NIEDRIG	PARTIE CENTRALE, VARIO IV, BA
14	UC-537	FODRING, VARIOFLEX IV	FOOT-RING, VARIOFLEX IV	FUSS-RING, VARIOFLEX IV	REPOSE-PIED VARIOFLEX IV
15	AC-761	KIT, RYGSTOP, VARIOF. III+IV	KIT, BACKREST STOP, VF III+IV	KIT, RÜCKENSTOPP, VF III + IV	KIT, BUTEE DOSSIER, VF III+I
16	MC-494	AKSEL. KONUSB/GASFJ. SÆDEKASSE	AXLE, GAS SPRING, SEAT	ACHSE, GASFEDER, SITZ	ARBRE, RESSORT GAZ, SIEGE
17	MB-017	MØTRIK, SUGESL. HOLDER, M10X1	NUT, SUCT.HOSE HOLDER, M10X1	SCHRAUBENM., SAUGSCHL. HALTER	ECROU, PORTE-TUYAUX ASPIR.
18	SD-456	RYGBØJLE, VARIO IV INCL. BØSNING	BACK SUPPORT, VARIOFLEX IV	BÜGEL, RÜCKLEHNE VARIOFLEX IV	ETRIER, DOSSIER, VARIOFLEX I
19	FA-287	STOPFJEDER (CLIPS) VARIO	STOP SPRING (CLIP) VARIO	STOPPFEDER (KLAMMER), VARIO	RESSORT BLOCAGE (CLIPS) VARI
20	UC-530	GLIDELEJE, GLY.PG 141610 F	JOURNAL BEARING	GLEITLAGER	PALIER A GLISSEMENT
21	FA-288	BØSNING FOR Ø8 AKSEL	SOCKET FOR Ø8 AXLE	BUCHSE FÜR Ø8 ACHSE	DOUILLE POUR AXE Ø8
22	MC-623	BESLAG, RYG, STOP, VARIO. III+IIII	FITTINGS, BACK. STOP, VARIO	BESCHLAG, RÜCKEN, STOPP, VARIO	GARNITURE, DOS, ARRET, VARIO
23	UC-511	SKRUE M 4 X 35 UH-SPAX	SCREW M 4 X 35	SCHRAUBE M 4 X 35	VIS M 4 X 35
24	UC-767	GASFJEDER 300 N SÆDE/RYG/95	GAS SPRING SEAT/BACK/95	GASFEDER SITZ/RÜCKEN/95	RESSORT A GAZ SIEGE/DOS 95
25	FA-538	LEJEAKSEL Ø 14 TVÆRMØTRIK	BEARING AXLE Ø 14 NUT	LAGERACHSE Ø 14 SCHRAUBENM.	ARBRE PALIER Ø 14 ECROU
26	UC-171	MØTRIK M 10 X 1	NUT M 10 X 1	SCHRAUBENMUTTER M 10 X 1	ECROU M 10 X 1

New type 2003 with steel cross - neue Typ 2003 mit Stahlkreuz

1 UC-898 Fodkryds, stål	Foot cross, steel	Fusskreutz, Stahl
2 UC-901 Centerpart, for stålkyds, høj	Center part for steel cross, high	Mittelteil für Stahlkreuz, hohes model
3 UC-899 Hjul for stålkyds	Wheel for steel cross	Rolle für Stahlkreuz
13 UC-900 Centerpart for stålkyds, std	Center part, steel cross, std	Mittelteil f. Stahlkreuz, Std.
14 UC-902 Fodring for stålkyds	Foot ring for steel cross	Fussring für Stahlkreuz



No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	UC-521	BAGPLADE RYG VARIOFLEX IV	REAR PLATE, BACKREST, VARIOF. IV	HINTERPLATTE, RÜCKLEHNE, VF IVV	PLAQUE ARR., DOSSIER VARIOF. I
2	UC-520	UNDERPLADE, SÆDE VARIO FLEX IV	UNDER PLATE, SEAT, VARIOFLEX IV	UNTERPLATTE, SITZ, VARIOF. IV	PLAQUE DESSOUS, SIEGE VF IV
3	UC-655	POLSTER, RYG, VARIO. IV WB-230	UPHOLST. BACK, VARIO. IV WB-230	POLSTER, RÜCK. VARIO. IV WB-230	REMBOUR. DOSSIER VF IV. WB-230
3	UC-652	POLSTER, RYG VARIO. IV WB-263	UPHOLST. BACK, VARIO. IV WB-263	POLSTER, RÜCK. VARIO. IV WB-263	REMBOUR. DOSSIER VF IV. WB-263
3	UC-653	POLSTER, RYG VARIO. IV WB-264	UPHOLST. BACK, VARIO. IV WB-264	POLSTER, RÜCK. VARIO. IV WB-264	REMBOUR. DOSSIER VF IV. WB-264
3	UC-654	POLSTER, RYG VARIO. IV WB-265	UPHOLST. BACK, VARIO. IV WB-265	POLSTER, RÜCK. VARIO. IV WB-265	REMBOUR. DOSSIER VF IV. WB-265
3	UC-656	POLSTER, RYG VARIO. IV WB-234	UPHOLST. BACK, VARIO. IV WB-234	POLSTER, RÜCK. VARIO. IV WB-234	REMBOUR. DOSSIER VF IV. WB-234
3	UC-793	POLSTER, RYG VARIO IV BA-210	UPHOLST. BACK, VARIO. IV BA-210	POLSTER, RÜCK. VARIO. IV BA-210	REMBOUR. DOSSIER VF IV. BA-210
3	UC-794	POLSTER, RYG VARIO IV BA-211	UPHOLST. BACK, VARIO. IV BA-211	POLSTER, RÜCK. VARIO. IV BA-211	REMBOUR. DOSSIER VF IV. BA-211
3	UC-795	POLST. RYG VARIO. IV BA-208	UPHOLST. BACK, VARIO. IV BA-208	POLSTER, RÜCK. VARIO. IV BA-208	REMBOUR. DOSSIER VF IV. BA-208
3	UC-796	POLSTER, RYG VARIO IV BA-209	UPHOLST. BACK, VARIO. IV BA-209	POLSTER, RÜCK. VARIO. IV BA-209	REMBOUR. DOSSIER VF IV. BA-209
4	UC-650	POLSTER, SÆDE VARIO. IV WB-230	UPHOLST. SEAT, VARIO. IV WB-230	POLSTER, SITZ. VARIO. IV WB-230	REMBOUR. SIEGE VARIO. IV WB-23
4	UC-647	POLSTER SÆDE VARIO. IV WB-263	UPHOLST. SEAT, VARIO. IV WB-263	POLSTER, SITZ. VARIO. IV WB-263	REMBOUR. SIEGE VARIO. IV WB-26
4	UC-648	POLSTER, SÆDE VARIOF. IV WB-264	UPHOLST. SEAT, VARIO. IV WB-264	POLSTER, SITZ. VARIO. IV WB-264	REMBOUR. SIEGE VARIO. IV WB-26
4	UC-649	POLSTER, SÆDE VARIOF. IV WB-265	UPHOLST. SEAT, VARIO. IV WB-265	POLSTER, SITZ. VARIO. IV WB-265	REMBOUR. SIEGE VARIO. IV WB-26
4	UC-651	POLSTER SÆDE VARIOF. IV WB-234	UPHOLST. SEAT, VARIO. IV WB-234	POLSTER, SITZ. VARIO. IV WB-234	REMBOUR. SIEGE VARIO. IV WB-23
4	UC-789	POLSTER, SÆDE VARIO IV BA-210	UPHOLST. SEAT, VARIO. IV BA-210	POLSTER, SITZ. VARIO. IV BA-210	REMBOUR. SIEGE VARIO. IV BA-21
4	UC-790	POLSTER, SÆDE VARIO IV BA-211	UPHOLST. SEAT, VARIO. IV BA-211	POLSTER, SITZ. VARIO. IV BA-211	REMBOUR. SIEGE VARIO. IV BA-21
4	UC-791	POLST, SÆDE VARIO IV BA-208	UPHOLST. SEAT, VARIO. IV BA-208	POLSTER, SITZ. VARIO. IV BA-208	REMBOUR. SIEGE VARIO. IV BA-20
4	UC-792	POLSTER, SÆDE VARIO IV BA-209	UPHOLST. SEAT, VARIO. IV BA-209	POLSTER, SITZ. VARIO. IV BA-209	REMBOUR. SIEGE VARIO. IV BA-20
5	AB-745	SKUMSTRIMMEL, SÆDE, VARIOFLEX	FOAM STRIP, SEAT, VARIOFLEX	SCHAUMSTOFFSTREIFEN, SITZ, VF	BANDE MOUSSE, SIEGE, VARIOFL




No	NUMBER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	AC-262	SVINGARM,BAKKE ERGOFLEX (PUR)*	SWIVEL ARM TABLE ERGOFLEX(PUR*	SCHWENKARM TISCH ERGOFLEX(PUR*	BRAS TOURNANT TABLE.ERGO(PUR
1	AC-640	SVINGARM,BAKKE ERGOFLEX (ALU)*	SWIVEL ARM, TABLE ERGOFLEX(ALU*	SCHWENKARM,TISCH ERGOFLEX(ALU*	BRAS TOURNANT,TABL.ERGO(ALU)
2	FA-318	SKÆRM, LYNKOBLINGER,9002 HØJRE	COVER, QUICK CONNECTIONS,9002	ABDECKUNG, SCHNELLKUPPL., 9002	CACHE,RACCORDES RAPIDES, 9002
2	FA-570	SKÆRM FOR LYNKOBLINGER VENSTRE	COVER, QUICK CONNECTIONS, LEFT	ABDECKUNG, SCHNELLKUPPL., LINKE	CACHE,RACCORDES RAPIDES, GAUCH
3	MR-216	DÆKNAP, HVID	COVER PLUG, WHITE	ABDECKSTÖPSEL, WEISS	BOUCHON BLANC
4	UB-056	BØSNING + MØTRIK, FOR SVINGARM	BUSH + NUT,	BUCHSE +SCHRAUBENMUTTER,	DOUILLE + ECROU,
5	SD-118	SKUFFEINDLÆG,ERGOFLEX NY MODEL	DRAWER INSERT, ERGOF.NEW MODEL	SCHUBL.EINSATZ,ERGOF.NEUES MOD	COMPART.ALVÉO.,ERGO.NOUV.MOD
6	FA-067	KAPPEGLASSÆT, A 5 STK.	SET OF STOPP. VIALS (5 VIALS)	KAPPENGLÄSERSATZ (5 GLÄSER)	JEU DE FLACONS A PHARM.(5 FL
6	FA-144	KAPPEGLAS, RØDT	STOPPERED VIAL,RED	KAPPENGLAS, ROT	FLACON A PHARMACIE ROUGE
6	FA-145	KAPPEGLAS, BLÅT	STOPPERED VIAL,BLUE	KAPPENGLAS, BLAU	FLACON A PHARMACIE BLEU
6	FA-146	KAPPEGLAS, GRØNT	STOPPERED VIAL, GREEN	KAPPENGLAS, GRÜN	FLACON A PHARMACIE VERT
6	FA-147	KAPPEGLAS, GULT	STOPPERED VIAL, YELLOW	KAPPENGLAS, GELB	FLACON A PHARMACIE JAUNE
6	FA-148	KAPPEGLAS, SORT	STOPPERED VIAL, BLACK	KAPPENGLAS, SCHWARZ	FLACON A PHARMACIE NOIR
7	FA-361	VATPELLETHOLDER	COTTON PELLET HOLDER	WATTEPELLETS BEHÄLTER	PORTE-BOULETTES DE COTON
8	FA-362	BORSTANDER, ERGOFLEX	BUR BLOCK, ERGOFLEX	BOHRERSTÄNDER, ERGOFLEX	PORTE-FRAISES, ERGOFLEX
9	FA-363	UNDERLAG, BORSTANDER	UNDERLAY, BUR BLOCK	UNTERLAGE, BOHRSTÄNDER	BASE, PORTE-FRAISES
10	FA-627	HOLDER, FLEXVIEW ERGOFLEX	HOLDER, FLEXVIEW ERGOFLEX	HALTER, FLEXVIEW ERGOFLEX	BRAS, FLEXVIEW ERGOFLEX
11	UA-568	PINOLSKRUE, M 5 X 6	POINTED SCREW M 5 X 6	REITSTOCKSCHRAUBE M 5 X 6	VIS POINTUE M 5 X 6
12	UC-008	BØSNING Ø 25 X 30	BUSH Ø 25 X 30	BUCHSE Ø 25 X 30	DOUILLE Ø 25 X 30
13	UR-060	FIBERPAKNING, 32X25,1X1, SORT	FIBRE GASKET, 32X25,1X1, BLACK	FASERDICHTUNG,32X25,1X1SCHWARZ	JOINT EN FIBRE,32X25,1X1,NOI
14	FA-355	DREJEOPHÆNG, FLEXVIEW	SWIVEL SUSPENSION, FLEXVIEW	DREHGELENK, FLEXVIEW	FIXATION PIVOTANTE, FLEXVIEW
15	UA-849	PINOLSKRUE M 8 M/KUGLELÅS	POINTED SCREW,M8+BALL BLOCKING	REITSTOCKSCHRAUBE M8+KUGELSCH.	VIS POINTUE M 8 + OBTURATION
16	YA-088	UNBRAKØNGLE M 2,5	ALLEN KEY M 2,5	INBUSS-SCHLÜSSEL M 2,5	CLE BTR M 2,5
17	AB-893	SPÆNDINGSUDTAG, FLEXVIEW, ERGOF.	POWER CONNECTOR, FLEXVIEW, ERGOF	SPANNUNG, FLEXVIEW, ERGOFLEX	SORTIE TENSION,FLEXVIEW,ERGO
18	AB-895	TRAFØ INKL. MONTRINGSPL., ERGOF	TRANSFORMER + MOUNTING PLATE	TRAFØ + MONTAGEPLATTE, ERGOFLEX	TRANSFORMATEUR + PLAQUE, ERGO
18	AB-217	TRAFØ, FLEX VIEW, TESTET, 100V.	TRANSFORMER, FLEXVIEW, 100V	NUR FUER JAPAN	ONLY FOR JAPAN
19	UA-820	MØTRIK M 5 Z	NUT M 5 Z	SCHRAUBENMUTTER M 5 Z	ECROU M 5 Z
20	UA-160	O-RING 14,3 X 2,4	O-RING 14.3 X 2.4	O-RING 14,3 X 2,4	JOINT TORIQUE 14,3 X 2,4
21	MR-209	O-RING 17,1 X 1,6	O-RING 17.1 X 1.6	O-RING 17,1 X 1,6	JOINT TORIQUE 17,1 X 1,6




No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
	YA-443	LAKPRØVE INTEGRAL	COLOUR SAMPLE,LACQUERED METAL	FARB PROBE, METALL	TENTIER, METAL LAQUE
	YA-444	SKAIPRØVE INTEGRAL	COLOUR SAMPLE SKAI	FARBPROBE SKAI	TEINTIER SKAI
	YA-970	LAKPRØVE INTEGRAL 96	COLOUR SAMPLE,LACQUER.METAL 96	FARB PROBE, METALL 96	TEINTIER, METAL LAQUE 96
	YA-971	SKAIPRØVE INTEGRAL 96	COLOUR SAMPLE SKAI 96	FARBPROBE, SKAI 96	TEINTIER SKAI 96
	YR-026	FEDTSPRØJTE DARINA R2 STOL	GREASE GUN DARINA R2 CHAIR	FETTSPRITZE DARINA R2 STUHL	GRAISSEUR DARINA R2 FAUTEUIL
	YR-032	ENGANGSSPRØJTE MED PARALIQ	SYRINGE WITH PARALIQ GREASE	SPRITZE MIT PARALIQ FETT	SERINGUE AVEC GRAISSE PARALI
1	YR-035	FLEX VAC CLEAN 12 FL. X 800 ML	FLEX VAC CLEAN, 12 X 800 ML	FLEX VAC CLEAN, 12 X 800 ML	FLEX VAC CLEAN, 12 X 800 ML
2	YR-002	SILICONEFEDT	SILICONE GREASE	SILIKONFETT	GRAISSE SILICONE
3	YR-001	FLEX MAKE-UP 150 ML	FLEX MAKE-UP 150 ML	FLEX MAKE-UP 150 ML	MAKE-UP FLEX 150 ML
5	YR-037	REPARATIONSLAK, RAL 9018	REPAIR PAINT, LIGHT GREY, 9018	REPARATURLACK, HELLGRAU, 9018	LAQUE POUR RETOUCHES,GRIS CL
5	YR-038	REPARATIONSLAK, BA-026	REPAIR PAINT, GREY	REPARATURLACK, GRAU	LAQUE POUR RETOUCHES, GRISE
5	YR-039	REPARATIONSLAK, BA-015	REPAIR PAINT, WHITE	REPARATURLACK, WEISS	LAQUE POUR RETOUCHES, BLANCH
5	YR-040	REPARATIONSLAK, BA-049	REPAIR PAINT, BLUE	REPARATURLACK, BLAU	LAQUE POUR RETOUCHES, BLEUE
5	YR-041	REPARATIONSLAK, BA-048	REPAIR PAINT, GREEN	REPARATURLACK, GRÜN	LAQUE POUR RETOUCHES, VERTE
5	YR-003	REPARATIONSLAK, BA-102	REPAIR PAINT, RED	REPARATURLACK, ROT	LAQUE POUR RETOUCHES, ROUGE
5	YR-055	REPARATIONSLAK, BA-206	REPAIR PAINT, BA-206	RAPARATURLACK, BA-206	LAQUE POUR RETOUSHES, BA-20
5	YR-056	REPARATIONSLAK, BA-207	REPAIR PAINT, BA-207	REPARATURLACK, BA-207	LAQUE POUR RETOUCHES, BA-207
5	YR-057	REPARATIONSLAK, BA-205	REPAIR PAINT, BA-205	REPARATURLACK, BA-205	LAQUE POUR RETOUCHES, BA-205
5	YR-058	REPARATIONSLAK, BA-204	REPAIR PAINT, BA-204	REPARATURLACK, BA-204	LAQUE POUR RETOUCHES, BA-204
5	YR-059	REPARATIONSLAK, BA-203	REPAIR PAINT, BA-203	REPARATURLACK, BA-203	LAQUE POUR RETOUCHES, BA-203
5	YR-060	REPARATIONSLAK, BA-202	REPAIR PAINT, BA-202	REPARATURLACK, BA-202	LAQUE POUR RETOUCHES, BA-202
5	YR-061	REPARATIONSLAK, CLASSIC	REPAIR PAINT, CLASSIC	REPARATURLACK, CLASSIC	LAQUE POUR RETOUCHES, CLASSI
6	SC-686	SPRAYOLIE, FLEXLUX H/V 12 FL.	SPRAY OIL,FLEXLUX (12 BOTTLES)	SPRAYÖL, FLEXLUX (12 FLASCHEN)	SPRAY LUBRIF.FLEXLUX(12 BOUT
7	YA-004	UNBRAKONØGLE M 4	ALLEN KEY M 4	INBUSS-SCHLÜSSEL M 4	CLE BTR M 4
7	YA-050	UNBRAKONØGLE M 3	ALLEN KEY M 3	INBUSS-SCHLÜSSEL M 3	CLE BTR M 3
7	YA-051	UNBRAKONØGLE M 5	ALLEN KEY M 5	INBUSS-SCHLÜSSEL M 5	CLE BTR M 5
7	YA-052	UNBRAKONØGLE M 6	ALLEN KEY M 6	INBUSS-SCHLÜSSEL M 6	CLE BTR M 6
7	YA-053	UNBRAKONØGLE M 8	ALLEN KEY M 8	INBUSS-SCHLÜSSEL M 8	CLE BTR M 8
7	YA-088	UNBRAKONØGLE M 2,5	ALLEN KEY M 2,5	INBUSS-SCHLÜSSEL M 2,5	CLE BTR M 2,5
7	YA-185	UNBRAKONØGLE 5/32"	ALLEN KEY WG 5/32"	INBUSS-SCHLÜSSEL WG 5/32"	CLE BTR WG 5/32"
8	ML-059	SIKRING FLINK 1,6AF	FUSE 1,6 A	SICHERUNG 1,6 A	FUSIBLE 1,6 A
8	ML-060	SIKRING TRÆG 1 A	FUSE T 1 A	SICHERUNG, T 1 A	FUSIBLE, T 1 A
8	ML-061	SIKRING TRÆG 1,6A	FUSE T 1,6 A	SICHERUNG T 1,6 A	FUSIBLE T 1,6 A
8	ML-062	SIKRING TRÆG 630 MA	FUSE 630 MA	SICHERUNG 630 MA	FUSIBLE 630 MA
8	ML-063	SIKRING TRÆG 2 A	FUSE T 2 A	SICHERUNG T 2 A	FUSIBLE T 2 A
8	ML-064	SIKRING TRÆG 3,15 A	FUSE 3,15 A	SICHERUNG 3,15 A	FUSIBLE 3,15 A
8	ML-065	SIKRING TRÆG 5 A	FUSE T 5 A	SICHERUNG T 5 A	FUSIBLE T 5 A
8	ML-548	SIKRING TRÆG 8 A	FUSE T 8 A	SICHERUNG T 8 A	FUSIBLE T 8 A
8	ML-066	SIKRING TRÆG 2,5 A	FUSE T 2,5 A	SICHERUNG T 2,5 A	FUSIBLE T 2,5 A
8	ML-067	SIKRING FLINK 2A	FUSE 2 A	SICHERUNG 2 A	FUSIBLE 2 A
8	ML-068	SIKRING TRÆG 4 A	FUSE T 4 A	SICHERUNG T 4 A	FUSIBLE T 4 A
8	ML-069	SIKRING TRÆG 6,3 A	FUSE T 6,3 A	SICHERUNG T 6,3 A	FUSIBLE T 6,3 A
8	ML-158	SIKRING TRÆG 32 MA	FUSE T 32 MA	SICHERUNG T 32 MA	FUSIBLE T 32 MA
8	ML-159	SIKRING FLINK 3,15A	FUSE F 3,15 A	SICHERUNG F 3,15 A	FUSIBLE F 3,15 A
8	MH-373	SIKRING TRÆG 10A	FUSE T 10 A	SICHERUNG T 10 A	FUSIBLE T 10 A
8	MH-427	SIKRING TRÆG 200 MA	FUSE T 200 MA	SICHERUNG T 200 MA	FUSIBLE T 200 MA
8	MH-428	SIKRING TRÆG 500 MA	FUSE T 500 MA	SICHERUNG T 500 MA	FUSIBLE T 500 MA
8	MH-430	SIKRING FLINK 6,3 A	FUSE F 6,3 A	SICHERUNG F 6,3 A	FUSIBLE F 6,3 A
8	MH-435	SIKRING TRÆG 315 MA	FUSE T 315 MA	SICHERUNG T 315 MA	FUSIBLE T 315 MA
9	ML-280	KABELBINDER 95 MM	CABLE BINDER 95 MM	KABELBINDER 95 MM	SERRE-CABLE 95 MM
9	ML-281	KABELBINDER 150 MM	CABLE BINDER 150 MM	KABELBINDER 150 MM	SERRE-CABLE 150 MM
9	ML-282	KABELBINDER 280 MM	CABLE BINDER 280 MM	KABELBINDER 280 MM	SERRE-CABLE 280 MM
10	MC-453	FASTNØGLE 13/22 MM	FORK SPANNER 13/22 MM	GABELSCHLÜSSEL 13/22 MM	CLE PLATE 13/22 MM
11	MC-500	FASTNØGLE NV 13	FORK SPANNER NV 13	GABELSCHLÜSSEL NV 13	CLE PLATE NV 13
12	FA-041	MANOMETER 4 HUL	MANOMETER 4-HOLE	MANOMETER 4-LOCH	MANOMÈTRE 4-VOIES

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
13	FA-400	FLOWMETER, KØLELUFT	FLOWMETER, COOLING AIR	MENGENMESSER, KÜHLLUFT	DEBITMETRE, AIR DE REFROIDIS
14	YR-031	TRI-FLOW SMØRREMIDDEL 100 ML	TRI-FLOW LUBRICANT 100 ML	TRI-FLOW SCHMIERMITTEL 100 ML	LUBRIFIANT TRI-FLOW 100 ML
15	UC-342	FEDT, PRESSOL 40 G. NR.:10164	GREASE, PRESSOL 40 G NR.:10164	FETT, PRESSOL 40 G NR.: 10164	GRAISSE PRESSOL 40 G NR.:101
16	AC-892	TESTFLOWMETER, VAND	TESTFLOWMETER, WATER	TESTMENGENMESSER, WASSER	TESTDEBITMETRE, EAU
17	SD-240	MANOMETER QF 4+4, INTEGRAL.	MANOMETER QF.4+4,INTEGRAL	MANOMETER QF.4+4,INTEGRAL	MANOMETRE QF.4+4,INTEGRAL
18	AC-188	SERVICEKIT ANTISTATISK	SERVICE KIT, ANTISTATIC	SERVICESATZ, ANTISTATISCH	SERVICE KIT ANTI-STATIQUE
19	UC-600	HÅNDELSREM KOMPLET FOR AC-188	HAND STRAP COMPLETE FOR AC-188	ARMBAND KOMPLETT FÜR AC-188	BRACELET COMPLET POUR AC-188
20	MR-114	SELVKLÆBENDE KABELAFLASTNING	CABLE RELIEF, ADHESIVE	KABELENTLASTUNG, SELBSTKLEBEND	SERRE-CÂBLE AUTOCOLLANT
21	MR-142	SELVKLÆBENDE KABELAFL. 9,5 MM	CABLE RELIEF, ADHESIVE 9.5 MM	KABELENTLAST,SELBSTKLEB. 9,5MM	SERRE-CABLE AUTOCOLLANT 9,5M
22	UC-158	KABELCLIPS 15 X 17 M/KLÆBEFOD	CABLE CLIPS 15X17, ADHESIVE	KABELCLIPS 15X17,SELBSTKLEBEND	CLIPS CÂBLE 15X17 AUTOCOLLAN
23	SD-391	NØGLE,OVERPART,REDUKTIONSVENT.	KEY, UPPER VALVE, REDUCTION V.	SCHLÜSSEL,VENTIL OBEN, REDUZ.	CLE,VANNE SUPÉRIEURE,DÉTENT.
24	SD-392	NØGLE,BUNDVENT.REDUKTIONSVENT.	KEY,BOTTOM VALVE, REDUCTION V.	SCHLÜSSEL,BODENVENTIL,REDUZIER	CLE,VANNE INFÉRIEURE,DETENT.
25	SD-393	NØGLE,BEHOLDER,LUFTFILTER	KEY, RECEPTACLE, AIR FILTER	SCHLÜSSEL, BEHÄLTER,LUFTFILTER	CLE, BIDON, FILTRE A AIR
MH-847	Kabletester til signalkabel	Cable tester for signal cable	Kabletester für Signalkabel	Adapteur, test cable signal	
SD-388	Nøgle til justering, nakkestøtte	Key for adjustment head rest	Schlüssel für Justierung Kopfstütze	Cle, tetiere	



No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	AD-108	PRINTSÆT, FLEXVIEW TYPE 287	PRINTED CIRCUIT, SET, TYPE 287	PLATINE, SATZ, TYP 287	CARTE, JEU, TYP 287
2	AD-105	PRINT, FLEXVIEW CONTROL, TY.287	PRINTED CIRC. "FV CONTROL" 287	PLATINE "FLEXVIEW CONTROL" 287	CARTE "FLEXVIEW CONTROL" 287
3	AD-106	PRINT, FLEXVIEW SWITCH PANEL	PR.CIRC"FLEXVIEW SWITCH PANEL"	PLATINE"FLEXVIEW SWITCH PANEL"	CARTE "FLEXVIEW SWITCH PANEL"
4	FA-392	RØRREFLEKTOR, FLEXVIEW	TUBE REFLECTOR, FLEXVIEW	REFLEKTOR, FLEXVIEW	RÉFLECTEUR, FLEXVIEW
5	FA-396	ENDEREFLEKTOR, FLEXVIEW	END REFLECTOR, FLEXVIEW	ENDREFLEKTOR, FLEXVIEW	RÉFLECTEUR D'EXTRÉMITÉ
6	MA-129	ALU-STAG, FLEXVIEW TYPE 287	SUPPORT PIN, ALUMINIUM, TYPE 287	STIFT, ALUMINIUM, TYP 287	GOUPILLE, ALUMINIUM, TYPE 28
7	MA-139	GUMMIBEN FLEXVIEW	RUBBER LEG FLEXVIEW	GUMMIFUSS FLEXVIEW	PIED EN CAOUTCHOUC FLEXVIEW
8	MA-236	ACRYLPLADE, FLEXVIEW TY.287	ACRYLIC PLATE, FLEXVIEW T. 287	AKRYLLPLATTE, FLEXVIEW TYP 287	PLAQUE ACRYLIQUE, FLEXVIEW 28
9	MH-107	KABEL, FLEXVIEW/KIRURGI VENTIL	CABEL, COILED, FLEXVIEW	SPIRALLEITUNG, FLEXVIEW	CORDON TÉLÉPHONE, FLEXVIEW
10	UA-706	PLADESKRUE 3,5 X 4,5 PHPX-B	SCREW 3.5 X 4.5 PHPX-B	SCHRAUBE 3,5 X 4,5 PHPX-B	VIS 3,5 X 4,5 PHPX-B
11	UA-728	SKRUE KB 3,0 X 6	SCREW KB 3.0 X 6	SCHRAUBE KB 3,0 X 6	VIS KB 3,0 X 6
12	UA-855	SKRUE M 5 X 16	SCREW M 5 X 16	SCHRAUBE M 5 X 16	VIS M 5 X 16
13	UC-018	PRINTKLEMME, FLEXVIEW	TERMINAL POINT, FLEXVIEW	KLEMME, FLEXVIEW	TERMINAL, FLEXVIEW
14	MK-135	LYSSTOFRØR TYND 8W/HVID	FLUORESCENT TUBE THIN 8W/WHITE	LEUCHTSTOFFRÖHRE DÜNN 8W/WEISS	TUBE LUMINESCENT MINCE 8W/BL
15	UC-057	FILMHOLDER, FLEXVIEW	FILM HOLDER, FLEXVIEW	FILMHALTER, FLEXVIEW	PORTE-FILM, FLEXVIEW
16	AD-109	PRINTSÆT, FLEXVIEW TYPE 210	PRINTED CIRCUIT, SET, TYPE 210	PLATINE, SATZ, TYP 210	CARTE, JEU, TYP 210
17	AD-103	PRINT, FLEXVIEW CONTROL, TY.210	PRINTED CIRC. "FC CONTROL" 210	PLATINE "FLEXVIEW CONTROL" 210	CARTE "FLEXVIEW CONTROL" 210
18	FA-319	RØRREFLEKTOR, FLEXVIEW	TUBE REFLECTOR, FLEXVIEW	REFLEKTOR, FLEXVIEW	RÉFLECTEUR, FLEXVIEW
19	FA-395	BÆRETAP, FLEXVIEW 50/210	PIN, FLEXVIEW 50/210	ZAPFEN, FLEXVIEW 50/210	TENON, FLEXVIEW 50/210
20	MA-083	ALU-STAG, FLEXVIEW TYPE 210	SUPPORT PIN, ALUMINIUM, TYPE 210	STIFT, ALUMINIUM, TYP 210	GOUPILLE, ALUMINIUM, TYPE 21
21	MA-235	AKRYLPLADE, FLEXVIEW 210	ACRYLIC PLATE, FLEXVIEW 210	AKRYLLPLATTE, FLEXVIEW 210	PLAQUE ACRYLIQUE, FLEXVIEW 21
22	MK-056	LYSSTOFRØR TYND 6W/33 (HVID)	FLUORESCENT THIN 6W/33, WHITE	LEUCHTSTOFFRÖHRE DÜNN 6W/33, W.	TUBE LUMINESCENT MINCE 6W/33
23	UA-611	SKRUE, GKN M 8 X 25	SCREW, GKN M 8 X 25	SCHRAUBE, GKN M 8 X 25	VIS, GKN M 8 X 25
24	AD-110	PRINTSÆT, FLEXVIEW TYPE 50	PRINTED CIRCUIT, SET, TYPE 50	PLATINE, SATZ, TYP 50	CARTE, JEU, TYP 50
25	AD-102	PRINT, FLEXVIEW CONTROL TY. 50	PRINTED CIRC. "FV CONTROL" 50	PLATINE "FLEXVIEW CONTROL" 50	CARTE "FLEXVIEW CONTROL" 50
26	AD-111	MB-256 FLEXVIEW SWITCH PANEL	PR.CIRC"FLEXVIEW SWITCH PANEL"	PLATINE"FLEXVIEW SWITCH PANEL"	CARTE "FLEXVIEW SWITCH PANEL"
27	MA-086	ALU-STAG, FLEXVIEW TYPE 50	SUPPORT PIN, ALUMINIUM, TYPE 50	STIFT, ALUMINIUM, TYP 50	GOUPILLE, ALUMINIUM, TYPE 50
28	MA-137	ACRYLPLADE, FLEXVIEW TY. 50	ACRYLIC PLATE, FLEXVIEW T. 50	AKRYLLPLATTE, FLEXVIEW TYP 50	PLAQUE ACRYLIQUE, FLEXVIEW 5
33	ML-164	TLF-STIK JACK 4POL	TELEPHONE POINT JACK 4-POLE	TELEPHONANSCHLUSS JACK 4-POLIG	PRISÉ POUR TEL. JACK 4 POLE
34	AB-667	HF TRANSFORMATOR, FLEXVIEW	TRANSFORMER, FLEXVIEW	TRANSFORMATOR, FLEXVIEW	FRANSFORMATEUR, FLEXVIEW
34	AE-217	TRAFO, FLEX VIEW, TESTET, 100V.	TRANSFORMER, FLEXVIEW, 100V	NUR FUER JAPAN	ONLY FOR JAPAN
37	FA-447	PRIMÆR ENDEBUND FLEXVIEW 9002	END COVER, FLEXVIEW, 9002	ENDABDECKUNG, FLEXVIEW, 9002	CARTER FLEXVIEW, 9002
38	FA-448	SEKUND. ENDEBUND FLEX VIEW 9002	END COVER, FLEXVIEW, 9002	ENDABDECKUNG, FLEXVIEW, 9002	CARTER FLEXVIEW, 9002
39	MA-107	PROFIL FLEXVIEW 9002	PROFILE FLEXVIEW 9002	PROFIL FLEXVIEW 9002	PROFIL FLEXVIEW 9002
40	MA-134	PROFIL FLEXVIEW 9002	PROFILE FLEXVIEW 9002	PROFIL FLEXVIEW 9002	PROFIL FLEXVIEW 9002




Flex Dental A/S
Tandrenser, Supprasson
Zahnreiniger, Supprasson

Scaler, Supprasson
Détartreur, Détartreur

1999-06-28
YA-001
0.6.26.B

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	AC-771	BALANCEOPHÆNG SUPPRAS. SCALER	BALANCED SUSPENSION, SUPRASSON	BALANCEAUFHÄNGUNG, SUPRASSON	POUET SUPRASSON
2	SD-387	STYRING SUPRASSON SP3055-B	CONTROL, SUPRASSON SP3055-B	STEUERUNG, SUPRASSON SP3055-B	COMMANDE SUPRASSON SP3055-B
3	AD-426	PRINT, SUPRASSON INTERFACE	PR. CIRCUIT, SUPRASS. INTERFACE	PLATINE, SUPRASSON INTERFACE	CARTE SUPRASSON INTERFACE
4	AC-441	KABEL, SUPRASSON POWER, FC	CABLE, SUPRASSON POWER, FC	KABEL, SUPRASSON LEISTUNG, FC	CABLE SUPRASSON PUISSANCE, F
5	AC-442	KABEL, SUPRASSON QUICK FLEX, FC	CABLE, SUPRASSON QUICKFLEX, FC	KABEL, SUPRASSON QUICKFLEX, FC	CABLE SUPRASSON QUICKFLEX, F
6	AC-443	KABEL, SUPRASSON CONTROL, FC	CABLE, SUPRASSON CONTROL, FC	KABEL, SUPRASSON STEUERUNG, FC	CABLE, COMMANDE SUPRASSON, F
7	SC-927	INSTRUMENT, SUPRASSON, NR. 1	INSTRUMENT, SUPRASSON, NO. 1	INSTRUMENT, SUPRASSON, NR. 1	INSTRUMENT, SUPRASSON, NO. 1
8	SC-928	INSTRUMENT, SUPRASSON, NR. 2	INSTRUMENT, SUPRASSON, NO. 2	INSTRUMENT, SUPRASSON, NR. 2	INSTRUMENT, SUPRASSON, NO 2
9	SD-018	INSTRUMENT, SUPRASSON, NR. 3	INSTRUMENT, SUPRASSON, NO. 3	INSTRUMENT, SUPRASSON, NR. 3	INSTRUMENT, SUPRASSON, NO 3
10	SD-022	VÆRKTØJ, INSTRUMENT, SUPRASSON	TOOL, INSTRUMENT, SUPRASSON	WERKZEUG, INSTRUMENT, SUPRASSON	OUTIL, INSTRUMENT, SUPRASSON
11	SD-021	HÅNDSTYKKE, SUPRASSON, F00225	HANDPIECE, SUPRASSON	HANDSTÜCK, SUPRASSON	PIÈCE À MAIN SUPRASSON
12	MC-860	BESLAG, STYRING, SCALER, SUPRA	COVER, OIL SUMP, SUPRASSON	BESCHLAG, ÖLWANNE, SUPRASSON	COUVERCLE, HUILE EXCEDENTAIR
13	BA-087	ADAPTOR, HYGIFLEXT, SUPRASSON	ADAPTOR, HYGIFLEXT, SUPRASSON	ADAPTOR, HYGIFLEXT, SUPRASSON	ADAPTATEUR, HYGIFLEXT, SUPRASS



No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
	SD-507	KIT, DÚRR FONTÆNEVENTIL	KIT,DÚRR SPITOONVALVE	KIT,DÚRR SPEIBECKENVENTIL	KIT,DÚRR VANNE DE CRACHOIR
	SD-509	KIT, DÚRR FONTÆNEVENT. CLASSIC	DÚRR SPITOONVALVE CLASSIC	DÚRR SPEIBECKENVENTIL CLASSIC	DÚRR VANNE DE CRACHOIR CLASS
1	UC-694	FONTÆNEVENTIL, DÚRR 7110	VALVE, CUSPIDOR, DÚRR 7110	VENTIL, SPEIFONTÆNE, DÚRR 7110	VANNE, CRACHOIR, DÚRR 7110
2	UC-695	SVØMMERHUS KOMPLET, DÚRR 7110	FLOAT CHAMBER, COMPLETE, DÚRR	SCHWIMMERGEHÄUSE,KOMPLETT,DÚRR	CHAMBRE FLOTTEUR COMPLET,DÚRR
3	UC-696	REED-KONTAKT, SVØM. DÚRR 7110	REED SWITCH, FLOAT, DÚRR 7110	REED KONTAKT, SCHWIMMER, DÚRR	CONTACT REED, FLOTTEUR DÚRR
4	UC-697	HOVEDPRINT, DÚRR 7110	MAIN PRINTED CIRCUIT,DÚRR 7110	HAUPTPLATINE, DÚRR 7110	PLATINE GENERALE DÚRR 7110
5	UC-698	PRINT,MOTOROVERVÅGN. DÚRR 7110	PR. CIRCUIT, MOTOR SUPERVISION	PLATINE, MOTORÜBERWACHUNG 7110	CARTE,SUPERVISION MOTEUR 711
6	UC-699	TASTATUR, DÚRR 7110	KEYBOARD, DÚRR 7110	TASTATUR, DÚRR 7110	JEU DE CONTACTS DÚRR 7110
7	UC-701	SL. OVERLØB SP.ENHED DÚRR 7110	OVERFLOW HOSE,FLUSHING UNIT	ÜBERLAUF SCHL.,SPÜLEINHEIT	DEBORDEMENT TUYAU,UNITE RINC
8	UC-702	MICROSWITCH,BEHOLDER,DÚRR 7110	MICROSWITCH, RECEPTACLE, DÚRR	MIKROSCHALTER, BEHÄLTER, DÚRR	MICRO-INTERRUPTEUR,BIDON,DÚRR
9	UC-703	KONTRAKLAP, AFLØB, DÚRR 7110	NON-RETURN FLAP, OUTLET, 7110	ABSPERRKLAPPE, ABFLUSS, 7110	CLAPET DE RETENUE, ECOULEMEN
10	UC-704	VINKEL,SVØMMERHUS,DÚRR 7110	ANGLE, FLOAT CHAMBER, DÚRR	WINKEL,SCHWIMMERGEHÄUSE, DÚRR	COUDE, CHAMBRE FLOTTEUR, DÚRR
11	UC-705	VINKEL, FONTÆNEVENTIL,DÚRR 7110	ANGLE, CUSPIDOR VALVE, DÚRR	WINKEL, SPEIFONTÆNEVENTIL	ANGLE, VANNE CRACHOIR, DÚRR
12	UC-706	VINKEL, SUG, DÚRR 7110	INKL.SI ANGLE,SUCTION,STRAINER INCL.	WINKELSTÜCK KOMPLETT + SIEB	ANGLE,ASPIRATION,FILTRE INCL
13	UC-686	SLANGE Ø4 LUFT, DÚRR 7110	TUBE Ø4 AIR, DÚRR 7110	SCHLAUCH Ø4 LUFT, DÚRR 7110	TUYAU Ø4, AIR, DÚRR 7110
14	UC-661	GUMMIVINKEL, SUG, DÚRR 7110	ANGLE,RUBBER,SUCTION,DÚRR 7110	WINKEL,GUMMI,ABSAUGER,DÚRR7110	COUDE CAOUTCHOUC,ASPIR.,DÚRR
15	UC-663	GUMMIMUFFE Ø20, DÚRR 7110	SLEEVE, RUBBER, Ø20, DÚRR 7110	MUFFE, GUMMI, Ø20, DÚRR 7110	MANCHON CAOUTC.,Ø20,DÚRR 711
16	UC-664	AMALGAMAFFALDS BEHOLDER, DÚRR	AMALGAM COLLECTOR, DÚRR 7110	AMALGAMBEHÄLTER, DÚRR 7110	RECEPTACLE AMALGAME DÚRR 711
17	ML-035	SIKRING, PLUG-IN 3,15AT Ø10	FUSE, PLUG-IN 3,15AT Ø 10	SICHERUNG	FUSIBLE
18	AE-081	STUDS, TILSL. SUG INKL.O-RINGE	CONNECTOR, SUCTION, + O-RINGS	STUTZEN, ABSAUGER, + O-RINGE	RACCORD, ASPIRATION, +JOINTS
19	AC-808	AMALGAMSEPARATOR INTEGRAL 92 *	AMALGAM SEPARATOR INT. 92 *	AMALGAMABSCHIEDER INT. 92 *	SEPARATEUR AMALGAME INT. 92
19	AE-226	AMALGAMSEPARATOR, 100 VOLT	AMALGAM SEPARATOR, 100V	NUR FÜER JAPAN	ONLY FOR JAPAN
20	AE-312	TILSLUTNING, FOR AFLØB	DRAIN CONNECTION	ABFLUSSANSCHLUSS	RACCORD D'ÉCOULEMENT
21	MC-900	BESLAG, TASTATUR, DÚRR	MOUNT.FITTING,DISPLAY,DÚRR	BESCHLAG,TASTATUR,DÚRR	SUPPORT,DISPLAY,DÚRR
22	MC-928	STUDS Ø20-Ø25 ,DÚRR SEPARATOR	CONNECTION Ø20-Ø25 DÚRR SEPAR.	VERBINDUNG Ø20-Ø25 DÚRR, SEP.	RACCORD Ø20-Ø25 DÚRR SEPAR.
23	MH-327	KLRK 2P RIACON 31049102	CONNECTOR 2P RIACON 0201.4901B	STECKER 2P RIACON 0201.4901B	PRISE 2P RIACON 0201.4901B
24	UA-114	SKRUE M 4,2 X 16 PH PLASTITE	SCREW M 4,2 X 16 PH PLASTITE	SCHRAUBE M 4,2 X 16 PLASTITE	VIS M 4,2 X 16 PH PLASTITE
25	UB-006	TANDSKIVE Ø 3 TYPE A	WASHER Ø 3	SCHIBE Ø 3	DISQUE Ø 3
26	UB-086	TANDSKIVE Ø 4 TYPE I	TOOTHED DISK Ø 4 TYPE I	ZAHNRADSCHIBE Ø 4 TYP I	DISQUE DENTÉ Ø 4 TYPE I
27	UC-155	SLANGE 6/4 BLÅ POLYURETHAN	HOSE 6/4 BLUE POLYURETHAN	SCHLAUCH,6/4,BLAU POLYURETHAN	TUYAUX,6/4,BLEU POLYURETHAN
28	UC-662	SPENDEBÅND Ø25	TIGHTENSTRAP Ø25	SCHLAUCHSCHELLE Ø25	RUBAN DE SERRAGE Ø25
29	UC-666	FINGERMØTRIK, M3 MESSING	MILLED NUT, M3, BRASS	RAENDELMUTTER, M3,MESSING	ECROU SPECIALE, M3, LETON
30	UC-687	T-STYKKE M. SNAP, 6 MM	T-CONNECTOR, 6 MM.	T-STUECK, 6 MM.	RACCORD EN T, 6 MM.
31	UC-688	REDUKTION M. SNAP, 6 MM - 4 MM	REDUCTION, 6 MM - 4 MM.	REDUZIERUNG, 6 MM - 4 MM.	REDUCTION, 6 MM - 4 MM.
32	UC-742	BLINDPROP Ø 4MM	COVER PLUG Ø 4 MM.	BLINDSTOPFEN Ø 4 MM.	BOUCHON Ø 4 MM
33	AC-870	DÚRR AMALGAM SEP.220V CABLE	DÚRR AMALGAM SEP.220V CABLE	DÚRR AMALGAMSEP. 220V KABEL	DÚRR SEP.D'AMALGAME 220V CÂB
34	AC-871	DÚRR AMAL.SEP.SIGNAL CABLE	DÚRR AMAL.SEP. SIGNAL CABLE	DÚRR AMAL.SEP. SIGNALKABEL	DÚRR SEP.D'AMALG. SIGNAL CÂB
35	AC-872	DÚRR AMALGAM SEP.24V CABLE	DÚRR AMALGAM SEP. 24V CABLE	DÚRR AMALGAMSEP. 24V KABEL	DÚRR SEP. D'AMALGAME 24V CÂB

AE-057 Manifold, Select Vac/Amalgam Ssep

Left-side bridge
Standard
3



Right-side bridge
Lefthand
2




1



4

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	MC-351	DÆKSEL, ARM TIL SKÆRM	COVER, SCREEN ARM	ABDECKUNG, SCHIRM ARM	OPERCULE, BRAS ECRAN
1	AE-452	DÆKSEL, ARM TIL SKÆRM, CLASSIC	COVER, SCREEN ARM	ABDECKUNG, SCHIRM ARM	OPERCILE, BASE ECRAN
2	AC-079	NTC/1 SKÆRM HØJRE.	* NTC/1 SCREEN, RIGHT	* NTC/1 SCHIRM, RECHTS	* ECRAN NTC/1 DROIT
2	AE-303	NTC/2 SKÆRM HØJRE,	* NTC II SCREEN, RIGHT	* NTC II SCHIRM, RECHTS	* ECRAN NTC II, DROIT
2	AE-737	NTC/2 SKÆRM HØJRE, CLASSIC	* NTC/2 SCREEN, RIGHT CLASSIC	* NTC/2 SCHIRM, RECHTS CLASSIC	* ECRAN NTC/2 DROIT CLASSIC
2	AE-305	FC/V SKÆRM HØJRE	* FC/V SCREEN, RIGHT	* FC/V SCHIRM, RECHTS	* ECRAN FC/V, DROIT
2	AE-738	FC/V SKÆRM HØJRE CLASSIC	* FC/V SCREEN, RIGHT CLASSIC	* FC/V SCHIRM, RECHTS CLASSIC	* ECRAN FC/V DROIT CLASSIC
3	AC-834	NTC/1 SKÆRM, VENSTRE.	* NTC-SCREEN, INTEGRAL, LEFT	* NTC-SCHIRM, INTEGRAL, LINKS	* ECRAN NTC, INTEGRAL, GAUCHE
3	AE-304	NTC/2 SKÆRM VENSTRE	* NTC/2 SCREEN, LEFT	* NTC/2 SCHIRM, LINKS	* ECRAN NTC/2 GAICHE
3	AE-739	NTC/2 SKÆRM VENSTRE CLASSIC	* NTC/2 SCREEN, LEFT CLASSIC	* NTC/2 SCHIRM, LINKS CLASSIC	* ECRAN NTC/2 GAICHE CLASSIC
3	AE-306	FC/V SKÆRM VENSTRE	* FC/V SCREEN, LEFT	* FC/V SCHIRM, LINKS	* ECRAN FC/V, GAUCHE
3	AE-740	FC/V SKÆRM VENSTRE CLASSIC	* FC/V SCREEN, LEFT CLASSIC	* FC/V SCHIRM, LINKS CLASSIC	* ECRAN FC/V GAICHE CLASSIC
4	AE-308	TRIMMEPINDESÆT	TRIMMING PIN SET	TRIMMSTIFTEN	TIGES DE RÉGLAGE



Patientstol, Type 93: Polster. Kunststofdele

Patient Chair, Type 93: Upholstery. Plast Parts

YA-001


Patientenstuhl, Typ 93: Polster. Kunststoffteile

Fauteuil, typ 93: Rembourrage. Pièces en matière synthétique

0.6.29.B


No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
WB-004	BA-137	TILLÆG, PLAST, BENSTØTTE	PROTECTION COVER OF PLASTIC	SCHUTZBEZUG AUS KUNSTSTOFF	PROTECTION EN PLASTIC
1	UC-227	POLST. NAKKE, STOL I+II WB-230	UPHOLSTERY, HEADREST, I+I WB-230	POLSTER, KOPFSTÜTZE, I+II WB-230	REMBOURRAGE, TETIERE, I+IIWB-2
1	UC-230	POLST. NAKKE, STOL I+II WB-234	UPHOLSTERY, HEADREST, I+IIWB-234	POLSTER, KOPFSTÜTZE, I+II WB-234	REMBOURRAGE, TETIERE, I+IIWB-2
1	UC-631	POLST. NAKKE, STOL I+II WB-264	UPHOLSTERY, HEADREST, I+IIWB-264	POLSTER, KOPFSTÜTZE, I+II WB-264	REMBOURRAGE, TETIERE, I+IIWB-2
1	UC-635	POLST. NAKKE, STOL I+II WB-263	UPHOLSTERY, HEADREST, I+IIWB-263	POLSTER, KOPFSTÜTZE, I+II WB-263	REMBOURRAGE, TETIERE, I+IIWB-2
1	UC-639	POLST. NAKKE, STOL I+II WB-265	UPHOLSTERY, HEADREST, I+IIWB-265	POLSTER, KOPFSTÜTZE, I+II WB-265	REMBOURRAGE, TETIERE, I+IIWB-2
1	UC-778	POLST. NAKKE, STOL I+II BA-210	UPHOLSTERY, HEADREST, I+IIBA-210	POLSTER, KOPFSTÜTZE, I+II BA-210	REMBOURRAGE, TETIERE, I+IIBA-2
1	UC-797	POLSTER, NAKKE, I+II BA-208	UPHOLSTERY, HEADREST, I+IIBA-208	POLSTER, KOPFSTÜTZE, I+II BA-208	REMBOURRAGE, TETIERE, I+IIBA-2
1	UC-815	POLSTER, NAKKE, I+II BA-211	UPHOLSTERY, HEADREST, I+IIBA-211	POLSTER, KOPFSTÜTZE, I+II BA-211	REMBOURRAGE, TETIERE, I+IIBA-2
1	UC-811	POLST. NAKKE, STOL I+II BA-209	UPHOLSTERY, HEADREST, I+IIBA-209	POLSTER, KOPFSTÜTZE, I+II BA-209	REMBOURRAGE, TETIERE, I+IIBA-2
2	UC-779	POLS. RYG, STOL II BA-210	UPHOLSTERY, BACK, II BA-210	POLSTER, RÜCKEN, II BA-210	REMBOUR.DOSSIER, II BA-210
2	UC-798	POLSTER, RYG, STOL II BA-208	UPHOLSTERY, BACK, II BA-208	POLSTER, RÜCKEN, II BA-208	REMBOUR.DOSSIER, II BA-208
2	UC-816	POLSTER, RYG, STOL II BA-211	UPHOLSTERY, BACK, II BA-211	POLSTER, RÜCKEN, II BA-211	REMBOUR.DOSSIER, II BA-211
2	UC-812	POLSTER, RYG, STOL II BA-209	UPHOLSTERY, BACK, II BA-209	POLSTER, RÜCKEN, II BA-209	REMBOUR.DOSSIER, II BA-209
2	UC-711	POLSTER, RYG, STOL II WB-230	UPHOLSTERY, BACK, II WB-230	POLSTER, RÜCKEN, II WB-230	REMBOUR.DOSSIER, II WB-230
2	UC-719	POLSTER, RYG, STOL II WB-234	UPHOLSTERY, BACK, II WB-234	POLSTER, RÜCKEN, II WB-234	REMBOUR.DOSSIER, II WB-234
2	UC-715	POLSTER, RYG, STOL II WB-264	UPHOLSTERY, BACK, II WB-264	POLSTER, RÜCKEN, II WB-264	REMBOUR.DOSSIER, II WB-264
2	UC-717	POLSTER RYG, STOL II WB-263	UPHOLSTERY, BACK, II WB-263	POLSTER, RÜCKEN, II WB-263	REMBOUR.DOSSIER, II WB-263
2	UC-713	POLSTER, RYG, STOL II WB-265	UPHOLSTERY, BACK, II WB-265	POLSTER, RÜCKEN, II WB-265	REMBOUR.DOSSIER, II WB-265
3	UC-780	POLS. SÆDE, STOL II BA-210	UPHOLSTERY, SEAT, II BA-210	POLSTER, SITZ, II BA-210	REMBOUR., SIEGE, II BA-210
3	UC-799	POLSTER, SÆDE, STOL II BA-208	UPHOLSTERY, SEAT, II BA-208	POLSTER, SITZ, II BA-208	REMBOUR., SIEGE, II BA-208
3	UC-817	POLSTER, SÆDE, STOL II BA-211	UPHOLSTERY, SEAT, II BA-211	POLSTER, SITZ, II BA-211	REMBOUR., SIEGE, II BA-211
3	UC-813	POLSTER, SÆDE, STOL II BA-209	UPHOLSTERY, SEAT, II BA-209	POLSTER, SITZ, II BA-209	REMBOUR., SIEGE, II BA-209
3	UC-712	POLSTER, SÆDE, STOL II WB-230	UPHOLSTERY, SEAT, II WB-230	POLSTER, SITZ, II WB-230	REMBOUR., SIEGE, II WB-230
3	UC-720	POLSTER, SÆDE, STOL II WB-234	UPHOLSTERY, SEAT, II WB-234	POLSTER, SITZ, II WB-234	REMBOUR., SIEGE, II WB-234
3	UC-716	POLSTER, SÆDE, STOL II WB-264	UPHOLSTERY, SEAT, II WB-264	POLSTER, SITZ, II WB-264	REMBOUR., SIEGE, II WB-264
3	UC-718	POLSTER, SÆDE, STOL II WB-263	UPHOLSTERY, SEAT, II WB-263	POLSTER, SITZ, II WB-263	REMBOUR., SIEGE, II WB-263
3	UC-714	POLSTER, SÆDE, STOL II WB-265	UPHOLSTERY, SEAT, II WB-265	POLSTER, SITZ, II WB-265	REMBOUR., SIEGE, II WB-265
4	UC-781	POLST. BEN, STOL I+II BA-210	UPHOLSTERY, LEGREST, I+IIBA-210	POLSTER, BEINSTÜTZE, I+IIBA-210	REMBOURRAGE, APPUIEJ. I+IIBA-2
4	UC-810	POLSTER, BEN, STOL I+II BA-208	UPHOLSTERY, LEGREST, I+IIBA-208	POLSTER, BEINSTÜTZE, I+IIBA-208	REMBOURRAGE, APPUIEJ. I+IIBA-2
4	UC-818	POLSTER, BEN, I+II BA-211	UPHOLSTERY, LEGREST, I+IIBA-211	POLSTER, BEINSTÜTZE, I+IIBA-211	REMBOURRAGE, APPUIEJ. I+IIBA-2
4	UC-814	POLSTER, BEN, STOL I+II BA-209	UPHOLSTERY, LEGREST, I+IIBA-209	POLSTER, BEINSTÜTZE, I+IIBA-209	REMBOURRAGE, APPUIEJ. I+IIBA-2
4	UC-273	POLSTER, BEN, STOL I+II WB-230	UPHOLSTERY, LEGREST, I+IIWB-230	POLSTER, BEINSTÜTZE, I+IIWB-230	REMBOURRAGE, APPUIEJ. I+IIWB-2
4	UC-276	POLSTER, BEN, STOL I+II WB-234	UPHOLSTERY, LEGREST, I+IIWB-234	POLSTER, BEINSTÜTZE, I+IIWB-234	REMBOURRAGE, APPUIEJ. I+IIWB-2
4	UC-633	POLSTER, BEN, STOL I+II WB-264	UPHOLSTERY, LEGREST, I+IIWB-264	POLSTER, BEINSTÜTZE, I+II WB-264	REMBOURRAGE, APPUIE. I+II WB-2
4	UC-637	POLSTER, BEN, STOL I+II WB-263	UPHOLSTERY, LEGREST, I+IIWB-263	POLSTER, BEINSTÜTZE, I+IIWB-263	REMBOURRAGE, APPUIE. I+II WB-2
4	UC-641	POLSTER, BEN, STOL I+II WB-265	UPHOLSTERY, LEGREST, I+II WB-265	POLSTER, BEINSTÜTZE, I+II WB-265	REMBOURRAGE, APPUIE. I+II WB-2
5	MC-875	SKØRTE, CENTERBOX, H, STOL II	COVERING, CENTRE BOX, RIGHT	ABDECKUNG, MITTELKASTEN, RECHTS	COUVERTURE, BOITE CENT., DROIT
5	MC-876	SKØRTE, CENTERBOX, V, STOL II	COVERING, CENTRE BOX, LEFT	ABDECKUNG, MITTELKASTEN, LINKS	COUVERTURE, BOITE CENT, GAUCH
6	MC-848	DAKSEL, STOL II HØJRE INTEG.	COVER, RIGHT II	ABDECKUNG, RECHTS II	CARTER, DROIT II
6	ME-142	DAKSEL, STOL II HØJRE CLASSIC	COVER, RIGHT II CLASSIC	ABDECKUNG, RECHTS II CLASSIC	CARTET, DROIT II CLASSIC
6	MC-847	DAKSEL, STOL II VENSTRE INTEGR.	COVER, LEFT II	ABDECKUNG, LINKS II	CARTET, GAUCHE II
6	ME-141	DAKSEL, STOL II VENSTRE CLASSIC	COVER, LEFT II CLASSIC	ABDECKUNG, LINKS II CLASSIC	CARTET, GAUCHE II CLASSIC
7	MC-870	PLUGS HØJRE ARMREST, STOL II	ARMREST CAP CHAIR II RIGHT	DECKSTÖPSEL, RECHTE ARMLEHNE	BOUCHON, ACCOUDOIR DROIT
7	MC-977	PLUGS VENSTRE ARMREST, STOL II	ARMREST CAP CHAIR II LEFT	DECKSTÖPSEL, LINKTE ARMLEHNE	BOUCHON, ACCOUDOIR GAUCHE
8	ME-131	DAKSEL, ØVERSTE STOL II CLASSIC	COVER, ARM, UPPER CLASSIC	ABDECKUNG, ARM, OBEN CLASSIC	CARTER, BRAS DESSUS CLASSIC
8	MC-712	DAKSEL, PARALLELARM, ØVERSTE II	COVER, ARM, UPPER	ABDECKUNG, ARM, OBEN	CARTER, BRAS DESSUS
9	ME-132	DAKSEL, NEDERST STOL II CLASSIC	COVER, ARM, LOWER CLASSIC	ABDECKUNG, ARM, UNTEN CLASSIC	CARTER, BRAS DESSOUS CLASSIC
9	MC-721	DAKSEL, PARALLELARM, NEDERSTE II	COVER, ARM, LOWER	ABDECKUNG, ARM, UNTEN	CARTER, BRAS DESSOUS
10	ME-133	DAKSEL, HØJRE, II CLASSIC INT.	COVER, RIGHT II CLASSIC	ABDECKUNG, RECHT II CLASSIC	CARTER, DROIT II CLASSIC
10	MC-714	DAKSEL, HØJRE, STOL II INTEG.	COVER, RIGHT II	ABDECKUNG, RECHT II	CARTER, DROIT II
10	ME-134	DAKSEL, VENST. II CLASSIC INT.	COVER, LEFT II CLASSIC	ABDECKUNG, LINKS II CLASSIC	CARTER, GAUCHE II CLASSIC
10	MC-715	DAKSEL, VENSTRE, STOL II INTEG	COVER, LEFT II	ABDECKUNG, LINKS II	CARTER, GAUCHE II
11	ME-135	FRONTDAKSEL, STOL II CLASSIC	FRONT PLATE II CLASSIC	FRONTPLATTE II CLASSIC	PANNEAU AVANT, II CLASSIC
11	MC-978	FRONTDAKSEL, AFBRYDER, STOL II	FRONT PLATE, CHAIR II	FRONTPLATTE, STUHL II	PANNEAU AVANT, FAUTEUIL II

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
12	ME-136	DÆKSEL, FODK.H.STOL II CLASSIC	COVER, RIGHTHANDED II CLASSIC	ABDECK.RECHTHÄNDER II CLASSIC	CARTER, DROIT II CLASSIC
12	MC-720	DÆKSEL, FODK.HØJRE, STOL II	COVER, RIGHTHANDED, II	ABDECK. RECHTHÄNDER II	CARTER, DROIT II
12	MC-980	DÆKSEL FODK.STOL II VENSTRE	COVER, LEFTHANDED, II	ABDECK.,LINKSHÄNDER II	CARTER, GAUCHER II
12	ME-137	DÆKSEL FODK.STOL II V. CLASSIC	COVER, LEFTHANDED, II CLASSIC	ABDECK.LINGSHÄNDER II CLASSIC	CARTER, GAUCHER II CLASSIC
13	ME-138	BASIS, RYG, STOL II CLASSIC	BASE,BACKREST, II CLASSIC	BASIS, RÜCKLEHNE, II CLASSIC	BASE, DOSSIER, II CLASSIC
13	UC-709	BASIS, RYG, STOL II	BASE, BACKREST, CHAIR II	BASIS, RÜCKLEHNE, STUHL II	BASE, DOSSIER,FAUTEUIL II
14	ME-139	BASIS, SÆDE, STOL II CLASSIC	BASE, SEAT, II CLASSIC	BASIS, SITZ, II CLASSIC	BASE, SIEGE, II CLASSIC
14	UC-710	BASIS, SÆDE, STOL II	BASE, SEAT, CHAIR II	BASIS, SITZ, STUHL II	BASE, SIEGE, FAUTEUIL II
15	ME-140	BASIS, BEN. STOL I-II CLASSIC	BASE, LEGREST, I-II CLASSIC	BASIS, BEINSTÜTZE I-II CLASSIC	BASE,APPUIE-JAMBES,I-II CLAS
15	UC-257	BASIS, BENSTØTTE, STOL I-II	BASE,LEGREST,PATIENT CHAIR II	BASIS,BEINSTÜTZE, STUHL II	BASE,APPUIE-JAMBES,FAUTEUIL
16	SD-423	PROGRAMKNAP FOR STOL I,II	PROGRAM BOTTON, CHAIR II	PROGRAMMIERKNOPF, STUHL II	BOUTON PROGRAMMES,FAUTEUIL I
17	AC-424	DÆKNAP STOLEBASIS INTEG. I,II	COVER, CHAIR BASE	ABDECKUNG, STUHLBASIS	CAPOT, BASE DU FAUTEUIL
18	AC-429	DÆKNAP, RAL 9018	COVER PLUG,9018,AC-336,AC-590	DECKSTÖPSEL,9018,AC-336,AC-590	BOUCHON, 9018, AC-336,AC-590
19	SD-428	BASIS PLADE INTEGR. II STOL	BASE BOARD, CHAIR II	BASISPLATTE, STUHL II	PLAQUE DE MONTAGE,FAUTEUIL I
20	SD-475	BØSNING FOR POLSTER STOL II	BUSHING FOR UPHOLSTERY, CH. II	BUCHSE, POLSTER, STUHL II	DOUILLE,POUR SELLERIE FAUT.I
21	SD-474	TAP FOR POLSTER STOL II	FIXING PIN M8 FOR UPHOLSTERY	ZAPPEN, POLSTER,STUHL II	AXE POUR SELLERIE,FAUTEUIL I
22	SD-476	TAP FOR POLSTER STOL II LANG	FIXING PIN FOR UPHOLSTERY, LONG	ZAPPEN, POLSTER, STUHL II, LANG	AXE POUR SELLERIE,FAUT.II,LO
23	MC-755	DÆKSEL ABC MEKANIK HØJRE	COVER, ABC MECHANICS, RIGHT	ABDECKUNG, ABC MECHANIK,RECHTS	CARTER, MECANIQUE ABC, DROIT
23	MC-756	DÆKSEL MEKANIK ABC VENSTRE	COVER, ABC MECHANICS, LEFT	ABDECKUNG, ABC MECHANIK, LINKS	CARTER, MECANIQUE ABC, GAUCH




Flex Dental A/S		1999-06-28
Patientstol,type 93: Mekaniske dele. Elektriske dele	Patient Chair,Type 93: Mechanical Parts. Electrical Parts	YA-001
Patientenstuhl,Typ 93: Mechanik. Elektronik	Fauteuil, type 93: Pièces mécaniques. Pièces électroniques	0.6.30.B

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZICHNUNG	DESCRIPTION
	1	MC-979 NAKKESTØTTE, INTEGRAL, HØJRE *	HEAD REST,INTEGRAL,RIGHT-HAND*	KOPFSTÜTZE,INTEGRAL,RECHTSH. *	TETIERE,INTEGRAL, DROITIER
	1	MC-791 NAKKESTØTTE, INTEGRAL VENSTRE*	HEADREST INTEGRAL, LEFT-HAND *	KOPFSTÜTZE INTEGRAL,LINKSHÄN.*	TETIERE,INTEGRAL, GAUCHER
1992-1998 1999 - new type	2	SD-412 PROGRAMVÆLGER STOL II	PEDAL, CHAIR II	PEDAL, STUHL II	PEDALE, FAUTEUIL II
	2	AB-272 PROGRAMVÆLGER II INCL. KABEL	PEDAL, CHAIR II, INCL. CABLE	PEDAL, STUHL II, MIT KABEL	PEDALE, FAUTEUIL II, AVEC CABL
	3	BA-136 TILLÆG, HØJRE ARMLÆN, STOL	ARM REST, RIGHT, CHAIR II	ARMLEHNE, RECHTS, STUHL II	ACCOUDOIR, DROIT, FAUTEUIL I
	3	BA-135 TILLÆG, VENSTRE ARMLÆN, STOL	LEFT ARMREST, PATIENT CHAIR/2	LINKE ARMLEHNE,PATIENTENST./2	ACCOUDOIR GAUCHE DU FAUTEUIL
	4	MC-873 GLIDER, BENSTØTTE, STOL II	SLIDING PART, FOOTREST, INT.	GLEITTEIL, FUSSLEHNE, INT.	PART DE GLISSEMENT, APPUI-PIE
	5	MC-983 GLIDESKINNE, BENSTØTTE STOL II	SLIDE GUIDEWAY OF LEGREST II	GLEITSCHIENE, FUSSLEHNE, INT.	TIGE DE GLISSEMENT, APPUI-PIE
	6	MH-242 HOVEDAFBRYDER, STOL/UNIT 220V	MAIN SWITCH, UNIT/CHAIR 220V	HAUPTSCHALTER, UNIT/STUHL 220V	INTERRUPTEUR, UNIT/FAUTEUIL22
	6	MH-745 HOVEDAFBRYDER, STOL/UNIT 100V	MAIN SWITCH, UNIT/CHAIR 100V	ONLY FOR JAPAN	ONLY FOR JAPAN
	7	SD-388 TAPNØGLE, JUST, NAKKESTØTTE	KEY, ADJUSTMENT HEADREST	SCHLÜSSEL,REGULIER. KOPFSTÜTZE	CLÉF, REGLAGE TÊTIÈRE
	11	SC-842 LØFTEMOTOR, FLEX CHAIR	* LIFTING MOTOR, FLEX CHAIR	* HUBMOTOR, FLEX STUHL	* MOTEUR DE LEVAGE, FAUT.FLEX
	11	SD-432 LØFTEMOTOR, INTEG/STOL II 100V	LIFTING MOTOR CHAIR II 100V	NUR FÜR JAPAN	ONLY FOR JAPAN
	12	SC-843 RYGMOTOR, FLEX CHAIR	* BACKREST MOTOR, FLEX CHAIR	* RÜCKLEHNENMOTOR, FLEX STUHL	* MOTEUR DU DOSSIER, FAUT.FLEX
	12	SD-433 RYGMOTOR, INTEG/STOL II 100V	BACKREST MOTOR, CHAIR II,100V	NUR FÜR JAPAN	ONLY FOR JAPAN
	13	SD-380 UDLØSERKNAP, NAKKESTØTTE	RELEASE NUT, HEAD REST	AUSLÖSEMUTTER, KOPFSTÜTZE	ECROU DE DECLENCHEMENT,TETIE
	14	SD-484 KAMSKIVE FOR MC846 STOL II	ADJUSTING CAM OF CHAIR II BACK	SCHBIBE,NOCKENW.MC-846,STUHLII	CAME POUR MC-846 FAUTEUIL II
	15	SD-416 UDLØSERARM, NAKKESTØTTE	RELEASE ARM, HEAD REST	AUSLÖSE ARM, KOPFSTÜTZE	DECLENCHEMENT BRAS, TETIERE
	16	SD-417 TRYKSTANG FOR LÅSEPAL STOLE	PRESSURE BAR, LOCKING PAWL	DRUCKSTANGE, SPERRHAKEN	BATON DE PRESSION, CLIQUET
	17	SD-418 FJEDER FOR LÅSEPAL STOL	SPRING, LOCKING PAWL,INTEGRAL	FEDER, SPERRHAKEN, INTEGRAL	RESSORT,CLIQUET FERMETURE
	18	SD-419 RØRSTIFT 3MM FOR LÅSEPAL STOLE	PIN 3 MM, LOCKING PAWL	STIPT 3 MM, SPERRHAKEN	GOUPILLE 3MM,CLIQUET FERMETU
	19	SD-420 LÅSEPAL NAKKESTØTTE STOL	LOCKING PAWL, INTEGRAL CHAIR	SPERRHAKEN, INTEGRAL STUHL	CLIQUET FERMETURE,FAUTEUIL
	20	SD-421 UNBRAKOSKRUB FOR NAKKEST. STOL	SCREW FOR HEAD REST	SCHRAUBE FÜR KOPFSTÜTZE	VIS, TETIERE
	21	SD-353 SPINDELMØTRIK, PLAST INT II	SPINDLE NUT, PLASTIC INT. II	SPINDELSCHR.MUTTER,PLAST,IN II	ECROU DE TRANSLATION, INT. I
	22	SD-480 PLASTKNAP FOR SD-423	PLASTIC BUTTON FOR SD-423	PLASTKNOPF FÜR SD-423	BOUTON PLASTIQUE,POUR SD-423
	23	ME-059 KNOP, PROGRAMVÆLGER STOL II	KNOB, PROGR.SELECTOR CHAIR II	KNOPF,PROGRAMMWÄHLER,STUHL II	BOUTON,SELECTEUR PROGR,FAUT
	24	AD-427 PRINT, STOL II	* PR.CIRCUIT,CHAIR CON.,INT./2 *	* PLATINE,CHAIR CONTROL,INT./2 *	* CARTE, CHAIR CONTROL,INT./2
	24	AD-436 PRINT,CHAIR CONTROL II 100V	PCB, CHAIR CONTROL II, 100V	NUR FÜR JAPAN	ONLY FOR JAPAN
	25	AD-428 PRINT,CHAIR CONNECTOR STOL II*	PR.CIRCUIT,CHAIR CONN.INT./2 *	* PLATINE,CHAIR CONN. INT./2 *	* CARTE,CHAIR CONNECTOR INT/2
	26	UC-480 STØFFILTER 6-A EVOK EBF6-6	NOISE FILTER	STÖRFILTER	FILTRE ANTIPARASITES
	36	UA-131 SKRUE, ARMLÆN INTEGRAL RUSTFRI	SCREW, ARM REST INTEGRAL	SCHRAUBE, ARMLEHNE INTEGRAL	VIS, ACCOUDOIR INTEGRAL
	37	MC-682 DÆKKAPPE, PROGRAMVÆLGER STOL	FOOTCONTROL CAP CHAIR II	ABDECKUNG, FUSSANLASSER, STUHL	CARTER, PEDALE, FAUTEUIL
	38	AB-238 E-PROM PROGRAMMERET, AD-427/436	E-PROM PROGRAMMDD, AD-427	E-PROM PROGRAMMIERT, AD-427	E-PROM PROGRAMME, AD-427
	39	MH-428 SIKRING TRÆG 500 MA	FUSE T 500 MA	SICHERUNG T 500 MA	FUSIBLE T 500 MA
	40	ML-069 SIKRING TRÆG 6,3 A	FUSE T 6.3 A	SICHERUNG T 6,3 A	FUSIBLE T 6,3 A
	41	SD-500 KIT TRÆKARM STOL II	KIT, ARMS.SCREWS, CHAIR II	KIT, ZUGARM+SCHRAUBEN,STUHL II	KIT,BRAS+VIS, FAUTEUIL II
	42	MC-724 BÆRESKINNE, PARALLELARMD,HØJRE	SUPPORTING BAR, ARM COVER,RIGHT	TRAGSCHIENE, ARMDECKPL,RECHTS	TIGE SUPPORT,CACHE BRAS,DROI
	43	MC-725 BÆRESKINNE, PARALLELARMD,VENST	SUPPORTING BAR, ARM COVER,LEFT	TRAGSCHIENE, ARMDECKEL,LINKS	TIGE SUPPORT,CACHE BRAS,GAUC
	44	SC-886 PLASTGREB, INTEG.NAKKESTØT.I.II	LOCKING HANDLE FOR HEADREST	VERSCHLUSSGRIFF FÜR KOPFSTÜTZE	POIGNEE FERMETURE POUR TETIE
	45	SD-423 PROGRAMKNAP FOR STOL I,II	PROGRAM BOTTON, CHAIR II	PROGRAMMIERKNOPF, STUHL II	BOUTON PROGRAMMES,FAUTEUIL I
	47	SD-490 VÆRKTØJ F/SD-412 PROGRAMVÆLGER	TOOL FOR SD-412, FOOTCONTROL	WERKZEUG,FÜR SD-412,FUSSPEDAL	OUTIL POUR SD-412, PEDAL
	48	SD-440 SIKKERHEDS SWITCHES STOL II	SAFETY SWITCH CHAIR II	SICHERHEITSSCHALTER STUHL II	INTERRUPTEUR DE SECURITE,F.I
	49	SD-450 SWITCH, RYGLÆN, INCL. KABEL	MICROSWITCH WITH CABLE	MIKROSCHALTER MIT KABELSATZ	INTERRUPTEUR AVEC CABLAGE
	50	MC-822 TRÆKARM, INTEGRALSTOL II	PAINTED FOOTREST MOV. LEVER II	PLEUELSTANGE, LACKIERT, II	BARRE DE COUPLAGE, II
	51	MC-736 SÆDELØFTBAKSEL INTEG.STOL II	LEVER, SEAT, CHAIR II	Achse SITZFLÄCHE, STUHL II	BARRE, SIEGE, FAUTEUIL II
	52	SC-849 MICROAFBRYDER, FLEX CHAIR	MICROSWITCH, FLEX CHAIR	MIKROSCHALTER, FLEX STUHL	MICRO-INTERRUPT.FAUTEUIL FLE
	53	SD-472 FJEDER FOR PROGAMVÆLGER TY II	FOOTCONTROL SPRING, INT. II	FEDER, PROGRAMMWÄHLSCHALTER	RESSORT,PEDALE, FAUTEUIL II
	54	MC-738 AKSEL Ø8X40, INTEGRALSTOL II	INTEGRAL CHAIR II PIN Ø8X40	ACHSE,Ø8X40,INT.STUHL II	AXE, Ø8X40,FAUTEUIL INT. II
	55	SD-473 RULLE FOR INTEGRAL STOL II	INTEGRAL CHAIR II ROLL	ROLLE, INT. STUHL II	ROULEMENT POUR FAUTEUIL II
	56	MC-846 AKSEL, CENTERBOX, SVBJST,	FOOTREST LEVER ROD/BAR, INT.II	ACHSE, INT.II	ARBRE, INT.II
	57	MC-737 LØFTEARM, BENSTØTTE, STOL II	ROD/BAR, FOOTREST LEVER, IN.II	HEBELARM, FUSS-STUETZE,STUHL II	BRAS DE LEVAGE, APPUI-PIEDS,
	SD-454	SPINDELMØTRIK, INTEGRAL II	COMPLETE BACKREST LEAKAGE II	SPINDELMUTTER STUHL II	ECROU DE TRANSLATION II
	SD-504	LEJEBOLT STOL II	SCREW TCCE 8 - 6X20 UMBRACO	LAGERBOLZEN TCCE 8 - 6X20	VIS AXIALE TCCE 8 - 6X20
	UC-437	Fiberskive	Fiber disc	Fiberscheibe	
	SD-413	Kondensator, Rygmotor	Condensator, bagrest motor	Kondensator, Rückenmotor	
	SD-431	Kondensator, Løftemotor	Condensator, lifting motor	Kondensator, Hubmotor	
	24 MC-898	Isoleringsplade print-427	Isolation plate PCB AD-427	Isolierungsplatte Platine AD-427	
	24 ME-047	Plastafskærmning print AD-427	Plastic cover PCB AD-427	Plastabdeckung Platine AD-427	




No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	PH-099	INSTRUMENT, UNIVERSAL	UNIVERSAL INSTRUMENT	UNIVERSALINSTRUMENT	INSERT UNIVERSEL
2	PH-105	INSTRUMENT, VENSTRE	INSTRUMENT, LEFT	INSTRUMENT, LINKS	INSERT, COURBÉ À GAUCHE
3	PH-106	INSTRUMENT, LIGE	INSTRUMENT, STRAIGHT	INSTRUMENT, GERADE	INSERT DROIT
4	PH-107	INSTRUMENT, HØJRE	INSTRUMENT, RIGHT	INSTRUMENT, RECHTS	INSERT, COURBÉ À DROITE
5	PH-117	INSTRUMENT, AMALGAM, HØJRE	INSTRUMENT, AMALGAM, RIGHT	INSTRUMENT, AMALGAM, RECHT	INSERT, AMALGAME, COURBÉ À DR
6	PH-118	INSTRUMENT, AMALGAM, VENSTRE	INSTRUMENT, AMALGAM, LEFT	INSTRUMENT, AMALGAM, LINKS	INSERT, AMALGAME, COURBÉ À G
7	PH-119	INSTRUMENT, PERIO, STANDARD	INSTRUMENT, PERIODONTAL	INSTRUMENT, PARODONTAL	INSERT PAROD., COURBÉ
9	SA-048	KRONEFJERNER	CROWN REMOVER	KRONENLÖSER	INSERT DESCHELLEUR
10	PH-123	INSTRUMENT, THIN LINE, LIGE	INSTRUMENT, THIN LINE STRAIGHT	INSTRUMENT, SCHLANK, GERADE	INSERT DROIT EXTRA FIN
11	PH-124	INSTRUMENT, THIN LINE, HØJRE	INSTRUMENT, THIN LINE, RIGHT	INSTRUMENT, SCHLANK, RECHTS	INSERT XTRA FIN COURBÉ A DRO
12	PH-125	INSTRUMENT, THIN LINE, VENSTRE	INSTRUMENT, THIN LINE, LEFT	INSTRUMENT, SCHLANK, LINKS	INSERT XTRA FIN COURBÉ A GAU
13	PH-126	INSTRUMENT, CEM	INSTRUMENT, CEM	ZEMENTIERUNGsinstrument	INSERT POUR LA SCHELLEMENT
14	PH-113	INSTRUMENT, ENDO, UDEN FILE	INSTRUMENT, ENDO, WITHOUT FILES	INSTRUMENT, ENDO, OHNE FEILEN	INSERT ENDO, SANS LIMES
15	PH-115	ENDOFIL, PH-113 (10 X 3 STK.)	ENDO-FILES, PH-113 (10 X 3)	ENDO-FEILE, PH-113 (10 X 3)	LIMES ENDO, PH-113 (10 X 3)
16	PH-116	INSTRUMENT SÆT, ENDO, INCL. FILE	INSTRUMENT, SET, ENDO WITH FILES	INSTR. SÆT, ENDO, MIT FEILEN	JEU D'INSERTS ENDO, AVEC LIME
17	PH-104	FERRITSTAV	FERRITE ROD	FERRITKERN	BÂTON EN FERRITE
18	PH-108	TANG, FERRIT SKIFT	FERRITE WRENCH	FERRITZANGE	PINCE POUR CHANGER LA FERRIT
19	PH-109	VÆRKTØJ, FERRIT	TOOL, FERRITE	WERKZEUG, FERRIT	TIGE POUR ECHANGE DES FERRIT
20	PH-110	HOLDER, INSTRUMENTER, TEFLON	HOLDER, INSTRUMENTS, TEFLON	INSTRUMENTENHALTER, TEFLON	CARQUOIS, INSERTS, TÉFLON
21	PH-122	STERIL.RØR, INSTR. TANDR. TY.O	STERILIZATION CAP TYPE O	STERILISATIONSKAPPE TYP O	PROTEGE-FERRITE TYPE O
22	SC-714	O-RING, FLEX/ODONTOSON, LILLE	O-RING, FLEX/ODONTOSON, SMALL	O-RING, FLEX/ODONTOSON, KLEIN	JOINT TORIQ.FLEX/ODONT., PETI
23	SC-715	O-RING, FLEX/ODONTOSON, STOR	O-RING, FLEX/ODONTOSON, BIG	O-RING, FLEX/ODONTOSON, GROSS	JOINT TORIQ.FLEX/ODONT., GRAN
24	AD-443	STYRING, TANDRENSER TYPE O *	CONTROL SCALER TYPE O	* STEUERUNG, ZEG TYP O	* GENERATEUR DETARTREUR TY O
25	AE-263	CODE CHIP, TANDRENSER TYPE O	CODE CHIP, SCALER TYPE O	CODE CHIP, ZEG TYP O	BLOC-CODE, DETARTREUR TYPE O
26	PH-127	LÆBEBESKYTTER, TANDR. TYPE O	LIP PROTECTION, SCALER TYPE O	LIPPENSCHUTZ, ZEG TYP O	PROTECTEUR LEVRE, DETARTREUR
27	UC-759	10 STK PLASTHÆTTE, CEM INSTR.	10 PLASTIC CAPS, CEM. INSTR.	10 PLASTIKKAPPEN, ZEMENT.INST.	10 CAPUCHON PLAST. INS. SCEL
28	ML-564	SIKRING FLINK 0.8A	FUSE 0.8 A	SICHERUNG 0,8 A	FUSIBLE 0,8 A
29	HT-115	HÅNDSTYKKE, TANDRENSER TYPE O*	HANDPIECE, SCALER TYPE O	* HANDSTÜCK, ZAHNREINIGER TY O *	* PIECE A MAIN DETARTREUR TY.O
30	SC-811	VÆRKTØJ, ENDOINSTRUMENT	TOOL, ENDO-INSTRUMENT	WERKZEUG, ENDO-INSTRUMENT	OUTIL, INSERT ENDO

29 208206 Gevinding for håndstykke Connector for hand piece Gewinde für Handstück



No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZICHTUNG	DESCRIPTION
1	AE-259	BLANDETANK, VAC ULTRA	MIXING TANK, VAC ULTRA	MISCHBEHÄLTER, VAC ULTRA	RECIPIENT A MELANGER, VACULTR
2	AE-260	TRAGT KOMPLET, VAC ULTRA	COMPLETE FUNNEL, VAC ULTRA	KOMPLETTER TRICHTER, VAC ULTRA	ENTONNOIR COMPLET, VAC ULTRA
3	AE-262	GENNEMFØRING M. DYSE OG RØR	FEED-THROUGH WITH NOZZLE+TUBE	DURCHFÜHRUNG MIT DÜSE UND ROHR	CARTER PASSE-TUYAUX + RACCOR
4	MB-280	FITTINGS, BANJO, 6/4 X M 5	FITTINGS, BANJO, 6/4 X M 5	FITTINGS, BANJO, 6/4 X M 5	RACCORD, BANJO, 6/4 X M 5
5	MC-231	OMLØBER Ø 43	UNION Ø 43	ÜBERWURFMUTTER Ø 43	ECROU Ø 43
6	UC-398	PAKNING, OMLØBER	GASKET, UNION	DICHTUNG, ÜBERWURFMUTTER	JOINT D'ECROU
7	AC-143	KABEL, HYGIFLEX VAC PCB	CABLE, HYGIFLEX VAC	KABEL, HYGIFLEX VAC	CABLE, HYGIFLEX VAC
8	AE-261	NIVEAUAFBRYDER, VAC ULTRA	LEVEL SENSOR, VAC ULTRA	FÜLLSTANDSMESSE, VAC ULTRA	CAPTEUR DE NIVEAU
9	AD-416	PRINT, HYGIFLEX VAC	PRINTED CIRCUIT "HYGIFLEX VAC"	PLATINE "HYGIFLEX VAC"	CARTE "HYGIFLEX VAC"
10	ML-183	MAGNETVENTIL, STANDER, FLANGE	SOLENOID VALVE, UNIT STAND	MAGNETVENTIL, STANDKÖRPER	ELECTROVANNE, BORNE
11	SD-292	ANKER, VAND, ML-183+AC-524	PLUNGER, WATER, ML-183+AC-524	PLUNGER, WASSER, ML-183+AC-524	PLONGEUR, EAU, ML-183+AC-524
12	SG-233	PAKNING, FOR ANKER SÆDS SD-292	GASKET, SEAT, VALVE SD-292	DICHTUNG, SITZ PLUNGER SD-292	JOINT, SIEGE, VANNE SD-292
13	AC-521	BESLAG, HYGIFLEXVAC	* FITTINGS, HYGIFLEX VAC, MEC.	* FITTINGS, HYGIFLEX VAC, MECH.	* GARNITURE, HYGIFLEX VAC, MEC.
14	FA-752	O-RING Ø 9.75 X 1.78 NBR 70	O-RING Ø 9.75 X 1.78 NBR 70	O-RING Ø 9.75 X 1.78 NBR 70	JOINT TORIQUE Ø9.75 X 1.78
15	MR-155	INDLØB TRAGT	HOSE FOR FUNNEL	SCHLAUCH, TRICHTER	



No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	SD-464	SPRAYDÆKSEL, TURBINE TYPE B	SPRAY COVER, TURBINE TYPE B	SPRAYPLATTE, TURBINE TYP B	CARTER TROUE, TURBINE TYPE B
2	SD-465	O-RING, SPRAYDÆKSEL, TYPE B	O-RING, SPRAY COVER, TYPE B	O-RING, SPRAYPLATTE, TYP B	JOINT TORIQUE, CARTER TROUE,
3	SD-466	MØTRIK, SPRAYDÆKSEL, TYPE B	NUT, SPRAY COVER, TYPE B	MUTTER, SPRAYPLATTE, TYP B	ECROU, CARTER TROUE, TYPE B
4	SA-983	VÆRKTØJ, SPRAYDÆKSELMØTRIK, B	TOOL, SPRAY COVER NUT, TYPE B	WERKZEUG, MUTTER SPRAYPLATTE, B	OUTIL, ECROU, CARTER TROUE,
5	SA-984	VÆRKTØJ, SPRAYDÆKSEL, TYPE B	TOOL, SPRAY COVER, TURBINE T.B	WERKZEUG, SPRAYPLATTE, TYP B	OUTIL, CARTER TROUE, TYPE B
6	SD-469	LYSLEDER, TURBINE TYPE B	OPTIC FIBRES, TURBINE TYPE B	LICHTLEITER, TURBINE TYP B	TIGE OPTIQUE, TURBINE TYPE B
7	SD-467	O-RING, SPRAYVAND/LUFT, TYPE B	O-RING, SPRAY WATER,AIR,TYPE B	O-RING, SPRAYWASSER/-LUFT, B	JOINT TORIQUE,BAU,AIR, TYPE
8	SD-468	O-RING, DRIVLUFT, TYPE B	O-RING, DRIVING AIR, TYPE B	O-RING, TREIBLUFT, TYP B	JOINT TOR.,L'AIR DE FONCT.,
9	HE-005	PÆRE, LUX MOTOR/TURBINE	BULB, LUX MOTOR/TURBINE	BIRNE, LUX MOTOR/TURBINE	AMPOULE, MOTEUR LUX/TURBINE
10	UC-736	BAGSTYKKE, TURBINE, TYPE B+P	QF PART, TURBINE, TYPE B+P	QF-SEKTION, TURBINE, TYP B+P	PARTIE ARRIERE,TURBINE, B+P
11	SA-985	KIT, SPRAYDÆKSEL, TURBINE B	KIT, SPRAY COVER, TURBINE B	KIT, SPRAYPLATTE, TURBINE B	KIT, CARTER TROUE, TURBINE B
12	SA-051	FEDTSPRØJTE, FLEXLUX TURBINE	GREASE GUN, FLEXLUX TURBINE	FETTSPRITZE, FLEXLUX TURBINE	GRAISSEUR, TURBINE FLEXLUX
13	SA-052	FEDT, 1 PATRON, FLEX TURBINE	GREASE, 1 CARTRIDGE, TURBINE	FETT, 1 PATRONE, FLEX TURBINE	GRAISSE, 1 CARTOUCHE, TURBIN
13	SD-318	FEDT 6 PATRON, FLEX TURBINE	GREASE, 6 CARTRIDGES, TURBINE	FETT, 6 PATRONEN, FLEX TURBINE	GRAISSE, 6 CARTOUCHES,TURBIN
14	SC-973	RENSENÅL, 1-HUL, BIEN AIR GUL	NOZZLE CLEANER,1-HOLE, YELLOW	REINIGUNGSNADEL,1-LOCH, GELB	AIGUILLE,1-VOIE,JAUNE,BIENAI
15	SD-429	FJEDER, TURBINE, TYPE P	SPRING, INTEGRAL TURBINE	FEDER, INTEGRAL TURBINE	RESSORT, TURBINE INTEGRAL
16	SC-740	O-RING, UNIFIX/QF, STOR	O-RING, UNIFIX COUPLING, BIG	O-RING, UNIFIX KUPPLUNG, ,GROSS	JOINT TORIQUE,RAC.UNIFIX,GRA
17	HT-070	QF KOBLING FOR 4+4 SLANGE	QUICKFLEX COUPLING, 4+4 TUBE	QUICKFLEX KUPPLUNG, 4+4 SCHL.	RACCORD QUICKFLEX, CORDON 4+

10 UC-721 Kontaktsæt for UC-736

Contact set for UC-736

Kontaktsatz für UC-736




Flex Dental A/S
Håndinstrumentbord, Holder
Handinstrumententische, Adapter

Hand instrument table, Adapter
Plateaux, Adaptateur

1999-06-28
YA-001
0.6.34.B

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	AC-066	ENKBELT BAKKE, INTEGRAL	HAND INSTR. TABLE, 1 TRAY, INT.	HANDINSTR.TISCH, 1 TRAY, INT.	TABLETTE, 1 BOÎTE NORMAL., IN
2	AC-067	DOBBLT BAKKE, INTEGRAL	HAND INSTR. TABLE, 2 TRAYS, INT.	HANDINSTR.TISCH, 2 TRAYS, INT.	TABLETTE, 2 BOÎTES NORMAL., I
3	AC-380	AFTAGELIG BAKKEHOLDER INTEGRAL	TRAY HOLDER, REMOVABLE	TABLETHALTER, ABNEHMBAR	PORTE-TABLETTE AMOVIBLE
4	AC-262	SVINGARM, BAKKE BRGOFLEX (PUR)*	SWIVEL ARM, TABLE BRGOFLEX (PUR)*	SCHWENKARM, TISCH BRGOFLEX (PUR)*	BRAS TOURNANT, TABL. BRGO (PUR)
5	AC-644	AFTAGELIG BAKKEH. ALU-B., INT	TRAY HOLDER, REMOVABLE, INT.	TABLETHALTER, ABNEHMBAR, INT.	PORTE-TABLETTE AMOVIBLE, INT
6	AC-640	SVINGARM, BAKKE BRGOFLEX (ALU)*	SWIVEL ARM, TABLE BRGOFLEX (ALU)*	SCHWENKARM, TISCH BRGOFLEX (ALU)*	BRAS TOURNANT, TABL. BRGO (ALU)
7	FA-069	BAKKE FOR 2 KASSETTER	HAND INSTR. TABLE, 2 NORM TRAYS	HANDINSTR.TISCH, 2 NORMTRAYS	TABLETTE POUR 2 BOÎTES NORMA
8	FA-070	ALU-BAKKE, 1 KASSETTE	HAND INSTR. TABLE, 1 NORM TRAY	HANDINSTR.TISCH, 1 NORMTRAY	TABLETTE POUR 1 BOÎTE NORMAL
9	UA-444	SKRUE, GKN M 4 X 10, BLANK	SCREW GKN M 4 X 10, PLAIN	SCHRAUBE GKN M 4 X 10, BLANK	VIS GKN M 4 X 10, BRUTE
10	AE-936	OVERDEL BAKKEHOLDER INTEGRAL	TOP PART TRAY HOLDER INTEGRAL	OBERTEIL TABLETHALTER INTEGR.	PARTIE SUP. PORTE-TABLETTE
11	AE-937	UNDERDEL BAKKEHOLDER ALU	BUTTOM PART TRAY HOLDER ALU	UNTERTEIL TABLETHALTER ALU	PARTIE INF. PORTE-TABLETTE A
12	MC-469	FINGERSKRUE FOR BAKKEOPHÆNG	THUMBSCREW, TRAY SUSPENSION	PLÜGELSCHRAUBE, TABLETTAUFHÄNG.	VIS TÊTE MOLETEE, FIXAT.PLATE
13	UB-056	BØSNING + MØTRIK, FOR SVINGARM	BUSH + NUT,	BUCHSE +SCHRAUBENMUTTER,	DOUILLE + ECROU,
14	AE-938	UNDERDEL BAKKEHOLDER PUR	BOTTOM PART TRAY HOLDER PUR	UNTERTEIL TABLETHALTER PUR	PARTIE INF. PORTE-TABLETTE P
15	AE-921	SVINGARM, BAKKE FLEXFORM (ALU)*	SWIVEL ARM, TABLE FLEXFORM ALU*	SCHWENKARM, TISCH FLEXFORM ALU*	BRAS TOURN. TABL. FLEXFORM ALU



Sugeslangeholder, fontæne. Lynkoblinger

Suction hose holder, cuspidor. Quick Connections


YA-001

Saugschlauchhalter, Fontæne. Schnellkupplungen

Porte-tuyaux d'aspiration, crachoir. Raccords rapides


0.6.35.B

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	AB-250	KABEL, SUGESL.HOLDER PÅ FONT.	CABLE, SUCTION AT THE CUSPIDOR	KABEL, ABSAUGER AN FONTÅNE	CABLE,ASPIR. SUR LE CRACHOIR
2	AD-403	PRINT, SUGESLANGEHOLDER	* PRINTED CIRC."SUCTION HOLDER"	* PLATINE "SUCTION HOLDER"	* CARTE "SUCTION HOLDER"
3	ME-037	DÆKPLADE, FONTÆNE M. SUG	PLATE, SUCTION AT THE CUSPIDOR	PLATTE, ABSAUGER AN FONTÅNE	PLAQUE,ASPIR. SUR LE CRACHOI
4	MB-038	REGULERINGSKØB, VAND FORNIKLET	REGULATING COCK, WATER	REGLERHAHN, WASSER	ROBINET DE RÉGLAGE, EAU
5	SC-252	SPINDEL,REGULERINGSKØB MB-038	SPINDLE,REGULATING COCK MB-038	SPINDEL, REGLERHAHN MB-038	ARBRE, ROBINET DE RÉGL. MB-0
6	SC-271	KNAP, VANDREGULERINGSKØB	PLUG	STÖPSEL	BOUCHON
7	HR-006	STIK, HAN, KØBLING 3/16"SLANGE	PLUG, MALE,CONNECT. 3/16" TUBE	STECKER,MÄNNL.,KUPPLUNG 3/16"S	PRISE MÂLE,RACCORD,3/16"CORD
8	HR-003	LYNKØBLING, SORT	QUICK CONNECTION, BLACK	SCHNELLKUPPLUNG, SCHWARZ	RACCORD RAPIDE, NOIR
8	HR-004	LYNKØBLING, GUL	QUICK CONNECTION, YELLOW	SCHNELLKUPPLUNG, GELB	RACCORD RAPIDE, JAUNE
8	HR-005	LYNKØBLING, RØD	QUICK CONNECTION, RED	SCHNELLKUPPLUNG, ROT	RACCORD RAPIDE, ROUGE
8	HR-007	LYNKØBLING, GRØN	QUICK CONNECTION, GREEN	SCHNELLKUPPLUNG, GRÜN	RACCORD RAPIDE, VERT
8	HR-041	LYNKØBLING, BLÅ	QUICK CONNECTION, BLUE	SCHNELLKUPPLUNG, BLAU	RACCORD RAPIDE, BLEU
9	MC-878	TOP,SUG,STOL	TOP, SUCTION, CHAIR	OBerteil, ABSAUGER, STUHL	PARTIE SUP.,ASPIR.,FAUTEUIL
10	MC-186	MODTAGER,HÅNDTAG, INSTRUMENTBRO	HOLDER, HANDLE, INSTR. BRIDGE	EINSTECKROHR,GRIF,INST.TRÄGER	TUYAU SUPP.,POIGNEE,PORTE-IN
11	MC-111	HÅNDTAG, BRO/LAMPE, 1 STK.	HANDLE, INSTRUMENT BRIDGE	HANDGRIF, INSTRUMENTENTRÄGER	POIGNEE, PORTE-INSTRUMENT
11	SD-399	HÅNDTAG, BRO/LAMPE, 4 STK.	HANDLES (4) BRIDGE/LAMP	HANDGRIFFE (4) TRÄGER/LAMPE	POIGNEES(4) PORTE-INSTR./ECL
11	AC-371	HÅNDTAG, BRO/LAMPE, 24 STK.	HANDLES (24), BRIDGE/LAMP	HANDGRIFFE (24), TRÄGER/LAMPE	POIGNEES(24)PORTE-INSTR./ECL
12	MC-683	HOLDER, SUGENIPLER, STOL	HOLDER, SUCTION NIPPLE, CHAIR	HALTER,SAUGSCHLAUCHANS.,STUHL	SUPPORT,EMBOUT ASPIR.FAUTEUI
13	MC-943	BESLAG, SLANGEH.,STOL,NCS 4502	FITTINGS,SUCTION,CHAIR,9018	BESCHLAG, SAUG, STUHL, 9018	FIXATION,ASPIR.,FAUTEUIL,901
14	UC-347	SILICONESLANGE Ø 21 X 16,8	SILICONE TUBE Ø 21 X 16,8	SILIKONSCHLAUCH Ø 21 X 16,8	TUYAU SILICONE Ø 21 X 16,8
15	MY-003	SPIRAL 40 CM, MESSING	SPIRAL 40 CM, BRASS	SPIRALE 40 CM, MESSING	SPIRALE 40 CM, LAITON
16	AB-162	KABEL,UNIT-CHAIR TIL SUG. II	CABLE,UNIT-CHAIR FOR SUCT.II	KABEL,UNIT-STUHL FÜR SAUG.II	CÂBLE,UNIT-FAUTEUIL P.ASPIR.



No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	SD-510	KAPPE, SPRØJTE S, MINILIGHT	COVER, SYRINGE S	DECKKAPPE, SPRITZE S	CARTER, SERINGUE S
2	SD-511	BALANCEOPHÆNG, SPRØJTE S 3F	* BALANCESUSPENSION S 3-WAY SYR.	BALANCEAUFHÄNG S 3-FUNKT. SPR.	FOUET, SERINGUE S 3-FONCTIONS
2	SD-512	BALANCEOPHÆNG, SPRØJTE S 6F	* BALANCESUSPENSION S 6-WAY SYR.	BALANCEAUFHÄNG S 6-FUNKT. SPR.	FOUET, SERINGUE S 6-FONCTIONS
4	SD-513	VENTILBLOK, SPRØJTE S 3F OG 6F	VALVE BLOCK, SYRINGE S 3F & 6F	VENTILBLOCK, SPRITZE S 3F & 6F	BLOCK VANNE, SERINGUE S 3F & 6F
5	SD-514	MELLEMSTYKKE, SPRØJTE S 6F	INTERMEDIATE PART, SYRINGE S 6F	ZWISCHENSTÜCK, SPRITZE S 6F	PRISE INTERMED., SERINGUE S 6
6	SD-516	OMLØBER, SPIDS, SPRØJTE S	UNION, TIP, SYRINGE S	ÜBERWURFM., SPITZE, SPRITZE S	ECROU, EMBOUT, SERINGUE S
7	SD-517	MICROSWITCHE VAND SPRØJTE S 6F	MICROSWITCH WATER SYRINGE S 6F	MIKROSCHALTER WASSER SPRITZE S	MICRO-INTERRUPT. EAU SERINGUE
8	SD-518	MICROSWITCHE MAIN SPRØJTE S 6F	MICROSWITCH MAIN SYRINGE S 6F	MIKROSCHALTER HAUPT SPRITZE S	MICRO-INTERRUPT. PRINC. SERINGUE
9	SD-519	FORENDE, SPRØJTE S 3F OG 6F	FRONT PART, SYRINGE S 3F & 6F	VORDERENDE, SPRITZE S 3F & 6F	PIECE AVANT, SERINGUE S 3F & 6F
10	SD-520	KNAPPER, SPRØJTE S 3F OG 6F	KNOBS, SYRINGE S 3F & 6F	KNÖPFE, SPRITZE S 3F & 6F	BOUTONS, SERINGUE S 3F & 6F
11	SD-521	FORBINDER, SPRØJTE S 3F OG 6F	CONNECTION, SYRINGE S 3F & 6F	VERBINDUNG, SPRITZE S 3F & 6F	RACCORD, SERINGUE S 3F & 6F
12	SD-522	O-RING, Ø4, SPRØJTE S 3F OG 6F	O-RING, Ø4, SYRINGE S 3F & 6F	O-RING, Ø4, SPRITZE S 3F & 6F	JOINT TORIQUE, Ø4, SERINGUE
13	SD-523	STUDS, SPRØJTE S 3F OG 6F	CONNECTOR, SYRINGE S 3F & 6F	STÜTZEN, SPRITZE S 3F & 6F	RACCORD, SERINGUE S 3F & 6F
14	SD-524	O-RING, Ø5, SPRØJTE S 3F OG 6F	O-RING, Ø5, SYRINGE S 3F & 6F	O-RING, Ø5, SPRITZE S 3F & 6F	JOINT TORIQUE, Ø5, SERINGUE
15	SD-525	FJEDER, SPRØJTE S 3F OG 6F	SPRING, SYRINGE S 3F & 6F	FEDER, SPRITZE S 3F & 6F	RESSORT, SERINGUE S 3F & 6F
16	SD-526	STIFT VAND, SPRØJTE S 3F OG 6F	PIN WATER, SYRINGE S 3F & 6F	STIFT WASSER, SPRITZE S 3F & 6F	BROCHE EAU, SERINGUE S 3F & 6F
17	SD-527	STIFT LUFT, SPRØJTE S 3F OG 6F	PIN AIR, SYRINGE S 3F & 6F	STIFT LUFT, SPRITZE S 3F & 6F	BROCHE AIR, SERINGUE S 3F & 6F
18	SD-528	O-RING, Ø3, SPRØJTE S 3F OG 6F	O-RING, Ø3, SYRINGE S 3F & 6F	O-RING, Ø3, SPRITZE S 3F & 6F	JOINT TORIQUE, Ø3, SERINGUE
19	SD-529	PLAST FOR SWITCH, SPRØJTE S 6F	PLASTIC F. SWITCH, SYRINGE S 6F	KUNSTSTOFF F. SCHALTER, SPRITZE	PLAST. P. INTERRUPTEUR, SERINGUE
20	SD-214	SPIDS, INTEGRAL SPRØJTE	TIP, SYRINGE	SPITZE, SPRITZE	EMBOUT, SERINGUE
21	BA-341	ADAPTOR, HYGIFLEX, SPRØJTE S	ADAPTER, HF THERMO SYRINGE S	ADAPTER, HF THERMO SPRITZE S	ADAPTATEUR, HF THERMO SERINGUE

6 SA-988 O-Ring for SD-516 O-Ring for SD-516 O-Ring für SD-516




Flex Dental A/S
Vandvarmer, vandglasfylder
Wasserwärmer, Glasfüller

Water heater, cup filler
Rècharffeur d'eau, jet de verre


1999-06-28
YA-001
0.6.37.B

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	AH-036	PRINT, "FLEX BOILER CONTROL"	PR. CIRC. "FLEX BOILER CONTROL"	PLATINE "FLEX BOILER CONTROL"	CARTE "FLEX BOILER CONTROL"
2	AC-836	VARMTVANDSBEHOLDER, LILLE	BOILER	WARMWASSERBOILER	CHAUFFER-EAU
3	ML-017	THERMOSTAT MED RESET	THERMOSTAT WITH RESET BUTTON	THERMOSTAT MIT RÜCKSCHALTUNG	THERMOSTAT+BOUTON DE "RESET"
4	ML-018	NTC THERMISTOR/TEMPERATURFØLER	NTC TEMPERATURE DETECTOR	NTC TEMPERATURFÜHLER	CAPTEUR DE TEMPERATURE NTC
5	FA-238	VARMEFOLIE	HEAT CONDUCTING FOIL	WÄRMELEITENDE FOLIE	FEUILLE CONDUCTRICE DE CHALE
6	AB-398	TEMPERATURFØLER	TEMPERATURE DETECTOR	TEMPERATURFÜHLER	SONDE DE TEMPÉRATURE

BC



No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	AE-047	VENTILBLOK BRO 2 OPHÆNG BC	VALVE BLOCK, 2 SUSPENSIONS, BC	VENTILBLOK, 2 AUFHÄNGUNGEN,BC	BLOC VANNE, 2 FOUETS, BC
2	AE-046	VENTILBLOK BRO 3 OPHÆNG BC	VALVE BLOCK, 3 SUSPENSIONS, BC	VENTILBLOK, 3 AUFHÄNGUNGEN,BC	BLOC VANNE, 3 FOUETS, BC
3	AE-031	BLOK MED KUGLER + NIPLER, BC	BLOCK WITH BALLS+NIPPLES, BC	BLOCK MIT KUGELN+ANSATZ, BC	BLOC AVEC BOULLES+EMBOUTS, B
4	AE-079	SELECTVENTIL, U/MAGNETV. BC	* VALVE BLOCK, NO SOLENOID, BC	* VENTILBLOK OHNE MAGNETV.,BC	* BLOCK VANNE SANS ELECTROV.BC
5	MC-946	TRÆKTRÅD FOR FC-BASIC	STRING, BC	ZUGDRAHT, BC	FIL, BC
6	MC-927	BØSNING TIL VANDREGULERING	BUSH, WATER REGULATION	BUCHSE, WASSERREGULIERUNG	DOUILLE, REGLAGE EAU
7	AE-033	PROP Ø5 MED O-RING	PLUG Ø5 WITH O-RING	STÖPSEL Ø5 MIT O-RING	BOUCHON Ø5 AVEC JOINT TORIQU
8	AE-045	PROP Ø4 MED O-RING	PLUG Ø4 WITH O-RING	STÖPSEL Ø4 MIT O-RING	BOUCHON Ø4 AVEC JOINT TORIQU
9	MC-925	REGULERINGSNØGLE FOR BC	REGULATING KEY, BC	EINSTELLSCHLÜSSEL, BC	CLE DE REGLAGE BC
10	MC-952	LÅSEBLIK, PROPPER,SELECTVENTIL	LOCKING PLATE,PLUGS,SELECT VAC	HALTEBLECH, STÖPSEL,SELECT VAC	PLAQUETTE VERROUILLAGE,SELEC
11	MC-948	BLÆNDEKSEL, OPHÆNG	COVER PLATE, INSTR. SUSPENSION	VERBLENDPLATTE,INST.AUFHÄNGUNG	CACHE FOUET
12	MC-949	LÅSEPLADE FOR BLÆNDEKSEL	LOCKING PLATE, COVER PLATE	HALTEBLECH, VERBLENDPLATTE	PLAQUETTE VERROUILLAGE, CACH
13	UC-683	PINOLSKRUE M 3 X 3	POINTED SCREW M 3 X 3	REITSTOCKSCHRAUBE M 3 X 3	VIS POINTUE M 3 X 3
14	AC-115	REGULERINGSVENTIL VAND/LUFT	REGULATING VALVE WATER/AIR	REGULIERVENTIL WASSER/LUFT	VANNE DE REGLAGE EAU/AIR
15	AC-164	MAGNETVENTIL,2 VEJS, M/STIK	SOLENOID VALVE, 2-WAY/PLUG	MAGNETVENTIL, 2-WEG/STECKER	ELEVTVANNE, 2-VOIES/PRISE
16	AC-165	MAGNETVENTIL,3 VEJS, M/STIK	SOLENOID VALVE,3-WAY/CONNECTOR	MAGNETVENTIL, 3-WEG/STECKER	ELEVTVANNE, 3-VOIES/PRISE
17	MR-150	PAKNING, SÆDE, SELECT VENTIL	GASKET, SEAT, SELECT VALVE	DICHTUNG, SITZ, SELECT VENTIL	JOINT, SIEGE, VALVE SELECT
18	SA-555	ANKER	PLUNGER	PLUNGER	PLONGEUR
19	SA-556	FJEDER	SPRING	FEDER	RESSORT
20	SD-308	O-RING	O-RING	O-RING	JOINT TORIQUE
21	SD-307	O-RING,12 X 1,5 B	O-RING	O-RING	JOINT TORIQUE
22	SD-303	O-RING, 12,5 X 1,1 N -O B75	O-RING	O-RING	JOINT TORIQUE
23	UC-041	O-RING 12 X 1,5 NB70	O-RING Ø 12 X 1.5	O-RING Ø 12 X 1,5	JOINT TORIQUE Ø 12 X 1,5
24	MR-227	O-RING Ø 3 X 1	O-RING Ø 3 X 1	O-RING Ø 3 X 1	JOINT TORIQUE Ø 3 X 1




Flex Dental A/S		1999-06-28
Sugeslangeholder, patientstol	Suction hose holder, patient chair	YA-001
Saugschlauchhalter, Patientenstuhl	Porte-tuyaux d'aspiration, fauteuil dentaire	0.6.39.B

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	MC-111	HÅNDTAG, BRO/LAMPE, 1 STK.	HANDLE, INSTRUMENT BRIDGE	HANDGRIFF, INSTRUMENTENTRÄGER	POIGNEE, PORTE-INSTRUMENT
1	SD-399	HÅNDTAG, BRO/LAMPE, 4 STK.	HANDLES (4) BRIDGE/LAMP	HANDGRIFFE (4) TRÄGER/LAMPE	POIGNEES (4) PORTE-INSTR./ECL
1	AC-371	HÅNDTAG, BRO/LAMPE, 24 STK.	HANDLES (24), BRIDGE/LAMP	HANDGRIFFE (24), TRÄGER/LAMPE	POIGNEES (24) PORTE-INSTR./ECL
2	MC-186	MODTAGER, HÅNDTAG, INSTRUMENTBRO	HOLDER, HANDLE, INSTR. BRIDGE	EINSTECKROHR, GRIFF, INST. TRÄGER	TUYAU SUPP., POIGNEE, PORTE-IN
3	UC-618	SKRUE M 8 X 10 JERN ELFORZ	SCREW M 8 X 10	SCHRAUBE M 8 X 10	VIS M 8 X 10
4	MC-683	HOLDER, SUGENIPLER, STOL	HOLDER, SUCTION NIPPLE, CHAIR	HALTER, SAUGSCHLAUCHANS., STUHL	SUPPORT, EMBOUT ASPIR. FAUTEUI
5	MC-943	BESLAG, SLANGEH., STOL, NCS 4502	FITTINGS, SUCTION, CHAIR, 9018	BESCHLAG, SAUG, STUHL, 9018	FIXATION, ASPIR., FAUTEUIL, 901
6	MC-878	TOP, SUG, STOL	TOP, SUCTION, CHAIR	OBERTAIL, ABSAUGER, STUHL	PARTIE SUP., ASPIR., FAUTEUIL
7	MC-933	COVER FOR SUG	COVER, SUCTION	ABDECKUNG, ABSAUGER	CARTER, ASPIRATION
8	AD-403	PRINT, SUGESLANGEHOLDER	* PRINTED CIRC. "SUCTION HOLDER" *	PLATINE "SUCTION HOLDER"	* CARTE "SUCTION HOLDER"
9	UC-347	SILICONESLANGE Ø 21 X 16,8	SILICONE TUBE Ø 21 X 16,8	SILIKONSCHLAUCH Ø 21 X 16,8	TUYAU SILICONE Ø 21 X 16,8
10	MY-003	SPIRAL 40 CM, MESSING	SPIRAL 40 CM, BRASS	SPIRALE 40 CM, MESSING	SPIRALE 40 CM, LAITON
11	AB-160	KABEL, SUGESLANGEHOLDER, STOL2	CABLE, SUCT. HOSE HOLDER, CHAIR2	KABEL, SAUGSCHL. HALTER, STUHL2	CÂBLE, PORTE-TUYAUX ASPI., FAU
12	AB-161	KABEL I STOL2 F/ SUGESL. HOLDER	CABLE I CHAIR2 F/ SUCT. HOSE H.	KABEL I STUHL2 F/ SAUGSL. HALTER	CÂBLE I FAUT. 2 P/ PORTE-TUY. A
13	AB-162	KABEL, UNIT-CHAIR TIL SUG. II	CABLE, UNIT-CHAIR FOR SUCT. II	KABEL, UNIT-STUHL FÜR SAUG. II	CÂBLE, UNIT-FAUTEUIL P. ASPIR.

14 AC-849 Connector

Cable connector

Kabelstecker




Flex Dental A/S
 Fontæneventil: Dúrr
 Speitontæneventil: Dúrr

Cuspidor valve: Dúrr
 Vanne crachoir: Dúrr

1999-06-28
 YA-001
 0.6.40.B

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	AE-957	T-STYKKE Ø32 M. OVERLØB	T-CONNECTOR Ø32 OVERFLOW	T-STÜCK Ø32, ÜBERFLUSS	CONNECTEUR T, Ø32 DEBORDEMENT
2	AE-958	SLANGE, FONTÆNEVENTIL-GULV	HOSE, CUSPIDOR VALVE-FLOOR	SCHLAUCH, SPEIFONTÄNENV.-BODEN	TUYAU, VANNE, CRACHOIR-
3	AE-966	ADAPTOR, Ø20 - Ø25, M. OVERLØB	ADAPTOR, Ø20-Ø25, OVERFLOW	ADAPTOR, Ø20 - Ø25, ÜBERLAUF	ADAPTEUR, DEBORDEMENT Ø20-Ø2
4	UC-822	ADAPTOR, Ø36-Ø32 M. UDTAG, SUG	PLAST FITTING COMP.WITH O-RING	PLAST FITTING COMP.WITH O-RING	RACCORD ASPI Ø36/32
5	UC-823	FONTÆNEVENTIL, DÚRR	* CUSPIDOR VALVE, DÚRR	SPEIFONTÄNENVENTIL, DÚRR	VANNE, CRACHOIR, DÚRR
6	UC-828	OVERGANG Ø32 - Ø36	ADAPTOR Ø32 - Ø36		PIÈCE INTERMÉDIARE Ø32 - Ø34
7	UC-838	PRINT, FONTÆNEVENTIL	PRINTED CIRC. CUSPIDOR VALVE	PLATINE, SPEIFONTÄNENVENTIL	CARTE "VANNE, CRACHOIR"
8	UC-839	TASTATUR, FONTÆNEVENTIL	KEYBOARD, CUSPIDOR VALVE	TASTATUR, SPEIFONTÄNENVENTIL	JEU DE CONTACTS, VANNE CRACHO




Flex Dental A/S
Spuleenhet: Type 98, Dürr
Spüleinheit: Typ 98, Dürr

Flushing unit: Type 98, Dürr
Unité de rinçage: Type 98, Dürr

1999-06-28
YA-001
0.6.41.B

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	AB-954	SLANGE, TILSLUTN. VÅDT SUG	HOSE, CONNECTION, WET SUCTION	SCHLAUCH, ANSCHLUSS, NASSABS. S.	CONNEXION TUYAU ASPI HUMIDE
2	UC-821	SPULEENHED TYPE 7100	FLUSHING UNIT, TYPE 7100	SPÜLEINHEIT, TYPE 7100	UNIT DE RINCAGE, TYPE 7100
3	UC-822	ADAPTOR, Ø36-Ø32 M. UDTAG, SUG	PLAST FITTING COMP. WITH O-RING	PLAST FITTING COMP. WITH O-RING	RACCORD ASPI Ø36/32
4	UC-835	MAGNETVENTIL, SPULEENHED, 7100	SOLENOID VALVE, FLUSHING U, 7100	MAGNETVENTIL, SPÜLEINHEIT, 7100	UNIT DE RINCAGE, VANNE CRA, 71
5	UC-836	PRINT, SPULEENHED, 7100	PRINTED C., FLUSHING UNIT, 7100	PLATINE, SPÜLEINHEIT, 7100	CARTE, UNIT DE RINCAGE, 7100
6	UC-837	FILTER, SPULEENHED, 7100	FILTER, FLUSHING UNIT, 7100	FILTER, SPÜLEINHEIT, 7100	FILTRE, UNIT DE RINCAGE, 710




No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	AD-462	PRINT INST.SELEK.INTEGRAL 5 *	PRINTED C, INST.SELEC.INTEG5 *	PLATINE, INST.SELEC.INTEG 5 *	CARTE, INSTR.SELEC.INTEGRAL5
2	AD-466	PRINT, INSTR. CONTR. INTEG/5 *	PR.CIRCUIT, INSTR.CONTROL TY 5	PLATINE, INSTR.CONTROL, TY 5	CARTE, INSTR.CONTROL TY 5
3	AD-468	PRINT LC ADAPT 1 INTEGRAL 5	PRINTED C, LC ADAPT 1 INTEGR 5	PLATINE, K.LICHT ADAPT 1 INT 5	CARTE, ADAPT.L.Á COMPO 1 INT
3	AD-471	PRINT LC ADAPT 2 INTEGRAL 5	PRINTED C, LC ADAPT 2 INTEGR 5	PLATINE, K.LICHT ADAPT 2, INT 5	CARTE, ADAPT L.Á COMPO 2, INT
4	AE-751	DÆKSEL INSTR.BRO TY5 RAL 9018	COVER, INSTR.BRIDGE, TYP.5 9018	ABDECK. INSTR. TRÄGER, TYP.5 9018	CARTER, PORTE-INSTR, TYPE5 90
4	AE-752	DÆKSEL INSTR.BRO TY5 CLASSIC	COVER, INSTR.BRIDGE, T5 CLASSIC	ABDECK. INSTR. TRÄGER, T5 CLASSIC	CARTER, PORTE-INSTRU, T5 CLASS
5	AB-753	BALANCEOPHÆNG, TURBINE MIDWEST*	BALANCED SUSPENSION, MID-WEST	BALANCEAUFHÄNGUNG, MID-WEST	FOUET, MID-WEST
6	AF-015	E-PROM PROGRAMMERET, NTC5, 440	E.PROM PROGRAMMED, NTC 5, 440	E.PROM PROGRAMMIERT, NTC 5, 440	E.PROM PROGRAMMÉE, NTC 5, 44
6	AF-017	E-PROM PROG.JAPAN NTC5, AD-446	E-PROM PROG.JAPAN NTC5, AD-446	E-PROM PROG.JAPAN NTC5, AD-446	E-PROM PROG.JAPAN NTC5, AD-4
6	AF-019	E-PROM PROGRAMM. FC-V5 AD-440			
6	AF-027	E-PROM PROGRAMM. FC-V5 AD-446	EPROM PROGRAMM. FC-V5, AD-446	EPROM PROGRAMM. FC-V5, AD-446	EPROM PROGRAMM. FC-V5, AD-446
6	AF-029	E-PROM PROGRAMM.BC/FC5 AD-440	EPROM PROGRAMM. BC/FC5 AD-440	EPROM PROGRAMM. BC/FC5 AD-440	EPROM PROGRAMM. BC/FC5 AD-44
6	AF-032	E-PROM PROG. JAPAN FC5, AD-446	EPROM PROG. JAPAN FC5, AD-446	EPROM PROG. JAPAN FC5, AD-446	EPROM PROG. JAPAN FC5, AD-44
7	ME-182	INSTRUMENTANLÆG, 5 INSTR.	INSTRUMENT PAD INTEGRAL TY 5	INSTRUMENTENABLAG INT.TY 5	REPOSE-INSTRUMENT INTEGRALTY
8	ME-184	GUMMILISTE, 5 INSTRUMENTER	RUBBER, 5 INSTRUMENT	GUMMILEISTE, 5 INSTRUMENTEN	LISTEAU, CAOUTCHOUC, 5 INSTR.
9	ME-183	BASIS HOLDER, INSTRUMENTOPHÆNG	BASE HOLDER, INSTRUMENT SUSPENS	HALTER, INSTRUMENT AUFHÄNGUNG	SUPPORT, SUSPENSION, INSTRUME
10	ME-188	AFDÆKNING, AD-462	COVERING, AD-464	ABDECKUNG, AD-464	COUVERTURE, AD-464
11	UC-859	SKRUE, PT UHX 3,5X9,5 Z POSI	SCREW, PT UHX 3,5X9,5 Z POSI	SCHRAUBE, PT UHX 3,5X9,5 Z POS	VIS, PT UHX 3,5 X 9,5 Z POSI

5 AE-245 Rulle m/fjeder Turbine
12 MH-862 Kodeblock, Scaler / 5

Rocker piece Turbine
Code block, scaler / 5

Kippteil Turbine
Kodeblock, ZEG / 5




Flex Dental A/S
Tandrenser: Amdent
Zahnreiniger: Amdent

Scaler: Amdent
Detartreur: Amdent

1999-06-28
YA-001
0.6.43.B

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	AD-467	PRINT SCALER ADAPT INTEGRAL 5	PR.CIRC, SCALER ADAPT, INTEG.5	PLATINE, ZEG ADAPTER, INTEGR.5	CARTE, ADAPT.DETARTR, INTEGR
2	AE-756	RULLE M.MAGN.+ FJEDER, SCALER	ROCKER PIECE, SCALER	KIPPTTEIL, SCALER	PIECE BASCUL., SCALER
3	AE-761	VENTIL, VAND, INTEGRAL/5 BC	WATER VALVE, INTEGRAL BC/5	WASSERVENTIL, INTEGRAL BC/5	VANNE EAU, INTEGRAL BC/5
4	AE-762	VENTIL, VAND INTEGRAL/5 FC/NTC	VALVE,WATER INTEGRAL/5 FC/NTC	VENTIL,WASSER INTEGRAL/5FC/NTC	VANNE, EAU, INTEGRAL/5PC/NTC
5	AE-763	KABEL, SCALERTILPAS. INTEGRAL5	CABLE, SCALER ADAPT. INTEGRAL5	KABEL, ZEG ADAPTER, INTEGRAL 5	CÂBLE, ADAPT.DETAR, INTEHRAL
6	AE-785	STYRING, TANDRENSER AMDENT *	CONTROL, SCALER AMDENT	STEUERUNG, ZEG AMDENT	GENERATEUR, DETARTREUR AMDEN
7	AE-786	KABEL TIL AMDENSCALER INTEGRA5	CABLE, AMDENT SCALER,INTEGRAL5	KABEL, AMDENT ZEG,INTEGRAL 5	CÂBLE, DETART.AMDENT,INTEGRA
8	AE-789	BALANCEOPHÆNG, AMDENT SCALER *	BALNCED SUSP. AMDENT SCALER	BALANCEAUFHÄNG.AMDENT SCALER	FOUET, AMDENT SCALER
9	AE-793	ADAPTOR, HF THERMO, AMDENT	ADAPTOR, HF THERMO, AMDENT	ADAPTER, HF THERMO, AMDENT	ADAPTATEUR, HF THERMO, AMDEN
10	HE-116	AMDENT SPIDS NR.31, PADDLE	AMDENT TIP NR. 31, PADDLE	AMDENT SPITZE NR. 31, PADDLE	AMDENT INSERT NR. 31, PADDLE
11	HE-117	AMDENT SPIDS NR.32, HOE	AMDENT TIP NR. 32, HOE	AMDENT SPRITZE NR. 32, HOE	AMDENT INSERT NR. 32, HOE
12	HE-118	AMDENT SPIDS NR.33, PROBE	AMDENT TIP NR. 33, PROBE	AMDENT SPITZE NR. 33, PROBE	AMDENT INSERT NR. 33, PROBE
13	HE-119	AMDENT SPIDS NR.31L, PROBE	AMDENT TIP NR 31L, PROBE	AMDENT SPITZE NR. 31L, PROBE	AMDENT INSERT NR. 31L. PROBE
14	HE-120	AMDENT SPIDS NR.33R, PROBE	AMDENT TIP NR. 33R, PROBE	AMDENT SPITZE NR. 33R, PROBE	AMDENT INSERT NR. 33R, PROBE
15	HE-121	AMDENT SPIDS NR.34 KRONEFJERNER	AMDENT TIP NR.34,CROWN REMOVER	AMDENT SPITZE NR.34,KRONENLÖSE	AMDENT INSERT NR.34,DESCELLE
16	HE-122	AMDENT SPIDS NR.34S KRONEFJERN	AMDENT TIP 34S, CROWN REMOVER	AMDENT SPITZE 34S,KRONENLÖSER	AMDENT INSERT 34S,DESCELLEUR
17	HE-123	AMDENT SPIDS NR.35, THIXO	AMDENT TIP NR. 35, THIXO	AMDENT SPITZE NR. 35 THIXO	AMDENT INSERT NR. 35, THIXO
18	HE-124	AMDENT SPIDS NR.35TA, THIXO	AMDENT TIP NR. 35TA, THIXO	AMDENT SPITZE NR. 35TA, THIXO	AMDENT INSERT NR. 35TA, THIX
19	HE-125	AMDENT SPIDS NR.35TB, THIXO	AMDENT TIP NR. 35TB, THIXO	AMDENT SPITZE NR. 35TB, THIXO	AMDENT INSERT NR. 35TB, THIX
20	HE-126	AMDENT SPIDS NR.35TC, THIXO	AMDENT TIP NR. 35TC, THIXO	AMDENT SPITZE NR. 35TC, THIXO	AMDENT INSERT NR. 35TC, THIX
21	HE-127	AMDENT SPIDS NR.37, THIN	AMDENT TIP NR. 37, THIN	AMDENT SPITZE NR. 37, THIN	AMDENT INSERT NR. 37, THIN
22	HE-128	AMDENT SPIDS NR.37L, THIN	AMDENT TIP NR. 37L, THIN	AMDENT SPITZE NR. 37L, THIN	AMDENT INSERT NR. 37L, THIN
23	HE-129	AMDENT SPIDS NR.37R, THIN	AMDENT TIP NR. 37R, THIN	AMDENT SPITZE NR. 37R, THIN	AMDENT INSERT NR. 37R, THIN
24	HE-130	AMDENT SPIDS NR.39, POWER	AMDENT TIP NR. 39, POWER	AMDENT SPITZE NR. 39, POWER	AMDENT INSERT NR. 39, POWER
25	HE-131	AMDENT O-RING FOR/GRØN 9,0X1,5	AMDENT O-RING F/GREEN, 9,0X1,5	AMDENT O-RING F/GRÜN, 9,0X1,5	AMDENT TORIQUE F/VERT, 9,0X1,
25	HE-132	AMDENT O-RING BAG/SORT10,1X1,6	AMDENT O-RING B/BLACK 10,1X1,6	AMDENT O-RING HINTER 10,1X1,6	AMDENT TORIQUE CAPOT 10,1X1,




Flex Dental A/S
 LC-lampe: Demetron
 Kompositlicht: Demetron

ComPOSITE Lamp: Demetron
 Lampe à composite: Demetron

1999-06-28
 YA-001
 0.6.44.B

No	NUMBER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	AE-760	STYRING LC-LAMPE, DEMETRON *	CONTROL LC, DEMETRON	STEUERUNG LC, DEMETRON	GENERATEUR LC, DEMETRON
2	AE-795	BALAN.OPHÆNG DEM KUN BAGKAPPE*	BALANCED SUSP.DEMETRON	BALANCEAUFHÄNGUNG, DEMETRON	FOUET, DEMETRON
3	SD-630	PÆRE, DEMETRON	BULB, DEMETRON	BIRNE, DEMETRON	AMPOULE, DEMETRON
4	SD-632	BAGPART, HUS, DEMETRON	REAR HOUSE, DEMETRON	HINTERGEHÄUSE, DEMETRON	CAPOT ARRIÈRE,LAMPE à DEMETR
5	SD-633	SKRUE FOR LC-LAMPE, DEMETRON	SCREW FOR LC, DEMETRON	SCHRAUBE, LC DEMETRON	VIS, LC DEMETRON
6	SD-634	FORSTYKKE , DEMETRON LC	FRONT COVER, DEMITRON LC	VORDERGEHÄUSE, DEMITRON LC	CAPOT AVANT L-A-C DEMETRON
7	SD-635	FRONT DEKSEL, DEMETRON	COVER, FRONT, DEMETRON	DECKEL, FRONT, DEMETRON	CARTER, FRONT, DEMETRON
8	SD-636	SKRUE FOR FRONT - DEMETRON	SCREW, FRONT, DEMETRON	SCHRAUBE, FRONT, DEMETRON	VIS, FRONT, DEMETRON
9	SD-637	FILTER, DEMETRON	FILTER, DEMETRON	FILTER, DEMETRON	FILTRE, DEMETRON
10	SD-638	VARMEFILTER, DEMETRON	FILTER, HEAT, DEMETRON	FILTER, HITZE,DEMETRON	FILTRE DE CHALEUR, DEMETRON
11	SD-639	O-RING, FRONT, DEMETRON	O-RING, FRONT, DEMETRON	O-RING, FRONT, DEMETRON	JOINT TORIQUE, FRONT, DEMETR
12	SD-640	BLÆSER KOMPLET, DEMETRON	FAN, COMPLETE, DEMETRON	VENTILATOR, KOMPLET, DEMETRON	VENTILATEUR, COMPLET, DEMETR
13	SD-641	LAMPESOKKET KOMPLET, DEMETRON	LAMP SOCKET COMPLETE, DEMETRON	LAMPESOKKET, DEMETRON	DOUILLE POUR LAMPE, DEMETRON
14	SD-642	SOKKET FOR LAMPE, DEMETRON	LEG, LAMP, DEMETRON	SOKKET, LAMP, DEMETRON	FIXEMENT, LAMP, DEMETRON
15	SD-643	PLADE FOR LAMPE, DEMETRON	COVER, LAMP, DEMETRON	PLATTE, LAMPE, DEMETRON	PLAQUE, LAMPE, DEMETRON
16	SD-644	PRINTKORT FOR LAMPE, DEMETRON	PRINTED CIRC. LC DEMETRON	PLATINE, LC DEMETRON	CARTE, LC DEMETRON
17	SD-645	FASTGØRELSESELEMENT, DEMETRON	CONNECTION PIECE, DEMETRON	ANSLUSSTEIL, DEMETRON	FIXATION POUR LAMPE DEMETRON
18	SD-676	LYSSTAV FOR DEMETRON LC-LAMPE	FIBRE ROD, DEMETRON LC-LAMPE	FIBERLICHTSTAB, DEMETRON LC	TIGE EN FIBRE, DEMETRON LC




Flex Dental A/S
 Motor: Type 98, MC3
 Motor: Typ 98, MC3

Motor: Type 98, MC3
 Moteur: Type 98, MC3

1999-06-28
 YA-001
 0.6.45.B


No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	AE-927	OPHÆNG, MOTOR, TYPE MC3	* MOTOR TUBE, TYPE MC3	MOTORSCHLAUCH, TYPE MC3	CORDON MOTEUR, TUPE MC3
2	HE-005	PÆRE, LUX MOTOR/TURBINE	BULB, LUX MOTOR/TURBINE	BIRNE, LUX MOTOR/TURBINE	AMPOULE, MOTEUR LUX/TURBINE
3	HE-138	BESKYTTELSKAPPE, MOTOR MC3	JACKET, MOTOR MC3	GEHÄUSE, MOTOR MC3	CARTER, MOTEUR MC3
4	HE-139	KOBLINGSRØR BIEN AIR MOTOR MC3	COOPLING TUBE,BAIR MOTOR,MC3	KUPPLUNGSROHR,BAIR MOTOR,MC3	TUYAU, RACCORDEMENT,BA MOT M
5	HE-140	O-RING, BIEN AIR MOTOR, MC3	O-RING, BIEN AIR MOTOR, MC3	O-RING, BIEN AIR MOTOR, MC3	JOINT TORIQUE, MOTEUR BA, MC
6	SD-629	KAPPE, MOTOR MC 3	COVER, MOTOR MC 3	KAPPE, MOTOR MC 3	CAPOT MOTEUR MC 3
7	HE-142	KUL, BIEN AIR MOTOR, MC3	CARBON, BIEN AIR MOTOR, MC3	KOHLE, BIEN AIR MOTOR, MC3	CHARBON, BIEN AIR MOTOR, MC3
8	HE-143	FJEDER, BIEN AIR MOTOR, MC3	SPRING, BIEN AIR MOTOR, MC3	FEDER, BIEN AIR MOTOR, MC3	RESSORT, BIEN AIR MOTOR, MC3
9	HE-144	KULHOLDER, BIEN AIR MOTOR, MC3	CARBON HOLDER, B AIR MOTOR,MC3	KOHLEHALTER,BIEN AIR MOTOR,MC3	SUPPORT,CHARBON,B.AIR MOTE,M
10	HE-145	PAKNING, BIEN AIR MOTOR, MC3	GASKET, BIEN AIR MOTOR, MC3	PACKUNG, BIEN AIR MOTOR, MC3	JOINT, MOTEUR BIEN AIR, MC3
11	HT-099	MOTOR, MC3, INCL. EMBALLAGE	* MOTOR, MC3	MOTOR, MC3	MOTEUR, MC3
12	SA-024	O-RING, BIEN AIR/FLEX MOTOR	O-RING, BIEN AIR/FLEX MOTOR	O-RING, BIEN AIR/FLEX MOTOR	JOINT TOR.,MOTEUR BIEN A./FL
1	AE-264	Rulle m/fjeder MC3	Rocker piece MC3	Kippteil MC3	



No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	AE-064	ADAPTOR, FODKONTAKT, BC	ADAPTOR, FOOT CONTROL, BC	ADAPTER, FUSSANLASSER, BC	ADAPTATEUR, PEDALE, BC
2	AE-757	AFBLÆNDING, JOYSTICK, FODK.	COVER, JOYSTICK, FOOTCONTROL	ABDECKUNG, JOYSTICK, FUSSANLASS.	COUVERTURE, JOYSTICK, PEDALE.
3	AE-758	SÆT, BØJLE, FODKONTAKT	KIT, HOOP, FOOTCONTROL TYPE II	FUSSANLASSERBÜGEL TYPE II	ETRIER, PEDALE TYPE II
4	AE-792	FODKONTAKT, PEDALARM MED STIFT	FOOT CONTROL, PEDAL ARM W PIN	FUSSANLASSER, PEDALARM M STIFT	PEDALE, LEVIER AVEC AXE
5	AE-969	PRINT INKL. BESLAG OG POTMETER	PRINTBOARD COMPL.FOOTCONTRL II	PLATINE, KOMPL.FUSSANLASSER II.	PLATINE COMPL. PEDALE II.
6	AE-974	JOYSTICK, FODKONTAKT	JOYSTICK, FOOTCONTROL II.	JOYSTICK, FUSSANLASSER II	JOYSTICK, PEDALE II
7	AE-975	FODKONTAKT, FC/NTC, TYPE II *	FOOTCONTROL, FC/NTC, TYPE II	FUSSANLASSER, FC/NTC, TYPE II	PEDAL, FC/NTC, TYPE II
8	AE-976	FODKONTAKT, BC, TYPE II *	FOOTCONTROL, BC, TYPE II	FUSSANLASSER, BC, TYPE II	PEDALE, BC, TYPE II
9	AE-977	KABEL FODKONTAKT	CABLE, FOOTCONTROL II	KABEL, FUSSANLASSER II	CABLE, PEDALE II
10	FA-752	O-RING Ø 9.75 X 1.78 NBR 70	O-RING Ø 9.75 X 1.78 NBR 70	O-RING Ø 9.75 X 1.78 NBR 70	JOINT TORIQUE Ø9.75 X 1.78
11	ME-161	FJEDER, PEDALARM, FODKONTAKT	SPRING, CRANK, FOOTCONTROL II	FEDER, PEDALARM, ANLASSER II	RESORT, BRAS DE PEDALE II
12	ME-162	BØSNING PEDALARM, FODKONTAKT	BUSH, CRANK, FOOTCONTROL II	BÜCHSE, PEDALARM, ANLASSER II	DOUILLE, BRAS DE PEDALE II
13	MV-018	GUMMI, FODKONTAKT	RUBBER, FOOTCONTROL II	GUMMIE, FUSSANLASSER II	CAOUTCHOUC, PEDALE II
14	MV-021	DÆKSEL, FODKONTAKT	COVER, FOOTCONTROL II	KAPPE, FUSSANLASSER II	CAPOT, PEDALE II
15	UA-058	SKRUE M 3 X 6 UH GEVINDF.	SCREW M 3 X 6	SCHRAUBE M 3 X 6	VIS M 3 X 6
16	UA-093	SCOTCH PRODUKTFØDDER	RUBBER FOOT, FOOTCONTROL II	GUMMIEFUSS, FUSSANLASSER II	PIED CAOUTCHOUC, PEDALE II
17	UA-588	SKRUE M 6 X 25 CH.INDV.S.Z.	SCREW, M 6 X 25	SCHRAUBE, M 6 X 25	VIS, M 6 X 25
18	UA-859	SKRUE, ECO-SYN(T) Ø3.5X25	SCREW, ECO-SYN, 3.5 X 25.	SCHRAUBE, ECO-SYN, 3.5 X 25	VIS, ECO-SYN, 3.5 X 25
19	UA-860	SKRUE M 4 X 20 UHZX	SCREW, M 4 X 20	SCHRAUBE, M 4 X 20	VIS, M 4 X 20

Tech News 16:

AE-899 Modifikationskit, horisontale bevæg. Modification kit, horisontal move. Modifikationskit, horizontale Bewegung
SD-691 Modifikationskits, vertikale bevæg. Modification kits, verticale move. Modifikationskits, Vertikale Bewegung



Flex Dental A/S
Tandrenser: Odontoson
Zahnreiniger: Odontoson

Scaler: Odontoson
Détartreur: Odontoson

1999-06-28
YA-001
0.6.47.B


No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	AD-467	PRINT SCALER ADAPT INTEGRAL 5	PR.CIRC, SCALER ADAPT, INTEG.5	PLATINE, ZEG ADAPTER, INTEGR.5	CARTE, ADAPT.DETARTR, INTEGR
2	AE-754	BALANCEPHÆNG ODONTOSONSCALER*	BALANCED SUSP. ODONTOSONSCALER	BALANCEAUFHÄNG, ODONTOSONSCALER	FOUET ODONTOSON SCALER
3	AE-756	RULLE M.MAGN.+ FJEDER, SCALER	ROCKER PIECE, SCALER	KIPPTTEIL, SCALER	PIECE BASCUL., SCALER
4	AE-761	VENTIL, VAND, INTEGRAL/5 BC	WATER VALVE, INTEGRAL BC/5	WASSERVENTIL, INTEGRAL BC/5	VANNE EAU, INTEGRAL BC/5
5	AE-762	VENTIL, VAND INTEGRAL/5 FC/NTC	VALVE, WATER INTEGRAL/5 FC/NTC	VENTIL, WASSER INTEGRAL/5 FC/NTC	VANNE, EAU, INTEGRAL/5 FC/NTC
6	AE-763	KABEL, SCALERTILPAS. INTEGRAL5	CABLE, SCALER ADAPT. INTEGRAL5	KABEL, ZEG ADAPTER, INTEGRAL 5	CÂBLE, ADAPT.DETAR, INTEHRAL
7	AE-779	ADAPTOR, HF THERMO, ODONTOSON	ADAPTOR, HF THERMO, ODONTOSON	ADAPTER, HF THERMO, ODONTOSON	ADAPTATEUR, HF THERMO, ODONTOS
8	AE-780	KABEL TIL GOOFSCALER INTEGRAL5	CABLE, GOOF SCALER, INTEGRAL 5	KABEL, GOOF SCALER, INTEGRAL 5	CÂBLE, GOOF SCALER, INTEGRAL
9	AE-782	STYRING, TANDRENSER GOOF	* CONTROL, SCALER GOOF	STRUERUNG, ZEG GOOF	GENERATEUR, DETARTREUR GOOF
10	HT-144	ODONTOSON TANDRENSER GRÅ U. QF			

Flex Dental A/S
Tandrenser: EMS
Zahnreiniger: EMS

Scaler: EMS
D\etartreur: EMS

1999-06-28
YA-001
0.6.48.B

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	AD-467	PRINT SCALER ADAPT INTEGRAL 5	PR.CIRC, SCALER ADAPT, INTEG.5	PLATINE, ZEG ADAPTER, INTEGR.5	CARTE, ADAPT.DETARTR, INTEGR
2	AE-755	BALANCEOPHÆNG EMS SCALER	* BALANCED SUSPENSION EMS SCALER	BALANCEAUFHÄNGUNG EMS SCALER	FOUET EMS SCALER
3	AE-756	RULLE M.MAGN.+ FJEDER, SCALER	ROCKER PIECE, SCALER	KIPPTAIL, SCALER	PIECE BASCUL., SCALER
4	AE-759	STYRING, TANDRENSER EMS	* CONTROL, SCALER EMS	STEUERUNG, ZEG EMS	GENERATEUR, DETARTREUR EMS
5	AE-761	VENTIL, VAND, INTEGRAL/5 BC	WATER VALVE, INTEGRAL BC/5	WASSERVENTIL, INTEGRAL BC/5	VANNE EAU, INTEGRAL BC/5
6	AE-762	VENTIL, VAND INTEGRAL/5 FC/NTC	VALVE, WATER INTEGRAL/5 FC/NTC	VENTIL, WASSER INTEGRAL/5FC/NTC	VANNE, EAU, INTEGRAL/5FC/NTC
7	AE-763	KABEL, SCALERTILPAS. INTEGRAL5	CABLE, SCALER ADAPT. INTEGRAL5	KABEL, ZEG ADAPTER, INTEGRAL 5	CÂBLE, ADAPT.DETAR, INTEHRAL
8	AE-778	ADAPTOR, HF THERMO, EMS	ADAPTOR, HF THERMO, EMS	ADAPTER, HF THERMO, EMS	ADAPTEUR, HF THERMO, EMS
9	AE-781	KABEL TIL EMSSCALER INTEGRAL5	CABLE, EMS SCALER, INTEGRAL 5	KABEL, EMS ZEG, INTEGRAL 5	CÂBLE, DETARTR. EMS, INTEGRAL
10	HE-133	SPIDS TYPE A, EMS PIEZON SCALER	TIP TYPE A, EMS PIEZON SCALER	SPITZE TYP A, EMS PIEZON SCALER	INSERT TYP A, EMS PIEZON SCAL
10	HE-134	SPIDS TYPE B, EMS PIEZON SCALER	TIP TYPE B, EMS PIEZON SCALER	SPITZE TYP B, EMS PIEZON SCALER	INSERT TYP B, EMS PIEZON SCAL
10	HE-135	SPIDS TYPE C, EMS PIEZON SCALER	TIP TYPE C, EMS PIEZON SCALER	SPITZE TYP C, EMS PIEZON SCALER	INSERT TYP C, EMS PIEZON SCAL
11	SD-663	HÅNDSTYKKE, EMS SCALER	* HANDPIECE, EMS SCALER	HANDSTÜCK, EMS SCALER	PIÈCE MAIN. EMS SCALER
12	SD-664	VÆRKTØJ, UDSKIFT AF SPIDS EMS	TOOL, CHANGE OF SCALER TIP, EMS	WERKZEUG, ZEG-SPITZE, EMS	OUTIL, CHANGEMENT INSERT, EMS
13	SD-693	O-Ring (10 Stk)	O-rings EMS (10 Pcs)	O-Ring EMS (10 Stck)	



Flex Dental A/S		1999-06-28
Videoskærm, 12" video	Monitor, 12" video	YA-001
Monitor, 12" Video	Moniteur, vidéo 12"	0.6.49.B

No	NUMMER	BESKRIVELSE	DESCRIPTION	BEZEICHNUNG	DESCRIPTION
1	AE-993	FÆRDIG STRØMFORSYNING FOR 12"	POWER SUPPLY FOR 12" SCREEN	TRANSFORMATOR FÜR 12" SCHIRM	ADAPTATEUR SECTEUR, ECRAN 12"
2	ZE-998	MONITOR 12", OMBYTNING	SCREEN 12", EXCHANGE	MONITOR 12", UMTAUSCH	ECRAN 12", ECHANGE
3	ZE-996	SKÆRM ARM MED KABLER, CLASSIC	SCREEN ARM WITH CABLE, CLASSIC	MONITOR ARM MIT KABEL, CLASSIC	BRAS ÉCRAN AVEC CABLE, CLASS
3	ZE-995	SKÆRM ARM MED KABLER, INTEGRAL	SCREEN ARM WITH CABLE, 9018	MONOTOR ARM MIT KABEL, 9018	BRAS ECRAN AVEC CÂBLE, 9018
4	ML-566	VGA SKÆRM FORLÆNGERKABEL, 5M	EXTENSION FOR VGA SCREEN CABLE	KABEL FÜR MONITOR, MULTIMED. 5M	EXTENSION CÂBLE ECRAN VGA

Monitor 15" - Type Mermaid

1 AF-052	Strømforsyning 15"	Power supply 15"	Trafo 15"
2 AF-919	Monitor 15"	Screen 15"	Flachbildschirm 15"
3 ZE-996	= AE-996		
3 AE-997	Kabel for arm med stik	Cable for arm with connector	Kabel für Arm mit Steckern
4 ZE-995	= AE-995		
5 AE-588	Dæksel, arm	Cover, arm	Deckel, arm
5 UA-076	Skrue for dæksel	Screw for cover	Schraube für Deckel

NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No

AA-030	1.8.32.B	40	AB-138	1.8.13.B	13	AB-197	1.8.33.B	1	AB-254	1.8.27.B	2
AA-112	1.8.32.B	40	AB-139	1.8.15.B	7	AB-197	1.8.35.B	6	AB-260	1.8.25.B	18
AA-469	1.8.25.B	36	AB-141	1.8.15.B	9	AB-199	1.8.33.B	12	AB-262	1.8.10.B	32
AB-003	1.8.05.B	11	AB-142	1.8.15.B	9	AB-200	1.8.34.B	1	AB-262	1.8.27.B	3
AB-007	1.8.11.B	10	AB-143	1.8.15.B	10	AB-201	1.8.34.B	3	AB-264	1.8.10.B	30
AB-009	1.8.15.B	12	AB-144	1.8.15.B	10	AB-202	1.8.30.B	1	AB-266	1.8.10.B	33
AB-009	1.8.16.B	11	AB-145	1.8.18.B	9	AB-202	1.8.34.B	9	AB-266	1.8.27.B	4
AB-012	1.8.25.B	14	AB-146	1.8.19.B	30	AB-203	1.8.34.B	19	AB-269	1.8.29.B	14
AB-017	1.8.24.B	7	AB-147	1.8.22.B	4	AB-204	1.8.34.B	20	AB-270	1.8.32.B	48
AB-017	1.8.25.B	12	AB-148	1.8.34.B	6	AB-205	1.8.34.B	21	AB-276	1.8.01.B	8
AB-018	1.8.25.B	15	AB-149	1.8.34.B	22	AB-206	1.8.35.B	1	AB-277	1.8.01.B	10
AB-027	1.8.10.B	14	AB-150	1.8.29.B	5	AB-207	1.8.35.B	8	AB-279	1.8.29.B	18
AB-029	1.8.10.B	14	AB-151	1.8.29.B	9	AB-208	1.8.35.B	9	AB-289	1.8.38.B	1
AB-040	1.8.38.B	1	AB-152	1.8.34.B	4	AB-209	1.8.35.B	11	AB-290	1.8.38.B	1
AB-048	1.8.12.B	16	AB-153	1.8.29.B	23	AB-210	1.8.39.B	2	AB-291	1.8.38.B	3
AB-048	1.8.13.B	20	AB-154	1.8.29.B	24	AB-211	1.8.39.B	3	AB-292	1.8.38.B	3
AB-049	1.8.16.B	10	AB-155	1.8.29.B	28	AB-212	1.8.39.B	4	AB-293	1.8.38.B	3
AB-053	1.8.32.B	37	AB-156	1.8.29.B	32	AB-213	1.8.39.B	5	AB-298	1.8.28.B	16
AB-075	1.8.10.B	14	AB-158	1.8.29.B	38	AB-214	1.8.39.B	6	AB-395	1.8.11.B	13
AB-077	1.8.12.B	17	AB-160	1.8.31.B	2	AB-215	1.8.39.B	9	AB-395	1.8.27.B	10
AB-078	1.8.19.B	26	AB-161	1.8.31.B	2	AB-216	1.8.37.B	3	AB-396	1.8.18.B	15
AB-078	1.8.20.B	26	AB-162	1.8.31.B	2	AB-217	1.8.39.B	10	AB-398	0.6.37.B	6
AB-100	3.2.02.B	1	AB-163	1.8.31.B	2	AB-218	1.8.39.B	11	AB-398	1.8.06.B	26
AB-104	3.2.02.B	8	AB-164	1.8.31.B	2	AB-219	1.8.37.B	9	AB-400	1.8.01.B	1
AB-105	1.8.05.B	16	AB-165	1.8.31.B	2	AB-219	1.8.40.B	2	AB-404	1.8.01.B	8
AB-106	1.8.08.B	5	AB-166	1.8.31.B	2	AB-220	1.8.37.B	10	AB-406	1.8.01.B	10
AB-107	1.8.08.B	8	AB-167	1.8.36.B	6	AB-220	1.8.40.B	3	AB-436	1.8.09.B	8
AB-108	1.8.08.B	9	AB-168	1.8.32.B	8	AB-222	1.8.41.B	1	AB-438	1.8.09.B	10
AB-109	1.8.08.B	17	AB-174	1.8.32.B	20	AB-223	1.8.41.B	2	AB-441	1.8.25.B	17
AB-111	1.8.09.B	8	AB-174	1.8.33.B	14	AB-224	1.8.41.B	3	AB-442	1.8.25.B	16
AB-112	1.8.09.B	10	AB-175	1.8.32.B	22	AB-225	1.8.41.B	7	AB-449	1.8.25.B	19
AB-113	1.8.09.B	4	AB-176	1.8.32.B	21	AB-226	1.8.41.B	7	AB-451	1.8.01.B	5
AB-114	1.8.09.B	5	AB-177	1.8.32.B	26	AB-228	1.8.19.B	20	AB-469	1.8.01.B	5
AB-120	1.8.11.B	13	AB-178	1.8.32.B	26	AB-229	1.8.28.B	6	AB-487	1.8.38.B	2
AB-121	1.8.11.B	36	AB-179	1.8.32.B	28	AB-230	1.8.28.B	8	AB-494	1.8.38.B	2
AB-122	1.8.11.B	37	AB-180	1.8.32.B	28	AB-231	1.8.37.B	8	AB-512	1.8.38.B	2
AB-123	1.8.11.B	10	AB-181	1.8.32.B	33	AB-231	1.8.40.B	1	AB-514	1.8.38.B	4
AB-124	1.8.11.B	11	AB-182	1.8.32.B	34	AB-232	1.8.40.B	12	AB-516	1.8.38.B	4
AB-125	1.8.11.B	24	AB-183	1.8.32.B	41	AB-233	1.8.40.B	13	AB-526	1.8.38.B	4
AB-126	1.8.11.B	26	AB-186	1.8.32.B	44	AB-234	1.8.40.B	18	AB-527	1.8.32.B	36
AB-128	1.8.11.B	29	AB-187	1.8.01.B	1	AB-235	1.8.40.B	18	AB-528	1.8.32.B	36
AB-129	1.8.11.B	32	AB-188	1.8.01.B	5	AB-236	1.8.40.B	18	AB-529	3.2.01.B	6
AB-130	1.8.11.B	33	AB-189	1.8.01.B	5	AB-237	1.8.40.B	18	AB-532	1.8.10.B	30
AB-131	1.8.12.B	11	AB-190	1.8.02.B	18	AB-238	1.8.29.B	2	AB-536	1.8.28.B	33
AB-132	1.8.12.B	11	AB-191	1.8.02.B	19	AB-239	1.8.29.B	4	AB-542	1.8.33.B	24
AB-133	1.8.12.B	12	AB-192	1.8.02.B	21	AB-240	1.8.29.B	22	AB-550	1.8.31.B	2
AB-134	1.8.12.B	12	AB-193	1.8.02.B	15	AB-242	1.8.33.B	24	AB-556	1.8.31.B	2
AB-135	1.8.12.B	13	AB-194	1.8.03.B	1	AB-244	1.8.27.B	1	AB-560	1.8.31.B	2
AB-136	1.8.12.B	13	AB-195	1.8.19.B	3	AB-246	1.8.29.B	41	AB-561	1.8.29.B	28
AB-137	1.8.13.B	12	AB-196	1.8.29.B	43	AB-253	1.8.29.B	42	AB-564	1.8.31.B	2

NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No
AB-565	1.8.29.B	23	AB-762	1.8.27.B	3	AB-906	1.8.01.B	10	AC-117	0.6.03.B	18
AB-568	1.8.31.B	2	AB-766	1.8.10.B	33	AB-908	1.8.38.B	20	AC-117	0.6.17.B	16
AB-569	1.8.29.B	5	AB-766	1.8.27.B	4	AB-913	1.8.18.B	9	AC-117	0.6.18.B	11
AB-572	1.8.31.B	2	AB-771	1.8.22.B	23	AB-914	3.2.02.B	8	AC-121	0.6.08.B	11
AB-575	1.8.29.B	43	AB-772	1.8.22.B	23	AB-936	1.8.09.B	8	AC-130	0.6.07.B	8
AB-576	1.8.31.B	2	AB-777	1.8.22.B		AB-938	1.8.09.B	10	AC-131	0.6.07.B	15
AB-585	1.8.41.B	1	AB-778	1.8.22.B		AB-944	1.8.38.B	21	AC-132	0.6.07.B	16
AB-590	1.8.41.B	7	AB-791	7.3.01.B	29	AB-955	1.8.39.B	2	AC-137	0.6.02.B	65
AB-599	1.8.40.B	18	AB-794	7.3.01.B	30	AB-956	1.8.39.B	3	AC-138	0.6.02.B	23
AB-612	1.8.32.B	22	AB-796	1.8.10.B	30	AB-957	1.8.39.B	4	AC-141	0.6.07.B	34
AB-617	1.8.24.B		AB-806	1.8.38.B	5	AB-958	1.8.39.B	6	AC-142	0.6.13.B	2
AB-621	1.8.24.B		AB-807	1.8.38.B	5	AB-959	1.8.39.B	9	AC-143	0.6.05.B	1
AB-622	1.8.24.B		AB-808	1.8.38.B	5	AB-960	1.8.39.B	10	AC-143	0.6.32.B	7
AB-623	1.8.24.B		AB-809	1.8.38.B	7	AB-962	1.8.22.B	18	AC-144	0.6.20.B	
AB-624	1.8.29.B	42	AB-810	1.8.38.B	7	AB-964	1.8.39.B	11	AC-146	0.6.04.B	41
AB-626	1.8.29.B	40	AB-811	1.8.38.B	7	AB-978	1.8.26.B	16	AC-147	0.6.09.B	11
AB-627	1.8.29.B	40	AB-812	1.8.38.B	6	AB-978	1.8.28.B	26	AC-152	0.6.04.B	24
AB-629	1.8.29.B	41	AB-813	1.8.38.B	6	AC-009	1.8.33.B	25	AC-153	0.6.04.B	4
AB-635	1.8.28.B	18	AB-814	1.8.38.B	6	AC-009	1.8.35.B	7	AC-154	0.6.04.B	5
AB-640	1.8.25.B	11	AB-815	1.8.38.B	8	AC-011	0.6.10.B	17	AC-159	0.6.04.B	27
AB-641	1.8.28.B	19	AB-816	1.8.38.B	8	AC-016	1.8.38.B	18	AC-160	0.6.02.B	1
AB-645	1.8.32.B	2	AB-817	1.8.38.B	8	AC-024	1.8.11.B	38	AC-162	0.6.17.B	24
AB-647	1.8.41.B	11	AB-819	1.8.01.B	8	AC-025	0.6.07.B	3	AC-164	0.6.10.B	9
AB-650	1.8.32.B	2	AB-825	1.8.32.B	36	AC-026	0.6.07.B	6	AC-164	0.6.38.B	15
AB-661	4.3.01.B	4	AB-829	1.8.29.B	5	AC-033	0.6.04.B	22	AC-165	0.6.10.B	19
AB-665	1.8.32.B	12	AB-831	1.8.29.B	28	AC-035	0.6.09.B	20	AC-165	0.6.38.B	16
AB-667	0.6.25.B	34	AB-832	1.8.29.B	23	AC-035	0.6.99.B		AC-166	0.6.02.B	4
AB-667	7.3.01.B	34	AB-833	1.8.29.B	9	AC-040	0.6.04.B	40	AC-167	0.6.02.B	5
AB-669	1.8.01.B	5	AB-834	1.8.29.B	42	AC-041	0.6.02.B	64	AC-168	0.6.02.B	6
AB-677	6.3.01.B	1	AB-835	1.8.29.B	40	AC-043	0.6.06.B	21	AC-170	0.6.18.B	3
AB-679	6.3.01.B		AB-836	1.8.29.B	41	AC-066	0.6.34.B	1	AC-175	0.6.08.B	3
AB-680	1.8.32.B	13	AB-839	1.8.29.B	43	AC-067	0.6.34.B	2	AC-176	0.6.13.B	1
AB-682	6.3.01.B	12	AB-843	1.8.29.B	45	AC-076	0.6.10.B	21	AC-179	0.6.13.B	19
AB-688	1.8.04.B	7	AB-853	1.8.41.B	2	AC-079	0.6.28.B	2	AC-181	0.6.13.B	5
AB-691	1.8.09.B	4	AB-857	1.8.41.B	3	AC-079	0.6.99.B		AC-182	0.6.05.B	2
AB-702	1.8.15.B	12	AB-858	1.8.41.B	1	AC-081	0.6.09.B	9	AC-183	0.6.05.B	11
AB-705	1.8.29.B	18	AB-871	1.8.40.B	13	AC-097	0.6.13.B	17	AC-187	0.6.09.B	21
AB-705	1.8.33.B	8	AB-875	1.8.10.B	14	AC-103	0.6.10.B	3	AC-188	0.6.24.B	18
AB-710	1.8.11.B	28	AB-878	1.8.16.B	10	AC-103	0.6.99.B		AC-209	0.6.01.B	21
AB-715	1.8.11.B	28	AB-879	1.8.10.B	30	AC-104	0.6.09.B	6	AC-213	0.6.01.B	18
AB-720	1.8.11.B	28	AB-885	1.8.04.B	1	AC-105	0.6.04.B	17	AC-214	0.6.01.B	18
AB-725	1.8.11.B	28	AB-893	0.6.23.B	17	AC-106	0.6.04.B	2	AC-215	0.6.01.B	18
AB-727	1.8.10.B	14	AB-893	4.3.01.B	18	AC-107	0.6.04.B	1	AC-216	0.6.01.B	17
AB-729	1.8.10.B	14	AB-895	0.6.23.B	18	AC-108	0.6.04.B	1	AC-217	0.6.01.B	17
AB-730	1.8.11.B	28	AB-895	4.3.01.B	19	AC-110	0.6.04.B	9	AC-218	0.6.01.B	17
AB-738	1.8.15.B	8	AB-900	1.8.01.B	1	AC-111	0.6.04.B	10	AC-219	0.6.01.B	21
AB-739	1.8.15.B	8	AB-901	1.8.40.B	12	AC-112	0.6.00.B	9	AC-220	0.6.01.B	21
AB-741	1.8.32.B	2	AB-902	1.8.41.B	7	AC-114	0.6.10.B	4	AC-224	0.6.07.B	18
AB-758	1.8.29.B	38	AB-903	1.8.41.B	7	AC-115	0.6.10.B	6	AC-225	0.6.07.B	18
AB-762	1.8.10.B	32	AB-904	1.8.05.B	10	AC-115	0.6.38.B	14	AC-226	0.6.07.B	18

NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No

AC-227	0.6.07.B	23	AC-369	0.6.04.B	15	AC-521	0.6.32.B	13	AC-795	0.6.07.B	8
AC-228	0.6.07.B	23	AC-370	0.6.04.B	3	AC-521	0.6.99.B		AC-796	0.6.07.B	15
AC-229	0.6.07.B	23	AC-371	0.6.09.B	4	AC-522	0.6.06.B	9	AC-808	0.6.27.B	19
AC-236	0.6.00.B	9	AC-371	0.6.17.B	1	AC-524	0.6.02.B	28	AC-808	0.6.99.B	
AC-237	0.6.00.B	9	AC-371	0.6.18.B	1	AC-525	0.6.02.B	29	AC-834	0.6.28.B	3
AC-238	0.6.00.B	9	AC-371	0.6.35.B	11	AC-527	0.6.08.B	9	AC-834	0.6.99.B	
AC-239	0.6.00.B	15	AC-371	0.6.39.B	1	AC-527	0.6.17.B	40	AC-836	0.6.37.B	2
AC-240	0.6.00.B	15	AC-378	0.6.13.B	20	AC-528	0.6.13.B	10	AC-870	0.6.27.B	33
AC-241	0.6.00.B	15	AC-379	0.6.13.B	21	AC-529	0.6.13.B	11	AC-871	0.6.27.B	34
AC-242	0.6.00.B	14	AC-380	0.6.34.B	3	AC-530	0.6.09.B	6	AC-872	0.6.27.B	35
AC-243	0.6.00.B	14	AC-384	0.6.12.B	7	AC-533	0.6.09.B	10	AC-884	0.6.09.B	
AC-244	0.6.00.B	14	AC-387	0.6.09.B	21	AC-541	0.6.17.B	17	AC-886	0.6.20.B	
AC-248	0.6.12.B	5	AC-388	0.6.17.B	14	AC-541	0.6.18.B	16	AC-892	0.6.24.B	16
AC-249	0.6.02.B	24	AC-401	0.6.12.B	4	AC-542	0.6.11.B	1	AC-893	0.6.09.B	20
AC-250	0.6.02.B	63	AC-402	0.6.12.B	6	AC-566	0.6.11.B	8	AC-894	0.6.01.B	18
AC-252	0.6.17.B	7	AC-404	0.6.18.B	4	AC-568	0.6.10.B	18	AC-895	0.6.01.B	17
AC-253	0.6.07.B	18	AC-405	0.6.11.B	8	AC-568	0.6.99.B		AC-897	0.6.07.B	23
AC-257	0.6.17.B	15	AC-415	0.6.07.B	35	AC-586	0.6.12.B	2	AC-899	0.6.01.B	21
AC-259	0.6.05.B	13	AC-416	0.6.18.B	10	AC-589	0.6.02.B	59	AC-901	0.6.07.B	18
AC-260	0.6.05.B	16	AC-420	0.6.04.B	42	AC-590	0.6.00.B	16	AC-902	0.6.01.B	4
AC-262	0.6.23.B	1	AC-421	0.6.18.B	10	AC-591	0.6.06.B	8	AC-903	0.6.00.B	14
AC-262	0.6.34.B	4	AC-423	0.6.04.B	39	AC-595	0.6.09.B	20	AC-904	0.6.00.B	15
AC-262	0.6.99.B		AC-424	0.6.20.B	21	AC-599	0.6.08.B	12	AC-938	0.6.00.B	9
AC-264	0.6.02.B	47	AC-424	0.6.29.B	17	AC-599	0.6.09.B	23	AC-965	0.6.12.B	21
AC-265	0.6.01.B	4	AC-429	0.6.00.B	1	AC-601	0.6.02.B	44	AC-967	0.6.09.B	9
AC-266	0.6.01.B	4	AC-429	0.6.03.B	10	AC-622	0.6.01.B	21	AC-985	0.6.17.B	7
AC-267	0.6.01.B	4	AC-429	0.6.29.B	18	AC-626	0.6.17.B	31	AD-001	1.8.03.B	17
AC-268	0.6.01.B	4	AC-432	0.6.11.B		AC-627	0.6.17.B	36	AD-002	1.8.03.B	12
AC-275	0.6.08.B	1	AC-433	0.6.11.B	6	AC-628	0.6.06.B	5	AD-005	1.8.12.B	2
AC-276	0.6.08.B	2	AC-436	0.6.12.B	20	AC-630	0.6.06.B	14	AD-005	1.8.13.B	2
AC-279	0.6.04.B	16	AC-441	0.6.26.B	4	AC-632	0.6.10.B	30	AD-005	1.8.15.B	1
AC-280	0.6.04.B	15	AC-442	0.6.26.B	5	AC-640	0.6.23.B	1	AD-005	1.8.16.B	1
AC-282	6.3.01.B	1	AC-443	0.6.26.B	6	AC-640	0.6.34.B	6	AD-005	1.8.19.B	25
AC-283	6.3.01.B	2	AC-447	0.6.12.B	1	AC-640	0.6.99.B		AD-005	1.8.25.B	21
AC-285	6.3.01.B	12	AC-451	0.6.12.B	2	AC-644	0.6.34.B	5	AD-005	1.8.26.B	15
AC-286	6.3.01.B	5	AC-457	0.6.13.B	6	AC-664	0.6.17.B	34	AD-006	1.8.19.B	5
AC-287	6.3.01.B		AC-467	0.6.00.B	3	AC-712	0.6.04.B	32	AD-006	1.8.21.B	2
AC-293	0.6.09.B	1	AC-468	0.6.00.B	4	AC-718	0.6.21.B	5	AD-006	1.8.24.B	6
AC-316	0.6.08.B	11	AC-469	0.6.00.B	7	AC-720	0.6.21.B	8	AD-006	1.8.25.B	13
AC-317	0.6.02.B	43	AC-469	0.6.12.B	18	AC-738	0.6.03.B	12	AD-007	1.8.15.B	2
AC-329	0.6.17.B	5	AC-479	0.6.01.B	19	AC-741	0.6.17.B	35	AD-008	1.8.12.B	1
AC-331	0.6.07.B	23	AC-480	0.6.01.B	19	AC-743	0.6.13.B	17	AD-008	1.8.26.B	14
AC-332	0.6.00.B	14	AC-481	0.6.01.B	19	AC-745	0.6.10.B	27	AD-009	1.8.20.B	7
AC-333	0.6.00.B	15	AC-482	0.6.01.B	19	AC-756	0.6.08.B		AD-012	1.8.04.B	6
AC-336	0.6.00.B	16	AC-483	0.6.01.B	20	AC-761	0.6.21.B	15	AD-015	1.8.05.B	1
AC-345	0.6.02.B	2	AC-484	0.6.02.B	22	AC-763	0.6.00.B	2	AD-018	1.8.28.B	1
AC-352	0.6.06.B	2	AC-517	0.6.07.B	29	AC-763	0.6.01.B	3	AD-019	1.8.13.B	1
AC-366	0.6.19.B	1	AC-518	0.6.07.B	30	AC-771	0.6.09.B	20	AD-025	1.8.05.B	2
AC-367	0.6.04.B	34	AC-520	0.6.07.B	32	AC-771	0.6.26.B	1	AD-027	1.8.18.B	14
AC-368	0.6.04.B	16	AC-521	0.6.05.B	10	AC-775	0.6.02.B	30	AD-028	1.8.22.B	1

NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No

AD-029	1.8.23.B	14	AD-134	1.8.12.B	19	AD-462	0.6.42.B	1	AE-238	0.6.30.B	38
AD-038	1.8.23.B	15	AD-134	1.8.13.B	25	AD-462	0.6.99.B		AE-239	0.6.03.B	13
AD-047	1.8.04.B	8	AD-134	1.8.28.B	39	AD-466	0.6.42.B	2	AE-239	0.6.06.B	12
AD-049	1.8.03.B	11	AD-135	1.8.03.B	12	AD-466	0.6.99.B		AE-243	0.6.06.B	21
AD-051	1.8.19.B	21	AD-136	1.8.03.B	12	AD-467	0.6.43.B	1	AE-245	0.6.09.B	26
AD-052	1.8.20.B	1	AD-228	1.8.22.B	21	AD-467	0.6.47.B	1	AE-245	0.6.16.B	21
AD-054	1.8.27.B	5	AD-402	0.6.09.B	2	AD-467	0.6.48.B	1	AE-246	0.6.08.B	11
AD-058	1.8.07.B	18	AD-402	0.6.99.B		AD-468	0.6.42.B	3	AE-250	0.6.35.B	1
AD-060	1.8.03.B	12	AD-403	0.6.09.B	12	AD-471	0.6.42.B	3	AE-259	0.6.32.B	1
AD-066	1.8.41.B	14	AD-403	0.6.35.B	2	AD-575	1.8.08.B	38	AE-260	0.6.32.B	2
AD-067	1.8.41.B	14	AD-403	0.6.39.B	8	AD-792	7.3.01.B	31	AE-261	0.6.32.B	8
AD-068	1.8.41.B	14	AD-403	0.6.99.B		AE-013	0.6.06.B	15	AE-262	0.6.32.B	3
AD-069	1.8.16.B	2	AD-404	0.6.09.B	8	AE-031	0.6.38.B	3	AE-263	0.6.09.B	31
AD-078	1.8.19.B	24	AD-404	0.6.99.B		AE-033	0.6.38.B	7	AE-263	0.6.31.B	25
AD-089	1.8.25.B	23	AD-408	0.6.09.B	7	AE-034	0.6.02.B	62	AE-264	0.6.09.B	26
AD-102	0.6.25.B	25	AD-408	0.6.99.B		AE-045	0.6.38.B	8	AE-264	0.6.16.B	21
AD-102	7.3.01.B	25	AD-409	0.6.07.B	25	AE-046	0.6.38.B	2	AE-265	0.6.10.B	15
AD-103	0.6.25.B	17	AD-409	0.6.99.B		AE-047	0.6.38.B	1	AE-271	0.6.09.B	20
AD-103	7.3.01.B	17	AD-412	0.6.03.B	6	AE-064	0.6.13.B	15	AE-271	0.6.99.B	
AD-104	1.8.18.B	7	AD-412	0.6.99.B		AE-064	0.6.46.B	1	AE-272	0.6.30.B	2
AD-105	0.6.25.B	2	AD-413	0.6.03.B	6	AE-071	0.6.03.B	12	AE-273	0.6.03.B	14
AD-105	7.3.01.B	2	AD-413	0.6.99.B		AE-073	0.6.03.B	18	AE-274	0.6.21.B	9
AD-106	0.6.25.B	3	AD-416	0.6.05.B	7	AE-073	0.6.17.B	16	AE-275	0.6.21.B	12
AD-106	7.3.01.B	3	AD-416	0.6.32.B	9	AE-073	0.6.18.B	11	AE-276	0.6.09.B	28
AD-108	0.6.25.B	1	AD-417	0.6.19.B	2	AE-075	0.6.03.B	2	AE-277	0.6.09.B	38
AD-108	7.3.01.B	1	AD-417	0.6.99.B		AE-079	0.6.38.B	4	AE-281	0.6.09.B	41
AD-109	0.6.25.B	16	AD-419	0.6.02.B	25	AE-079	0.6.99.B		AE-287	0.6.03.B	2
AD-109	7.3.01.B	16	AD-419	0.6.99.B		AE-081	0.6.27.B	18	AE-303	0.6.28.B	2
AD-110	0.6.25.B	24	AD-422	0.6.11.B	5	AE-103	0.6.02.B		AE-303	0.6.99.B	
AD-110	7.3.01.B	24	AD-423	0.6.09.B	3	AE-105	0.6.01.B	19	AE-304	0.6.28.B	3
AD-111	0.6.25.B	26	AD-423	0.6.99.B		AE-107	0.6.03.B	9	AE-304	0.6.99.B	
AD-111	7.3.01.B	26	AD-426	0.6.26.B	3	AE-108	0.6.03.B	9	AE-305	0.6.28.B	2
AD-116	1.8.09.B	17	AD-427	0.6.30.B	24	AE-110	0.6.13.B	17	AE-305	0.6.99.B	
AD-121	1.8.03.B	19	AD-427	0.6.99.B		AE-111	0.6.13.B	5	AE-306	0.6.28.B	3
AD-121	1.8.05.B	22	AD-428	0.6.30.B	25	AE-124	0.6.04.B	18	AE-306	0.6.99.B	
AD-127	1.8.04.B	2	AD-428	0.6.99.B		AE-129	0.6.13.B	4	AE-308	0.6.28.B	4
AD-128	1.8.28.B	23	AD-436	0.6.30.B	24	AE-129	0.6.99.B		AE-312	0.6.27.B	20
AD-129	1.8.19.B	23	AD-439	0.6.03.B	6	AE-160	0.6.39.B	11	AE-322	0.6.08.B	8
AD-129	1.8.28.B	24	AD-443	0.6.31.B	24	AE-161	0.6.39.B	12	AE-331	0.6.01.B	18
AD-130	1.8.16.B	12	AD-443	0.6.99.B		AE-162	0.6.35.B	16	AE-332	0.6.01.B	18
AD-130	1.8.28.B	35	AD-443	0.6.99.B		AE-162	0.6.39.B	13	AE-333	0.6.01.B	18
AD-131	1.8.22.B	2	AD-444	0.6.09.B	8	AE-165	0.6.13.B	27	AE-334	0.6.01.B	18
AD-131	1.8.28.B	36	AD-444	0.6.99.B		AE-166	0.6.13.B		AE-335	0.6.01.B	18
AD-132	1.8.19.B	23	AD-446	0.6.03.B	17	AE-194	0.6.02.B		AE-336	0.6.01.B	18
AD-132	1.8.28.B	37	AD-447	0.6.03.B	4	AE-203	0.6.03.B	9	AE-337	0.6.01.B	17
AD-133	1.8.12.B	18	AD-447	0.6.03.B	17	AE-217	0.6.23.B	18	AE-338	0.6.01.B	17
AD-133	1.8.13.B	24	AD-447	0.6.99.B		AE-217	0.6.25.B	34	AE-339	0.6.01.B	17
AD-133	1.8.15.B	14	AD-447	0.6.99.B		AE-226	0.6.27.B	19	AE-340	0.6.01.B	17
AD-133	1.8.16.B	13	AD-448	0.6.02.B	67	AE-232	0.6.09.B	14	AE-341	0.6.01.B	17
AD-133	1.8.28.B	38	AD-448	0.6.03.B	15	AE-237	0.6.03.B	9	AE-342	0.6.01.B	17

NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No

AE-349	0.6.07.B	23	AE-658	0.6.00.B	14	AE-779	0.6.47.B	7	AE-981	0.6.17.B	5
AE-350	0.6.07.B	23	AE-659	0.6.00.B	14	AE-780	0.6.47.B	8	AE-982	0.6.17.B	5
AE-351	0.6.07.B	23	AE-660	0.6.00.B	14	AE-781	0.6.48.B	9	AE-993	0.6.49.B	1
AE-352	0.6.07.B	23	AE-661	0.6.00.B	14	AE-782	0.6.47.B	9	AE-995	0.6.99.B	
AE-353	0.6.07.B	23	AE-662	0.6.00.B	15	AE-782	0.6.99.B		AE-996	0.6.99.B	
AE-354	0.6.07.B	23	AE-663	0.6.00.B	15	AE-785	0.6.43.B	6	AE-998	0.6.99.B	
AE-361	0.6.01.B	21	AE-664	0.6.00.B	15	AE-785	0.6.99.B		AF-015	0.6.42.B	6
AE-362	0.6.01.B	21	AE-665	0.6.00.B	15	AE-786	0.6.43.B	7	AF-017	0.6.42.B	6
AE-363	0.6.01.B	21	AE-666	0.6.00.B	15	AE-789	0.6.43.B	8	AF-019	0.6.42.B	6
AE-364	0.6.01.B	21	AE-667	0.6.00.B	15	AE-789	0.6.99.B		AF-027	0.6.42.B	6
AE-365	0.6.01.B	21	AE-674	0.6.01.B	16	AE-790	0.6.02.B		AF-029	0.6.42.B	6
AE-366	0.6.01.B	21	AE-681	0.6.03.B	9	AE-792	0.6.46.B	4	AF-032	0.6.42.B	6
AE-373	0.6.07.B	18	AE-694	0.6.17.B	34	AE-793	0.6.43.B	9	AG-100	1.8.05.B	20
AE-374	0.6.07.B	18	AE-702	0.6.00.B	16	AE-795	0.6.44.B	2	AG-102	1.8.05.B	9
AE-375	0.6.07.B	18	AE-737	0.6.28.B	2	AE-795	0.6.99.B		AG-107	1.8.05.B	9
AE-376	0.6.07.B	18	AE-737	0.6.99.B		AE-908	0.6.03.B	9	AG-301	1.8.05.B	3
AE-377	0.6.07.B	18	AE-738	0.6.28.B	2	AE-909	0.6.03.B	9	AG-552	1.8.33.B	19
AE-378	0.6.07.B	18	AE-739	0.6.28.B	3	AE-915	0.6.09.B	40	AG-552	1.8.35.B	21
AE-379	0.6.01.B	4	AE-739	0.6.99.B		AE-915	0.6.16.B		AG-553	1.8.33.B	19
AE-380	0.6.01.B	4	AE-740	0.6.28.B	3	AE-921	0.6.34.B	15	AG-553	1.8.35.B	21
AE-381	0.6.01.B	4	AE-740	0.6.99.B		AE-921	0.6.99.B		AG-555	1.8.33.B	18
AE-382	0.6.01.B	4	AE-745	0.6.22.B	5	AE-927	0.6.45.B	1	AG-556	1.8.33.B	18
AE-383	0.6.01.B	4	AE-751	0.6.42.B	4	AE-927	0.6.99.B		AH-006	1.8.03.B	14
AE-384	0.6.01.B	4	AE-752	0.6.42.B	4	AE-936	0.6.34.B	10	AH-007	1.8.08.B	2
AE-441	0.6.00.B	8	AE-753	0.6.09.B	20	AE-937	0.6.34.B	11	AH-008	1.8.09.B	17
AE-442	0.6.01.B	5	AE-753	0.6.42.B	5	AE-938	0.6.34.B	14	AH-010	1.8.12.B	6
AE-443	0.6.01.B	12	AE-753	0.6.99.B		AE-939	0.6.01.B	19	AH-010	1.8.13.B	4
AE-452	0.6.09.B	42	AE-754	0.6.47.B	2	AE-940	0.6.01.B	19	AH-010	1.8.15.B	4
AE-452	0.6.28.B	1	AE-754	0.6.99.B		AE-941	0.6.01.B	19	AH-010	1.8.16.B	4
AE-453	0.6.17.B	15	AE-755	0.6.48.B	2	AE-942	0.6.01.B	19	AH-011	1.8.12.B	7
AE-454	0.6.17.B	14	AE-755	0.6.99.B		AE-943	0.6.01.B	19	AH-011	1.8.13.B	5
AE-461	0.6.00.B	9	AE-756	0.6.43.B	2	AE-944	0.6.08.B	1	AH-011	1.8.15.B	5
AE-462	0.6.00.B	9	AE-756	0.6.47.B	3	AE-945	0.6.08.B	2	AH-011	1.8.16.B	6
AE-463	0.6.00.B	9	AE-756	0.6.48.B	3	AE-946	0.6.09.B	6	AH-013	1.8.12.B	5
AE-464	0.6.00.B	9	AE-757	0.6.46.B	2	AE-947	0.6.17.B	7	AH-013	1.8.13.B	7
AE-465	0.6.00.B	9	AE-758	0.6.46.B	3	AE-948	0.6.17.B	7	AH-013	1.8.15.B	6
AE-466	0.6.00.B	9	AE-759	0.6.48.B	4	AE-949	0.6.01.B	19	AH-013	1.8.16.B	7
AE-581	0.6.18.B	10	AE-759	0.6.99.B		AE-954	0.6.41.B	1	AH-014	1.8.18.B	1
AE-582	0.6.18.B	10	AE-760	0.6.44.B	1	AE-957	0.6.40.B	1	AH-015	1.8.18.B	7
AE-583	0.6.17.B	5	AE-760	0.6.99.B		AE-958	0.6.40.B	2	AH-016	1.8.20.B	3
AE-584	0.6.17.B	36	AE-761	0.6.43.B	3	AE-966	0.6.40.B	3	AH-017	1.8.20.B	4
AE-585	0.6.17.B	31	AE-761	0.6.47.B	4	AE-969	0.6.46.B	5	AH-018	1.8.20.B	4
AE-586	0.6.17.B	35	AE-761	0.6.48.B	5	AE-972	0.6.10.B	8	AH-019	1.8.21.B	6
AE-600	0.6.18.B	4	AE-762	0.6.43.B	4	AE-973	0.6.10.B	14	AH-020	1.8.21.B	6
AE-610	0.6.00.B	2	AE-762	0.6.47.B	5	AE-974	0.6.46.B	6	AH-021	1.8.22.B	5
AE-610	0.6.01.B	3	AE-762	0.6.48.B	6	AE-975	0.6.46.B	7	AH-022	1.8.23.B	12
AE-615	0.6.09.B	43	AE-763	0.6.43.B	5	AE-975	0.6.99.B		AH-023	1.8.29.B	31
AE-616	0.6.00.B	13	AE-763	0.6.47.B	6	AE-976	0.6.46.B	8	AH-024	1.8.31.B	10
AE-656	0.6.00.B	14	AE-763	0.6.48.B	7	AE-976	0.6.99.B		AH-025	1.8.03.B	10
AE-657	0.6.00.B	14	AE-778	0.6.48.B	8	AE-977	0.6.46.B	9	AH-026	1.8.19.B	1

NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No

AH-027	1.8.21.B	1	AK-063	0.6.17.B	6	FA-064	1.8.22.B	7	FA-145	4.3.01.B	7
AH-029	1.8.33.B	7	AK-063	1.8.01.B	15	FA-065	1.8.22.B	9	FA-146	0.6.23.B	6
AH-030	1.8.35.B	12	AK-064	0.6.17.B	6	FA-067	0.6.23.B	6	FA-146	1.8.40.B	4
AH-031	1.8.09.B	17	AK-064	1.8.01.B	15	FA-067	1.8.40.B	5	FA-146	4.3.01.B	7
AH-036	0.6.37.B	1	AK-065	0.6.17.B	6	FA-069	0.6.34.B	7	FA-147	0.6.23.B	6
AH-036	1.8.06.B	27	AK-065	1.8.01.B	15	FA-069	1.8.37.B	4	FA-147	1.8.40.B	4
AH-036	1.8.07.B	19	AK-301	1.8.01.B	14	FA-069	4.3.01.B	21	FA-147	4.3.01.B	7
AH-041	1.8.12.B	4	BA-061	0.6.11.B		FA-070	0.6.34.B	8	FA-148	0.6.23.B	6
AH-041	1.8.13.B	6	BA-062	0.6.11.B		FA-070	1.8.37.B	5	FA-148	1.8.40.B	4
AH-041	1.8.15.B	13	BA-082	0.6.06.B	20	FA-073	1.8.23.B	1	FA-148	4.3.01.B	7
AH-041	1.8.16.B	5	BA-083	0.6.06.B	19	FA-074	1.8.23.B	13	FA-152	3.2.02.B	6
AH-044	1.8.25.B	22	BA-085	0.6.06.B	17	FA-083	1.8.30.B	2	FA-155	1.8.28.B	32
AH-119	1.8.21.B	6	BA-086	0.6.06.B	18	FA-084	1.8.30.B	6	FA-156	1.8.18.B	10
AH-121	6.3.01.B	6	BA-087	0.6.06.B	16	FA-085	1.8.30.B	7	FA-157	1.8.39.B	12
AH-209	1.8.09.B	18	BA-087	0.6.26.B	13	FA-086	1.8.30.B	11	FA-158	1.8.39.B	13
AH-209	1.8.29.B	49	BA-135	0.6.30.B	3	FA-087	1.8.30.B	9	FA-159	1.8.39.B	14
AH-214	1.8.18.B	1	BA-136	0.6.30.B	3	FA-088	1.8.30.B	8	FA-160	1.8.32.B	3
AH-242	1.8.20.B	21	BA-137	0.6.20.B		FA-091	1.8.10.B	5	FA-161	1.8.32.B	11
AH-314	1.8.18.B	1	BA-137	0.6.29.B		FA-096	1.8.32.B	1	FA-162	1.8.31.B	13
AK-001	1.8.01.B	12	BA-341	0.6.36.B	21	FA-099	1.8.29.B	34	FA-166	1.8.21.B	12
AK-032	1.8.01.B	12	FA-010	1.8.01.B	18	FA-103	1.8.33.B	5	FA-167	1.8.25.B	20
AK-033	1.8.01.B	13	FA-019	1.8.29.B	8	FA-104	1.8.33.B	3	FA-168	1.8.27.B	11
AK-034	1.8.01.B	14	FA-022	1.8.40.B	10	FA-106	1.8.40.B	19	FA-172	1.8.34.B	24
AK-035	1.8.01.B	12	FA-023	1.8.40.B	11	FA-107	1.8.05.B	15	FA-174	1.8.27.B	12
AK-037	1.8.01.B	14	FA-032	1.8.05.B	21	FA-109	1.8.33.B	2	FA-175	1.8.27.B	13
AK-049	0.6.17.B	6	FA-038	1.8.07.B	13	FA-109	1.8.35.B	22	FA-176	1.8.11.B	29
AK-049	1.8.01.B	15	FA-040	1.8.12.B		FA-111	1.8.28.B	9	FA-176	1.8.27.B	14
AK-050	0.6.17.B	6	FA-040	1.8.13.B		FA-113	1.8.28.B	20	FA-179	1.8.27.B	6
AK-050	1.8.01.B	15	FA-040	1.8.14.B		FA-118	1.8.31.B	1	FA-195	1.8.05.B	18
AK-050	3.2.02.B	12	FA-041	0.6.24.B	12	FA-119	1.8.31.B	7	FA-197	1.8.27.B	7
AK-052	0.6.17.B	6	FA-041	1.8.12.B		FA-120	1.8.28.B	10	FA-198	1.8.20.B	6
AK-052	1.8.01.B	15	FA-041	1.8.13.B		FA-121	1.8.28.B	11	FA-201	1.8.40.B	17
AK-053	0.6.17.B	6	FA-041	1.8.14.B		FA-124	1.8.40.B	19	FA-203	1.8.05.B	17
AK-053	1.8.01.B	15	FA-043	1.8.12.B	9	FA-125	1.8.28.B	10	FA-223	1.8.25.B	27
AK-054	0.6.17.B	6	FA-043	1.8.13.B	22	FA-127	1.8.28.B	10	FA-224	1.8.25.B	28
AK-054	1.8.01.B	15	FA-044	1.8.18.B	2	FA-129	1.8.28.B	10	FA-225	1.8.25.B	29
AK-055	0.6.17.B	6	FA-045	1.8.18.B	3	FA-133	1.8.29.B	33	FA-230	1.8.25.B	30
AK-055	1.8.01.B	15	FA-046	1.8.18.B	4	FA-134	1.8.28.B	11	FA-230	1.8.42.B	2
AK-056	0.6.17.B	6	FA-047	1.8.18.B	5	FA-135	1.8.32.B	27	FA-231	1.8.25.B	31
AK-056	1.8.01.B	15	FA-048	1.8.18.B	6	FA-136	1.8.32.B	29	FA-238	0.6.37.B	5
AK-057	0.6.17.B	6	FA-050	1.8.42.B	1	FA-138	1.8.32.B	46	FA-238	1.8.06.B	24
AK-057	1.8.01.B	15	FA-051	1.8.18.B	12	FA-140	1.8.33.B	17	FA-247	1.8.29.B	35
AK-058	0.6.17.B	6	FA-052	1.8.18.B	12	FA-141	1.8.37.B	7	FA-256	1.8.07.B	20
AK-058	1.8.01.B	15	FA-054	1.8.21.B	3	FA-142	1.8.37.B	6	FA-257	1.8.29.B	47
AK-059	0.6.17.B	6	FA-055	1.8.21.B	4	FA-143	1.8.39.B	1	FA-270	6.3.01.B	5
AK-059	1.8.01.B	15	FA-056	1.8.21.B	5	FA-144	0.6.23.B	6	FA-274	6.3.01.B	4
AK-060	0.6.17.B	6	FA-057	1.8.21.B	8	FA-144	1.8.40.B	4	FA-277	6.3.01.B	7
AK-060	1.8.01.B	15	FA-058	1.8.21.B	9	FA-144	4.3.01.B	7	FA-287	0.6.21.B	19
AK-061	0.6.17.B	6	FA-059	1.8.21.B	10	FA-145	0.6.23.B	6	FA-288	0.6.21.B	21
AK-061	1.8.01.B	15	FA-063	1.8.21.B	11	FA-145	1.8.40.B	4	FA-293	1.8.28.B	10

NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No

FA-294	1.8.28.B	11	FA-544	1.8.20.B	22	FH-109	1.8.19.B	19	FS-152	1.8.24.B	
FA-297	1.8.01.B	18	FA-546	1.8.29.B	33	FH-110	0.6.31.B	20	FS-262	1.8.24.B	
FA-299	1.8.28.B	10	FA-551	1.8.27.B	8	FH-110	1.8.19.B	17	HA-003	1.8.31.B	3
FA-308	6.3.01.B	8	FA-557	1.8.29.B	35	FH-113	0.6.31.B	14	HA-006	1.8.13.B	27
FA-309	6.3.01.B	9	FA-558	1.8.40.B	6	FH-113	1.8.19.B	33	HA-006	1.8.16.B	15
FA-318	0.6.23.B	2	FA-560	1.8.40.B	6	FH-115	0.6.31.B	15	HA-006	1.8.28.B	4
FA-318	4.3.01.B	1	FA-569	1.8.18.B	8	FH-115	1.8.19.B	34	HA-007	1.8.01.B	4
FA-319	0.6.25.B	18	FA-570	0.6.23.B	2	FH-116	0.6.31.B	16	HA-009	1.8.41.B	4
FA-319	7.3.01.B	18	FA-572	1.8.39.B	1	FH-116	1.8.19.B	33	HA-010	1.8.05.B	6
FA-336	1.8.28.B	17	FA-605	1.8.28.B	14	FH-117	0.6.31.B	5	HA-011	1.8.29.B	20
FA-337	1.8.29.B	33	FA-614	6.3.03.B	4	FH-117	1.8.19.B	14	HA-011	1.8.34.B	7
FA-344	7.3.01.B	35	FA-624	1.8.28.B	34	FH-118	0.6.31.B	6	HA-012	1.8.29.B	1
FA-346	7.3.01.B	36	FA-627	0.6.23.B	10	FH-118	1.8.19.B	15	HA-013	1.8.06.B	5
FA-355	0.6.23.B	14	FA-627	4.3.01.B	11	FH-119	0.6.31.B	7	HA-013	1.8.42.B	3
FA-355	4.3.01.B	15	FA-638	1.8.22.B	22	FH-119	1.8.19.B	13	HA-016	1.8.01.B	9
FA-361	0.6.23.B	7	FA-639	1.8.22.B	10	FH-120	1.8.19.B	10	HA-018	1.8.01.B	7
FA-361	4.3.01.B	8	FA-640	1.8.22.B	11	FH-122	0.6.31.B	21	HA-019	1.8.01.B	3
FA-362	0.6.23.B	8	FA-642	1.8.08.B	1	FH-123	0.6.31.B	10	HA-022	3.2.02.B	9
FA-362	4.3.01.B	9	FA-643	1.8.08.B	1	FH-124	0.6.31.B	11	HA-023	1.8.01.B	13
FA-363	0.6.23.B	9	FA-652	1.8.25.B	38	FH-125	0.6.31.B	12	HA-028	1.8.01.B	17
FA-363	4.3.01.B	10	FA-653	1.8.25.B	37	FH-126	0.6.31.B	13	HA-037	3.2.02.B	10
FA-367	0.6.21.B	4	FA-665	1.8.08.B	1	FH-127	0.6.31.B	26	HA-038	1.8.01.B	17
FA-391	1.8.27.B	19	FA-742	1.8.08.B	1	FR-300	1.8.10.B	19	HA-039	1.8.01.B	4
FA-392	0.6.25.B	4	FA-752	0.6.32.B	14	FR-301	1.8.10.B	19	HA-044	1.8.31.B	5
FA-392	7.3.01.B	4	FA-752	0.6.46.B	10	FR-302	1.8.10.B	19	HA-045	1.8.01.B	3
FA-395	0.6.25.B	19	FB-001	1.8.10.B	16	FR-320	1.8.10.B	21	HA-054	1.8.01.B	13
FA-395	7.3.01.B	19	FB-101	1.8.10.B	16	FR-321	1.8.10.B	21	HA-057	1.8.09.B	1
FA-396	0.6.25.B	5	FC-005	0.6.07.B	26	FR-322	1.8.10.B	21	HA-058	1.8.09.B	1
FA-396	7.3.01.B	5	FE-001	1.8.24.B		FR-350	1.8.10.B	18	HA-059	1.8.09.B	2
FA-400	0.6.24.B	13	FE-002	1.8.24.B		FR-351	1.8.10.B	18	HA-063	1.8.09.B	1
FA-402	1.8.09.B	16	FE-003	1.8.24.B		FR-352	1.8.10.B	18	HA-067	1.8.09.B	2
FA-403	1.8.09.B	16	FE-004	1.8.24.B		FR-370	1.8.10.B	20	HA-068	1.8.09.B	2
FA-408	6.3.01.B	10	FE-101	1.8.24.B		FR-371	1.8.10.B	20	HA-074	1.8.09.B	3
FA-409	6.3.01.B	11	FE-103	1.8.24.B		FR-372	1.8.10.B	20	HA-083	1.8.02.B	1
FA-418	1.8.31.B	1	FH-004	1.8.13.B	26	FS-003	1.8.15.B	11	HA-091	1.8.29.B	7
FA-419	1.8.29.B	8	FH-004	1.8.16.B	14	FS-004	1.8.16.B	8	HA-095	1.8.09.B	3
FA-421	1.8.31.B	7	FH-004	1.8.28.B	3	FS-005	1.8.13.B	8	HA-103	6.3.03.B	6
FA-425	1.8.32.B	1	FH-099	0.6.31.B	1	FS-025	1.8.15.B	11	HA-104	6.3.03.B	7
FA-430	1.8.32.B	1	FH-101	1.8.19.B	6	FS-035	1.8.16.B	8	HA-105	6.3.03.B	5
FA-445	1.8.32.B	14	FH-104	0.6.31.B	17	FS-057	1.8.19.B	32	HA-107	6.3.03.B	1
FA-446	1.8.32.B	15	FH-104	1.8.19.B	12	FS-067	1.8.20.B	19	HA-108	6.3.03.B	2
FA-447	0.6.25.B	37	FH-105	0.6.31.B	2	FS-070	1.8.19.B	30	HA-109	6.3.03.B	3
FA-447	7.3.01.B	37	FH-105	1.8.19.B	10	FS-071	1.8.19.B	31	HA-115	1.8.11.B	23
FA-448	0.6.25.B	38	FH-106	0.6.31.B	3	FS-075	1.8.13.B	8	HA-123	1.8.09.B	15
FA-448	7.3.01.B	38	FH-106	1.8.19.B	11	FS-079	1.8.25.B	11	HA-125	1.8.41.B	6
FA-455	1.8.21.B	7	FH-107	0.6.31.B	4	FS-080	1.8.27.B	9	HA-125	6.3.01.B	3
FA-491	1.8.10.B	5	FH-107	1.8.19.B	13	FS-085	1.8.19.B	32	HA-138	1.8.41.B	12
FA-538	0.6.21.B	25	FH-108	0.6.31.B	18	FS-086	1.8.16.B	9	HA-139	1.8.09.B	15
FA-542	1.8.01.B	18	FH-108	1.8.19.B	18	FS-087	1.8.16.B	9	HA-140	1.8.01.B	17
FA-543	1.8.20.B	23	FH-109	0.6.31.B	19	FS-090	1.8.27.B	9	HA-142	1.8.41.B	4

NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No

HA-144	1.8.38.B	14	HE-005	0.6.45.B	2	HE-125	0.6.43.B	19	HR-034	1.8.10.B	27
HA-147	0.6.08.B	4	HE-005	1.8.17.B	8	HE-126	0.6.43.B	20	HR-035	1.8.10.B	27
HA-150	1.8.29.B	46	HE-024	1.8.24.B	1	HE-127	0.6.43.B	21	HR-041	0.6.35.B	8
HA-153	1.8.01.B	17	HE-025	1.8.24.B	4	HE-128	0.6.43.B	22	HR-041	1.8.42.B	56
HA-160	0.6.02.B	26	HE-026	1.8.24.B	5	HE-129	0.6.43.B	23	HR-042	1.8.10.B	27
HA-160	0.6.99.B		HE-030	1.8.24.B	2	HE-130	0.6.43.B	24	HR-045	1.8.10.B	27
HA-163	5.3.01.B	22	HE-041	1.8.24.B	8	HE-131	0.6.43.B	25	HR-050	1.8.10.B	15
HA-169	1.8.09.B	15	HE-042	1.8.24.B	9	HE-132	0.6.43.B	25	HR-051	1.8.09.B	
HA-173	1.8.09.B	15	HE-043	1.8.24.B	10	HE-133	0.6.48.B	10	HR-051	1.8.34.B	10
HA-175	1.8.09.B	3	HE-044	1.8.24.B	11	HE-134	0.6.48.B	10	HR-051	1.8.41.B	10
HA-179	1.8.39.B	20	HE-045	1.8.24.B	12	HE-135	0.6.48.B	10	HR-052	3.2.01.B	10
HA-200	1.8.10.B	10	HE-046	1.8.24.B	13	HE-138	0.6.45.B	3	HR-055	1.8.25.B	4
HA-202	1.8.10.B	26	HE-047	1.8.24.B	14	HE-139	0.6.45.B	4	HR-110	1.8.37.B	2
HA-211	1.8.01.B	17	HE-048	1.8.24.B	15	HE-140	0.6.45.B	5	HR-150	1.8.10.B	15
HA-215	1.8.41.B	9	HE-049	1.8.24.B	16	HE-142	0.6.45.B	7	HR-209	1.8.25.B	9
HA-250	1.8.29.B	46	HE-050	1.8.24.B	17	HE-143	0.6.45.B	8	HR-210	1.8.10.B	5
HA-300	1.8.10.B	10	HE-052	1.8.24.B	19	HE-144	0.6.45.B	9	HR-211	1.8.10.B	36
HA-301	1.8.10.B	11	HE-053	1.8.24.B	20	HE-145	0.6.45.B	10	HR-212	1.8.10.B	36
HA-302	1.8.10.B	26	HE-054	1.8.24.B	21	HH-110	1.8.03.B	3	HR-250	1.8.10.B	5
HA-304	0.6.17.B	4	HE-056	1.8.24.B	22	HH-111	1.8.03.B	3	HR-257	1.8.10.B	22
HA-339	1.8.28.B	27	HE-058	1.8.24.B	23	HH-115	1.8.05.B	19	HR-310	1.8.10.B	5
HA-341	1.8.41.B	4	HE-059	1.8.24.B	3	HH-151	1.8.05.B	13	HR-311	1.8.10.B	36
HA-349	1.8.08.B	37	HE-060	1.8.24.B	24	HH-152	1.8.05.B	13	HR-312	1.8.10.B	36
HA-402	1.8.01.B	4	HE-061	1.8.24.B	25	HH-153	1.8.05.B	14	HR-357	1.8.10.B	22
HA-403	1.8.01.B	3	HE-065	1.8.24.B	26	HH-154	1.8.05.B	14	HT-005	1.8.07.B	12
HA-420	1.8.05.B	7	HE-066	1.8.20.B	16	HR-003	0.6.35.B	8	HT-006	3.2.02.B	16
HA-447	1.8.05.B	6	HE-067	1.8.20.B	17	HR-003	1.8.42.B	56	HT-009	1.8.30.B	3
HA-461	1.8.40.B	18	HE-068	1.8.20.B	18	HR-004	0.6.35.B	8	HT-010	1.8.30.B	15
HA-462	1.8.28.B	18	HE-100	1.8.26.B	12	HR-004	1.8.42.B	56	HT-012	1.8.05.B	12
HB-020	1.8.34.B	11	HE-101	1.8.26.B	13	HR-005	0.6.35.B	8	HT-017	1.8.11.B	14
HB-030	1.8.34.B	12	HE-102	1.8.26.B	11	HR-005	1.8.42.B	56	HT-019	1.8.42.B	6
HB-101	1.8.10.B	25	HE-103	1.8.26.B	2	HR-006	0.6.35.B	7	HT-020	1.8.06.B	21
HB-102	1.8.10.B	25	HE-104	1.8.26.B	3	HR-006	1.8.42.B	61	HT-021	1.8.07.B	12
HB-103	1.8.10.B	25	HE-105	1.8.26.B	6	HR-007	0.6.35.B	8	HT-022	1.8.17.B	1
HB-104	1.8.10.B	25	HE-106	1.8.26.B	7	HR-007	1.8.42.B	56	HT-023	1.8.17.B	1
HB-105	1.8.10.B	25	HE-107	1.8.19.B	29	HR-008	1.8.10.B	5	HT-034	1.8.17.B	1
HB-201	1.8.10.B	25	HE-107	1.8.26.B	8	HR-009	1.8.25.B	9	HT-035	1.8.17.B	1
HB-202	1.8.10.B	25	HE-108	1.8.26.B	9	HR-013	1.8.10.B	27	HT-036	1.8.17.B	1
HB-203	1.8.10.B	25	HE-109	1.8.26.B	10	HR-014	1.8.10.B	27	HT-037	1.8.17.B	1
HB-204	1.8.10.B	25	HE-110	1.8.26.B	1	HR-015	1.8.10.B	27	HT-038	1.8.17.B	1
HB-205	1.8.10.B	25	HE-111	1.8.26.B	4	HR-016	1.8.10.B	27	HT-039	1.8.17.B	1
HB-501	1.8.34.B	8	HE-116	0.6.43.B	10	HR-017	1.8.10.B	27	HT-040	1.8.14.B	1
HE-001	1.8.13.B	15	HE-117	0.6.43.B	11	HR-019	1.8.10.B	27	HT-042	1.8.17.B	1
HE-002	1.8.13.B	14	HE-118	0.6.43.B	12	HR-020	1.8.24.B	27	HT-055	1.8.14.B	2
HE-003	1.8.13.B	17	HE-119	0.6.43.B	13	HR-021	1.8.10.B	27	HT-057	1.8.06.B	28
HE-003	1.8.17.B	8	HE-120	0.6.43.B	14	HR-022	1.8.10.B	27	HT-060	1.8.14.B	3
HE-004	1.8.13.B	16	HE-121	0.6.43.B	15	HR-023	1.8.10.B	27	HT-061	1.8.14.B	3
HE-004	1.8.14.B	11	HE-122	0.6.43.B	16	HR-024	3.2.01.B	7	HT-062	1.8.14.B	4
HE-005	0.6.14.B	4	HE-123	0.6.43.B	17	HR-032	1.8.10.B	27	HT-065	0.6.14.B	8
HE-005	0.6.33.B	9	HE-124	0.6.43.B	18	HR-033	1.8.10.B	27	HT-065	0.6.99.B	

Indeks over alle Flex reservedelskataloger

Index für alle Flex Ersatzteilkataloge

Index for all spare part catalogues

Index de tous catalogues Flex des pièces détachées

Indice per tutti i cataloghi dei pezzi di ricambio

NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No

HT-070	0.6.14.B	20	MA-127	6.3.02.B	3	MB-161	1.8.42.B	23	MC-144	0.6.00.B	8
HT-070	0.6.33.B	17	MA-128	6.3.02.B	4	MB-163	1.8.42.B	24	MC-161	0.6.06.B	1
HT-092	0.6.14.B	23	MA-129	0.6.25.B	6	MB-164	1.8.42.B	25	MC-165	0.6.06.B	11
HT-092	0.6.99.B		MA-129	7.3.01.B	6	MB-166	1.8.42.B	27	MC-176	0.6.09.B	32
HT-097	0.6.14.B	22	MA-132	6.3.01.B	2	MB-180	1.8.30.B	17	MC-176	0.6.10.B	28
HT-097	0.6.99.B		MA-134	0.6.25.B	40	MB-181	1.8.32.B	50	MC-177	0.6.13.B	3
HT-099	0.6.45.B	11	MA-134	7.3.01.B	40	MB-200	1.8.08.B	18	MC-186	0.6.09.B	5
HT-099	0.6.99.B		MA-137	0.6.25.B	28	MB-201	1.8.08.B	19	MC-186	0.6.17.B	2
HT-115	0.6.31.B	29	MA-137	7.3.01.B	28	MB-254	1.8.11.B	11	MC-186	0.6.18.B	19
HT-115	0.6.99.B		MA-139	0.6.25.B	7	MB-255	1.8.11.B	16	MC-186	0.6.35.B	10
HT-144	0.6.47.B	10	MA-139	7.3.01.B	7	MB-262	1.8.11.B	22	MC-186	0.6.39.B	2
MA-001	1.8.29.B	10	MA-149	6.3.02.B	2	MB-280	0.6.05.B	12	MC-188	0.6.04.B	11
MA-003	1.8.30.B	4	MA-150	6.3.02.B	1	MB-280	0.6.32.B	4	MC-190	0.6.04.B	12
MA-005	1.8.39.B	15	MA-152	6.3.02.B	2	MB-280	1.8.11.B	1	MC-191	0.6.14.B	6
MA-007	1.8.29.B	30	MA-153	6.3.02.B	1	MB-280	1.8.42.B	28	MC-195	0.6.00.B	10
MA-008	1.8.07.B	11	MA-154	6.3.02.B	2	MB-281	1.8.11.B	2	MC-196	0.6.04.B	13
MA-020	1.8.38.B	17	MA-155	6.3.02.B	1	MB-281	1.8.42.B	28	MC-197	0.6.04.B	14
MA-032	1.8.08.B	25	MA-235	0.6.25.B	21	MB-283	1.8.11.B	6	MC-200	0.6.01.B	22
MA-035	1.8.39.B	16	MA-235	7.3.01.B	21	MB-284	0.6.02.B	50	MC-205	0.6.17.B	13
MA-042	1.8.41.B	13	MA-236	0.6.25.B	8	MB-284	1.8.11.B	7	MC-205	0.6.18.B	6
MA-046	1.8.11.B	39	MA-236	7.3.01.B	8	MB-288	1.8.11.B	3	MC-223	0.6.04.B	38
MA-046	1.8.19.B	27	MA-299	1.8.33.B	22	MB-321	1.8.42.B	29	MC-231	0.6.05.B	5
MA-046	1.8.42.B	7	MA-354	0.6.11.B		MB-322	1.8.02.B	10	MC-231	0.6.32.B	5
MA-072	1.8.32.B	16	MB-006	1.8.15.B	15	MB-322	1.8.42.B		MC-245	0.6.07.B	11
MA-076	6.3.02.B	1	MB-007	1.8.15.B	16	MB-329	0.6.10.B	24	MC-248	0.6.04.B	35
MA-077	6.3.02.B	1	MB-038	0.6.35.B	4	MB-382	0.6.10.B	13	MC-249	0.6.17.B	27
MA-079	6.3.02.B	1	MB-038	1.8.42.B	63	MC-028	0.6.08.B	8	MC-261	0.6.09.B	45
MA-083	0.6.25.B	20	MB-040	1.8.06.B	15	MC-033	0.6.08.B		MC-263	0.6.04.B	26
MA-083	7.3.01.B	20	MB-041	1.8.06.B	17	MC-047	0.6.08.B	5	MC-264	0.6.01.B	10
MA-085	1.8.28.B	12	MB-100	1.8.07.B	8	MC-063	0.6.00.B	5	MC-265	0.6.01.B	11
MA-086	0.6.25.B	27	MB-101	1.8.07.B	6	MC-063	0.6.12.B	19	MC-266	0.6.09.B	46
MA-086	7.3.01.B	27	MB-109	1.8.42.B	8	MC-070	0.6.17.B	25	MC-274	0.6.00.B	18
MA-087	6.3.02.B	1	MB-110	1.8.42.B	9	MC-080	0.6.04.B	25	MC-274	0.6.01.B	23
MA-088	1.8.11.B	19	MB-111	1.8.42.B	10	MC-101	0.6.07.B	19	MC-286	0.6.02.B	51
MA-089	1.8.11.B	18	MB-112	1.8.42.B	11	MC-102	0.6.07.B	20	MC-346	0.6.13.B	8
MA-091	6.3.02.B	2	MB-113	1.8.42.B	12	MC-105	0.6.04.B	19	MC-347	0.6.13.B	9
MA-092	1.8.34.B	2	MB-114	1.8.30.B	14	MC-108	0.6.07.B	24	MC-351	0.6.09.B	42
MA-093	1.8.10.B	41	MB-115	1.8.42.B	13	MC-109	0.6.09.B	27	MC-351	0.6.28.B	1
MA-093	1.8.35.B	14	MB-116	1.8.42.B	13	MC-111	0.6.09.B	4	MC-352	0.6.09.B	22
MA-094	6.3.02.B	2	MB-118	1.8.42.B	14	MC-111	0.6.17.B	1	MC-359	0.6.13.B	12
MA-096	6.3.02.B	2	MB-119	1.8.42.B	15	MC-111	0.6.18.B	1	MC-361	0.6.13.B	14
MA-097	6.3.02.B	2	MB-125	1.8.42.B	16	MC-111	0.6.35.B	11	MC-363	0.6.13.B	13
MA-107	0.6.25.B	39	MB-126	1.8.42.B	17	MC-111	0.6.39.B	1	MC-367	0.6.13.B	16
MA-107	7.3.01.B	39	MB-127	0.6.06.B	6	MC-112	0.6.05.B	6	MC-370	0.6.13.B	5
MA-118	6.3.02.B	1	MB-127	1.8.42.B	18	MC-128	0.6.09.B	13	MC-377	0.6.05.B	14
MA-120	6.3.02.B	1	MB-128	0.6.06.B	7	MC-131	0.6.07.B	37	MC-381	0.6.02.B	53
MA-122	6.3.02.B	1	MB-128	1.8.42.B	19	MC-133	0.6.01.B	18	MC-390	1.8.02.B	16
MA-124	6.3.02.B	2	MB-129	1.8.42.B	20	MC-134	0.6.01.B	17	MC-390	1.8.08.B	36
MA-125	6.3.02.B	2	MB-159	1.8.42.B	21	MC-140	0.6.01.B	5	MC-397	0.6.08.B	10
MA-126	6.3.02.B	2	MB-160	1.8.42.B	22	MC-141	0.6.01.B	12	MC-436	0.6.07.B	5

NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No

MC-453	0.6.24.B	10	MC-791	0.6.20.B	18	ME-131	0.6.29.B	8	MH-745	0.6.03.B	11
MC-454	0.6.03.B	1	MC-791	0.6.30.B	1	ME-132	0.6.29.B	9	MH-745	0.6.19.B	14
MC-457	0.6.09.B	25	MC-791	0.6.99.B		ME-133	0.6.29.B	10	MH-745	0.6.30.B	6
MC-457	0.6.16.B	22	MC-804	0.6.09.B	17	ME-134	0.6.29.B	10	MK-008	1.8.19.B	22
MC-458	0.6.09.B	15	MC-805	0.6.09.B	18	ME-135	0.6.29.B	11	MK-010	7.3.01.B	32
MC-458	0.6.16.B	24	MC-822	0.6.30.B	50	ME-136	0.6.29.B	12	MK-056	0.6.25.B	22
MC-459	0.6.09.B	39	MC-846	0.6.30.B	56	ME-137	0.6.29.B	12	MK-056	7.3.01.B	22
MC-459	0.6.10.B	29	MC-847	0.6.29.B	6	ME-138	0.6.29.B	13	MK-106	1.8.10.B	17
MC-465	0.6.07.B	10	MC-848	0.6.29.B	6	ME-139	0.6.29.B	14	MK-135	0.6.25.B	14
MC-469	0.6.34.B	12	MC-859	0.6.01.B	7	ME-140	0.6.29.B	15	MK-135	7.3.01.B	14
MC-494	0.6.21.B	16	MC-860	0.6.26.B	12	ME-141	0.6.29.B	6	MK-254	1.8.22.B	16
MC-496	0.6.21.B	6	MC-870	0.6.29.B	7	ME-142	0.6.29.B	6	MK-300	1.8.03.B	13
MC-500	0.6.24.B	11	MC-873	0.6.30.B	4	ME-161	0.6.46.B	11	MK-306	1.8.10.B	17
MC-525	0.6.12.B	16	MC-875	0.6.29.B	5	ME-162	0.6.46.B	12	MK-350	1.8.03.B	2
MC-540	0.6.01.B	1	MC-876	0.6.29.B	5	ME-182	0.6.42.B	7	MK-350	1.8.04.B	3
MC-540	0.6.04.B	30	MC-878	0.6.35.B	9	ME-183	0.6.42.B	9	MK-500	1.8.03.B	18
MC-551	0.6.11.B	3	MC-878	0.6.39.B	6	ME-184	0.6.42.B	8	MK-549	1.8.03.B	18
MC-562	0.6.17.B	18	MC-900	0.6.27.B	21	ME-188	0.6.42.B	10	ML-005	1.8.28.B	2
MC-562	0.6.18.B	17	MC-925	0.6.38.B	9	MH-107	0.6.11.B	4	ML-006	1.8.28.B	2
MC-563	0.6.12.B	17	MC-927	0.6.38.B	6	MH-107	0.6.25.B	9	ML-007	1.8.28.B	2
MC-584	0.6.02.B	54	MC-928	0.6.27.B	22	MH-107	7.3.01.B	9	ML-008	0.6.09.B	16
MC-590	0.6.00.B	6	MC-933	0.6.39.B	7	MH-215	0.6.03.B	13	ML-008	1.8.29.B	39
MC-590	0.6.12.B	3	MC-938	0.6.07.B	38	MH-215	0.6.06.B	12	ML-017	0.6.37.B	3
MC-605	0.6.10.B	23	MC-943	0.6.35.B	13	MH-242	0.6.02.B	45	ML-017	1.8.06.B	25
MC-623	0.6.21.B	22	MC-943	0.6.39.B	5	MH-242	0.6.03.B	11	ML-018	0.6.37.B	4
MC-666	0.6.07.B	18	MC-946	0.6.38.B	5	MH-242	0.6.19.B	14	ML-024	1.8.44.B	1
MC-673	0.6.21.B	7	MC-948	0.6.09.B	43	MH-242	0.6.30.B	6	ML-030	1.8.08.B	3
MC-678	0.6.13.B	22	MC-948	0.6.38.B	11	MH-327	0.6.03.B	8	ML-030	1.8.10.B	4
MC-682	0.6.30.B	37	MC-949	0.6.09.B	44	MH-327	0.6.27.B	23	ML-030	1.8.35.B	13
MC-683	0.6.35.B	12	MC-949	0.6.38.B	12	MH-339	0.6.18.B	21	ML-035	0.6.12.B	22
MC-683	0.6.39.B	4	MC-952	0.6.38.B	10	MH-341	0.6.18.B	20	ML-035	0.6.27.B	17
MC-712	0.6.29.B	8	MC-977	0.6.29.B	7	MH-373	0.6.24.B	8	ML-044	3.2.01.B	
MC-714	0.6.29.B	10	MC-978	0.6.29.B	11	MH-427	0.6.03.B	5	ML-048	1.8.19.B	2
MC-715	0.6.29.B	10	MC-979	0.6.19.B	6	MH-427	0.6.24.B	8	ML-051	1.8.11.B	35
MC-720	0.6.29.B	12	MC-979	0.6.20.B	18	MH-428	0.6.19.B	29	ML-053	1.8.19.B	4
MC-721	0.6.29.B	9	MC-979	0.6.30.B	1	MH-428	0.6.24.B	8	ML-056	1.8.06.B	20
MC-724	0.6.30.B	42	MC-979	0.6.99.B		MH-428	0.6.30.B	39	ML-058	1.8.07.B	17
MC-725	0.6.30.B	43	MC-980	0.6.29.B	12	MH-429	0.6.03.B	5	ML-058	1.8.32.B	6
MC-736	0.6.30.B	51	MC-983	0.6.30.B	5	MH-430	0.6.09.B	30	ML-058	1.8.36.B	5
MC-737	0.6.30.B	57	MC-984	0.6.13.B	23	MH-430	0.6.24.B	8	ML-059	0.6.24.B	8
MC-738	0.6.30.B	54	MC-987	0.6.15.B	4	MH-432	0.6.18.B	7	ML-060	0.6.24.B	8
MC-748	0.6.01.B	2	MC-988	0.6.13.B	28	MH-433	0.6.18.B	7	ML-060	1.8.03.B	6
MC-748	0.6.17.B	19	ME-010	0.6.01.B	13	MH-435	0.6.03.B	5	ML-061	0.6.09.B	29
MC-748	0.6.18.B	12	ME-014	0.6.03.B	16	MH-435	0.6.24.B	8	ML-061	0.6.24.B	8
MC-755	0.6.29.B	23	ME-015	0.6.06.B	1	MH-438	0.6.03.B	7	ML-062	0.6.24.B	8
MC-756	0.6.29.B	23	ME-017	0.6.21.B	17	MH-439	0.6.03.B	3	ML-063	0.6.24.B	8
MC-758	0.6.17.B	28	ME-023	0.6.09.B	35	MH-535	0.6.03.B	13	ML-063	1.8.03.B	6
MC-778	0.6.01.B	6	ME-037	0.6.07.B	24	MH-539	0.6.12.B	14	ML-064	0.6.24.B	8
MC-779	0.6.01.B	16	ME-037	0.6.35.B	3	MH-720	1.8.08.B	4	ML-064	1.8.03.B	6
MC-791	0.6.19.B	6	ME-059	0.6.30.B	23	MH-745	0.6.02.B	45	ML-065	0.6.03.B	5

NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No

ML-065	0.6.24.B	8	ML-412	1.8.44.B	22	MM-252	1.8.23.B	3	MR-110	1.8.08.B	45
ML-065	3.2.02.B	18	ML-413	1.8.44.B	22	MM-253	1.8.23.B	4	MR-112	1.8.02.B	9
ML-066	0.6.03.B	5	ML-414	1.8.44.B	22	MM-254	1.8.23.B	5	MR-113	1.8.02.B	8
ML-066	0.6.24.B	8	ML-415	1.8.44.B	20	MM-255	1.8.23.B	6	MR-114	0.6.24.B	20
ML-066	1.8.03.B	6	ML-416	1.8.44.B	22	MM-256	1.8.23.B	7	MR-115	0.6.18.B	5
ML-067	0.6.24.B	8	ML-417	1.8.44.B	20	MM-257	1.8.23.B	8	MR-115	3.2.02.B	11
ML-068	0.6.03.B	5	ML-418	1.8.44.B	20	MM-258	1.8.23.B	9	MR-117	1.8.27.B	15
ML-068	0.6.24.B	8	ML-421	1.8.44.B	18	MM-259	1.8.23.B	10	MR-122	1.8.06.B	1
ML-068	1.8.03.B	6	ML-425	1.8.44.B	19	MM-260	1.8.23.B	11	MR-129	1.8.27.B	16
ML-069	0.6.03.B	5	ML-427	1.8.44.B	20	MR-002	1.8.02.B	20	MR-131	1.8.06.B	11
ML-069	0.6.19.B	30	ML-429	1.8.44.B	20	MR-005	1.8.09.B	7	MR-132	1.8.06.B	8
ML-069	0.6.24.B	8	ML-430	1.8.44.B	20	MR-006	1.8.09.B	9	MR-133	1.8.06.B	9
ML-069	0.6.30.B	40	ML-431	1.8.44.B	20	MR-009	1.8.09.B	14	MR-134	1.8.06.B	12
ML-069	1.8.03.B	6	ML-434	1.8.03.B	9	MR-009	1.8.29.B	48	MR-135	0.6.10.B	22
ML-070	1.8.07.B	16	ML-447	1.8.08.B	40	MR-009	1.8.32.B	30	MR-135	1.8.06.B	6
ML-100	1.8.44.B	2	ML-448	1.8.08.B	39	MR-012	1.8.12.B	3	MR-136	1.8.06.B	10
ML-103	1.8.03.B	5	ML-449	1.8.08.B	41	MR-012	1.8.13.B	3	MR-137	1.8.06.B	7
ML-110	1.8.03.B	4	ML-457	1.8.32.B	47	MR-012	1.8.15.B	3	MR-139	1.8.06.B	3
ML-158	0.6.24.B	8	ML-457	1.8.33.B	10	MR-012	1.8.16.B	3	MR-140	1.8.06.B	4
ML-159	0.6.24.B	8	ML-461	1.8.44.B	15	MR-016	1.8.40.B	7	MR-141	1.8.06.B	2
ML-164	0.6.25.B	33	ML-462	1.8.29.B	51	MR-020	1.8.32.B	4	MR-142	0.6.24.B	21
ML-164	7.3.01.B	33	ML-462	1.8.44.B	17	MR-030	1.8.29.B	19	MR-143	1.8.06.B	1
ML-183	0.6.02.B	27	ML-463	1.8.44.B	14	MR-030	1.8.32.B	32	MR-150	0.6.10.B	20
ML-183	0.6.05.B	8	ML-464	1.8.44.B	16	MR-031	1.8.33.B	21	MR-150	0.6.38.B	17
ML-183	0.6.32.B	10	ML-467	1.8.29.B	50	MR-032	1.8.29.B	19	MR-203	0.6.10.B	16
ML-204	1.8.44.B	3	ML-467	1.8.44.B	10	MR-034	1.8.25.B	32	MR-203	1.8.11.B	25
ML-215	1.8.44.B	8	ML-468	1.8.44.B	11	MR-039	1.8.36.B	1	MR-206	1.8.06.B	19
ML-217	1.8.44.B	7	ML-469	1.8.03.B	8	MR-040	1.8.36.B	2	MR-206	1.8.32.B	23
ML-218	1.8.44.B	6	ML-469	1.8.44.B	12	MR-041	0.6.02.B	49	MR-207	1.8.29.B	21
ML-250	1.8.28.B	29	ML-470	1.8.44.B	9	MR-043	1.8.25.B	33	MR-207	1.8.34.B	5
ML-280	0.6.24.B	9	ML-471	1.8.44.B	13	MR-044	1.8.25.B	34	MR-208	1.8.10.B	40
ML-280	1.8.44.B	23	ML-488	1.8.42.B	30	MR-045	1.8.25.B	35	MR-208	1.8.11.B	12
ML-281	0.6.24.B	9	ML-489	1.8.42.B	31	MR-047	1.8.11.B	5	MR-209	0.6.07.B	17
ML-281	1.8.44.B	24	ML-513	1.8.28.B	30	MR-049	1.8.09.B	7	MR-209	0.6.23.B	21
ML-282	0.6.24.B	9	ML-548	0.6.24.B	8	MR-050	1.8.11.B	4	MR-209	1.8.29.B	16
ML-282	1.8.44.B	25	ML-564	0.6.31.B	28	MR-051	1.8.11.B	8	MR-209	1.8.30.B	12
ML-290	1.8.44.B	4	ML-566	0.6.49.B	4	MR-052	1.8.09.B	9	MR-209	1.8.32.B	35
ML-306	1.8.28.B	31	MM-016	1.8.22.B	3	MR-061	1.8.32.B	17	MR-209	1.8.35.B	20
ML-350	1.8.03.B	16	MM-016	1.8.28.B	5	MR-075	0.6.04.B	3	MR-210	0.6.04.B	29
ML-352	1.8.03.B	16	MM-044	1.8.22.B	15	MR-075	0.6.07.B	28	MR-210	1.8.01.B	16
ML-372	1.8.08.B	46	MM-100	1.8.03.B	7	MR-075	1.8.32.B	18	MR-210	1.8.35.B	18
ML-400	1.8.44.B	20	MM-120	1.8.04.B	4	MR-076	0.6.07.B	21	MR-211	1.8.30.B	13
ML-401	1.8.44.B	20	MM-128	1.8.04.B	4	MR-077	0.6.07.B	4	MR-212	1.8.37.B	12
ML-403	1.8.44.B	20	MM-212	1.8.22.B	20	MR-101	0.6.13.B	7	MR-214	0.6.07.B	7
ML-404	1.8.44.B	20	MM-220	1.8.06.B	22	MR-101	1.8.08.B	12	MR-215	1.8.08.B	28
ML-406	1.8.44.B	20	MM-226	1.8.03.B	7	MR-102	1.8.08.B	13	MR-216	0.6.23.B	3
ML-408	1.8.44.B	20	MM-228	1.8.04.B	4	MR-103	1.8.08.B	14	MR-216	1.8.40.B	7
ML-409	1.8.44.B	22	MM-230	1.8.07.B	14	MR-104	1.8.08.B	15	MR-216	4.3.01.B	2
ML-410	1.8.44.B	22	MM-231	1.8.07.B	15	MR-105	1.8.08.B	16	MR-217	0.6.07.B	22
ML-411	1.8.44.B	22	MM-251	1.8.23.B	2	MR-109	1.8.10.B	42	MR-218	0.6.07.B	1

NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No

MR-220	1.8.07.B	5	MS-222	1.8.13.B	12	SA-046	1.8.25.B	3	SC-301	5.3.03.B	2
MR-221	1.8.06.B	23	MS-305	1.8.43.B	12	SA-048	0.6.31.B	9	SC-302	5.3.03.B	9
MR-226	0.6.07.B	2	MS-306	1.8.43.B	13	SA-048	1.8.19.B	16	SC-303	5.3.03.B	3
MR-227	0.6.10.B	7	MS-307	1.8.43.B	10	SA-049	1.8.17.B	7	SC-304	5.3.03.B	5
MR-227	0.6.38.B	24	MS-308	1.8.43.B	11	SA-050	1.8.17.B	7	SC-305	5.3.03.B	6
MR-250	1.8.07.B	1	MS-309	1.8.43.B	14	SA-051	0.6.33.B	12	SC-306	5.3.03.B	11
MR-255	1.8.07.B	2	MS-312	1.8.43.B	15	SA-051	1.8.14.B	8	SC-307	5.3.03.B	13
MR-256	1.8.07.B	2	MS-313	1.8.43.B	16	SA-052	0.6.33.B	13	SC-308	5.3.03.B	14
MR-257	1.8.32.B	5	MS-320	1.8.33.B	6	SA-052	1.8.14.B	9	SC-309	5.3.03.B	15
MR-257	1.8.36.B	4	MS-320	1.8.35.B	23	SA-053	1.8.13.B	9	SC-310	5.3.03.B	16
MR-260	1.8.07.B	9	MS-320	1.8.43.B	8	SA-053	1.8.14.B	7	SC-311	5.3.03.B	27
MR-265	1.8.07.B	10	MS-321	1.8.43.B	9	SA-054	1.8.25.B	10	SC-312	5.3.03.B	30
MR-266	1.8.07.B	7	MS-340	1.8.43.B	7	SA-055	1.8.17.B	7	SC-313	5.3.03.B	31
MR-270	1.8.07.B	3	MS-341	0.6.04.B	15	SA-056	1.8.17.B	7	SC-314	5.3.03.B	32
MR-302	1.8.08.B	10	MS-347	1.8.32.B	42	SA-407	1.8.05.B	25	SC-315	5.3.03.B	26
MR-303	1.8.08.B	11	MS-347	1.8.43.B	6	SA-434	1.8.05.B	26	SC-316	5.3.03.B	40
MR-309	1.8.10.B	42	MS-348	0.6.04.B	16	SA-448	1.8.35.B	5	SC-317	5.3.03.B	41
MR-310	0.6.17.B	22	MS-350	1.8.43.B	4	SA-555	0.6.10.B	1	SC-318	5.3.03.B	28
MR-311	1.8.35.B	2	MS-352	1.8.43.B	3	SA-555	0.6.38.B	18	SC-319	5.3.03.B	29
MR-312	1.8.33.B	16	MS-355	1.8.31.B	4	SA-556	0.6.10.B	2	SC-321	5.3.03.B	20
MR-312	1.8.35.B	3	MS-355	1.8.43.B	1	SA-556	0.6.38.B	19	SC-322	5.3.03.B	33
MR-403	1.8.08.B	7	MS-356	1.8.43.B	2	SA-569	1.8.43.B		SC-323	5.3.04.B	43
MR-404	1.8.08.B	29	MS-358	1.8.43.B	5	SA-583	1.8.11.B	40	SC-324	5.3.04.B	62
MR-407	1.8.09.B	13	MV-018	0.6.46.B	13	SA-583	1.8.19.B	28	SC-325	5.3.04.B	63
MR-415	1.8.09.B	13	MV-021	0.6.46.B	14	SA-890	1.8.06.B	29	SC-326	5.3.04.B	86
MR-425	1.8.32.B	24	MY-001	1.8.05.B	8	SA-891	1.8.22.B	8	SC-327	5.3.04.B	46
MR-426	1.8.32.B	25	MY-002	1.8.29.B	44	SA-983	0.6.33.B	4	SC-328	5.3.04.B	47
MR-430	1.8.35.B	17	MY-003	0.6.35.B	15	SA-984	0.6.33.B	5	SC-329	5.3.04.B	49
MR-431	1.8.35.B	19	MY-003	0.6.39.B	10	SA-985	0.6.33.B	11	SC-330	5.3.04.B	50
MR-456	1.8.05.B	5	MY-003	1.8.32.B	43	SC-089	1.8.13.B	10	SC-331	5.3.04.B	51
MR-465	1.8.05.B	5	MY-004	1.8.32.B	43	SC-089	1.8.14.B	5	SC-332	5.3.04.B	56
MR-466	1.8.01.B	2	MY-005	1.8.29.B	44	SC-183	1.8.10.B	27	SC-333	5.3.04.B	57
MR-485	1.8.01.B	2	MY-006	1.8.29.B	44	SC-184	1.8.10.B	27	SC-334	5.3.04.B	58
MR-487	1.8.31.B	20	MY-007	1.8.29.B	44	SC-185	1.8.10.B	27	SC-335	5.3.05.B	108
MR-500	1.8.11.B	27	MY-030	1.8.12.B	10	SC-186	1.8.10.B	27	SC-336	5.3.05.B	107
MR-505	1.8.32.B	9	MY-030	1.8.13.B	23	SC-249	1.8.22.B	14	SC-337	5.3.05.B	115
MR-580	1.8.11.B	31	SA-002	1.8.17.B	3	SC-251	1.8.42.B	32	SC-339	5.3.05.B	104
MR-600	1.8.32.B	10	SA-012	1.8.17.B	2	SC-252	0.6.35.B	5	SC-341	5.3.05.B	102
MR-601	1.8.32.B	10	SA-024	0.6.14.B	3	SC-252	1.8.42.B	64	SC-342	5.3.05.B	96
MR-603	1.8.33.B	13	SA-024	0.6.45.B	12	SC-262	1.8.38.B	22	SC-343	5.3.05.B	88
MR-603	1.8.35.B	4	SA-024	1.8.17.B	5	SC-271	0.6.35.B	6	SC-344	5.3.05.B	101
MS-001	0.6.04.B	21	SA-030	1.8.17.B	7	SC-271	1.8.42.B	70	SC-345	5.3.05.B	90
MS-008	1.8.30.B	18	SA-035	0.6.14.B	18	SC-274	1.8.19.B	31	SC-346	5.3.05.B	93
MS-009	0.6.05.B	9	SA-036	0.6.14.B	17	SC-275	1.8.19.B	30	SC-347	5.3.05.B	100
MS-012	1.8.05.B	24	SA-038	1.8.17.B	7	SC-276	1.8.25.B	5	SC-348	5.3.05.B	94
MS-022	1.8.13.B	12	SA-042	0.6.14.B	7	SC-277	1.8.25.B	6	SC-349	5.3.05.B	122
MS-023	1.8.13.B	13	SA-042	1.8.17.B	4	SC-278	1.8.25.B	8	SC-350	5.3.05.B	121
MS-025	1.8.27.B	23	SA-043	1.8.25.B	1	SC-279	1.8.19.B	35	SC-351	5.3.05.B	125
MS-030	1.8.12.B	11	SA-044	1.8.25.B	7	SC-281	3.2.02.B	20	SC-352	5.3.05.B	128
MS-040	1.8.13.B	13	SA-045	1.8.25.B	2	SC-300	5.3.03.B	7	SC-353	5.3.05.B	129

NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No

SC-354	5.3.05.B	135	SC-406	5.3.03.B	4	SC-646	1.8.19.B	36	SC-871	1.8.13.B	18
SC-355	5.3.05.B	127	SC-407	5.3.03.B	10	SC-647	1.8.19.B	37	SC-872	1.8.12.B	14
SC-356	5.3.05.B	147	SC-408	5.3.03.B	12	SC-686	0.6.24.B	6	SC-872	1.8.13.B	19
SC-357	5.3.05.B	148	SC-409	5.3.03.B	23	SC-688	0.6.12.B	8	SC-883	1.8.29.B	53
SC-358	5.3.05.B	116	SC-410	5.3.03.B	24	SC-696	1.8.19.B	7	SC-886	0.6.19.B	39
SC-359	5.3.04.B	71	SC-411	5.3.03.B	25	SC-714	0.6.31.B	22	SC-886	0.6.30.B	44
SC-360	5.3.04.B	68	SC-412	5.3.03.B	34	SC-714	1.8.19.B	8	SC-890	0.6.17.B	41
SC-361	5.3.04.B	81	SC-413	5.3.03.B	35	SC-715	0.6.31.B	23	SC-895	0.6.17.B	38
SC-362	5.3.04.B	72	SC-414	5.3.03.B	36	SC-715	1.8.19.B	9	SC-895	3.2.01.B	13
SC-363	5.3.04.B	80	SC-415	5.3.03.B	37	SC-733	1.8.17.B	6	SC-896	6.3.01.B	
SC-364	5.3.04.B	78	SC-416	5.3.03.B	38	SC-740	0.6.14.B	1	SC-901	1.8.42.B	65
SC-365	5.3.06.B	151	SC-417	5.3.03.B	39	SC-740	0.6.33.B	16	SC-903	1.8.44.B	5
SC-366	5.3.06.B	152	SC-418	5.3.03.B	42	SC-740	1.8.13.B	11	SC-905	1.8.32.B	19
SC-368	5.3.06.B	155	SC-421	5.3.03.B	22	SC-740	1.8.14.B	6	SC-909	3.2.01.B	32
SC-369	5.3.06.B	156	SC-422	5.3.03.B	21	SC-748	0.6.14.B	13	SC-910	3.2.01.B	33
SC-370	5.3.06.B	157	SC-423	5.3.04.B	44	SC-750	0.6.14.B	25	SC-912	3.2.01.B	34
SC-371	5.3.06.B	158	SC-424	5.3.04.B	64	SC-811	0.6.31.B	30	SC-913	3.2.01.B	35
SC-372	5.3.05.B	140	SC-425	5.3.04.B	65	SC-811	1.8.19.B	40	SC-914	3.2.01.B	36
SC-373	5.3.05.B	134	SC-426	5.3.04.B	85	SC-812	1.8.42.B	39	SC-915	3.2.01.B	37
SC-374	5.3.05.B	142	SC-427	5.3.04.B	45	SC-813	1.8.42.B	40	SC-917	0.6.17.B	29
SC-375	5.3.05.B	143	SC-433	5.3.05.B	110	SC-817	0.6.14.B	26	SC-918	1.8.19.B	38
SC-376	5.3.05.B	144	SC-434	5.3.04.B	69	SC-821	0.6.14.B	24	SC-919	0.6.17.B	42
SC-377	5.3.07.B	186	SC-435	5.3.04.B	70	SC-839	5.3.01.B	3	SC-920	0.6.17.B	43
SC-378	5.3.08.B	190	SC-436	5.3.04.B	66	SC-840	5.3.01.B	5	SC-923	1.8.42.B	62
SC-379	5.3.08.B	191	SC-437	5.3.04.B	67	SC-841	0.6.19.B	7	SC-927	0.6.26.B	7
SC-380	5.3.08.B	192	SC-438	5.3.04.B	76	SC-841	5.3.01.B	6	SC-927	1.8.20.B	8
SC-381	5.3.08.B	193	SC-439	5.3.04.B	77	SC-842	0.6.19.B	3	SC-928	0.6.26.B	8
SC-382	5.3.08.B	194	SC-440	5.3.04.B	79	SC-842	0.6.30.B	11	SC-928	1.8.20.B	9
SC-383	5.3.08.B	195	SC-441	5.3.04.B	82	SC-842	0.6.99.B		SC-935	1.8.25.B	25
SC-384	5.3.08.B	196	SC-442	5.3.06.B	150	SC-842	5.3.01.B	1	SC-936	1.8.27.B	20
SC-385	5.3.08.B	197	SC-443	5.3.06.B	149	SC-843	0.6.19.B	4	SC-942	5.3.01.B	14
SC-386	5.3.08.B	198	SC-445	5.3.06.B	153	SC-843	0.6.30.B	12	SC-942	5.3.02.B	10
SC-387	5.3.08.B	199	SC-446	5.3.06.B	154	SC-843	0.6.99.B		SC-943	5.3.01.B	15
SC-388	5.3.08.B	200	SC-447	5.3.05.B	124	SC-843	5.3.01.B	2	SC-944	0.6.19.B	5
SC-389	5.3.08.B	201	SC-448	5.3.05.B	123	SC-845	0.6.19.B	10	SC-944	5.3.01.B	7
SC-390	5.3.08.B	202	SC-450	5.3.04.B	60	SC-845	5.3.01.B	8	SC-944	5.3.01.B	16
SC-391	5.3.08.B	203	SC-451	5.3.04.B	61	SC-846	5.3.01.B	9	SC-945	5.3.01.B	17
SC-392	5.3.08.B	204	SC-452	5.3.05.B	131	SC-847	5.3.01.B	10	SC-946	0.6.17.B	30
SC-393	5.3.08.B	206	SC-453	5.3.04.B	59	SC-848	5.3.01.B	11	SC-948	5.3.02.B	7
SC-394	5.3.08.B	207	SC-453	5.3.06.B	159	SC-849	0.6.19.B	9	SC-949	5.3.02.B	7
SC-395	5.3.08.B	208	SC-454	5.3.07.B	187	SC-849	0.6.30.B	52	SC-950	5.3.02.B	8
SC-396	5.3.08.B	209	SC-455	5.3.04.B	84	SC-849	5.3.01.B	12	SC-951	5.3.02.B	8
SC-397	5.3.08.B	210	SC-485	1.8.05.B	23	SC-850	5.3.02.B	1	SC-952	5.3.01.B	18
SC-398	5.3.08.B	211	SC-536	1.8.26.B	5	SC-851	5.3.02.B	2	SC-953	5.3.01.B	19
SC-399	5.3.08.B	213	SC-588	1.8.42.B	33	SC-852	5.3.01.B	13	SC-957	5.3.01.B	3
SC-400	5.3.08.B	205	SC-589	1.8.42.B	34	SC-853	5.3.02.B	3	SC-960	1.8.27.B	22
SC-401	5.3.08.B	212	SC-590	1.8.42.B	35	SC-854	5.3.02.B	5	SC-961	6.3.01.B	
SC-402	5.3.08.B	219	SC-591	1.8.42.B	36	SC-855	5.3.02.B	4	SC-962	6.3.01.B	
SC-403	5.3.08.B	218	SC-592	1.8.42.B	37	SC-856	5.3.02.B	6	SC-963	6.3.01.B	3
SC-404	5.3.03.B	1	SC-593	1.8.42.B	38	SC-871	1.8.12.B	15	SC-970	1.8.37.B	11

NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No

SC-973	0.6.33.B	14	SD-122	5.3.01.B	19	SD-188	5.3.02.B	5	SD-291	0.6.05.B	18
SC-973	1.8.14.B	10	SD-125	0.6.17.B	8	SD-189	5.3.02.B	4	SD-292	0.6.02.B	31
SC-980	6.3.03.B	8	SD-126	0.6.17.B	11	SD-190	5.3.02.B	4	SD-292	0.6.32.B	11
SC-981	1.8.18.B	13	SD-127	0.6.17.B	9	SD-191	5.3.02.B	4	SD-293	0.6.02.B	32
SC-983	1.8.42.B	45	SD-130	0.6.17.B	10	SD-192	5.3.02.B	4	SD-294	0.6.02.B	33
SC-984	1.8.42.B	50	SD-131	3.2.01.B	38	SD-193	5.3.02.B	4	SD-297	0.6.02.B	34
SC-991	1.8.25.B	26	SD-136	0.6.02.B	48	SD-194	5.3.02.B	4	SD-298	0.6.02.B	35
SC-994	1.8.42.B	55	SD-153	1.8.42.B	11	SD-195	5.3.02.B	4	SD-299	0.6.02.B	36
SC-998	0.6.14.B	15	SD-157	6.3.01.B	13	SD-196	5.3.02.B	4	SD-300	0.6.02.B	37
SD-005	0.6.14.B	5	SD-158	6.3.01.B	14	SD-197	5.3.02.B	4	SD-301	0.6.02.B	38
SD-005	1.8.17.B	9	SD-159	1.8.09.B	21	SD-205	0.6.15.B	14	SD-302	0.6.02.B	39
SD-018	0.6.26.B	9	SD-160	1.8.09.B	21	SD-210	0.6.14.B	12	SD-303	0.6.10.B	10
SD-018	1.8.20.B	10	SD-161	1.8.09.B	21	SD-213	0.6.16.B	1	SD-303	0.6.38.B	22
SD-019	1.8.20.B	11	SD-162	1.8.09.B	22	SD-213	0.6.99.B		SD-307	0.6.10.B	11
SD-020	1.8.20.B	11	SD-163	1.8.09.B	22	SD-214	0.6.16.B	3	SD-307	0.6.38.B	21
SD-021	0.6.26.B	11	SD-164	1.8.09.B	22	SD-214	0.6.36.B	20	SD-308	0.6.10.B	12
SD-021	1.8.20.B	13	SD-165	1.8.12.B	4	SD-215	0.6.16.B	4	SD-308	0.6.38.B	20
SD-022	0.6.26.B	10	SD-165	1.8.13.B	6	SD-216	0.6.14.B	10	SD-312	0.6.14.B	27
SD-022	1.8.20.B	14	SD-165	1.8.15.B	13	SD-217	0.6.15.B	2	SD-313	0.6.14.B	28
SD-023	1.8.20.B	15	SD-165	1.8.16.B	5	SD-218	0.6.15.B	3	SD-318	0.6.33.B	13
SD-024	1.8.20.B	19	SD-166	1.8.08.B	3	SD-220	0.6.15.B	7	SD-320	0.6.04.B	3
SD-025	1.8.20.B	24	SD-166	1.8.10.B	4	SD-221	0.6.15.B	8	SD-322	0.6.16.B	5
SD-026	1.8.20.B	25	SD-166	1.8.35.B	13	SD-222	0.6.15.B	9	SD-325	0.6.16.B	6
SD-028	1.8.26.B	17	SD-167	1.8.12.B	4	SD-223	0.6.15.B	10	SD-326	0.6.16.B	7
SD-029	1.8.19.B	39	SD-167	1.8.13.B	6	SD-226	0.6.16.B	2	SD-327	0.6.16.B	8
SD-030	1.8.23.B	1	SD-167	1.8.15.B	13	SD-226	0.6.99.B		SD-328	0.6.16.B	9
SD-031	1.8.23.B	17	SD-167	1.8.16.B	5	SD-227	0.6.15.B	1	SD-329	0.6.16.B	10
SD-032	1.8.23.B	18	SD-168	1.8.12.B	4	SD-228	0.6.09.B	20	SD-330	0.6.16.B	11
SD-034	5.3.01.B	20	SD-168	1.8.13.B	6	SD-228	0.6.16.B	19	SD-331	0.6.16.B	12
SD-035	5.3.01.B	21	SD-168	1.8.15.B	13	SD-228	0.6.99.B		SD-332	0.6.16.B	13
SD-036	1.8.27.B	24	SD-168	1.8.16.B	5	SD-229	0.6.09.B	20	SD-333	0.6.16.B	14
SD-037	1.8.27.B	25	SD-169	5.3.02.B	12	SD-229	0.6.16.B	19	SD-334	0.6.16.B	15
SD-038	1.8.27.B	26	SD-170	5.3.02.B	13	SD-229	0.6.99.B		SD-335	0.6.16.B	16
SD-039	1.8.27.B	21	SD-171	5.3.02.B	3	SD-232	0.6.02.B	12	SD-336	0.6.16.B	17
SD-040	1.8.33.B	15	SD-172	5.3.02.B	3	SD-232	0.6.06.B	10	SD-337	0.6.16.B	18
SD-041	1.8.30.B	10	SD-173	5.3.02.B	3	SD-233	0.6.02.B	15	SD-339	0.6.14.B	9
SD-042	1.8.29.B	52	SD-174	5.3.02.B	3	SD-235	0.6.02.B	61	SD-347	0.6.04.B	32
SD-043	1.8.29.B	52	SD-175	5.3.02.B	3	SD-236	0.6.02.B	21	SD-348	0.6.04.B	33
SD-105	1.8.25.B	39	SD-176	5.3.02.B	3	SD-237	0.6.02.B	18	SD-353	0.6.30.B	21
SD-108	5.3.01.B	23	SD-177	5.3.02.B	3	SD-240	0.6.24.B	17	SD-374	0.6.19.B	
SD-111	5.3.02.B		SD-178	5.3.02.B	3	SD-245	0.6.02.B	13	SD-380	0.6.19.B	12
SD-112	1.8.34.B		SD-179	5.3.02.B	3	SD-247	0.6.02.B	16	SD-380	0.6.30.B	13
SD-113	1.8.34.B		SD-180	5.3.02.B	5	SD-249	0.6.02.B	20	SD-381	0.6.19.B	40
SD-115	5.3.02.B	1	SD-181	5.3.02.B	5	SD-255	0.6.02.B	7	SD-382	0.6.02.B	46
SD-116	5.3.02.B	2	SD-182	5.3.02.B	5	SD-256	0.6.02.B	8	SD-384	5.3.01.B	
SD-117	5.3.02.B	11	SD-183	5.3.02.B	5	SD-257	0.6.02.B	9	SD-387	0.6.26.B	2
SD-118	0.6.23.B	5	SD-184	5.3.02.B	5	SD-258	0.6.02.B	10	SD-388	0.6.19.B	15
SD-118	4.3.01.B	6	SD-185	5.3.02.B	5	SD-259	0.6.02.B	11	SD-388	0.6.30.B	7
SD-119	4.3.01.B	6	SD-186	5.3.02.B	5	SD-276	0.6.14.B	2	SD-391	0.6.24.B	23
SD-121	5.3.01.B	18	SD-187	5.3.02.B	5	SD-291	0.6.02.B	57	SD-392	0.6.24.B	24

NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No
SD-393	0.6.24.B	25	SD-450	0.6.30.B	49	SD-525	0.6.36.B	15	UA-046	1.8.01.B	11
SD-394	0.6.19.B	16	SD-451	0.6.19.B	34	SD-526	0.6.36.B	16	UA-048	1.8.29.B	3
SD-395	0.6.19.B	17	SD-452	0.6.19.B	35	SD-527	0.6.36.B	17	UA-049	1.8.01.B	11
SD-396	0.6.19.B	18	SD-456	0.6.21.B	18	SD-528	0.6.36.B	18	UA-049	3.2.02.B	13
SD-397	0.6.17.B	8	SD-464	0.6.33.B	1	SD-529	0.6.36.B	19	UA-058	0.6.46.B	15
SD-399	0.6.09.B	4	SD-465	0.6.33.B	2	SD-532	0.6.08.B		UA-060	0.6.08.B	6
SD-399	0.6.17.B	1	SD-466	0.6.33.B	3	SD-532	0.6.99.B		UA-067	0.6.08.B	7
SD-399	0.6.18.B	1	SD-467	0.6.33.B	7	SD-533	0.6.17.B	32	UA-091	1.8.33.B	20
SD-399	0.6.35.B	11	SD-468	0.6.33.B	8	SD-629	0.6.45.B	6	UA-092	1.8.10.B	12
SD-399	0.6.39.B	1	SD-469	0.6.33.B	6	SD-630	0.6.44.B	3	UA-093	0.6.46.B	16
SD-400	0.6.04.B	32	SD-470	0.6.11.B	9	SD-632	0.6.44.B	4	UA-114	0.6.27.B	24
SD-401	0.6.07.B	31	SD-471	0.6.19.B	37	SD-633	0.6.44.B	5	UA-118	0.6.00.B	11
SD-402	0.6.12.B	15	SD-472	0.6.30.B	53	SD-634	0.6.44.B	6	UA-120	5.3.02.B	9
SD-403	0.6.19.B	20	SD-473	0.6.30.B	55	SD-635	0.6.44.B	7	UA-128	1.8.09.B	11
SD-410	0.6.15.B	17	SD-474	0.6.29.B	21	SD-636	0.6.44.B	8	UA-131	0.6.30.B	36
SD-411	0.6.15.B	19	SD-475	0.6.29.B	20	SD-637	0.6.44.B	9	UA-136	1.8.08.B	32
SD-411	0.6.99.B		SD-476	0.6.29.B	22	SD-638	0.6.44.B	10	UA-149	1.8.32.B	45
SD-412	0.6.30.B	2	SD-480	0.6.19.B	31	SD-639	0.6.44.B	11	UA-149	1.8.33.B	9
SD-414	0.6.08.B		SD-480	0.6.30.B	22	SD-640	0.6.44.B	12	UA-150	1.8.11.B	21
SD-414	0.6.99.B		SD-481	0.6.19.B	38	SD-641	0.6.44.B	13	UA-151	1.8.22.B	8
SD-415	0.6.15.B	20	SD-484	0.6.30.B	14	SD-642	0.6.44.B	14	UA-160	0.6.23.B	20
SD-415	0.6.99.B		SD-485	0.6.19.B	28	SD-643	0.6.44.B	15	UA-210	0.6.18.B	13
SD-416	0.6.19.B	21	SD-488	0.6.02.B	56	SD-644	0.6.44.B	16	UA-210	3.2.02.B	17
SD-416	0.6.30.B	15	SD-489	0.6.16.B	4	SD-645	0.6.44.B	17	UA-215	1.8.10.B	2
SD-417	0.6.19.B	22	SD-490	0.6.30.B	47	SD-663	0.6.48.B	11	UA-224	1.8.08.B	34
SD-417	0.6.30.B	16	SD-492	0.6.02.B	17	SD-663	0.6.99.B		UA-224	1.8.41.B	5
SD-418	0.6.19.B	23	SD-493	0.6.19.B	32	SD-664	0.6.48.B	12	UA-243	1.8.06.B	16
SD-418	0.6.30.B	17	SD-494	0.6.13.B	18	SD-676	0.6.44.B	18	UA-243	1.8.29.B	11
SD-419	0.6.19.B	24	SD-495	0.6.09.B	33	SG-233	0.6.02.B	55	UA-243	1.8.35.B	10
SD-419	0.6.30.B	18	SD-496	0.6.17.B	32	SG-233	0.6.32.B	12	UA-244	1.8.33.B	23
SD-420	0.6.19.B	25	SD-500	0.6.30.B	41	UA-002	1.8.30.B	16	UA-246	1.8.01.B	6
SD-420	0.6.30.B	19	SD-507	0.6.27.B		UA-004	1.8.08.B	26	UA-246	1.8.29.B	29
SD-421	0.6.19.B	26	SD-508	0.6.17.B	33	UA-005	1.8.10.B	3	UA-246	1.8.32.B	51
SD-421	0.6.30.B	20	SD-509	0.6.27.B		UA-007	1.8.10.B	3	UA-246	1.8.36.B	3
SD-422	0.6.17.B	33	SD-510	0.6.36.B	1	UA-009	1.8.10.B	31	UA-247	1.8.01.B	6
SD-423	0.6.19.B	13	SD-511	0.6.36.B	2	UA-012	0.6.13.B	24	UA-247	1.8.41.B	8
SD-423	0.6.29.B	16	SD-511	0.6.99.B		UA-018	1.8.10.B	38	UA-248	1.8.32.B	7
SD-423	0.6.30.B	45	SD-512	0.6.36.B	2	UA-019	1.8.28.B	25	UA-249	0.6.13.B	25
SD-426	0.6.15.B	16	SD-512	0.6.99.B		UA-022	1.8.22.B	17	UA-253	1.8.08.B	33
SD-427	0.6.20.B	20	SD-513	0.6.36.B	4	UA-024	1.8.29.B	6	UA-266	1.8.02.B	17
SD-428	0.6.29.B	19	SD-514	0.6.36.B	5	UA-025	0.6.02.B	42	UA-266	1.8.40.B	14
SD-429	0.6.33.B	15	SD-516	0.6.36.B	6	UA-025	0.6.18.B	18	UA-280	1.8.29.B	26
SD-432	0.6.30.B	11	SD-517	0.6.36.B	7	UA-025	1.8.02.B	4	UA-348	1.8.28.B	28
SD-433	0.6.30.B	12	SD-518	0.6.36.B	8	UA-025	3.2.02.B	5	UA-410	3.2.01.B	39
SD-436	0.6.19.B	19	SD-519	0.6.36.B	9	UA-027	1.8.11.B	17	UA-423	1.8.08.B	31
SD-437	0.6.19.B	27	SD-520	0.6.36.B	10	UA-027	1.8.31.B	9	UA-423	1.8.29.B	12
SD-439	0.6.14.B	11	SD-521	0.6.36.B	11	UA-028	1.8.18.B	11	UA-423	1.8.31.B	6
SD-440	0.6.30.B	48	SD-522	0.6.36.B	12	UA-030	1.8.11.B	20	UA-425	1.8.31.B	15
SD-446	0.6.04.B	37	SD-523	0.6.36.B	13	UA-035	1.8.39.B	19	UA-444	0.6.34.B	9
SD-449	0.6.19.B	36	SD-524	0.6.36.B	14	UA-038	1.8.27.B	17	UA-444	1.8.10.B	8

NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No

UA-452	1.8.29.B	3	UA-728	7.3.01.B	11	UC-041	0.6.38.B	23	UC-216	0.6.09.B	34
UA-452	1.8.34.B	23	UA-750	1.8.10.B	1	UC-049	0.6.06.B	3	UC-217	0.6.11.B	2
UA-455	1.8.28.B	21	UA-751	1.8.11.B	30	UC-051	0.6.07.B	9	UC-218	0.6.11.B	7
UA-470	1.8.09.B	12	UA-755	3.2.01.B	16	UC-053	0.6.07.B	12	UC-227	0.6.20.B	6
UA-472	1.8.10.B	29	UA-800	1.8.28.B	40	UC-057	0.6.25.B	15	UC-227	0.6.29.B	1
UA-526	1.8.08.B	30	UA-803	1.8.28.B	7	UC-057	7.3.01.B	15	UC-229	0.6.20.B	6
UA-544	1.8.08.B	6	UA-810	1.8.28.B	15	UC-058	0.6.07.B	13	UC-230	0.6.20.B	6
UA-544	1.8.10.B	8	UA-810	1.8.31.B	11	UC-059	0.6.07.B	14	UC-230	0.6.29.B	1
UA-544	1.8.31.B	8	UA-815	1.8.28.B	13	UC-059	0.6.07.B	27	UC-231	0.6.20.B	6
UA-547	3.2.01.B	8	UA-815	1.8.29.B	36	UC-062	0.6.04.B	36	UC-232	0.6.20.B	7
UA-548	1.8.40.B	16	UA-820	0.6.23.B	19	UC-063	0.6.04.B	7	UC-234	0.6.20.B	7
UA-561	1.8.27.B	18	UA-820	4.3.01.B	20	UC-064	0.6.04.B	8	UC-235	0.6.20.B	7
UA-564	1.8.32.B	31	UA-830	1.8.02.B	6	UC-068	0.6.07.B	36	UC-236	0.6.20.B	7
UA-566	1.8.40.B	8	UA-840	1.8.10.B	28	UC-101	0.6.18.B	8	UC-237	0.6.20.B	8
UA-567	1.8.31.B	14	UA-846	3.2.01.B	9	UC-105	0.6.06.B	4	UC-239	0.6.20.B	8
UA-567	1.8.40.B	9	UA-849	0.6.23.B	15	UC-110	1.8.08.B	44	UC-240	0.6.20.B	8
UA-568	0.6.23.B	11	UA-849	4.3.01.B	16	UC-120	0.6.02.B	60	UC-241	0.6.20.B	8
UA-568	4.3.01.B	12	UA-855	0.6.25.B	12	UC-121	1.8.06.B	14	UC-242	0.6.20.B	1
UA-571	1.8.06.B	13	UA-855	7.3.01.B	12	UC-124	0.6.06.B	8	UC-250	1.8.07.B	4
UA-572	1.8.10.B	29	UA-859	0.6.46.B	18	UC-126	0.6.02.B	58	UC-252	0.6.20.B	11
UA-579	1.8.38.B	16	UA-860	0.6.46.B	19	UC-126	0.6.04.B	23	UC-253	0.6.20.B	10
UA-579	1.8.39.B	18	UB-006	0.6.27.B	25	UC-126	0.6.05.B	17	UC-254	0.6.20.B	12
UA-584	1.8.38.B	11	UB-054	1.8.02.B	5	UC-127	0.6.04.B	31	UC-255	0.6.19.B	
UA-586	1.8.09.B	11	UB-056	0.6.23.B	4	UC-132	0.6.18.B	15	UC-255	0.6.20.B	13
UA-587	1.8.29.B	25	UB-056	0.6.34.B	13	UC-132	1.8.02.B	7	UC-256	0.6.20.B	14
UA-588	0.6.46.B	17	UB-056	4.3.01.B	5	UC-132	3.2.02.B	3	UC-257	0.6.20.B	2
UA-591	1.8.38.B	12	UB-057	1.8.32.B	49	UC-133	0.6.18.B	9	UC-257	0.6.29.B	15
UA-602	1.8.08.B	23	UB-060	1.8.32.B	39	UC-149	0.6.13.B	26	UC-261	1.8.08.B	20
UA-604	1.8.31.B	16	UB-061	1.8.32.B	39	UC-155	0.6.27.B	27	UC-269	0.6.20.B	15
UA-604	1.8.34.B	13	UB-062	1.8.32.B	39	UC-158	0.6.24.B	22	UC-270	0.6.20.B	16
UA-605	1.8.31.B	17	UB-070	0.6.18.B	14	UC-159	0.6.02.B	3	UC-271	0.6.19.B	11
UA-608	1.8.34.B	14	UB-070	1.8.02.B	3	UC-162	0.6.10.B	5	UC-272	0.6.20.B	17
UA-611	0.6.25.B	23	UB-070	1.8.39.B	21	UC-163	6.3.02.B	3	UC-273	0.6.20.B	19
UA-611	7.3.01.B	23	UB-070	3.2.02.B	4	UC-164	6.3.02.B	4	UC-273	0.6.29.B	4
UA-612	3.2.02.B	7	UB-085	1.8.08.B	24	UC-171	0.6.21.B	26	UC-275	0.6.20.B	19
UA-614	1.8.08.B	22	UB-086	0.6.27.B	26	UC-186	6.3.02.B	2	UC-276	0.6.20.B	19
UA-615	1.8.38.B	10	UB-086	1.8.29.B	37	UC-190	0.6.18.B	22	UC-276	0.6.29.B	4
UA-618	1.8.38.B	19	UB-092	1.8.38.B	9	UC-200	1.8.08.B	27	UC-277	0.6.20.B	19
UA-626	1.8.08.B	21	UB-094	3.2.02.B	14	UC-201	6.3.02.B	1	UC-278	0.6.15.B	11
UA-627	1.8.31.B	18	UC-007	1.8.10.B	13	UC-204	1.8.10.B	9	UC-279	0.6.12.B	9
UA-628	1.8.31.B	19	UC-008	0.6.23.B	12	UC-205	1.8.11.B	34	UC-281	0.6.12.B	10
UA-706	0.6.25.B	10	UC-008	4.3.01.B	13	UC-206	0.6.01.B	14	UC-282	0.6.12.B	11
UA-706	1.8.28.B	22	UC-010	1.8.34.B	15	UC-207	0.6.01.B	15	UC-283	0.6.12.B	12
UA-706	7.3.01.B	10	UC-011	1.8.34.B	16	UC-208	0.6.00.B	12	UC-289	0.6.05.B	3
UA-707	0.6.06.B	22	UC-018	0.6.25.B	13	UC-209	0.6.01.B	8	UC-291	0.6.20.B	4
UA-707	0.6.07.B	39	UC-018	7.3.01.B	13	UC-210	0.6.00.B	13	UC-293	0.6.20.B	5
UA-707	1.8.30.B	5	UC-032	0.6.02.B	14	UC-211	1.8.10.B	23	UC-294	0.6.20.B	3
UA-710	1.8.23.B	16	UC-033	0.6.02.B	19	UC-212	1.8.10.B	23	UC-298	0.6.08.B	
UA-718	1.8.02.B	2	UC-037	0.6.04.B	28	UC-213	1.8.10.B	23	UC-300	1.8.08.B	42
UA-728	0.6.25.B	11	UC-040	0.6.04.B	6	UC-214	0.6.01.B	9	UC-301	1.8.35.B	15

NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No

UC-304	1.8.08.B	43	UC-600	0.6.24.B	19	UC-703	0.6.27.B	9	UC-821	0.6.41.B	2
UC-305	1.8.10.B	7	UC-602	0.6.02.B	41	UC-704	0.6.27.B	10	UC-822	0.6.40.B	4
UC-321	0.6.15.B	5	UC-606	0.6.21.B	3	UC-705	0.6.27.B	11	UC-822	0.6.41.B	3
UC-325	6.3.01.B	3	UC-618	0.6.39.B	3	UC-706	0.6.27.B	12	UC-823	0.6.40.B	5
UC-331	0.6.17.B	3	UC-631	0.6.20.B	6	UC-707	0.6.05.B	15	UC-823	0.6.99.B	
UC-331	0.6.18.B	2	UC-631	0.6.29.B	1	UC-709	0.6.29.B	13	UC-828	0.6.40.B	6
UC-339	6.3.01.B	7	UC-632	0.6.20.B	7	UC-710	0.6.29.B	14	UC-835	0.6.41.B	4
UC-342	0.6.24.B	15	UC-633	0.6.20.B	19	UC-711	0.6.29.B	2	UC-836	0.6.41.B	5
UC-347	0.6.20.B		UC-633	0.6.29.B	4	UC-712	0.6.29.B	3	UC-837	0.6.41.B	6
UC-347	0.6.35.B	14	UC-634	0.6.20.B	8	UC-713	0.6.29.B	2	UC-838	0.6.40.B	7
UC-347	0.6.39.B	9	UC-635	0.6.20.B	6	UC-714	0.6.29.B	3	UC-839	0.6.40.B	8
UC-355	0.6.00.B	17	UC-635	0.6.29.B	1	UC-715	0.6.29.B	2	UC-843	0.6.09.B	19
UC-365	0.6.17.B	47	UC-636	0.6.20.B	7	UC-716	0.6.29.B	3	UC-844	0.6.09.B	24
UC-374	1.8.11.B	14	UC-637	0.6.20.B	19	UC-717	0.6.29.B	2	UC-859	0.6.42.B	11
UC-377	1.8.10.B	39	UC-637	0.6.29.B	4	UC-718	0.6.29.B	3	UD-013	1.8.11.B	15
UC-378	1.8.29.B	13	UC-638	0.6.20.B	8	UC-719	0.6.29.B	2	UD-013	1.8.25.B	24
UC-379	3.2.02.B	15	UC-639	0.6.20.B	6	UC-720	0.6.29.B	3	UR-003	1.8.12.B	8
UC-381	1.8.38.B	15	UC-639	0.6.29.B	1	UC-726	0.6.16.B	23	UR-003	1.8.13.B	21
UC-381	1.8.39.B	17	UC-640	0.6.20.B	7	UC-732	0.6.15.B	15	UR-051	1.8.22.B	19
UC-398	0.6.05.B	4	UC-641	0.6.20.B	19	UC-736	0.6.33.B	10	UR-053	1.8.34.B	18
UC-398	0.6.32.B	6	UC-641	0.6.29.B	4	UC-742	0.6.27.B	32	UR-055	1.8.09.B	19
UC-404	0.6.21.B	10	UC-642	0.6.20.B	8	UC-759	0.6.31.B	27	UR-060	0.6.23.B	13
UC-405	1.8.09.B	6	UC-647	0.6.22.B	4	UC-762	0.6.02.B	66	UR-060	4.3.01.B	14
UC-405	1.8.11.B	9	UC-648	0.6.22.B	4	UC-762	0.6.06.B	13	UR-062	3.2.02.B	19
UC-406	0.6.21.B	11	UC-649	0.6.22.B	4	UC-767	0.6.21.B	24	UR-063	1.8.34.B	17
UC-429	0.6.15.B	12	UC-650	0.6.22.B	4	UC-778	0.6.29.B	1	UR-064	1.8.38.B	13
UC-430	0.6.15.B	13	UC-651	0.6.22.B	4	UC-779	0.6.29.B	2	UR-064	1.8.39.B	7
UC-433	0.6.10.B	26	UC-652	0.6.22.B	3	UC-780	0.6.29.B	3	UR-065	3.2.01.B	22
UC-434	0.6.02.B	52	UC-653	0.6.22.B	3	UC-781	0.6.29.B	4	UR-109	1.8.29.B	27
UC-480	0.6.30.B	26	UC-654	0.6.22.B	3	UC-789	0.6.22.B	4	UR-110	1.8.33.B	4
UC-500	0.6.12.B	13	UC-655	0.6.22.B	3	UC-790	0.6.22.B	4	UR-115	1.8.29.B	17
UC-505	0.6.08.B	13	UC-656	0.6.22.B	3	UC-791	0.6.22.B	4	UR-150	1.8.03.B	15
UC-511	0.6.21.B	23	UC-661	0.6.27.B	14	UC-792	0.6.22.B	4	UR-150	1.8.04.B	5
UC-512	0.6.15.B	18	UC-662	0.6.27.B	28	UC-793	0.6.22.B	3	UR-591	1.8.39.B	8
UC-516	0.6.21.B	1	UC-663	0.6.27.B	15	UC-794	0.6.22.B	3	VB-001	4.3.02.B	1
UC-517	0.6.21.B	13	UC-664	0.6.27.B	16	UC-795	0.6.22.B	3	VB-002	4.3.02.B	2
UC-519	0.6.21.B	2	UC-665	0.6.15.B	21	UC-796	0.6.22.B	3	VB-003	4.3.02.B	3
UC-520	0.6.22.B	2	UC-666	0.6.27.B	29	UC-797	0.6.29.B	1	VB-004	4.3.02.B	4
UC-521	0.6.22.B	1	UC-683	0.6.38.B	13	UC-798	0.6.29.B	2	VB-005	4.3.02.B	5
UC-526	0.6.02.B	40	UC-686	0.6.27.B	13	UC-799	0.6.29.B	3	VB-006	4.3.02.B	6
UC-530	0.6.21.B	20	UC-687	0.6.27.B	30	UC-807	0.6.10.B	25	VB-007	4.3.02.B	7
UC-532	0.6.09.B	36	UC-688	0.6.27.B	31	UC-810	0.6.29.B	4	VB-008	4.3.02.B	8
UC-533	0.6.09.B	37	UC-694	0.6.27.B	1	UC-811	0.6.29.B	1	VB-009	4.3.02.B	9
UC-537	0.6.21.B	14	UC-695	0.6.27.B	2	UC-812	0.6.29.B	2	VB-010	4.3.02.B	10
UC-538	0.6.07.B	33	UC-696	0.6.27.B	3	UC-813	0.6.29.B	3	VB-011	4.3.02.B	11
UC-558	0.6.19.B		UC-697	0.6.27.B	4	UC-814	0.6.29.B	4	VB-012	4.3.02.B	12
UC-563	0.6.19.B	33	UC-698	0.6.27.B	5	UC-815	0.6.29.B	1	VB-013	4.3.02.B	13
UC-564	0.6.19.B		UC-699	0.6.27.B	6	UC-816	0.6.29.B	2	VB-014	4.3.02.B	14
UC-579	0.6.17.B	23	UC-701	0.6.27.B	7	UC-817	0.6.29.B	3	VB-015	4.3.02.B	15
UC-582	0.6.04.B	20	UC-702	0.6.27.B	8	UC-818	0.6.29.B	4	VB-016	4.3.02.B	16

NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No	NUMMER	SIDE	No
VB-017	4.3.02.B	17	WH-003	3.2.01.B	27	YR-037	0.6.24.B	5			
VB-018	4.3.02.B	18	WH-004	0.6.15.B	6	YR-038	0.6.24.B	5			
VB-019	4.3.02.B	19	WH-004	1.8.22.B	12	YR-039	0.6.24.B	5			
VB-020	4.3.02.B	20	WH-005	1.8.22.B	6	YR-040	0.6.24.B	5			
VB-021	4.3.02.B	21	WR-001	0.6.17.B	20	YR-041	0.6.24.B	5			
VB-022	4.3.02.B	22	WR-001	3.2.01.B	20	YR-055	0.6.24.B	5			
VB-023	4.3.02.B	23	WR-003	3.2.01.B	15	YR-056	0.6.24.B	5			
VB-024	4.3.02.B	24	WR-005	0.6.17.B	39	YR-057	0.6.24.B	5			
VB-025	4.3.02.B	25	WR-005	3.2.01.B	26	YR-058	0.6.24.B	5			
VB-026	4.3.02.B	26	WR-006	3.2.01.B	19	YR-059	0.6.24.B	5			
VB-027	4.3.02.B	27	WY-001	3.2.01.B	28	YR-060	0.6.24.B	5			
VB-028	4.3.02.B	28	XA-010	1.8.06.B	18	YR-061	0.6.24.B	5			
VB-029	4.3.02.B	29	YA-002	1.8.45.B	1	ZE-995	0.6.49.B	3			
VB-030	4.3.02.B	30	YA-004	0.6.24.B	7	ZE-996	0.6.49.B	3			
VB-031	4.3.02.B	31	YA-004	3.2.01.B	29	ZE-998	0.6.49.B	2			
VB-032	4.3.02.B	32	YA-050	0.6.24.B	7						
VB-033	4.3.02.B	33	YA-050	1.8.02.B	11						
VB-034	4.3.02.B	34	YA-051	0.6.24.B	7						
VB-035	4.3.02.B	35	YA-051	1.8.02.B	12						
VD-120	4.3.01.B	3	YA-052	0.6.24.B	7						
WE-004	0.6.17.B	45	YA-052	1.8.02.B	13						
WE-004	3.2.01.B	17	YA-053	0.6.24.B	7						
WE-005	3.2.01.B	1	YA-053	1.8.02.B	14						
WE-006	3.2.01.B	2	YA-079	1.8.45.B	4						
WE-007	3.2.01.B	2	YA-080	1.8.45.B	5						
WE-009	0.6.17.B	21	YA-088	0.6.23.B	16						
WE-009	3.2.01.B	21	YA-088	0.6.24.B	7						
WE-013	3.2.01.B	24	YA-088	4.3.01.B	17						
WE-014	0.6.17.B	37	YA-111	1.8.02.B	22						
WE-015	3.2.01.B	23	YA-111	1.8.09.B	20						
WE-016	0.6.17.B	46	YA-185	0.6.24.B	7						
WE-018	3.2.02.B	10	YA-443	0.6.24.B							
WE-020	3.2.01.B	5	YA-444	0.6.24.B							
WE-024	3.2.02.B	2	YA-970	0.6.24.B							
WE-026	3.2.01.B	4	YA-971	0.6.24.B							
WE-028	3.2.01.B	3	YR-001	0.6.24.B	3						
WE-031	3.2.01.B	2	YR-001	1.8.45.B	2						
WE-038	3.2.02.B	10	YR-002	0.6.24.B	2						
WE-040	3.2.02.B	16	YR-002	1.8.45.B	3						
WE-045	3.2.02.B	9	YR-003	0.6.24.B	5						
WE-046	3.2.02.B	9	YR-005	1.8.45.B	6						
WE-047	3.2.02.B	16	YR-006	1.8.45.B	6						
WE-099	3.2.01.B	30	YR-007	1.8.45.B	6						
WE-100	0.6.17.B	26	YR-008	1.8.45.B	6						
WE-100	3.2.01.B	31	YR-026	0.6.24.B							
WE-101	3.2.01.B	12	YR-028	1.8.45.B	6						
WH-001	0.6.17.B	12	YR-029	1.8.45.B	6						
WH-001	3.2.01.B	14	YR-031	0.6.24.B	14						
WH-002	3.2.01.B	18	YR-032	0.6.24.B							
WH-003	0.6.17.B	44	YR-035	0.6.24.B	1						